


PTC System Controls and Diagnostics



Pyramid Technical Consultants

Specifications (continued)

Example output peak shape from scintillation detector, +/-2 V vertical scale (60 keV X-rays, Cyberstar C17LA lanthanum chloride scintillation detector, non-inverting, 50 ohm termination)



PicoScope 6 19 Facoustants The Life Views Measurements Jools Hely 10 In dia dig 9 G (50 m/dw = 1) + 1 MS (1) 4 12 bits A (1 s20m/dw = 1) C = (5 (4 00 m/dw = 1)) + 1 MS (1) 4 12 bits A (1 s20m/dw = 1) C = (5 (4 00 m/dw = 1)) + 1 MS (1) 4 12 bits A (1 s20m/dw = 1) C = (5 (4 00 m/dw = 1)) + 1 MS (1) 4 12 bits A (1 s20m/dw = 1) + 1 C = (1) 5 (4 00 m/dw = 1) + 1 MS (1) 4 12 bits A (1) 4 12 bits (* * * 1of1 * Ø * x1 * * * * * * @ @ @ @ pico 16.0 - X ■ -26.85 uV -9.852 mV 9.825 mV 6 21 MN W -3.0 16.0 0.01 200.0 50.0 100.0 150.0 250.0 300.0 350.0 -50.0 Stopped 🕨 📕 Trigger Single 🖛 📌 🗛 🐨 🎜 📢 -1.638 mV 🔹 4 20 % 🔸 75 4 🔭 Messurements 🗊 💷 💷 Rulers 🛄 Notes 🖼 **PTC System Controls and Diagnostics**

Example corresponding direct detector signal, 20 mV vertical scale (50 ohm termination)



Pyramid Technical Consultants


CP15AL

Specifications (continued)

The CP15AL is designed with a compromise of bandwidth and gain to provide optimum pulse counting performance for fast photomultiplier-based detectors such as LaBr₃ and LaCl₃ scintillation detectors. However a pulse height analyzer that can handle the fast pulses from the C15AL can produce spectra for diagnostic purposes and photon energy analysis.

CP15AL


Pulse height spectrum example (Am-241 spectrum, Cyberstar C17LA lanthanum chloride scintillation detector)




	CP15A
Specifications (contin	ued)
Power input	+12V (+2/-4 V) DC, 70 mA and -12V (+4/-2 V), 70 mA
Displays	Power on LED
Controls	DC offset adjust 25-turn trimpot Gain adjust 25-turn trimpot (0 - 100%)
Case	Machined AI alloy case Ni plated with tapped hole pattern for mounting.
Dimensions	79 x 35 x 24 mm (see drawing)
Weight	95 g (3.4 oz)
Operating environment	10 to 35 C (15 to 25 C recommended to reduce drift and offset) , < 70% humidity, non-condensing, vibration < 0.1g all axes (1 to 100 Hz)
Shipping and storage environment	-10 to 50 C, < 80% humidity, non-condensing, vibration < 0.5 g all axes, 1 to 100 Hz
Connectors	
Signal input	Lemo coaxial size 00 receptacle (NIM-CAMAC CD/N)
Power	Lemo 0B four-pin female (mating cable connector FGG.0B.304) 1: -12 V 4: +12 V
	(view on CP15 connector / solder side of mating cable connector) 2: AGnd 3: AGnd (cable scrn)
Signal output	Lemo coaxial size 00 receptacle (NIM-CAMAC CD/N)
Ordering information	
_	Pre-amplifier suitable for LaBr $_3$ and LaCl $_3$ scintillator / photomultipliers.
Related products and	accessories
	Fast fixed gain pre-amplifier (suitable for LaCl ₃ , LaBr ₃ and plastic scintillation detectors)
	Very fast fixed gain pre-amplifier (suitable for plastic scintillation detectors, channel electron multipliers)
	Four channel pulse detector controller (fast window discriminators and scalers, supplies HV bias and pre-amp power).
PSU1212	+24 VDC to +/-12 VDC power converter
PSU24-25-1	+24 VDC power supply, universal, IEC C8 line connection

PTC System Controls and Diagnostics







Pyramid Technical Consultants