



BPM16-38

Beam Position Monitor

User Manual



Pyramid Technical Consultants, Inc.

135 Beaver Street, Suite 102, Waltham, MA 02452 USA

US: TEL: (781) 402-1700 ♦ FAX: (781) 402-1750 ♦ EMAIL: SUPPORT@PTCUSA.COM

Europe: TEL: +44 1273 492001

1 Contents

1	CONTENTS.....	2
2	FIGURES.....	4
3	SAFETY INFORMATION.....	5
4	MODELS.....	7
5	SCOPE OF SUPPLY.....	8
6	OPTIONAL ITEMS.....	9
6.1	READOUT.....	9
6.2	SIGNAL CABLES AND CABLE ACCESSORIES.....	9
7	INTENDED USE AND KEY FEATURES.....	10
7.1	INTENDED USE.....	10
7.2	KEY FEATURES.....	10
8	SPECIFICATION.....	11
9	INSTALLATION.....	15
9.1	PREPARATION AND HANDLING.....	15
9.2	MOUNTING.....	15
9.2.1	<i>Coordinate conventions.....</i>	<i>16</i>
9.2.2	<i>Beam direction.....</i>	<i>17</i>
9.2.3	<i>Rotation orientation.....</i>	<i>18</i>
9.3	CABLING AND SERVICES.....	19
9.3.1	<i>Connections.....</i>	<i>19</i>
9.3.2	<i>Cable and pipe arrangement.....</i>	<i>21</i>
9.3.3	<i>Alignment and connection to vacuum beamline.....</i>	<i>22</i>
9.4	VACUUM FORCE.....	23
10	AN OVERVIEW OF THE BPM16-38.....	24
10.1	IONIZATION CHAMBERS.....	24
10.1.1	<i>Signal formation.....</i>	<i>24</i>
10.1.2	<i>Gain calibration.....</i>	<i>25</i>
10.1.3	<i>Recombination.....</i>	<i>27</i>
10.2	POSITION READOUT.....	28
10.2.1	<i>Beam centroid.....</i>	<i>28</i>
10.2.2	<i>Beam width.....</i>	<i>31</i>
10.3	BPM16-38 MECHANICAL ASSEMBLY.....	31
10.3.1	<i>Ionization chamber electrodes.....</i>	<i>31</i>
10.3.2	<i>Gas pocket.....</i>	<i>32</i>
10.3.3	<i>Vacuum housing.....</i>	<i>32</i>
10.3.4	<i>Actuator system.....</i>	<i>34</i>
11	CONNECTORS.....	35
11.1	ELECTRICAL.....	35

11.1.1	Vertical sense signal output	35
11.1.2	Horizontal sense signal output	36
11.1.3	High voltage input	37
11.1.4	Actuator control	37
11.1.5	Redundant limit switch outputs.....	37
11.2	PNEUMATIC PRESSURE.....	38
11.3	FILL GAS	38
12	MAINTENANCE.....	39
12.1	PNEUMATIC CYLINDER LUBRICATION.....	39
12.1.1	Procedure	39
12.2	VACUUM BELLOWS REPLACEMENT	40
12.2.1	Procedure	41
12.3	MICROSWITCH ADJUSTMENT	41
12.3.1	Procedure	41
12.4	SPARE PARTS	41
13	FAULT-FINDING	43
14	RETURNS PROCEDURE.....	45
15	SUPPORT.....	46
16	DISPOSAL.....	47
17	REVISION HISTORY.....	48

2 Figures

Figure 1. BPM layout; sensor out of beam position (bellows extended). Dimensions mm.....	13
Figure 2. BPM layout continued; sensor out of beam position (bellows extended)). Dimensions mm.....	14
Figure 3. BPM mounting face.....	15
Figure 4. Ion optics standard coordinate convention (Transport).....	16
Figure 5. BPM16-38 internal coordinate system	17
Figure 6. BPM strip numbering relative to conventional beamline coordinates for beam entering via signal connector face (left) and via pneumatic connection face (right).....	18
Figure 7. Schematic installation.....	20
Figure 8. Fixed and moving cable connections.....	21
Figure 9. High flex cable connecting the actuator control connector to the solenoid valve.....	22
Figure 10. Ionization chamber signal formation.....	24
Figure 11. Approximate gain curve for the BPM16-38, protons in air at SATP	26
Figure 12. Measured chamber gains and corresponding centroid and width determinations, as a function of beam current of 228 MeV protons.....	27
Figure 13. Test result extract showing residual offsets of BPM sensor	28
Figure 14. Beam position in BPM sensor coordinates (blue).....	29
Figure 15. Displacement of BPM body coordinates(red) from BPM sensor coordinates (blue).....	30
Figure 16. Change to beamline coordinates (green).....	30
Figure 17. Ionization chamber electrode assembly.....	31
Figure 18. Electrode assembly mounted on spindle	32
Figure 19. Gas pocket and foil windows	32
Figure 20. Vacuum housing with sensor extended (in beam position), engaged on alignment pins.....	33
Figure 21. Vacuum housing with sensor retracted (out of beam position)	33
Figure 22. Complete BPM assembly with sensor extended (left) and retracted (right).....	34
Figure 23. Cylinder lubrication – step A.....	40
Figure 24. Cylinder lubrication – step B.....	40
Figure 25. Microswitch adjustment – step A. Upper switch (right) and lower switch (left).....	41

3 Safety Information

This unit is designed for compliance with harmonized electrical safety standard EN61010-1:2000. It must be used in accordance with its specifications and operating instructions. Operators of the unit are expected to be qualified personnel who are aware of electrical safety issues. The customer's Responsible Body, as defined in the standard, must ensure that operators are provided with the appropriate equipment and training.

The unit is designed to make measurements in **Measurement Category I** as defined in the standard.



CAUTION. Entrapment hazard.

The BPM is moved using pneumatic pressure during normal operation. There is a risk of fingers being trapped in pinch points. Never place body parts near the device when pneumatic pressure is present.



CAUTION. High voltage.

High voltage must be provided to this device for correct operation.

The high voltage is not exposed in the correctly assembled unit. A voltage of up to + or – 1000 V DC at 1.0 mA maximum is supplied to the BPM via the SHV connector. The hazardous live voltages delivered to the BPM are not accessible under the definitions of EN61010 but may nevertheless give a noticeable shock if misuse were to lead you to come into contact with them. The user must therefore exercise appropriate caution when servicing the device and when connecting cables. Power should be turned off before making any connections.

The body of the BPM should be grounded via its connection to the customer's beamline and/or mounting. If this is not the case, or if local regulations require, a safety ground must be securely connected to the ground lug on the case.



CAUTION. Radiation.

After use in a high-energy particle accelerator beamline, the BPM may become activated. Do not

work on the device, or move the device from a controlled area until it has been surveyed and declared safe by a qualified radiation supervisor.

Only Service Personnel, as defined in EN61010-1, should attempt to work on the disassembled unit, and then only under specific instruction from Pyramid Technical Consultants, Inc.

Some of the following symbols may be displayed on the unit, and have the indicated meanings.



Direct current



Earth (ground) terminal



Protective conductor terminal



Frame or chassis terminal



Equipotentiality



Supply ON



Supply OFF



CAUTION – RISK OF ELECTRIC SHOCK



CAUTION – RISK OF DANGER – REFER TO MANUAL



CAUTION – ENTRAPMENT HAZARD

4 Models

BPM16-38	Beam position monitor with sixteen sensing strips and 38mm sensing length in each axis.
-CF	With fixed tapped CF (Conflat™) 2 ¾ inch / DN40 flanges
-KF	With NW KF40 flanges
-FLO	With fittings for flow gas operation

Example:

BPM16-38-CF-FLO BPM with CF flanges and in/out flow gas fittings.

5 Scope of Supply

BPM16-38 model as specified in your order.

Four pin Weidmuller mating connector for redundant limit switch connection.

USB memory stick containing:

- BPM16-38 User manual

- Test data

Optional items as specified in your order.

6 Optional Items

6.1 Readout

I3200-XP10 32-channel electrometer with +1 kV bias supply and actuator control.

I3200-XP20 32-channel electrometer with +2 kV bias supply and actuator control.

6.2 Signal cables and cable accessories

CAB-D25M-5-D25F Cable, multiway low-noise, DSub 25 pin male to DSub 25 way female, 5 m.

CAB-D25M-10-D25F Cable, multiway low-noise, DSub 25 pin male to DSub 25 way female, 10 m.

CAB-SHV-5-SHV Cable, coaxial HV, SHV to SHV, 5 m.

CAB-SHV-10-SHV Cable, coaxial HV, SHV to SHV, 10 m.

CAB-D9M-5-D9F Cable, multiway, DSub 9 pin male to DSub 9 way female, 5 m.

CAB-D9M-10-D9F Cable, multiway, DSub 9 pin male to DSub 9 way female, 10 m.

7 Intended Use and Key Features

7.1 Intended Use

The BPM16-38 is intended to provide position and shape readout of high-energy ion beams, nominally proton beams in the energy range 30 to 350 MeV. It will form a part of a complete beam diagnostics suite in a particle accelerator system. The BPM performs a semi-destructive measurement of the beam. Although the beam passes through the device, the amount of scattering it receives is generally too great for the beam downstream to be used for other purposes. The BPM16-38 therefore includes a pneumatic actuator system to move the sensor element completely out of the beam path once the measurement has been made.

The BPM16-38 is intended to be incorporated into a high vacuum beamline, and is leak tested in the final stage of manufacture and test. It is not fully-UHV compatible because it includes elastomer seals, and thus cannot be baked to high temperatures. Moderate heating of the metal body to 70C maximum is acceptable to accelerate outgassing. The BPM may also be operated in atmosphere.

The operating environment should be clean and free of vibration. Users should be familiar with vacuum technology and low current measurement.

7.2 Key Features

- Low insertion length and choice of vacuum connectors.
- Small gap ionization chamber with 16 by 16 sensing strips.
- 38 mm by 38 mm sensitive area.
- Operation with atmospheric air or flow-through gas.
- Bias voltage up to 1 kV.
- High-vacuum compatible.
- Robust low-maintenance IC housing with bonded thin stainless steel windows.
- Integrated pneumatic actuation with dual limit-switch end-point sensing.
- Precision guides for in-beam location.
- 100000 cycle nominal bellows life.
- Compatible with the I3200 electrometer.

8 Specification

<i>Beam compatibility</i>		
	Species	Protons, deuterons, helium ions, carbon ions
	Energy range	30 MeV/nucleon to 500 500 MeV / nucleon
	Beam current density	10 pA cm ⁻² to 20 nA cm ⁻²
<i>Sensor</i>		
	Type	Parallel plate dual ionization chamber with multistrip cathodes
	HV bias range	500 V to 2000 V
	Sensitive area	38 mm by 38 mm Vacuum pipe internal diameter will define sensitive area.
	Sensitive volume	6mm anode – cathode spacing, each axis.
	Readout strip geometry	16 strips, equal width 2.38 mm on 2.534 mm pitch.
	Readout strip orientation	Horizontal sensing strips = sensing orthogonal to the axis along the centre of the vacuum bellows. Vertical sensing strips = sensing in the axis along the centre of the vacuum bellows.
<i>Vacuum</i>		
	Vacuum regime	High vacuum (1 e-8 mbar or greater).
	Bakeout	Maximum bakeout 70C (with screening of the electrical components and forced gas cooling of the sensor volume)
	Vacuum materials	Stainless steel, aluminium alloy, Viton O ring seals
	Bellows	Edge-welded stainless steel, rated lifetime > 100,000 cycles
	Beam path windows	50 µm stainless steel foil, diffusion bonded.
<i>Beampath materials</i>		
	Layers	1: 50 µm stainless steel foil entrance window 2: 10.8 mm fill gas 3: 152 µm FR4 fiberglass epoxy with 17.3 µm gold-flashed copper cathode strips 4: 6 mm fill gas (active volume) 5: 152 µm FR4 fiberglass epoxy with 17.3 µm gold-flashed copper anode planes each side 6: 6 mm fill gas (active volume) 7: 152 µm FR4 fiberglass epoxy with 17.3 µm gold-flashed copper cathode strips

		8: 10.8 mm fill gas 9: 50 μ m stainless steel foil exit window
<i>Actuator</i>		
	Type	Pneumatic, dual acting cylinder with DC solenoid control.
	Travel	58 mm
	Solenoid valve	Coil 24 VDC, 50 mA. DC voltage applied holds the sensor in the beam path.
	Pneumatic pressure	75 psi minimum, 110 psi maximum, clean dry air (CDA) or nitrogen. Maximum cylinder life is obtained with lubricated CDA.
	Limit sensing	Dual redundant microswitch sensing of fully in and fully out positions.
	Position reproducibility	Sensor positioning in beam reproduces to ± 0.1 mm in actuator motion direction, ± 0.25 mm in transverse dimension, referenced to BPM reference mounting surface.
	Position accuracy	Sensor positioning in beam within ± 0.25 mm of BPM body mechanical centre, in both axes, relative to BPM referenced to BPM reference mounting surface.
<i>Mechanical</i>		
	Insertion length	152.4 mm (6") flange face to flange face with CF flanges.
	Orientation	Operable in any orientation. Vertical orientation (actuator above beampipe) recommended for best position reproducibility.
	Weight	12 kg (26.4 lb)
	Operating environment	Clean and dust-free 5 to 35 C (15 to 25 C recommended), < 70% humidity, non-condensing Vibration < 0.1g all axes, 0.1 to 1000 Hz.
	Shipping and storage environment	Vacuum flanges should be sealed off and gas ports should be plugged. -10 to 50 C < 80% humidity, non-condensing Vibration < 1g all axes, 0.1 to 1000 Hz.
	Dimensions	See figures 1 and 2 for dimensions




Figure 1. BPM layout; sensor out of beam position (bellows extended). Dimensions mm.




Figure 2. BPM layout continued; sensor out of beam position (bellows extended)). Dimensions mm.

9 Installation

9.1 Preparation and handling

The BPM16-38 is shipped as a vacuum-tested clean assembly. Do not remove protective covers until the unit can be mounted in its operating position in a clean, dry area. Use gloves if you might touch any vacuum surfaces. Take care to protect knife edges on CF flanges, or sealing surfaces on KF flanges.

Similarly, any preparation or service work must be carried out in a clean, dry environment.



CAUTION. Entrapment hazard.

Once pneumatic pressure has been connected, there is a risk of unexpected movement of the actuator. Exercise due care.

9.2 Mounting

The BPM16-38 should be mounted via the reference face on its base, as this is where the sensor position is referenced. A pattern of four M6 threaded holes is provided on a precision machined surface. The weight of the device should be taken via this mounting, not via the beamline flanges. The mounting must be aligned accurately with the beamline ion optical elements, rigid and not subject to vibration. Depending upon your application, you may require that the mount provides position adjustment in the two transverse axes, so that the BPM can be adjusted relative optical elements in the beamline, to remove physically any offsets.




Figure 3. BPM mounting face.

9.2.1 Coordinate conventions

In the conventional beamline optics convention, z is always directed along the beam forward direction (thus tangent to beam trajectory when in a dipole field), the horizontal (X axis) is defined as orthogonal to z and outwards along the radius of curvature of local dipole magnets, and vertical (Y axis) being the orthogonal direction to Z and X (and thus the non-dispersion direction of dipole magnets). However your facility may adopt a different convention – check carefully to ensure you know how the BPM orientation you use relates to your coordinate system.




Figure 4. Ion optics standard coordinate convention (Transport).

The BPM16-38 has a natural internal XY coordinate system arising from the way the cathode strips are numbered, and this in turn is based on the assumption that an I3200 electrometer is used to read out the device using pin to pin cabling, with channels 1-16 used for the vertical sensing and 17-32 for the horizontal. Figure 5 shows this internal coordinate system.




Figure 5. BPM16-38 internal coordinate system

9.2.2 Beam direction

The BPM can be inserted into the beamline facing in either direction. The choice may be constrained by physical interferences with other beamline components.

The sensor chamber is symmetric, the only difference being which sensing axis the beam sees first as it passes through, and whether the positive X direction for conventional beam optics coordinates has increasing or decreasing strip numbering. Referring to figure 5, you can see that the internal coordinate system has the same sense as the conventional beamline optics system if the BPM is mounted with the pneumatic solenoid upstream. If you mount in the other direction, with the signal connectors upstream, you must remember to make the necessary inversion when converting from BPM coordinates to conventional beamline coordinates. As an example, figure 6 shows how the BPM strips are arranged relative to the conventional beamline coordinate system if the beam enters the BPM faces shown in figure 5.




Figure 6. BPM strip numbering relative to conventional beamline coordinates for beam entering via signal connector face (left) and via pneumatic connection face (right)

With high energy beams, the order in which the beam passes the vertical and horizontal halves of the sensor will make no difference to the reading. At and below the low end of the specification proton energy range, you may start to see the beam width apparently greater than it actually is in the downstream chamber, due to scattering in the entrance window and upstream chamber. This will have negligible effect on the centroid computation. If you use low energy X-rays to test the BPM on the bench, then you will see less integrated signal from the downstream half of the sensor, due to absorption.

9.2.3 Rotation orientation

The BPM will operate in any rotation orientation (rotation about the beam axis), to facilitate use on beamline gantry systems. The position reproducibility and accuracy are guaranteed for mounting with the actuator vertical in the absolute sense, either above or below the beamline. Note that the terms “horizontal-sensing” and “vertical sensing” lose their absolute meaning if the beamline or the BPM is rotated. Usually we refer to horizontal and vertical relative to the nearest dipole magnets (figure 4).

Beamline optics that use dipole, quadrupole and sextupole elements in standard relative orientations preserve the independence of the beam envelope trajectory in X and Y. This means we can look at the progression of X profiles along a sequence of BPMs independently of the Y profiles. If a BPM needs to be mounted at some arbitrary rotation angle to the local beamline coordinate system, then readout becomes complicated, and the data handling and display must manage this coordinate transformation. Rotation of the BPM by 90 degrees exchanges X and Y. Rotation by 180 degrees swaps the directions of X and Y. However a sensor mounted at an arbitrary angle mixes X and Y in the data, which complicates interpretation, and is therefore not recommended.

9.3 Cabling and services

9.3.1 Connections

The following cable connections are required:

<i>Function</i>	<i>Connector on the BPM16-38</i>	<i>Location on the BPM16-38</i>
Vertical axis readout	Dsub 25-way male	Moving upper body
Horizontal axis readout	Dsub 25-way male	Moving upper body
Bias voltage input	SHV receptacle	Moving upper body
Actuator control and limit switch readback	Dsub 9-way male	Static lower body
Redundant limit switch readback (optional)	Weidmuller 3.81mm 4-way terminal header male	Static lower body

The following service connections are required:

<i>Function</i>	<i>Connector on the BPM16-38</i>	<i>Location on the BPM16-38</i>
Pneumatic pressure in	¼" tubing	Moving upper body
Fill gas in (optional)	Threaded hole for 1/8" NPT fitting	Moving upper body
Fill gas out (optional)	Threaded hole for 1/8" NPT fitting	Moving upper body

Figure 7 shows a schematic installation using the Pyramid I3200 electrometer for readout. The I3200 is well-suited as it provides HV bias and actuator control as well as ionization chamber readout.




Figure 7. Schematic installation.

The location of the I3200 must reconcile the conflicting requirements to keep the signal cables short for best noise performance, but keep the electronics out of the radiation field. A maximum cable length of 10 m is recommended. Longer cables will still function, but signal to noise performance will degrade, and the capacitive load may introduce systematic calibration offsets. Refer to the I3200 user manual for further details on cable lengths. Low-noise signal cables must be used to minimize noise due to cable movement and vibration. These are available to special order from Pyramid Technical Consultants, Inc.

The actuator control link to the I3200 provides the switched 24 VDC for solenoid actuation, and reads back the “in” and “out” limit switches. A separate redundant set of limit switches allows the same information to be sent to an independent system, for example for safety interlocking. The BPM16-38 operates the same whether or not this optional connection is made.

The BPM16-38 can be supplied with pipe fittings to suit 1/8” OD pipe for flow gas feed and return as an option (-FLO option). Flow gas is not necessary for the sensor to function, but it may be preferred in some installations where the atmosphere is humid, or where large atmospheric pressure excursions are common.

9.3.2 Cable and pipe arrangement

Give careful attention to the routing of cables and pipe that connect to the BPM. All the connections that make to the moving upper body should be routed so that the movement of the cable or pipe when the actuator moves is allowed for, without any strain being placed on the connections. A strain-relief bracket is provided for the signal and HV cables, and the pneumatic feed pipe may also be routed via this point.

Do not tie the moving (signal and high voltage) and non-moving (actuator control) cables together within 1 metre (3') of the BPM. There should be no tendency for the cable to the non-moving part to move during actuation. Failure to observe this precaution may result in failure of the actuator connection due to flexing.




Figure 8. Fixed and moving cable connections.

The flexible cable connection between the actuator DSub connector on the static part of the BPM, and the solenoid valve must also be allowed to flex freely without risk of interference with other cables or services. Pay particular attention to this if the BPM is rotated out of the vertical, or would be rotated by a beamline gantry rotation.



Figure 9. High flex cable connecting the actuator control connector to the solenoid valve.

9.3.3 Alignment and connection to vacuum beamline

Each BPM is factory-tested for sensor positioning accuracy relative to the reference mounting surface, and the test data includes values for the residual offset in horizontal (X) and vertical (Y) axes. The errors will be less than 0.25 mm in X and Y. When you install the BPM you may wish to compensate these small residual errors if an adjustable mount permits it, or simply include the offsets in the software calibration of the sensor response.

The BPM must be surveyed and aligned relative to the ion optical components of the beamline (the magnets), not to the vacuum pipe. Generally the BPM16-38 must therefore be connected to the beam pipe via flexible sections on either side. The exception would be in beamlines where vacuum chambers are designed to be precision aligned to magnet pole pieces.

See section 10.2.1 for more information on how to read out beam position with appropriate corrections for residual offsets.

9.4 Vacuum force

When the BPM16-38 is under vacuum, atmospheric pressure will try to contract the bellows and thus move the sensor into the beam path. This force is overcome by the pneumatic pressure in normal operation. The default position for the actuator (solenoid un-energized) is for the sensor to be clear of the beam, so you should be aware that this condition will not be met if pneumatic pressure is lost.

10 An Overview of the BPM16-38

10.1 Ionization chambers

10.1.1 Signal formation

High energy ions pass through matter with relatively small lateral scattering and energy loss, but nevertheless leave a trail of ionization behind. In the ionization chamber the free charge that is created is separated by the applied bias voltage, with the positive air ions moving to the cathodes, which are grounded at the virtual earths of the readout preamplifiers, and the electrons (or negative ions formed by electron capture) moving to the anode. The resulting small current is measured by the readout electronics. The current from an individual ion is too small to measure, but for beam currents of a few 10's of pA or more, the aggregate current can be measured by sensitive electrometer electronics.




Figure 10. Ionization chamber signal formation.

The ionization chamber uses a parallel plate geometry to provide uniform gain over its active area. At higher applied bias voltages and with field intensifying geometry such as thin wires or points, the chamber would start operate in the proportional regime, where the signal is increased greatly by electron avalanching. This regime is generally avoided in high energy ion beamlines, because it is less stable, and the chamber is more prone to degradation from beam exposure.

Although the readout strips on the cathodes are separated by small gaps, all the charge arriving at the cathodes is routed by the electric field onto the strips. Therefore you can safely assume that the effective strip width is the same as the strip pitch, and the conversion to physical units is given by multiplying a result in strips by the strip pitch.

10.1.2 Gain calibration

The energy to produce an ion-electron pair is almost constant for a given target species, and very small compared to the energy carried by the ion, hence the minimal effect on the beam.

Air	34 eV
Nitrogen gas	36.4 eV
Oxygen gas	32.2 eV
Argon gas	26.3 eV
Helium gas	42.7 eV

The amount of ionization per ion in the gas filling of the BPM ionization chamber depends upon the gas composition (essentially constant if atmospheric air is used), the distance travelled through the gas, the pressure of the gas, and to a small extent, the energy of the ion. If the factors mentioned above are fixed by control or by calibration, for example against a Faraday cup collector, then the effective gain of the ionization chamber, (= chamber signal / ion beam current) is known and the BPM can be used to give a reasonable indication of beam current.

The approximate gain curve for the BPM16-38 with 6 mm electrode spacing, as a function of beam energy for protons in air at standard ambient temperature and pressure, is as shown in figure 10 below. A value of 75 was measured for the BPM16-38 for 228 MeV protons, with atmospheric air filling.



CAUTION

If the BPM16-38 is used for critical dosimetry applications, then you must use accurate gain values referenced to traceable standards, and regularly validated.

If the BPM is only being used to measure the position and shape of the beam, then we don't need to know the gain accurately, only that it is consistent across the chamber.





Figure 11. Approximate gain curve for the BPM16-38, protons in air at SATP.

Since the gain at any beam energy is a function of the gas density in the electrode gaps, it has to be corrected for variation in pressure and temperature, relative to the gain at some reference pressure and temperature where it is known. For example, if you know the gain at standard atmospheric temperature and pressure (SATP; Temperature_{SATP} = 298.15 K, Pressure_{SATP} = 100000 Pa), then the actual gain for another air pressure and temperature is given by

$$\text{Gain}_{\text{ACTUAL}} = \text{Gain}_{\text{SATP}} * 1 / [(\text{Pressure}_{\text{SATP}} / \text{Pressure}_{\text{ACTUAL}}) * (\text{Temperature}_{\text{ACTUAL}} / \text{Temperature}_{\text{SATP}})]$$

Temperatures must be in Kelvin, pressures can be in any convenient absolute unit.

10.1.3 Recombination

The upper beam current measurement limit of the ionization chamber is set by recombination of the ions and electrons before they can be collected on the electrodes, which reduces the measured current. Recombination is a function of the local beam current density, the electric field strength, the gas composition and pressure. It is mitigated in the BPM16-38 by the use of a small anode-cathode gap and a relatively high bias voltage.

At low current densities, corresponding to ion beam currents up to around 10 nA, recombination is negligible and can be safely ignored. Even at higher currents, the effect on the measurement of the beam centroid and width is comparatively little affected. The BPM16-38 has been tested with 1 kV bias at beam currents up to 35 nA and shown negligible difference in the computed beam parameters.




Figure 12. Measured chamber gains and corresponding centroid and width determinations, as a function of beam current of 228 MeV protons.

10.2 Position readout

The strip pitch of the BPM16-38, $S = 2.534$ mm, is relatively small compared to the beam width for the intended application, so that you will see signal on three strips at least. You can then use peak fitting or centre of mass calculation to determine the position of the peak to much less than one strip width, typically 10% of the strip width or less for normal beam currents and noise levels.

A centre of mass calculation (CoG) is simple to calculate and makes no assumptions about the shape of the peak. However it is unreliable if the whole peak is not included, or if excess background noise or offset is included. A fit to a Gaussian is often the best solution. Pyramid real time controller products include fast algorithms for peak fitting and centroid finding.

10.2.1 Beam centroid

The sensor strip geometry is controlled to high accuracy by the use of precision PCB manufacturing techniques. The centre of the pattern is between the eighth and ninth strips in each axis, with the strips numbered from 1 to 16 vertical, and 17 to 32 horizontal. Thus in a perfect system with no offsets, a perfectly centred beam would give a centroid reading of 24.500 strips in X and 8.5 in Y, which you would typically translate as a physical position of (0.000, 0.000) mm.

If you determine a peak position, (P_X, P_Y) expressed in strips, then the position in physical units in the BPM coordinate system (figure 5) is given by

$$X_S = (P_X - 24.5) * S$$

$$Y_S = (P_Y - 8.5) * S$$

To translate to the beamline coordinate system we must allow for any rotations of the BPM relative to the beamline coordinate system due to the way it is mounted, and also for any residual offsets between the BPM sensor coordinate system and the beamline coordinate system. The residual offsets between the sensor and the body of the BPM are determined during final factory test, and are given in the test results. The offsets are given in the BPM coordinate system, in microns.

Alignment Check	Measured Value	Results
Position 1 X	8.57	Nominal channel value 8.5
Position 1 Y	24.48	Nominal channel value 24.5
Position 2 X	8.52	Nominal channel value 8.5
Position 2 Y	24.52	Nominal channel value 24.5
dX (<250 microns)	114.21	PASS
dY (<250 microns)	0	PASS

Figure 13. Test result extract showing residual offsets of BPM sensor

Offsets of the BPM body relative to its ideal position in the beamline coordinate system will be determined during beamline survey.

The following example, with grossly exaggerated offsets, shows how a measured X position in sensor coordinates is translated to the beamline coordinate system. We'll take the more complicated case where the beamline X axis direction is opposite to the BPM X axis direction.

Firstly, the raw centroid, X_s , is found in the sensor coordinate system from the cathode strip signals, using a suitable centroid calculation.




Figure 14. Beam position in BPM sensor coordinates (blue)

Next we recognize that the beam position relative to BPM body must be adjusted by the residual offset, ΔX_0 , measured in the factory. This measurement tells us where the body is relative to the sensor coordinate system. It is simplest to express this as the displacement of the sensor centre relative to the body, which is the negation of the value dX given in the test data, and in the sensor coordinate system.




Figure 15. Displacement of BPM body coordinates (red) from BPM sensor coordinates (blue)

Finally we change to beamline coordinates. In this example this involves a change in the X axis direction, and a further measured residual offset, ΔX_1 , of the BPM body relative to the beamline, stated in the beamline coordinate system.




Figure 16. Change to beamline coordinates (green)

Thus the beam position in beamline coordinates is

$$X_B = \lambda(X_S - \Delta X_0) + \Delta X_1$$

where $\lambda=1$ if the sensor X axis is in the same direction as the beamline axis, and $\lambda=-1$ if it is in the opposed direction.

An alternative approach better suited to a final installation, is to assign absolute positions to the strips in beamline coordinates, with all scaling and offsets included, and store this in a software data file linked to the particular BPM. The beam centroid calculation will then return directly a value in beamline coordinates.

10.2.2 Beam width

Conversion of computed beam width to physical units is simply a matter of multiplying the width computed in strips to mm by multiplying by the strip pitch.

10.3 BPM16-38 Mechanical Assembly

10.3.1 Ionization chamber electrodes

The heart of the BPM16-38 is the dual ionization chamber electrode stack. This is manufactured from special thin FR4 PCB material using standard PCB manufacturing techniques, which give extremely high geometric accuracy. Cross hairs on the back sides of the cathodes are centred on the electrode strip patterns, allowing external optical alignment.




Figure 17. Ionization chamber electrode assembly.

The anode PCB in the middle is copper clad on both sides, to form the anode for both halves of the chamber. Note that the strips on the cathodes are orthogonal to the sensing direction (vertical strips sense in the horizontal axis, for example).

The electrode assembly is mounted on a hollow bore spindle that provides a route for the signal ribbon cables, HV bias coaxial cable and flow gas tube.




Figure 18. Electrode assembly mounted on spindle.

10.3.2 Gas pocket

A gas and vacuum-tight cap fits to the end of the spindle (shown transparent in the following figure).




Figure 19. Gas pocket and foil windows.

Thin stainless steel foil windows allow the passage of the beam. They are diffusion-bonded to the housing, and should require no maintenance for the life of the BPM. It is possible to assemble the BPM with the gas pocket removed if required for alignment checking. The crosshairs on the backs of the cathodes are then visible.

10.3.3 Vacuum housing

The gas pocket moves within the vacuum housing. At the end of travel, in the extended position, it engages on precision alignment pins to bring it into close alignment to the vacuum housing for

the X axis. The housing includes the precision mounting face. Stub tubes with alternative vacuum flanges can be mounted to suit the beamline standard.




Figure 20. Vacuum housing with sensor extended (in beam position), engaged on alignment pins




Figure 21. Vacuum housing with sensor retracted (out of beam position)

10.3.4 Actuator system

The vacuum housing is stationary, and the pneumatic actuator bears on it to move the sensor, gas pocket, spindle and the whole upper part of the assembly relative to it. An edge-welded bellows maintains vacuum while allowing the motion. Sensor connections pass up through the bellows and are terminated on bulkhead connectors. Dual microswitches make when the movement is within 4 mm of the end of travel. In the extended position, the travel limit is set by an adjustable hard stop, pre-set in manufacture, which defines the position accuracy in the Y axis.





Figure 22. Complete BPM assembly with sensor extended (left) and retracted (right)

11 Connectors

11.1 Electrical


11.1.1 Vertical sense signal output

DSub 25 pin male.



(External view on connector / solder side of mating plug)

View of cathode strips looking from the signal connectors side. Cathode strips are on the front side of the cathode from this view.




1	Strip 2	14	Strip 1
2	Strip 3	15	n/c
3	Strip 4	16	AGnd
4	Strip 5	17	AGnd
5	Strip 6	18	AGnd
6	Strip 7	19	AGnd
7	Strip 8	20	AGnd
8	Strip 9	21	AGnd
9	Strip 10	22	AGnd
10	Strip 11	23	AGnd
11	Strip 12	24	Strip 16
12	Strip 13	25	Strip 15
13	Strip 14		

The connector is compatible with a pin to pin connection to an I3200 electrometer. Strip 1 is connected to channel 1 and so on up to strip 16 being connected to channel 16.


11.1.2 Horizontal sense signal output

DSub 25 pin male.



(External view on connector / solder side of mating plug)

View of cathode strips looking from the signal connectors side. Cathode strips are on the back side of the cathode from this view.




1	Strip 18	14	Strip 17
2	Strip 19	15	n/c
3	Strip 20	16	AGnd
4	Strip 21	17	AGnd
5	Strip 22	18	AGnd
6	Strip 23	19	AGnd
7	Strip 24	20	AGnd
8	Strip 25	21	AGnd
9	Strip 26	22	AGnd
10	Strip 27	23	AGnd
11	Strip 28	24	Strip 32
12	Strip 29	25	Strip 31
13	Strip 30		

The connector is compatible with a pin to pin connection to an I3200 electrometer. Strip 17 is connected to channel 17 and so on up to strip 32 being connected to channel 32.


11.1.3 High voltage input

SHV receptacle. To mate with standard SHV connector.



11.1.4 Actuator control

Dsub 9 pin male.



(External view on connector / solder side of mating plug)


1	Solenoid power (switched +24 VDC in)	6	Limit switch common
2	Solenoid return	7	Limit switch in
3	n/c	8	n/c
4	n/c	9	n/c
5	Limit switch out		

The connector is compatible with a pin to pin connection to the actuator control port on an I3200 electrometer. When the actuator is retracted, pin 5 is connected to pin 6. When the actuator is extended (in beam), pin 7 is connected to pin 6.

11.1.5 Redundant limit switch outputs


Weidmuller 4 position 3.81 mm.

Two potential-free contact pairs. “Out” contacts close when the switch senses the actuator is retracted (sensor out of beam), “In” contacts close when the switch senses the actuator is extended (sensor in beam).



11.2 Pneumatic pressure

Push-fit connection (SMC “One-touch” / McMaster Carr “push to connect”) for ¼” od plastic hose.



11.3 Fill gas

Gas flow fittings are optional for customers who do not wish to operate the ionization chamber sensor in atmospheric air. 1/8” NPT threaded holes are provided in the upper body of the BPM for this purpose. The BPM16-38 can be ordered with the –FLO option which includes suitable fittings.

The feed tube which passes down to the ionization chamber volume should be connected on the inside to the fitting for the incoming gas. The outgoing gas is simply expelled from the internal gas space of the BPM. The flow should be regulated, and typical rates are very low, ≤ 1 sccm.

12 Maintenance

The BPM16-38 is designed for minimal maintenance. There are two recommended preventative maintenance procedures, and one adjustment that may be necessary.



CAUTION. Radiation.

Do not work in the beamline area, or on the BPM, until a survey has been completed by a qualified radiation supervisor and the radiation is known to be at acceptable levels.



CAUTION. Entrapment hazard.

Do not work on the BPM until pneumatic pressure has been removed.

12.1 Pneumatic cylinder lubrication

The cylinder should be lubricated with several drops of 10-weight engine oil every 500 hours of continuous operation. At 4 seconds typical per in and out movement, this corresponds to 450,000 cycles.

The procedure can be carried out with the BPM on the beamline.

12.1.1 Procedure

- 1) Remove pneumatic pressure and make safe using the working procedures for your site.
- 2) Disconnect the gas feed to the top of the cylinder.




Figure 23. Cylinder lubrication – step A.

- 3) Unscrew the pressure reducing valve, taking care not to change the setting.
- 4) Drop lubricant into the cylinder bore.




Figure 24. Cylinder lubrication – step B.

- 5) Reconnect the valve and pipe. Restore the pneumatic pressure and cycle the actuator several times to confirm correct operation.

12.2 Vacuum bellows replacement

The vacuum bellows are rated for 100,000 cycles of operation. You may wish to replace them as a preventative measure at or around this stage. The BPM must be removed from the beamline and partially disassembled to replace the bellows.

12.2.1 Procedure

The assembly requires precision alignment following replacement of the vacuum bellows. We therefore recommend that the BPM is returned to Pyramid for bellows replacement and refurbishment. Full details of the procedure will be supplied on request to users who have received appropriate training.

12.3 Microswitch adjustment

The microswitches should close when the actuator moves to within 2-4 mm of the hard end stop. You may need to adjust this if the switches are replaced, or if the setup is lost due to other work of accident. The BPM should not be under vacuum, so that the motion can be made manually.

12.3.1 Procedure

1) Move the BPM by hand to the fully retracted and fully extended positions. Ensure the microswitch actuation block contacts the switch rollers as shown.




Figure 25. Microswitch adjustment – step A. Upper switch (right) and lower switch (left)

- 2) Half-tighten the retaining screws, and move between the two positions to confirm that both switches at each end close between 4 and 2 mm of the hard stop. Adjust the positions as need to achieve this. Use a DVM to measure switch continuity.
- 3) When correct and consistent, tighten the screws and do a final check.

12.4 Spare parts

Item	Part number	Supplier
Vacuum bellows	Pyramid 11204461	Metal Flex Nevada Inc.
Pneumatic cylinder	Pyramid 165692	Bimba
Solenoid valve assy	SY5120-5LN-N7T-F2	SMC
Solenoid drive cable		Pyramid Technical Cons.
Microswitch	SS-01GL2D	Omron
Signal and bias PCB	9000-11112	Pyramid Technical Cons.
Actuator control PCB	9000-11128	Pyramid Technical Cons.

V sense cathode	9000-11078	Pyramid Technical Cons.
H sense cathode	9000-11076	Pyramid Technical Cons.
Anode	9000-11080	Pyramid Technical Cons.
CF flange assy	17219330	Pyramid Technical Cons.
Compressed gas tube ¼"	5181K232	McMaster-Carr
Ribbon cable assy	HF18U-12-ND	Digikey.
Anode cable assy		Pyramid Technical Cons.
Gas pocket assy	11204490	Pyramid Technical Cons.

13 Fault-finding

The BPM16-38 is designed to give you trouble-free service. We expect that a simple replacement policy will be followed for any units that fail on a beamline, and that failed units will be returned to Pyramid for refurbishment. However the following fault-finding is provided to help decide whether a BPM should be exchanged, and to guide repairs for customers who do not have a service arrangement.

<i>Symptom</i>	<i>Possible Cause</i>	<i>Confirmation</i>	<i>Solution</i>
Actuator does not move	Insufficient pneumatic pressure	Check pressure	Use correct pressure (75 to 110 psi).
	Mechanism is stuck.	Remove pneumatic pressure and attempt to move manually.	Maintain mechanism as necessary.
	Control cable disconnected	Check cable	Ensure cables are connected and secure.
	Cables or pipes are snagging	Visual inspection.	Ensure all incoming cables and services allow free actuator movement in all operating orientations.
	Solenoid failure	24 V across the solenoid? Solenoid magnetic field detectable?	Replace solenoid as necessary.
	Flow restrictors set incorrectly.	Check incoming gas pressure.	Adjust restrictors as necessary.
Limit switches not reading as expected.	Switches out of position.	Visual inspection.	Adjust and secure microswitches so that they close within 2 to 4mm of the end of travel.
	Switches faulty	Disconnect actuator and check switch function with a DVM.	Replace switch as necessary, re-set position.
Error in position readout	Mechanism not reaching end of travel and engaging with alignment pins.	Visual inspection – has actuator reached the hard stops?	Maintain mechanism as necessary.
	BPM is installed in an	Visual inspection;	Ensure calibration

	orientation that invalidates calibration values.	check calibration values being used.	values are correct.
	X and Y axis cables swapped.	Check cable	Ensure cables are connected and secure.
Small or no signal	HV bias not reaching the BPM.	Check HV supply. Check cable integrity.	Correct as necessary. If no external cause of missing HV, remove BPM and check internal connection.
Unstable signal	HV arcing due to high humidity.	Reduce HV and see if instability reduces. Improve air conditioning.	Improve air conditioning. Consider using flow gas.
	Triboelectric noise.	Check signal again when no vibration or movement.	Use low-noise cables. Do not use data from immediately after the actuator is moved.
Bad vacuum	Leak on seal	Leak check	Repair as necessary
	Leak on bellows	Aggravated by actuator movement. Leak check.	Repair as necessary.
	Leak on diffusion bond	Leak check (feed helium into gas pocket).	Replace gas pocket as necessary.

14 Returns procedure

Damaged or faulty units cannot be returned unless a Returns Material Authorization (RMA) number has been issued by Pyramid Technical Consultants, Inc. If you need to return a unit, contact Pyramid Technical Consultants at support@ptcusa.com, stating

- model
- serial number
- nature of fault



CAUTION. Radiation.

The unit cannot be shipped until it is certified to be below legal limits for radiation, and that it is clear of any chemical contamination.

An RMA will be issued, including details of which service center to return the unit to. We recommend that you retain the original shipping carton to use for returns.

15 Support

Manual and other documentation updates are available for download from the Pyramid Technical Consultants website at www.ptcusa.com. Technical support is available by email from support@ptcusa.com. Please provide the model number and serial number of your unit, plus relevant details of your application.

16 Disposal

We hope that the BPM16-38 gives you long and reliable service. The BPM16-38 is manufactured to be compliance with the European Union RoHS Directive 2002/95/EC, and as such should not present any health hazard, once any activation has decayed.



CAUTION. Radiation.

The BPM must not be released from a radiation controlled area until it has been surveyed and declared safe by a qualified Radiation Supervisor.

When your BPM16-38 has reached the end of its working life, you must dispose of it in accordance with local regulations in force. If you are disposing of the product in the European Union, this includes compliance with the Waste Electrical and Electronic Equipment Directive (WEEE) 2002/96/EC. Please contact Pyramid Technical Consultants, Inc. for instructions when you wish to dispose of the device.

17 Revision History

The release date of a Pyramid Technical Consultants, Inc. user manual can be determined from the document file name, where it is encoded yymmdd. For example, B10_UM_080105 would be a B10 manual released on 5 January 2008.

<i>Version</i>	<i>Changes</i>
BPM16-38_UM_120420	First general release
BPM16-38_UM_130521	Added information on gain curves and temperature and pressure correction.
BPM16-38_UM_131023	Added information and picture for cable strain relief Corrected gain temperature and pressure correction equation.
BPM16-38_UM_250722	Updated address and removed outdated references