CM100 Control and Charge Monitoring Unit

for Particle Therapy Control Rooms

User Manual





Pyramid Technical Consultants, Inc.

1050 Waltham Street Suite 200, Lexington, MA 02421 USA

US: TEL: (781) 402-1700 FAX: (781) 402-1750 EMAIL: <u>SUPPORT@PTCUSA.COM</u> Europe: TEL: +44 1273 492001

1 Table of Contents

1	ТАВ	LE OF CONTENTS	2	
2	FIGURES			
3	SAF	SAFETY INFORMATION		
	3.1	Standards	6	
	3.2	Power and Grounding	6	
	3.3	INTERNAL BATTERY	6	
	3.4	SAFETY CONSIDERATIONS FOR THE INTENDED APPLICATION	6	
	3.4.1	Interlock relays	6	
	3.4.2	Pulse Counting Cable Disconnect	6	
	3.4.3	Watchaog	/	
	5.5	JTWIDOLS	/	
4	MO	DELS	8	
5	500		٩	
5	300			
6	ΟΡΤ	IONAL ITEMS AND RELATED PRODUCTS	. 10	
	6.1	ELECTROMETERS	10	
	6.2	POWER SUPPLIES	10	
	6.3	BATTERY PACK	10	
	6.4	SIGNAL CABLES AND CABLE ACCESSORIES	10	
	6.5	DATA CABLES	10	
7	INTE	ENDED USE AND KEY FEATURES	. 11	
	7.1	INTENDED USE	11	
	7.2	Key Features	11	
8	SPE	CIFICATION	. 12	
9	INST		. 19	
5	0 1	Mounting	10	
	9.1		19	
	9.2	GROUNDING AND POWER SUPPLY	. 19	
	9.3	BATTERY	20	
	9.4 0 E		ZU	
	9.5 0.6		ZI	
	ס.צ ס.ד		ZI	
_	9.7		∠⊥	
1	D HOV	V THE CM100 WORKS AND HOW TO USE IT - AN OVERVIEW	. 22	
	10.1	TOUCH SCREEN INTERFACE	22	

10.2	INDEPENDENT TERMINATION SYSTEM AND PULSE COUNTING	22
10.3	AUDIBLE INDICATOR	24
10.4	Кеүѕwітсн	24
10.5	Start and Pause Buttons	24
10.6	EMERGENCY STOP BUTTON	25
10.7	SYNCHRONIZATION INTERFACE	25
10.8	SAFETY INTERLOCKS	25
10.0	ΒΑΤΤΕΡΥ ΒΑΓΥΙΙΟ	23
10.5		27
11 UPF	DATING THE CM100	27
11 012		20
11.1	Updating the Firmware	28
11.1.	1 Changing the SD card.	28
11.1.	.2 Copying program files	28
11.2	SYSTEM.XML	29
11.2.	1 Editing the File	29
11.2	.2 Attributes that can be Changed	29
11.3	CHANGING IP ADDRESS	29
11.4	Serial Number	30
11.5	LOOPBACK ADDRESS	30
12 STA	TE MACHINES	31
12.1	Operational State Machine	31
12.2	Εινοτιοναί State Machine	32
12 1165		24
15 036		54
13.1	STATUS LEDS	34
13.1.	.1 AC	34
13.1.	2 CMS	34
13.1.	.3 Cable	35
13.2	SYSTEM STATE	35
13.3	System Messages	35
13.4	Date/Time	35
13.5	Номе Виттол	35
13.6	Setup Button	35
13.6	1 Network	36
13.6.	2 About	37
13.6	.3 Screen settings	37
13.6	.4 Audio settings	38
13.6	.5 Connection settings	39
13.7	AUDIO VOLUME	39
13.8	Session Information	40
13.8	1 Count total	40
13.8	2 Unfiltered count	40
13.8	3 Count threshold	40
13.8	.4 Session ID	40
13.9	PULSE RATE	40

1	3.10	SPOT COUNT	40	
1	3.11	Error Messages	40	
14	EXA	MPLE USE CASE	43	
		0		
1	.4.1	UFF	44	
1	.4.2	STANDBY	44	
1	.4.3	PREP	45	
1	.4.4	I REAT	46	
1	.4.5	BEAM ON	47	
1	.4.6		48	
1	.4.7	PAUSING AN ONGOING TREATMENT	48	
1	.4.8	EMO Pressed	48	
1	.4.9	MOVING THE KEYSWITCH OFF THE TREAT POSITION DURING BEAM ON STATE	49	
15	CON	INECTORS	50	
1	.5.1	Rear Panel	50	
	15.1.	1 Ground lug	.50	
	15.1.	2 Power input, 24 VDC	.50	
	15.1.	3 Data links	.50	
	15.1.	4 Q Pulse inhibit	.51	
	15.1.	5 Pulse in	.51	
1	15.1.		.51	
Ŧ	.5.Z 15.2	1 Switches	52	
	15.2.	2 Battery	.53	
10				
10	CON		22	
1	.6.1	EPICS PROCESS VARIABLES	55	
17	FAU	ILT FINDING	61	
10			C2	
10	IVIAI		03	
1	.8.1	OPERATION CHECK	63	
1	.8.2	CHECKING THE BATTERY	63	
1	.8.3	BATTERY PACK REPLACEMENT	63	
19	RET	URNS PROCEDURE	64	
••			<u> </u>	
20	SUP	PUKI	65	
21	DISF	POSAL	66	
	REVISION HISTORY			
KEV	KEVISION HISTORY			

2 Figures

Figure 1 - CM100 top and side views. Dimensions mm	16
Figure 2 - CM100 top and rear views. Dimensions mm	17
Figure 3 - CM100 side view. Dimensions mm	18
Figure 4 - CM100 connections	20
Figure 5 - CM100 Block Diagram	22
Figure 6 - Pulse inhibit logic	23
Figure 7 - Keyswitch	24
Figure 8 - Start and pause buttons	25
Figure 9 - Emergency off button	25
Figure 10 - Safety interlocks – rear panel connectors	26
Figure 11 - CM100 state machine	31
Figure 12 - CM100 operational states	32
Figure 13 - Operation state conditions	32
Figure 14 - Functional states	32
Figure 15 - Functional state conditions	33
Figure 16 - Home screen	34
Figure 17 –Setup menu	36
Figure 18 - Network settings page	36
Figure 19 - About screen	37
Figure 20 - Screen settings	38
Figure 21 - Audio information	38
Figure 22 - Connection settings	39
Figure 23 - Audio volume	39
Figure 24 - Error message display and CLEAR button	42
Figure 25 - Sample sequence for CM100 use	43
Figure 26 - CM100 in the Off State	44
Figure 27 - CM100 in the Standby state	45
Figure 28 - CM100 in Prep state	45
Figure 29 - CM100 in Prep state ready to transition to Treat state	46
Figure 30 - CM100 in Ready state	47
Figure 31 - CM100 in Beam On state	47
Figure 32 - CM100 when the treatment is completed	48
Figure 33 - CM100 Emergency Off	49
Figure 34 - CM100 rear panel	50

3 Safety Information

3.1 Standards

This unit is designed for compliance with harmonized electrical safety standard EN61010-1:2000. It must be used in accordance with its specifications and operating instructions. Operators of the unit are expected to be qualified personnel who are aware of electrical safety issues. The customer's Responsible Body, as defined in the standard, must ensure that operators are provided with the appropriate equipment and training.

3.2 Power and Grounding

The unit is designed to operate from +24VDC power, with a typical current requirement of 1000 mA. A suitably rated power supply module is supplied as standard and is strongly recommended. Customers who make their own 24 V power provision should ensure that the supply cannot source more than 4200 mA, that it has the correct connector type, that it is rated for use in the environment, that it provides overcurrent and short circuit shutdown and that it has all necessary regulatory approvals.

A safety ground must be securely connected to the ground lug on the case.

3.3 Internal Battery

The unit has an internal battery that allows the unit to operate for up to 20 minutes after the power is removed. The battery is a lithium iron phosphate (LiFePO4) type, selected for its inherent safety, reliability, and low toxicity.

3.4 Safety considerations for the intended application

The CM100 is a component of an overall safety system for a proton therapy system. By itself, it cannot ensure safe delivery of dose. It must be connected to relevant signals and interlock systems, the configuration must be subjected to risk analysis and the connections must be checked and tested. The following are safety features of the CM100 that must be functional and used correctly. If there is any reason to suspect the features are not working as described, the proton therapy system must not be used for medical treatments until the problem is identified and resolved.

3.4.1 Interlock relays

The unit is equipped with two safety rated relays that can be used to interlock the beam. An Emergency Off Button (EMO) is provided that connects to dual parallel safety relays.

3.4.2 Pulse Counting Cable Disconnect

The TTL and fiber signals are set to be normally high. Disconnecting the cable is automatically detected by the system, opening the safety relays and placing the unit in the error state.

3.4.3 Watchdog

The unit is equipped with a watchdog timer, maintained by the ARM/Linux + PRU programs. The watchdog triggers within 2ms if the program is compromised, opening both safety relays and resetting the processor.

3.5 Symbols

Some of the following symbols may be displayed on the unit, and have the indicated meanings.

	Direct current
Ţ	Earth (ground) terminal
	Protective conductor terminal
\downarrow	Frame or chassis terminal
\triangleleft	Equipotentiality
	Supply ON
\bigcirc	Supply OFF
A	CAUTION – RISK OF ELECTRIC SHOCK

CAUTION - RISK OF DANGER - REFER TO MANUAL

4 Models

CM100Control and charge monitoring unit for particle therapy. De build is TTL on BNC for pulses and gate.	
-QN	Configured to count fast -ve pulses on BNC input.
-Q0	Configured to count optical pulses on ST fiber input.
-GO	Configured to respond to optical gate signals on ST fiber input.
-PM	6U panel mounting kit.

Example:

CM100-QO CM100 configured to count fast -ve pulses on ST fiber input.

5 Scope of Supply

CM100 model as specified in your order.

PSU24-100-1R 24 VDC 100-watt power supply with Redel locking connector, rated for medical use.

USB memory stick containing: CM100 User manual Test data

OEM customers will receive only components relevant to their application.

6 Optional Items and related products

6.1 Electrometers

1128 electrometer, single high current channel and 128 strip channel current readout electronics.

F100 electrometer, single high current channel current readout electronics.

XF03 accessory board for use with the F100 electrometer in proton therapy systems.

6.2 Power Supplies

PSU24-100-1R +24 VDC 40W PSU (universal voltage input, plug receptacle for standard IEC C14 three-pin socket) with output lead terminated in two-pin Redel PAG connector.

6.3 Battery Pack

BAT24-CM100 +24 V battery pack for CM100.

6.4 Signal Cables and Cable Accessories

CAB-BNC-xxx-BNC cable BNC terminated both ends, RG-58.

CAB-ST-xxHCS-ST Fiber-optic cable pair 200 μ m silica fiber ST terminated with color-coded sleeves, xx feet long.

6.5 Data Cables

CAB-RJ45-xxx-RJ45 Ethernet CAT5.

CAB-ST-xxHCS-ST Fiber-optic cable pair 200 μ m silica fiber ST terminated with color-coded sleeves, xx feet long.

7 Intended Use and Key Features

7.1 Intended Use

The CM100 was designed to meet relevant requirements of IEC 60601-2-64:2014 when installed in a proton therapy treatment control room and connected correctly to other subsystems. Its primary purpose is to provide a redundant means of recording total delivered charge, and making that data available in the event of a power outage. The CM100 accomplishes this though a battery-backed display and nonvolatile storage. In addition, the device includes the capability of counting individual spots so that it can satisfy the needs of PBS in spot-scanning mode when attached to compatible dosimetry systems.

The charge collection function is combined with a flexible set of switches that can be configured to satisfy IEC 60601-2-64:2014 with respect to hard-wired and relay-driven interfaces to arm, start, pause and stop the treatment process. Additional functions are available through the graphic touch-screen, allowing the device to be updated as new requirements arise.

The CM100 also provides an audible indicator that is proportional to the dose rate. This indicator is internal to the device but can also be transmitted to external speakers.

7.2 Key Features

Designed to meet the requirements of IEC 60601-2-64:2014.

Non-volatile charge recorder (> 30 minutes duration).

Key switch access for preparing and enabling irradiation.

Illuminated physical switches to start and pause irradiation.

Audible signal for dose being delivered.

Latching emergency stop button with direct connection to rear panel connector.

Relay to stop irradiation if total dose exceeds target plus allowed tolerance.

7" color LCD touch screen for real-time display of dose delivered, dose target and system state. Built-in self-diagnostics.

Automatic irradiation termination if internal failure detected.

Checks that critical cable connections are made.

Ethernet connectivity to compatible dosimetry and session management systems.

Maintains a record of pencil beam spot number when used with compatible dosimetry systems.

8 Specification

Ch	Charge recording			
	Signal type	Selectable from:		
		- TTL pulses		
		- Fiber optic 640 nm pulses		
		The dosimetry system that the CM100 works with defines the		
		quantum of dose that corresponds to one pulse.		
	Maximum pulse rate	2 MHz.		
	Dose per pulse	Configurable from 1e-5 to 1.0 MU recorded per pulse received.		
	Counting scheme	Upwards from zero to set dose.		
	Count limit	Automatic interlock if count reaches configurable limit, typically 120% of set dose.		
	Counter depth	64 bits		
	Data retention and	Data recorded to flash memory card.		
display CM100 display and function main minutes after power loss using bu pack.		CM100 display and function maintained for greater than 20 minutes after power loss using built-in re-chargeable battery pack.		
	Gating	Two inhibit inputs to allow charge monitoring into a gated counter channel to be temporarily disabled, for example if the beam is known to be absent by other means and any apparent dose would be spurious.		
		An example is to prevent counting apparent dose from an imaging X-ray shot when the particle beam is known to be blocked by a beam stop.		
		A parallel non-gated charge monitoring counter remains active.		
		Gate inputs purchase option select from:		
		- TTL pulses		
		- Fiber-optic 640 nm pulses		
		Dose pulse input Ungated counter		
		Gate inputs A Gated counter		
Au	Idio			
	Output	Internal speaker, adjustable volume (muting not allowed).		
Line out jack.		Line out jack.		

PTC2-2120709629-4089

	Sound	Selectable "tick" or "beep". Audio frequency of beep option selectable.	
	Scaling	Configurable pulse rate to audio tick rate mapping from 1e6 to 1e2 pulses per tick.	
	Source	Ungated channel.	
Ве	am disable		
	Pause	Safety-rated relay (Tyco SR4) with mechanically-guided contacts and sensing of welded contacts.	
		Normally open contacts. Relay closed if Max MU threshold is not exceeded, CM100 is in READY or BEAM ON state with no errors.	
	Overflow	Safety-rated relay (Tyco SR4) with mechanically-guided contacts and sensing of welded contacts.	
		Normally open contacts. Relay closed if Max MU threshold is not exceeded, CM100 is in READY or BEAM ON state with no errors.	
Pro	ocessor		
	Туре	TI Sitara AM335x (ARM Cortex A8) 1 GHz with dual PRU.	
	Operating system	Debian Linux.	
	Watchdog	Relays open (hardware action) and forced processor reset if watchdog is not tickled every millisecond.	
Self-test			
	Automated self-test (POST)	Automated tests of relay function, RAM and flash memory, battery function, Ethernet connection.	
		System operation prohibited if POST fails.	
	Other tests	Tests with user prompts for button function, emergency off, audio function, key switch function, touchscreen function, fiber- optic signalling.	
Со	nnectivity		
	Ethernet	Ethernet 10/100/1000 Mbps. Auto MDIX.	
		Embedded EPICS channel access server allows client software to monitor and control device function.	
	USB	USB port for device setup and diagnostics (qualified service technician access).	
		Connecting to the USB port creates a virtual network to a host system and appears at static IP address 192.168.7.2.	
	Dosimetry system	Direct fiber-optic communications channel to compatible pencil beam scanning scan and dose control systems to allow CM100 to	

		record spot number and system state.
		The connection is not mandatory; the CM100 will function in respects aside from spot number tracking if it is not present.
	Other devices	Fiber-optic communications to devices including magnetic field monitors (H20 device), power supply interfaces (M10, M40 devices) in normal or snooping mode for optional recording of data associated with beam scanning.
Ро	wer	
	Power input	24 V (+/- 2V) DC, 1000 mA typ, 1500 mA max.
	Battery backup	Operation of device including user interface continues for not less than 20 minutes if power is removed.
		Alert if battery pack is missing, low performance or not in good condition.
		Factory replaceable LiFePO4 battery assembly.
Са	se	
	Configuration	Desk-mounting console. See figures for dimensions.
	Protection rating	IP34 (proof against splashed liquid).
	Weight	4.5 kg (10.0 lb)
En	vironment	
	Intended location	Particle therapy treatment rooms (one CM100 per room).
	Operating environment	10 to 35 C (15 to 25 C recommended) , < 70% humidity, non- condensing, vibration < 0.2g all axes (1 to 100 Hz).
	Shipping and storage environment	-10 to 50 C, < 80% humidity, non-condensing, vibration < 1g all axes, 1 to 100 Hz.
Со	ntrols	
	Keyswitch	Three-position switch with key retention. Positions correspond to stages in an irradiation as defined in IEC 60601-2-64: Off / Prepare / Treat.
		Backlit labels indicate PREP and TREAT conditions.
		Key can only be removed in OFF position.
	Emergency Stop	Locking push switch with visual indication of actuation.
	Start	Pushbutton with green illuminated bezel indicating availability to start or resume irradiation.
	Pause	Pushbutton with blue illuminated bezel indicating availability to pause irradiation.

CM100 User Manual

User interface	1024 by 768 backlit color LCD capacitive touch screen, 7" (17.8 cm) diagonal.
Access rights	Clinical controls or service controls only enabled when authorization codes are received from a host system via Ethernet command.
Processor reset	Rear panel push-button.





Figure 1 - CM100 top and side views. Dimensions mm.

CM100 User Manual



Figure 2 - CM100 top and rear views. Dimensions mm.

CM100_UM_180228



Figure 3 - CM100 side view. Dimensions mm.

9 Installation

9.1 Mounting

The CM100 is intended for table top use in the treatment control room of a particle therapy facility. The unit should be placed on a level surface in an area easily accessible by the therapist. The unit ideally should be located close to the computer screens (clinical user interface) that are used to set up and monitor the treatment.

The unit is designed to prevent accidental movement via a set of rubber feet on the bottom, and is also heavy enough to prevent movement under normal use.

The mounting position should allow sufficient access to connectors and cable bend radii. 60 mm minimum clearance is recommended at the front and back of the device.

Best performance will be achieved if the CM100 is in a temperature-controlled environment. No forced-air cooling is required, but free convection should be allowed around the back and sides of the case.

An optional 19" rack mount configuration is available for the CM100. Details on the rack mount configuration are available on request.

9.2 Grounding and Power Supply

A grounding connection is not required for the CM100.

+24 VDC power should be provided from a suitably-rated power supply with the following minimum performance:

Output voltage	+24 +/- 0.5 VDC
Output current	1000 mA minimum, 4200 mA maximum
Ripple and noise	< 1% pk-pk, 1 Hz to 1 MHz
Regulation	< +/- 5%

The CM100 includes an internal automatically re-setting PTC fuse rated at 1.1 A that protects the internal circuitry. The battery charging circuit is not fused, so the external 24 V power supply must provide an overcurrent shutdown facility. The PSU24-100-1R power supply provides this protection.

9.3 Battery

The CM100 is shipped with the battery disconnected. Use of the battery is optional if you have connected the CM100 to a secure uninterruptable power supply (UPS).

If you wish to use the internal battery, unscrew the four screws on the bottom panel of the unit.



Attach the three Molex connectors connecting the battery to the PCB. Allow the battery to charge up for four hours, then check the battery as in the section "Checking the Battery".

9.4 Signal Source

The diagram below shows an example CM100 installation in schematic form.



The primary function of the CM100 is to record dose from one of the electrometers monitoring the treatment via an integral plane on an ionization chamber. In the diagram the signals are transmitted as pulses via a fiber-optic cable connected to the redundant readout electronics. On some models this is connection is made with a BNC cable to a separate dedicated connector on the unit. If a BNC cable is used, it should be 50 ohm terminated to prevent signal reflections.

A separate input from one of the scan/dose real time controllers specifies whether the pulses are to be ignored, in other words, are not considered as valid particles to be measured. In the diagram this connection is made via a fiber-optic cable connected to one of the scan/dose real time controllers.

9.5 Interlocks

Hardware relay signals generated by the CM100 (potential-free contacts or 24 V logic) are routed directly to one or more beam abort systems. The decisions to open these interlocks are made by the CM100 depending upon the state of the key-switch, buttons, device state, and beam overflow detection logic.

9.6 External Audio

A speaker jack allows the CM100 to be connected to an external speaker. The speaker must be self-powered. Note that the internal speaker remains active.

9.7 Client

The CM100 connects to a client system via a standard CAT5 ethernet cable. Data is exchanged with the client via the EPICS protocol.

10 How the CM100 Works and How to Use it - An Overview

The CM100 is specially designed for particle therapy systems, providing a physical interface to the therapist with the treatment system. The CM100 was specially designed to meet certain requirements of IEC 60601-2-64:2014.



The block diagram below shows the various functions of the CM100.



Figure 5 - CM100 Block Diagram

10.1 Touch Screen Interface

A 7" LCD color touch screen is used for the real-time display of dose delivered, dose target, spot number, and the system state. Certain setups can be modified via the touch screen interface.

10.2 Independent Termination System and Pulse Counting

The primary function of the CM100 is to provide an independent termination system for a particle therapy system in the case of a detected over-dose condition (see IEC 60601-2-64:2014 201.10.2.101.3.1.9 Independent termination system). This is accomplished by monitoring the applied dose and opening a hardware interlock in the event of a detected overdose.

The CM100 has an input (Q Pulse In) that counts raw pulses coming from a compatible electrometer such as the I128 or F100, where each pulse represents a known amount of dose. The pulses must be positive going minimum 1usec wide, maximum rate 2MHz. A cable-disconnect when using the BNC option will be automatically detected.

Typical particle therapy systems operate in one of two modes, Monitor Unit (MU) and Gigaproton (Gp). In MU mode a pulse would represent a fixed amount of charge as created in the ionization chamber (Coulombs), and in Gp mode a fixed number of beam particles. These values may also be corrected for temperature and pressure effects. The CM100 in any event is not concerned with the details of what the pulses mean; it merely counts what is sent to it.

The IEC standard requires that the device "only count when LIGHT ION BEAM enters the RADIATION HEAD". For example, the imaging process requires X-rays applied to the patient, which can be picked up by the ionization chamber. Magnetic effects from changing the scan magnet and dipole fields can also influence the readings. The CM100 provides an "inhibit" input that will inhibit the counting process for treatment particles in these cases. The CM100 accomplishes this via an "inhibit" input (Q Pulse Inhibit) that can be used to maintain a separate count of qualified pulses only. The diagram below shows this logic.



Figure 6 - Pulse inhibit logic

Two inhibit signals can be input to the CM100, and if either one is high then the beam is not inhibited.

The system thus maintains two counters: the total counts received by the unit since the last reset (unfiltered counter), and the counts allowed to pass through the inhibit filter (filtered counter). The filtered counter is used to see if the dose has overflowed. The unfiltered value is used to produce the audio output.

The termination logic starts by loading a target count, which is done automatically by the client via the EPICS interface. The target represents the <u>filtered</u> count above which the CM100 will open its "Overflow" relay to interrupt the beam and halt the treatment. This relay is normally closed, and is only opened in case of overflow or if the unit's watchdog is activated.

The filtered and unfiltered counts are published on the EPICS interface if the client wishes to display these values.

10.3 Audible Indicator

An audible "beep" is emitted by the CM100 proportional to the incoming <u>unfiltered</u> count received by the device. The beep volume can be modified in the unit's setup screens, but cannot be muted. The beep can be transmitted to an external, self-powered speaker via the external speaker jack on the unit.

10.4 Keyswitch

A key-switch is provided on the CM100 that can be used to interact with the workflow of the therapy session. The switch has three positions, OFF, PREP, and TREAT. The switch position impacts the state of the relays on the unit and is also reported via the EPICS interface.



Figure 7 - Keyswitch

A typical application would use the OFF position when the room is not active or the patient is being setup; the PREP position to load the treatment map from the planning system and to send the CM100 the overflow count; and the TREAT position to arm the system and allow the treatment to start.

The key itself provides a level of authorization to the system. The key can only be removed in the OFF state.

10.5 Start and Pause Buttons

Two buttons are provided that allow the therapist to start an irradiation or pause it. The buttons are surrounded by a colored ring – if the ring is green, then the button can be pressed; if the ring is red its function is suppressed.

CM100 User Manual



Figure 8 - Start and pause buttons

10.6 Emergency Stop Button

A latching emergency off button (EMO) is provided to allow the therapist to interrupt the treatment. The button connects directly to the EMO connections on the back of the unit. These contacts can be connected to a means of stopping the beam production.



Figure 9 - Emergency off button

10.7 Synchronization Interface

The CM100 is capable of snooping on the Pyramid inter-controller synchronization message to determine the present spot number. This is accomplished by a fiber input located on the rear connector labeled Data Link 1. Depending upon the arrangement, the signal may want to be retransmitted out the corresponding output fiber. The spot number is displayed on the touch screen interface.

Two other fiber pairs are also located on the back, Data Link 2 and Data Link 3. These presently are not used by the system.

10.8 Safety Interlocks

The back of the unit has several safety interlocks that can be used to interface to the overall therapy system.

CM100 User Manual



Figure 10 - Safety interlocks – rear panel connectors

The state of each of these relays depends upon various hardware and software states of the CM100. The states of SR4 safety relays are monitored internally to identify any malfunction such as welded contacts.

The following table is a summary of each of the interlocks and the conditions that determine the state of the contacts on the rear panel connectors:

Interlock	Hardware conditions	Software conditions
E-STOP, primary (direct connection to switch)	Open EMO switch depressed.	
E-STOP, second(ary) (SR4 safety relay)	Open when watchdog activated.	
Overflow relay (SR4 safety relay)	Open when watchdog activated.	Normally closed by software; opened when gated counts detected greater than configured overflow setting.
Pause relay (SR4 safety relay)	Pressed state transferred to software.	Normally opened by software; closed when in treatment state and pause button press detected.
Keyswitch OFF	Closed when keyswitch in the OFF position.	
Keyswitch PREP	Closed when keyswitch in the PREP position.	
Keyswitch TREAT	Closed when keyswitch in the TREAT position.	

10.9 Battery Backup

The unit is equipped with a lithium iron phosphate (LiFePO4) battery This type of battery was selected due to its inherent safety, reliability, and low toxicity.

The battery can provide power to the unit for at least twenty minutes in the event of a power interruption. In this event the therapist can manually record the total delivered dose and, where available, the spot number progress of the interrupted treatment.

The battery may degrade over time and should be checked periodically to determine that it is remains capable of powering the unit for at least twenty minutes. The recommended interval for this test is not longer than monthly.

10.10 Embedded Software

The CM100 runs the Linux operating system for the Beagle Bone Black board. The software can be updated remotely across the Ethernet interface.

11 Updating the CM100

The CM100 maintains its program and various options in files on its memory card. These files should be changed only by qualified service personnel. The procedure is listed here in case a change should be necessary, but it is strongly recommended that you contact Pyramid before proceeding with any change.



If the CM100 is in use in a medical therapy system, then following any change to setup, it is mandatory that the system should be tested and re-certified as fit for its purpose by a qualified person.

11.1 Updating the Firmware

The firmware can be updated via two methods.

11.1.1 Changing the SD card.

The SD flash memory card can be changed, although this requires opening the unit by removal of the top panel. This will allow the operating system and running program to be updated. When this is done, it is important to preserve the System.xml file from the previous card and copy it to the appropriate location (see below). It is recommended that this procedure only be carried out at Pyramid.

11.1.2 Copying program files

The device firmware applications, including ARM/Linux service and UI, and PRUSS, can be remotely updated over Ethernet. To do this, each of the files below needs to be copied to the indicated location:

File name	Location on CM100	Description
cm100_pru0.out	opt\pyramid\cm100\bin	PRUSS 0 code, responsible for counting
cm100_pru1.out	opt\pyramid\cm100\bin	PRUSS 1 code, responsible for spot detection
cm100.exe	opt\pyramid\cm100\bin	IG2 executable running in ARM, responsible for main logic, EPICS.
All files	opt\pyramid\cm100\Debug	GUI

11.2 System.xml

The System.xml file contains important setups relating to the performance of the CM100. This file should generally be changed under the guidance of Pyramid and we recommend contacting us before proceeding. Making incorrect entries can materially affect the performance of the unit.

11.2.1 Editing the File

To edit the file, perform the following steps:

- 1. The CM100 must be powered and on the network.
- 2. Find the **IP address** on the CM100 under **Setting->Network Settings**. Here we will assume the IP address is **192.168.100.123**
- 3. Access the file using the Windows Explorer at address \\192.168.100.123\bbbfilesystem\opt\pyramid\cm100\config.
- 4. Edit the file as needed and save back to the same location.
- 5. Restart CM100 using reset button on back of device

11.2.2 Attributes that can be Changed

The following attributes only should be changed:

Node	Attribute	Description
cm100state	userunits	MU – monitor unit mode.
		Gp – Gigaproton mode.
	spotcheck	true – enable spot display and sync checking during
		treatment
		false - disable
pulsecounter	pulsespertick	The number of unfiltered pulses that corresponds to
		a single audio beep.
	tonetickfreq	The frequency of the beep sound in Hz.
	prescaler	Relates the number of physical pulses to a single
		MU or Gp on the display.

11.3 Changing IP Address

The IP address can be changed by editing the ethernet communications file **/etc/network/interfaces**. In this file, you can change the primary network interface called **eth0**.

A static or dynamic IP can be set by using one of the following two lines:

iface eth0 inet dhcp iface eth0 inet static

// will give it a DHCP address.
// will give it a static address.

If using a static IP, you must follow that line with the following lines:

```
address 192.168.100.123
netmask 255.255.255.0
gateway 129.168.100.1
```

To activate the new settings, SSH into the device and run

ifdown eth0 ifup eth0

Otherwise, restart the CM100 using reset button on the back of device.

11.4 Serial Number

The CM100 serial number is kept in file **\opt\pyramid\cm100\Debug\Serial.txt**. Do not modify this file.

11.5 Loopback Address

The CM100 loopback address is kept in file **\opt\pyramid\cm100\Debug\IPAdress.txt**. Do not modify this file.

12 State Machines

The CM100 maintains two state machines that are used to control the device process workflow, an operational and a functional state. These states are reflected on the touchscreen interface and are also transmitted via the EPICS interface.



Figure 11 - CM100 state machine

12.1 Operational State Machine

The operational state machine reflects the operating state of the device.

State	Description
Running	The CM100 is running normally with no errors.
Error	The unit is in the error state, which occurs due to the following conditions:
	- Max MU/Gp threshold exceeded
	- Emergency Stop button latched
	- Watchdog timeout in CM100 processor

CM100 User Manual

- Pulse counter fail in BEAM ON state
- Max MU threshold invalid in BEAM ON state
- Pulse input cable-disconnect detected in BEAM ON state (BNC only)
- Fiber-optic communication error detected in BEAM ON state
 Host system communications synchronization heartbeat not present in BEAM ON state (if enabled)
- Spot sequence error detected (spot option must be active)

Figure 12 - CM100 operational states

The operational state is determined by several inputs, including the hardware itself and state information received over the EPICS interface. The following table summarizes the state requirements:

State	Key switch	Start button	Pause button	CMS heart beat	Ready permit	Clinical auth.	Max count and session ID	Overflow relay state	Pause relay state
Running	Any	Any	Any	Any	Any	Any	Any	Closed	Closed
Error	Any	Any	Any	Any	Not granted	Any	Not set	Open	Open

Figure 13 - Operation state conditions

12.2 Functional State Machine

Following is a list of CM100 functional states:

State	Description
Off	Keyswitch in OFF position.
Standby	Keyswitch in PREP position. Waiting for user authorization.
Preparatory	Keyswitch in PREP position. External particle therapy system and user preparing to treat.
Ready	Keyswitch in TREAT position. Waiting user start button press.
Beam On	Keyswitch in TREAT position. System is executing a treatment.

Figure 14 - Functional states

The functional state is determined by several inputs, including the hardware itself and state information received over the EPICS interface. The following table summarizes the state requirements:

CM100 User Manual

State	Key switch	Start button	Pause button	CMS heart beat	Ready permit	Clinical auth.	Max count and session ID	Overflow relay state	Pause relay state
Off	OFF	Not latched	Not latched	CMS heart beat	Not granted	Not granted	Not set	Closed	Closed
Standby	PREP	Not latched	Not latched	Any	Any	Not granted	Not set	Closed	Closed
Prep	PREP	Not latched	Not latched	Any	Any	Granted	Any	Closed	Closed
Ready	TREAT	Not latched	Any	Active	Granted	Granted	Set	Closed	Closed
Beam On	TREAT	Latched	Not latched	Active	Granted	Granted	Set	Closed	Closed if coming from Prep state. Open if coming from BeamOn state.

Figure 15 - Functional state conditions

13 User Interface

The CM100 has a 7" color touchscreen that provides a graphical user interface for the therapist. The touch screen allows the therapist to navigate screens and set certain parameters. The home screen after boot-up is shown below.



Figure 16 - Home screen

13.1 Status LEDs

13.1.1 **AC**

The AC LED shows the state of the AC power. When green, power is plugged in, when red power is disconnected.

13.1.2 **CMS**

The *CMS* LED is lit when the CMS is properly loaded and running. This module is necessary when the CM100 is operated in conjunction with the Pyramid Scan/Dose system.

13.1.3 **Cable**

The *Cable* LED is lit when the pulse counting cable is properly connected. The device can detect cable disconnect as one of its safety features. The cable must be connected for normal operation.

13.2 System State

This section displays the system state, which can take on the following values:

State	Description
Off	Keyswitch in OFF position.
Standby	Keyswitch in PREP position.
Prep	Keyswitch in PREP position.
Ready	Keyswitch in TREAT position.
Beam On	Keyswitch in TREAT position, treatment is active.

13.3 System Messages

The CM100 will display helpful messages and errors in this location from time to time.

13.4 Date/Time

This section displays the system date and time, and is only displayed if enabled in the Settings setup screen.

13.5 Home Button

The *Home* button returns the software to the home page, as shown above. This can be pressed while on any other screen to return home.

13.6 Setup Button

The *Setup* button gives access to the setup menu, shown below. Some parameters can be changed on these setups, others can only be viewed.

CM100 User Mo	anual		
AC CMS CABLE	Off Standby Prep	Ready Beam On Syst	em in OFF State
	Network Settings	About	Screen Settings
	-	í	
	Audio Info	Connection Info	
((()	(((

Figure 17 –Setup menu

Each button on the screen can be pressed to view the indicated setup or information page.

13.6.1 Network

The *Network Settings* button allows the unit's IP address, subnet mask, and default gateway to be displayed. These settings can only be changed by modifying the unit's XML file, described in a separate section of this manual.

O AC	О смs	CABLE	Off ▲	Standby	Prep Read	y Beam On	System in OFF State
		Netwoi	rk Set	ttings			
Î		IP Ac	dress	:			
C			192	168	100	15	
		Subr	net Ma	sk:			
			255	255	255	0	
		Defa	ult Ga	teway:			
	4		192	168	100	1	
))						



13.6.2 About

The *About* page give information about the version that is installed in the device.

O AC	О смs	CABLE	Off	Standby	Prep	Ready	Beam On	System in OFF State	
10		About	t						
1		Ver	sion:						
C			СМ100	1.2.0	0.0				
199			PRU0:	1.0.0	0.0				
			PRU1:	1.1.0	0.0				
			Hardwa	are: 2.0.0	0.0				
			IG2:	2.6.8	3.0				
))								

Figure 19 - About screen

The versions displayed are:

- CM100 Main software version for the device
- PRU0 Sub-processor used to handle counts and scaling
- PRU1 Sub-processor used to handle spot count information from synchronization message.
- Hardware the version of the hardware (PCB).
- IG2 library used to transmit data via EPICS.

13.6.3 Screen settings

The Screen Settings page modifies the visual behavior of the screens.

CM100 Use	er Manual							
	Off Off	Standby	Prep	Ready I	Beam On	System in OFF State	2/3/18	
AC CMS C	CABLE						15:34	
S	creen Sett	ings						
	Show Date/	Time:)~				
∎))								

Figure 20 - Screen settings

The following options are available:

• Show Date/Time – if enabled displays the date and time at the top right of the screen.

13.6.4 Audio settings

The *Audio Settings* page allows settings relating to the audio output to be modified and displayed.

AC	О смs	CABLE	Off ▲	Standby	Prep	Ready	Beam Or	n S	System in	OFF State	2/3/18 15:34
10	1	Audio	Infor	mation							
		1	Pulses p	oer Tick:	20			Pres	caler: 0.00	01	
			Tone Tio	k Freque 20(ncy: 00			Tick S	Source: ered		
)										

Figure 21 - Audio information

The following options are available (all read-only):

• Pulses per Tick – the number of detected unfiltered pulses that cause a single beep to be output.

PTC2-2120709629-4089

- Tone Tick Frequency The frequency of the tone emitted by the CM100 in Hz.
- Prescaler Relates the number of physical pulses to a single MU or Gp on the display.
- Tick Source The input associated with the tick, filtered or unfiltered.

13.6.5 Connection settings

When the *Connection Settings* button is pressed, the settings of the dip switches located on the bottom of the unit are displayed (see section 15.2.1). The settings of fiber vs. BNC for the pulse counting and inhibit A and B connections, and the disconnect method and mode are displayed.

AC CMS	CABLE Off Standby	Prep Ready Be	am On System in OF	F State 2/3/18 15:35
	Connection Informa	tion		
	Disconnect Mod	le:	Inhibit A:	
Q	Fiber		Fiber	
	Disconnect Met	hod:	Inhibit B:	
	Hardware Duty	Cycle	Fiber	
(()	Pulse Mode: Fiber			

Figure 22 - Connection settings

13.7 Audio Volume

When the *Audio Volume* button is pressed, a popup will appear that allows the volume to be increased or decreased. The sound cannot be completed muted, however.



Figure 23 - Audio volume

13.8 Session Information

Specific information relating to the session is included in this section.

13.8.1 Count total

The *Count Total* field displays the filtered count in Gp or MU., that is, the count after filtering according to the pulse inhibit gate.

13.8.2 Unfiltered count

The *Unfiltered Count* field displays the unfiltered count in Gp or MU., that is, the count received by the device irrespective of the pule inhibit gate.

13.8.3 Count threshold

The *Count Threshold* field is the threshold of filtered counts above which the CM100 will open the overflow interlock and signal an error. This value is set by the client software over the EPICS interface.

13.8.4 Session ID

The *Session ID* field is a unique identifier associated with a treatment. This value is set by the client software over the EPICS interface.

13.9 Pulse Rate

The *Pulse Rate* section displays the rate of dose delivery. A meter shows this graphically, as well as a numeric value.

13.10 Spot Count

The *Spot Count* field is displayed only if Spot checking is established in the setup. This value will be the spot number known across all energy layers.

13.11 Error Messages

The following error message can be displayed by the CM100:

Code	Text	Description
300	MU(Gp) threshold violation	Detected filtered pulses exceeds the MU or Gp threshold transmitted by the client.
301	EMO button pressed	The EMO button was depressed.
302	PRU ACK is NOT active	The PRU sub-processor is not active (fatal error).

303	PRU heartbeat not active	The PRU sub-processor heartbeat is not found (fatal error).
304	Counter hardware failure	The counter hardware has failed (fatal error).
305	Corrupt MU threshold	The MU threshold has been found to be corrupted.
306	Pulse cable disconnected	The pulse counter cable has been detected as disconnected.
307	CMS communication failure	The CMS client program communication has failed.
308	Keyswitch moved during BeamOn state	The keyswitch was moved during the treatment state, causing the session to terminate.
309	Timeout enabling PRU threshold check	An error was detected while setting the max threshold check.
310	Fiber spot error	Spot display has been enabled, and during the treatment the spot synchronization message was not detected as necessary.
311	Power disconnected	Power was disconnected while the unit was in the Beam On state.
400	<pre>processOpState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
401	<pre>processIntFuncState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
402	<pre>processOnState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
403	<pre>processPrepState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
404	<pre>processTreatState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
405	<pre>setOpState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
406	<pre>setIntFuncState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
407	<pre>setOnState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.

408	<pre>setPrepState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
409	<pre>setTreatState() unknown state</pre>	This error should not normally occur. Please contact Pyramid if one of these errors occurs.
500	keyswitch not in prep state when required	The keyswitch should be in the Prep state.
501	control system paused, waiting for resume	This error is normally not propagated to the GUI.
502	waiting on control system ready permit	This may occur when the Ready Permit is revoked by the client.

The error messages will be displayed with the error number across the entire portion of the screen. The user will need to move the keyswitch back to the PREP position and press the CLEAR button to continue.



Figure 24 - Error message display and CLEAR button

14 Example Use Case

Specific use of this device in the client system is flexible depending upon the configuration. The sequence diagram below shows a typical use case for the device.





PTC2-2120709629-4089

The sequence shown above shows a typical sequence that is used demonstrating interaction between the therapist, the CM100, and the treatment control system. The specific details of this sequence will be determined by the therapy system manufacturer.

14.1 Off

The system boots up in the *Off* state, waiting for action from the therapist and the client system.



Figure 26 - CM100 in the Off State

14.2 Standby

The first step required to prepare for treatment requires that the therapist (insert and) switch the key to the PREP position. This transitions the CM100 to the *Standby* state, ready to accept authorization from the client over the EPICS interface. A reminder message is displayed in the message area.



Figure 27 - CM100 in the Standby state

14.3 Prep

Once the authorization message is sent from the client, the system transitions to the *Prep* state. The client now must load the following information over the EPICS interface to the CM100:

- Session ID a unique numeric identifier associated with the treatment.
- Count Threshold total count (in MU or Gp) above which the overflow relay will open during the treatment, should the detected filtered count exceed this value.
- Ready Permit standard permit to authorize treatment.



Figure 28 - CM100 in Prep state

The message area shows which of these items have been received by lighting up the item. In the above example the Session ID and Count Threshold (Max) have been set and the system is still awaiting the Ready permit. The ready permit should only be transmitted after the session ID and max threshold have been sent, as well as any other conditions deemed necessary by the client system.

When all three parameters have been received, the CM100 displays the message "Turn Keyswitch to 'TREAT' informing the therapist of the next step. The TREAT label on the key switch will also be blinking at this point.



Figure 29 - CM100 in Prep state ready to transition to Treat state

14.4 Treat

Once the keyswitch is turned to the TREAT position, the CM100 enters the *Treat* state.

CM100 Us	er Manual					
AC CMS	Off Standby	Prep	Ready Beam On	Pres	ss 'START' When Read	y 2/3/18 15:44
	Count Total			-	500 400	600
		0	.000	GP	300	700
	Unfiltered Count		0.000	GP	100	900
	100).(000	GP	kCounts Per	Second
(()	Session ID		555			

Figure 30 - CM100 in Ready state

The ring around the Start button should be flashing green at this point, signaling to the therapist that the treatment can be started by pressing this button.

14.5 Beam On

After the Start button is pressed, the CM100 enters the *Beam On* state. The ring around the Pause button should flash blue, indicating that the treatment can be paused at any time. When beam is present the CM100 should begin to read pulses, displaying them in the Count Total and Unfiltered Count fields (note that these fields are always active, even when not in the Treat state).

	Off Standby Prep Ready Beam On Treatment In Progress	2/3/18 15:49
	Count Total	600
	99.000 GP 200	700
	Unfiltered Count 99.000 GP	900
	100.000 GP	econd O
(()	Session ID 555	

Figure 31 - CM100 in Beam On state

14.6 Treatment Complete

The client system signals the CM100 when the treatment is completed over the EPICS interface. When this occurs, the system goes back to the *Prep* state, retaining all information settings from the session.



Figure 32 - CM100 when the treatment is completed

After turning the keyswitch back to PREP, the ready permit is automatically revoked by the CM100. The client system should enforce at this point the conditions needed to start a new session, including sending a new session ID and Max threshold.

14.7 Pausing an Ongoing Treatment

The therapist can press the Pause button during the Beam On state at any time, pausing the treatment. In this case the Start button will flash green, indicating that the treatment can be restarted.

14.8 EMO Pressed

A special case occurs when the EMO button is pressed on the unit. When pressed, an error message will be displayed, and a CLEAR button appears. The system opens the EMO relays and signals the error to the client.

At this point the only available action is to move the keyswitch back to the PREP or OFF position, unlatch the EMO, press CLEAR on the touchscreen, and restart the entire setup process from the Prep state. The user authorization does not need to be retransmitted by the

client. If the system was delivering beam, then the session is terminated and cannot be restarted.



Figure 33 - CM100 Emergency Off

14.9 Moving the Keyswitch off the TREAT position during Beam On State

During the Beam On state, until the client notifies the CM100 that the treatment is finished, the therapist should not normally turn the keyswitch off the TREAT state. Doing so will cause the error "Keyswitch moved during BeamOn state" to be displayed, revoking the permit, session ID, and Max threshold. The therapist must turn the key back to PREP and press the CLEAR button at the top right of the touchscreen.

15 Connectors

15.1 Rear Panel

All connections for the CM100 are located on the rear panel, shown in the drawing below.



Figure 34 - CM100 rear panel

15.1.1 Ground lug

M4 threaded stud. To mate with M4 ring lug.

15.1.2 Power input, 24 VDC

Two-pin Redel PXG.M0.2GG.NG female. To mate with Redel PAG.M0.2 type or PFG.M0.2 type free plugs. Suitably terminated 24 V power supplies and leads are available from Pyramid Technical Consultants, Inc.



15.1.3 Data links

Six HFBR ST bayonets suitable for 1 mm plastic or 200 μ m silica fiber. 664 nm (visible red) light. Dark casing = receiver, light casing = transmitter.



Data Link 1 is used by the CM100 to intercept synchronization messages from the Pyramid Scan Dose system. This provides the spot information and a means to determine that the system is active. Note that it is not necessary to connect the output fiber if the synchronization signal is repeated with an X14 device.

Data Link 2 and Data Link 3 are provided for future enhancements.

15.1.4 **Q Pulse inhibit**

Three HFBR ST bayonets suitable for 200 μ m core step index silica fiber. 664 nm (visible red) light. Dark casing = receiver, light casing = transmitter.

Two Lemo coax 00 jack (female). Input termination 2.5 kohm.



The pulse inhibit is implemented via fiber (inputs A and B) or COAX (inputs A and B), depending upon the setting of the configuration switches located on the bottom of the unit. Only one input is required to inhibit the pulse counting. The Sync connection is not used.

15.1.5 Pulse in

One HFBR ST bayonet suitable for 200 μ m core step index silica fiber. 664 nm (visible red) light. Dark casing = receiver, light casing = transmitter.

One Lemo coax 00 jack (female). Input termination 2.5 kohm.



The pulse input is implemented via fiber or COAX input, depending upon the setting of the configuration switches located on the bottom of the unit.

15.1.6 Interlocks



Emergency stop

Phoenix 1827910 7-position 3.81 mm. Locking mating connector Phoenix 1851287 (supplied). Emergency off switch contacts. Opto-coupled digital input.

CM100 User Manual

1	n/c	5	Switch contact
2	Redundant contact	6	Opto input anode (1 k Ω)
3	Redundant contact	7	Opto input cathode (1 k Ω)
4	Switch contact		I

<u>Relay</u>

Phoenix 1827923 8-position 3.81 mm. Locking mating connector Phoenix 1851290 (supplied). Contact closure indication of pause state and counter exceeded limit state.

1	Closed in position 1	7	Closed in position 2
2	Position 1 common	8	Position 2 common
3	Open in position 1	9	Open in position 2
4	Closed in position 3	10	24 V rtn
5	Position 3 common	11	24 V rtn
6	Open in position 3	12	+24 VDC out

<u>Keyswitch</u>

Phoenix 1787056 12-position 3.5 mm. Locking mating connector Phoenix 1790331 (supplied). Changeover (SPDT) indication of keyswitch state to remote devices.

Position 1: Off Position 2: Prepare Position 3: Treat (direct connection)

1	Closed in position 1	7	Closed in position 2
2	Position 1 common	8	Position 2 common
3	Open in position 1	9	Open in position 2
4	Closed in position 3	10	24 V rtn
5	Position 3 common	11	24 V rtn
6	Open in position 3	12	+24 VDC out

15.2 Bottom Panel

The bottom of the unit has a removable hatch that gives access to various configurations and connects the battery.

15.2.1 Switches

Dip switches used for configuration.



Altering switch settings will affect the operation of the device. Changes should only be made by qualified personnel. If the CM100 is in use in a medical therapy system, then following any change to setup, it is mandatory that the system should be tested and re-certified as fit for its purpose by a qualified person.



1	INHIBIT A	Selects the Pulse Inhibit input between FIBER or BNC.
2	INHIBIT B	Selects the Pulse Inhibit input between FIBER or BNC.
3	PULSE TYPE	Selects the pulse cable type between FIBER or BNC.
4	PULSE TYPE	Must always be the same setting as 3 above.
5	DISCONNECT	Always set this to DUTY CYCLE.

15.2.2 Battery

The battery is connected to the bottom hatch. It connects to the PCB via three two pin Molex connectors. Each connector is unique, connecting the battery's power, thermistor, and indicator connections to the unit.

CM100 User Manual



16 Communications Interface

The CM100 connects to the client via an Ethernet interface. The device acts as an EPICS server, and all information with client applications is exchanged using this protocol. For more information on EPICS please go to the EPICS website at <u>http://www.aps.anl.gov/epics/index.php</u>.

16.1 EPICS Process Variables

The PV names will start with the CM100 unit's serial number so that they are unique. The following parameters are supported. Network Configuration

The CM100 uses EPICS communication over standard local area network hardware. Addressing is using the IP4 standard, and it supports static and dynamic (DHCP) address assignment. The device can be configured via the PTC DiagnosticG2, via the serial interface, or by your own host software using the appropriate procedure calls.

Most control and data acquisition systems are set up with fixed addresses assigned by the network administrator. It is also typical to isolate such networks from the internet to prevent unauthorized access, and to allow operation without firewalls which can disrupt communications.

In order for the host computer and the CM100 to communicate, they must be within the same subnet. It is typical to limit a local network to 256 addresses by setting the IP4 subnet mask to 255.255.255.0. Then the CM100 and the host must have the first three bytes of their addresses common, and must differ in the last byte. For example, the host could be 192.168.100.11 and the I128 at 192.168.100.20. The last byte must also not conflict with any other devices on the same subnet. Addresses with last byte 0 and 255 are reserved for special functions in TCP/IP. See the section on "File Setups" for information on changing the IP address of the unit.

PV Name	Description	Direction	Туре
[CM100serial#]_version	Version number, in the form: [byte0].[byte1].[byte2].[byte3] To be read-only for all clients.	R	Integer32
[CM100serial#]_arm_heartbeat	Heartbeat counter incremented by master CM100 service process running on ARM. This value should increment at a high rate and will overflow. To be read-only for all clients.	R	Integer32
[CM100serial#]_cm100_heartbeat	Heartbeat counter incremented by the master CM100 state machine process.	R	Integer32
[CM100serial#]_pru0_heartbeat	Heartbeat counter incremented by slave CM100 service process running on PRU0. This value should	R	Integer32

Following is a list of all available process variables:

	increment at a high rate and will overflow. To be read-only for all clients.		
[CM100serial#]_pru1_heartbeat	Heartbeat counter incremented by slave CM100 service process running on PRU1. This value should increment at a high rate and will overflow. To be read-only for all clients.	R	Integer32
[CM100serial#]_pause	State of physical PAUSE button. To be read/write for local CM100 UI, read-only for client applications.	R	Boolean
[CM100serial#]_start	State of physical START button. To be read/write for local CM100 UI, read-only for client application.	R	Boolean
[CM100serial#]_emo	State of physical EMO button. To be read-only for all clients.	R	Boolean
[CM100serial#]_pulse_discon_hw	Cable disconnect state, as detected by hardware.	R	Boolean
[CM100serial#]_ready_permit	Permit issued by user that allows the CM100 to transition out of the PREPARATORY state. 0x55555555 = PERMIT GRANTED other = PERMIT NOT GRANTED Other informational codes (ALL PERMIT NOT GRANTED): 0x1111111 = Control system process paused. 0x22222222 = Control system process finished.	RW	Integer32
[CM100serial#]_session_id	User-defined session ID. -1 = INVALID other = VALID Note: Changing this value will reset the treatment MU/GP integrator: [CM100serial#]_treat_count_mu	RW	Integer32
[CM100serial#]_cumulative_spot	Cumulative spot number, read from ScanDose fiber connection. -902 = Fiber timeout.	R	Integer32

PTC2-2120709629-4089

	-901 = Fiber message corrupt.		
	-900 = Fiber system uninitialized.		
	other = Valid spot number.		
[CM100serial#]_spot_count_enabled	1: Spot count check enabled	R	Boolean
	0: Spot count check disabled		
[CM100serial#]_count_reset	Reset count command. Calling this command will reset the unfiltered and filtered pulse count integrators to 0. To be written by client application or local UI. Argument ignored.	RW	Boolean
[CM100serial#]_key_position	State of keyswitch position readback. To be read-only for all clients.	R	Integer32
[CM100serial#]_relay1	Relay 1 command	RW	Boolean
[CM100serial#]_relay2	Relay 2 command	RW	Boolean
[CM100serial#]_count_unfiltered	Snapshot of unfiltered pulse count integrator, lower 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_unfiltered_ovf	Snapshot of unfiltered pulse count integrator, upper 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_filtered_qep	Snapshot of filtered pulse count integrator (QEP), lower 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_filtered_ovf_qep	Snapshot of filtered pulse count integrator (QEP), upper 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_filtered_dm	Snapshot of filtered pulse count integrator (DMTimer), lower 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_filtered_ovf_dm	Snapshot of filtered pulse count integrator (DMTimer), upper 32 bits. To be read-only for all clients.	R	Integer32
[CM100serial#]_count_rate	Snapshot of filtered/unfiltered (configurable) pulse count rate. To be read-only for all clients.	R	Integer32
[CM100serial#]_prescaler	Prescaler value, in units of MU/pulse.	R	Float
[CM100serial#]_treat_count_mu	Snapshot of treatment dose integrator, in units of MU/GP.	R	Float

[CM100serial#]_treat_count_mu_unfiltered	Snapshot of treatment dose integrator, in units of MU/GP, for unfiltered count	R	Float
[CM100serial#]_treat_count_max	Maximum MU/GP threshold for the irradiation (typically 120% of expected delivery).	RW	Float
[CM100serial#]_count_user_units	User Units of display counts. Will be either MU or GP	R	String
[CM100serial#]_audio_ppt	Number of pulses in the audio pulse integrator required to output an audio tick. The audio pulse integrator is reset after this threshold is reached. This parameter controls the audio tick rate to pulse rate relationship.	RW	Integer32
[CM100serial#]_audio_spt	Number of audio samples output per tick. This number is a factor in volume and tick sound.	RW	Integer32
[CM100serial#]_audio_tone_tick	Audio tick sound selection. 0 = Simple tick 1 = Sine wave tick	RW	Boolean
[CM100serial#]_audio_tone_tick_freq	Frequency of audio sine wave tick. This is only applicable if audio_tone_tick == 1.	RW	Integer32
[CM100serial#]_audio_int_period	Polling period for the audio process, in ms. The process will examine the audio pulse integrator at this period.	RW	Integer32
[CM100serial#]_audio_volume	Volume level (0-100%)	RW	Integer32
[CM100serial#]_func_state	Functional state of CM100. See state machine and sequence diagrams for more details. To be read-only for all clients. OFF = 0 STANDBY = 1 PREP = 2 READY = 3 BEAM ON = 4	R	Integer32
[CM100serial#]_op_state	Operational state of CM100. See state machine and sequence	R	Integer32

PTC2-2120709629-4089

	diagrams for more details. To be read-only for all clients. RUNNING = 0 SERVICE = 1 ERROR = 2		
[CM100serial#]_error	Array of error information: [0] = Error Code [1] = 1: Clearble error 0: Blocking error	R	Float[3]
[CM100serial#]_op_clear_error	Clear Operational Error command. Argument ignored.	RW	Boolean
[CM100serial#]_prepared_conditions	Boolean value indicating if the Prepared Conditions have been met and the user can turn the keyswitch to Treat	R	Boolean
[CM100serial#]_auth_code	Authorization code from external control system to CM100. Integer bitmask format. See User Authorization for details. BIT0 = CLINICAL AUTH BIT1 = SERVICE AUTH	RW	Integer32
[CM100serial#]_control_system_heartbeat	Heartbeat counter incremented by master CMS service process to monitor connection between DCS and CM100. This value should increment at a high rate and will overflow. To be written by CMS.	RW	Integer32
[CM100serial#]_ control_system_connected	Boolean value indicating state of heartbeat connection to external control system. TRUE = CONNECTED FALSE = NOT CONNECTED	R	Boolean
[CM100serial#]_ip_address	Integer representation of IP address	RW	Integer32
[CM100serial#]_subnet_mask	Integer representation of subnet mask	RW	Integer32
[CM100serial#]_default_gateway	Integer representation of default gateway	RW	Integer32

PTC2-2120709629-4089

[CM100serial#]_ethernet_set	Push down Ethernet settings command. Calling this command will push down the Ethernet's IP address, static/dynamic, subnet mask, and default gateway. To be written by client application or local UI. 1 = Set Ethernet Settings 0 = Cancel	RW	Boolean
[CM100serial#]_power_disconnected	1= AC power connected 0 = AC power not connected	R	Boolean

17 Fault Finding

Symptom	Possible Cause	Confirmation	Solution
Device will not boot or communicate	Damage to CM100.		Contact Pyramid Technical Consultants, Inc.
	Failed firmware update.		Contact Pyramid Technical Consultants, Inc.
	Network IP addresses between CM100 and client not compatible.	Check network addresses and IP masks are compatible.	Change addresses or configuration as needed.
	Duplicate Ethernet address of CM100.		Change CM100 address.
Pulses do not count	Bad pulse input cable.	Check "Cable LED" on GUI indicating that the cable is connected properly.	Replace cable.
	Broken fiber optic receiver used for incoming pulses.	Examine if there is any physical damage to receiver or if it is loose on the PCB.	Send unit back for repair.
	Broken fiber optic transmitter on client system used for pulse inhibit.	Look at pulse inhibit source and confirm it is stuck in the lit state.	Fix pulse inhibit transmitter.
	Damage to CM100.		Contact Pyramid Technical Consultants, Inc.
Pulse count is not accurate	Improperly terminated BNC cable.	Verify cable termination.	Terminate cable with 50 ohm.

PTC2-2120709629-4089

CM100 User Manual

	Interference from noise sources.	Investigate noise sources.	Eliminate noise source or reposition pulse input count cable.
	Client system is operating in different mode from CM100 display (MU/Gp).	Examine client system mode vs. label on CM100.	Set label on CM100 display to match client system mode.
	Damage to CM100.		Contact Pyramid Technical Consultants, Inc.
Pulse counts seen when no beam present	Electrometer connected to the CM100 has background offset.	Check raw data from electrometer.	Apply offset correction and/or locate and remove source of offset.
	Electrometer connected to the CM100 is detecting noise.	Check raw data from electrometer.	Remove source of noise.
System does not start when pressing the "Start" button.	Key switch not placed in "TREAT" position.	Examine key position.	Follow workflow procedures.
	Client has not properly loaded parameters and authorization due to network or other client issues.	Check network connection and client software.	Fix network connection or resolve client software issues.
	Damage to CM100.		Contact Pyramid Technical Consultants, Inc.
System does not pause when pressing the "Pause" button during a treatment.	Damage to CM100.		Contact Pyramid Technical Consultants, Inc.

18 Maintenance

There are no user-serviceable parts inside the CM100.

18.1 Operation check

A regular quality assurance check should be made to check that the following function as expected:

- Audio tone is clearly audible in the working environment
- Interlock relay (beam is stopped) opens on count overflow
- Disconnected pulse signal cable is detected by the unit
- Pause button functions (beam is stopped)
- Turning keyswitch from Treat back to Prep stops beam delivery
- EMO button functions (primary contacts open, connected systems are disabled)

The checks will be embedded in machine QA protocols for a particular facility.

18.2 Checking the Battery

To perform the battery check, remove the power cord to the CM100 with the unit powered on. Observe after 30 minutes that the system behaves normally. Reconnect the power cord.

18.3 Battery Pack Replacement

To replace the battery pack, unscrew the four screws on the bottom panel of the unit.



Remove the three Molex connectors connecting the battery to the PCB. Attach the new battery by reattaching the Molex connectors. Allow the battery to charge up for four hours, then check the battery as in the section "Checking the Battery".

It is recommended to replace the battery pack every three years or less.

19 Returns Procedure

Damaged or faulty units cannot be returned unless a Returns Material Authorization (RMA) number has been issued by Pyramid Technical Consultants, Inc. If you need to return a unit, contact Pyramid Technical Consultants at support@ptcusa.com, stating

- model
- serial number
- nature of fault

An RMA will be issued, including details of which service center to return the unit to.

20 Support

Manual and software driver updates are available for download from the Pyramid Technical Consultants website at <u>www.ptcusa.com</u>. A secondary site can be found at <u>www.ptceurope.com</u>. Technical support is available by email from support@ptcusa.com. Please provide the model number and serial number of your unit, plus relevant details of your application.

21 Disposal

We hope that the CM100 gives you long and reliable service. The CM100 is manufactured to be compliance with the European Union RoHS Directive 2002/95/EC, and as such should not present any health hazard. Nevertheless, when your CM100 has reached the end of its working life, you must dispose of it in accordance with local regulations in force. If you are disposing of the product in the European Union, this includes compliance with the Waste Electrical and Electronic Equipment Directive (WEEE) 2002/96/EC. Please contact Pyramid Technical Consultants, Inc. for instructions when you wish to dispose of the device.

Revision History

The release date of a Pyramid Technical Consultants, Inc. user manual can be determined from the document file name, where it is encoded YYMMDD. For example, CM100_UM_080105 would be a CM100 manual released on 5 January 2008.

Version	Changes
CM100_UM_180228	First general release.

DOCUMENT APPROVAL

This document has been reviewed and approved by the following individuals:

W. Nett/Author	 J. Iken
X J. Gordon	C. Pendleton
X J. Nett	P. Boisseau