

Evaluation of activity and function before and immediately after the provision of a microprocessor knee in individuals with transfemoral amputation

Prosthetics and Orthotics International, 49(4), 368–377. (2025).
DOI: 10.1097/PXR.0000000000000449. [OpenAccess](#)

Products

C-Leg, Kenevo

Major Findings

With microprocessor-controlled knees (MPKs) (C-Leg, Kenevo, Rheo, Plie) compared to non-MPKs:

- **MPK provision improved functional test outcomes and self-reported scores compared to nMPK use:**
- 33.3% higher activities-specific balance confidence (ABC, $p<0.05$)
 - 15.4% higher reintegration of normal living index (RNLI, $p<0.05$)
 - 15.3% higher walking performance in clinical tests (2MWT, $p<0.05$)
 - 30.4% higher patient-reported prosthesis function (PEQ total, $p<0.05$)




Figure 1. Mean values of self-reported outcomes (ABC, RNLI, PEQ total) and 2MWT. *indicates a statistically significant difference ($p<0.05$).

→ **Physical activity measured by number of steps, number of sit to stand events, active, upright and sedentary time did not exhibit statistically significant changes within the 4 weeks of the study**

Population	Subjects:	14 (12 male, 2 female)
	Previous nMPK:	OB: 3R60, 3R80, 3R95 Others: Mauch, Endolite PSPC, NOP4/NOP5, NOFM2, ESK, Smart IP, Paso
	Provided MPK:	C-Leg (9), Kenevo (1), Plie (3), Rheo (1)
	Amputation causes:	trauma (3), peripheral vascular disease (3), tumor (3), congenital (2), pain-related (2), infection (1)
	Mean age:	57 years (range: 41-78 years)
	Mean time since amputation:	22 years (range: 1-44 years)
	MFCL:	not reported

Study Design

Observational pre-post cohort study:



The individuals were assessed before and after a 4-week microprocessor-controlled knee (MPK) trial following routine clinical provision. The first data collection was performed with their previous nMPK knee. The data collection included self-reported scores and functional tests. Also, on seven consecutive days following both days of data collection, activity monitoring was performed.

Results

Functions and Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic Energy Consumption	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health Economics

Category	Outcomes	Results for MPKs vs nMPKs *	Sig. ^a					
Level Walking	Activities-Specific Balance Confidence Scale (ABC)	Mean self-reported balance confidence (ABC scores) significantly improved after MPK provision. 10 out of 14 patients reported an at least 20% higher score after MPK provision.	++					
		<table><tr><th>nMPK</th><th>MPK</th><th>Difference</th></tr><tr><td>60.1</td><td>82.05</td><td>33.3%</td></tr></table>	nMPK	MPK	Difference	60.1	82.05	33.3%
	nMPK	MPK	Difference					
	60.1	82.05	33.3%					
		Higher ABC scores with MPK correlated with - a higher distance in the 2MWT (r=0.62). - a lower time in the TUG (r=-0.57). - a higher number of sit to stand events (during actimetry) (r=0.69). Higher ABC scores with nMPK correlated with higher sedentary times (r=0.55).	++					
2 Minute Walk Test (2MWT) (m)	MPK use significantly increased short-distance walking capacity, measured with 2MWT. The improvement in distance was found in 12 out of 14 patients.	++						
	<table><tr><th>nMPK</th><th>MPK</th><th>Difference</th></tr><tr><td>96.6 ± 32.1</td><td>112.3 ± 34.4</td><td>15.3%</td></tr></table>	nMPK	MPK	Difference	96.6 ± 32.1	112.3 ± 34.4	15.3%	
nMPK	MPK	Difference						
96.6 ± 32.1	112.3 ± 34.4	15.3%						
	Higher 2MWT results during the use of MPKs correlated with - a higher number of steps (during actimetry) (r=0.58).	++						
Timed Up and Go Test (TUG) (s)	Similar values were measured at both time points	0						
	<table><tr><th>nMPK</th><th>MPK</th><th>Difference</th></tr><tr><td>29.8 ± 11.1</td><td>28.9 ± 15.9</td><td>8.7%</td></tr></table>	nMPK	MPK	Difference	29.8 ± 11.1	28.9 ± 15.9	8.7%	
nMPK	MPK	Difference						
29.8 ± 11.1	28.9 ± 15.9	8.7%						

Category	Outcomes	Results for MPKs vs nMPKs *	Sig. ^a						
		Lower TUG time during the use of MPKs and nMPKs correlated with - A higher number of steps (during actimetry) (MPK: r=-0.77; nMPK: r=-0.77). - A higher 2MWT (MPK: r=-0.83; nMPK: r=-0.59). - A higher upright time (MPK: r=-0.63; nMPK: r=-0.68). For nMPK use a higher TUG time was also correlated with lower active time (r=-0.69)	++						
Activity, Mobility, Activities of Daily Living (ADLs)	Actimetry	There were no significant changes in the activity-based parameters: number of steps, active time, upright time, sedentary time, number of sit to stand events (high degree of intersubject variability, limited mean difference between time points).	0						
Preference, Satisfaction, Quality of Life (QoL)	Reintegration to normal living index (RNLI)	MPK provision significantly enhanced perceived participation in daily and social life. For the majority the improvement was >10%. <table><tr><th>nMPK</th><th>MPK</th><th>Difference</th></tr><tr><td>77.4</td><td>89.6</td><td>15.4%</td></tr></table> A higher RNLI was correlated with a lower number of steps (r=-0.61) and lower active time (r=-0.72) with nMPK.	nMPK	MPK	Difference	77.4	89.6	15.4%	++
nMPK	MPK	Difference							
77.4	89.6	15.4%							
	Prosthesis Evaluation Questionnaire (PEQ)	Users perceived substantial and significant improvements in prosthesis performance and usability. <table><tr><th>nMPK</th><th>MPK</th><th>Difference</th></tr><tr><td>555.3</td><td>732.8</td><td>30.4%</td></tr></table> A higher PEQ score was correlated with higher sedentary time with nMPK (r=0.57)	nMPK	MPK	Difference	555.3	732.8	30.4%	++
nMPK	MPK	Difference							
555.3	732.8	30.4%							

^a No difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.). significance set at $p<0.05$; trends set at $0.1>p>0.05$.

* Results are displayed as mean or mean $\pm$ standard deviation if not specified otherwise.

Pearson correlation coefficients (r):

direction - Positive (0 and +1), negative (0 and -1), no correlation ($r=0$);

magnitude - Weak (<0.3), Moderate ($0.3-0.5$), Strong (>0.5).

Author's Conclusion

„Provision of an MPK improved self-reported mobility, function, and prosthetic satisfaction. This highlights an improved outcome after MPK provision within a relatively short period of follow-up observed in this study. There was a large degree of intra- and intersubject variability in the activity-based parameters monitored over a 7-day period, pre- and post-MPK provision. Future studies are needed to explore the role of activity monitoring and step count in assessing outcomes in individuals with lower limb loss. This should include patient-perceived benefit and an observed change in function.“ (Caggiari, 2025)

Author's Affiliation

¹ Skin Sensing Research Group, School of Health Sciences, University of Southampton, Southampton, United Kingdom

² The Dorset Prosthetic Centre, Royal Bournemouth Hospital, Bournemouth, United Kingdom

³ Portsmouth Enablement Centre, Portsmouth Hospitals NHS Trust, St Marys Community Health Campus, Portsmouth, United Kingdom

⁴ People Powered Prosthetics Research Group, University of Southampton, Southampton, United Kingdom

⁵ Bioengineering Science Research Group, Faculty of Engineering and Physical Sciences, University of Southampton, Southampton, United Kingdom

©2026, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.