

Agilium Freestep

Biomechanics – Gait analysis

Major Findings

With Agilium Freestep:

The knee adduction angular impulse is significantly reduced

→ Agilium Freestep in different adjustments: significant ($p < 0.05$) decreased impulse in all conditions; by 5% in neutral position, by 3% in varus position, by 7% in valgus position (Fantini-Pagani et al. 2013)

The knee adduction moment (KAM, fist peak) is significantly reduced

→ Patients: the KAM was reduced by 22% (Schmalz et al. 2011)

→ Subjects: the KAM was reduced by 11% to 20%


(Drewitz et al. 2017, Fantini-Pagani et al. 2013, Schmalz et al. 2011)

→ Agilium Freestep + 10mm shoe wedge: significant reduced KAM with both, medial and lateral shoe wedge (Schmalz et al. 2006)

→ Agilium Freestep + Insole: the KAM was reduced by 23% (Drewitz et al. 2017)

→ Agilium Freestep in different adjustments: significant ($p < 0.05$) decreased KAM in all conditions; by 11% in neutral position, by 8% in varus position, by 12% in valgus position (Fantini-Pagani et al. 2013)

Mean Maximum knee adduction moment during walking



++ = significant changes between the conditions with / without Agilium Freestep for healthy subjects and patients with knee OA ($p \leq 0.01$) [Schmalz et al. 2011]

Clinical Relevance

Osteoarthritis (OA) is a frequent cause of pain and disability in elderly people worldwide (Hinman et al. 2005) and its prevalence is expected to increase as the population ages (Ramsey et al. 2009).

The cause of OA is a non-physiological load distribution due to kinematic changes in the knee joint. A varus deformity of the leg axis leads to a long lever arm between the knee joint center and the vector of the ground reaction force, inducing a high external adduction moment. In the long term, these mechanisms can lead to damaged cartilage (Andriacci et al. 2004, Felson et al. 2009, Vincent et al. 2012).

Simply put, this means that the forces acting at the joint are high when the knee lever arm in the frontal plane is long due to varus malalignment. The length of the knee lever arm is correlated with a high knee adduction moment, and this adduction moment has been identified as the mechanism primarily responsible for a high compressive load in the knee (Ramsey et al. 2009).

The Agilium Freestep is used in unicompartmental OA. It is an AFO (ankle-foot-orthosis) that may reduce the eversion of the foot. Biomechanical studies showed that this shifts the vector of the ground reaction force closer to the center of the knee joint. As a result, the external adduction moment is reduced (Schmalz et al. 2006, Schmalz et al. 2011, Fantini-Pagani et al. 2013).

Summary

Four biomechanical studies have evaluated the effectiveness of the Agilium Freestep:

Schmalz et al. (2006) observed a significant reduction in the length of the knee lever arm to the vector of the ground reaction force (GRF) in the frontal plane and a significantly decreased knee adduction moment with Agilium Freestep in combination with shoe wedges.

Schmalz et al. (2011) investigated the Agilium Freestep without shoe wedges. Also here, a significantly reduced length of the knee lever arm to the GRF in the frontal plane and a significantly decreased knee adduction moment with Agilium Freestep could be shown.

Fantini-Pagani et al. (2013) conducted their study with subjects that showed a tendency towards a knee varus alignment. The former results could be confirmed. Significant reductions were seen in the knee lever arm, the knee adduction angular impulse and the knee adduction moment.

Drewitz et al. (2017) compared the effect of the Agilium Freestep to a combined use of Agilium Freestep with an insole as well as insoles without orthoses. The results showed that the orthosis reduced the knee adduction moment by 20% compared to no intervention. This effect could also be observed with the combination of orthosis and insole (reduction of 23%). An insole without orthosis increased the knee adduction moment by 1% compared to no intervention.

Thus, Agilium Freestep can alter the load distribution within the knee joint and thereby relieve the affected compartment. The use of this AFO, designed to offload the medial or lateral knee compartment, represents an alternative for conservative treatment of knee OA.

References of summarized studies

Drewitz, H., Schmalz, T., Wille, N. (2017). Konservative Behandlung der Varusgonarthrose mit einer Unterschenkelorthese. Conservative Treatment of Varus Osteoarthritis of the Knee with an Ankle-Foot-Orthosis. *OT: Orthopädie Technik*, 4: 36-41.

Fantini-Pagani, C. H., Willwacher, S., Benker, R., Brüggemann, G.-P. (2013). Effect of an ankle-foot orthosis on knee joint mechanics: A novel conservative treatment for knee osteoarthritis. *Prosthetics & Orthotics International*, *38*(6): 481-491.

Schmalz, T., Blumentritt, S., Drewitz, H. (2011). Die Nutzung von Unterschenkelorthesen im Rahmen der konservativen Behandlung der Gonarthrose. The application of orthoses for the lower leg in conservative treatment of gonarthrosis. *MOT: Medizinisch Orthopädische Technik*, *5*: 68-78.

Schmalz, T., Blumentritt, S., Drewitz, H., Freslier, M. (2006). The influence of sole wedges on frontal plane knee kinetics, in isolation and in combination with representative rigid and semi-rigid ankle-foot-orthoses. *Clinical Biomechanics (Bristol, Avon)*, *21*(6): 631-639. DOI: 10.1016/j.clinbiomech.2006.02.004.

Other References

Andriacchi, T. P., Mündermann, A., Lane Smith, R., Alexander, E. J., Dyrby, C. O., & Koo, S. (2004). A Framework for the in Vivo Pathomechanics of Osteoarthritis at the Knee. *Annals of Biomedical Engineering*, *32*(3): 447–457.

Felson, D. T. (2009). Developments in the clinical understanding of osteoarthritis. *Arthritis Research & Therapy*, *11*(1): 203-2013. DOI: 10.1186/ar2531.

Hinman, R. S., Crossley, K. M., McConnell, J., & Bennell, K. L. (2003). Efficacy of knee tape in the management of osteoarthritis of the knee: blinded randomized controlled trial. *BMJ: British Medical Journal*, *327*(7407): 135-140. DOI: 10.1136/bmj.327.7407.135.

O'Reilly, S. C., Jones, A., Muir, K. R., & Doherty, M. (1998). Quadriceps weakness in knee osteoarthritis: the effect on pain and disability. *Annals of the Rheumatic Diseases*, *57*(10): 588–594.

Ramsey, D. K., & Russell, M. E. (2009). Unloader braces for medial compartment knee osteoarthritis: implications on mediating progression. *Sports Health*, *8*: 416-426.

Vincent, K. R., Conrad, B. P., Fregly, B. J., & Vincent, H. K. (2012). The pathophysiology of osteoarthritis: a mechanical perspective on the knee joint. *PM&R: Physical Medicine and Rehabilitation*, *4*(5): S3-9. DOI: 10.1016/j.pmrj.2012.01.020.

© 2017, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.