
Reference

Safran Möller¹, David Rusaw¹, Kerstin Hagberg² and Nerrolyn Ramstrand¹

¹ School of Health and Welfare, Jonkoping University, Jonkoping, Sweden.

² Advanced Reconstruction of Extremities, Sahlgrenska University Hospital, Sweden and Department of Orthopaedics, Institute of Clinical Sciences; Sahlgrenska Academy, University of Gothenburg, Sweden.

Reduced cortical brain activity with the use of microprocessor-controlled prosthetic knees during walking

Prosthetics and Orthotics International (2019), Vol. 43(3): 257-265.

Products

Microprocessor-controlled prosthetic knee joints (MPK)^a vs. Non-microprocessor-controlled prosthetics knee joints (non-MPK)^b vs. Healthy controls

^a C-Leg (Ottobock), Genium (Ottobock), Rheo Knee (Ossur), Very Good Knee (Orthomobility)

^b 3R-80 (Ottobock), Mauch (Ossur), Total Knee (Ossur), Ultimate Knee (Ortho Europe)

Major Findings

→ **Less cortical brain activity while walking with an MPK compared to a non-MPK**

→ **Higher cortical brain activity while walking with a non-MPK compared to healthy controls**

→ **No difference in cortical brain activity while walking with an MPK compared to healthy controls**

→ **Prostheses users walked slower, used more steps on a 10 meter walking test than healthy individuals**

→ **Prostheses users covered less distances on a 6-minutes walk test than healthy individual**

Population

Subjects:	29 subjects (15 MPK-users, 14 non-MPK-users) 16 healthy controls
Previous prosthesis:	not reported
Amputation causes:	11 Tumour, 13 Trauma, 4 Vascular, 1 Infection
Mean age:	MPK-users: 50.1 ± 10.8 years Non-MPK-users: 50.8 ± 14.9 years Controls: 47 years (no individual data reported)
Mean time since amputation:	MPK-users: 16.1 ± 14.7 years Non-MPK-users: 19.1 ± 12.9 years
Prosthetic use score:	MPK-users: 84.9 ± 17.2 points Non-MPK-users: 69.4 ± 35.1 points




Figure 1, Experimental protocol

Participants attended a single session and were able to use their normal walking aid for the session. The experimental protocol was constructed as follows:

Prior to testing, participants were fitted with the fNIRS system (functional near-infrared spectroscopy; a wireless, portable, continuous wave system to record haemodynamic responses). fNIRS systems record the relative changes in concentration of oxygenated (oxyHb) and deoxygenated (de-oxyHb) haemoglobin, which results from neural activation. An increase in oxyHb concentration and a simultaneous decrease in de-oxyHb concentration would typically be observed as a response to regional brain activation.

Prior to the first testing, baseline measurements were collected for 30 seconds with the participants sitting quietly, with closed eyes. Following the baseline measures, participants then stood for another 30 seconds to allow to stabilization of fNIRS signals before given the signal to commence testing. Each walking trial was repeated for 4 times. A quiet rest of 30 seconds was given between each round.

Participants walked back and forth at a self-selected walking speed over a distance of 14 meters. Measurements of temporospatial data was recorded on the first 10 meters. After the fNIR probes were removed, they also fulfilled a 6 minutes walk test (on a 20 meter track rather than a 30 meter track due to space constraints).

Results

Functions and Activities					Participation			Environment	
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health, Economics
Category	Outcomes		Results for MPKs vs. Non-MPKs vs. Control				Sig.*		
Cognitive Demand	Haemodynamic response		Significantly increased oxyHb and decreased de-oxyHb in both prefrontal and motor cortices in individuals walking with a non-MPK when compared to controls (p<0.0167).						
			Significantly lower oxyHb concentrations in the prefrontal and motor cortex in individuals walking with an MPK compared to those walking with a non-MPK (p<0.0167).				--		

Category	Outcomes	Results for MPKs vs. Non-MPKs vs. Control	Sig.*
		No differences could be detected in de-oxyHb when comparing individuals walking with an MPK and non-MPK.	++
		No significant differences were seen in oxyHb between MPK users and controls.	0
		A significant decrease on de-oxyHb concentrations in the motor cortex for controls when compared to the MPK group could be detected ($p < 0.0167$).	0
		Significant differences could be shown for the time taken on the 10 meters walking distance when comparing controls to prosthesis users (non-MPK vs. controls $p < 0.001$; MPK vs. control $p < 0.001$).	--
Level Walking	Temporospatial	Significant differences could be shown for the number of steps taken on the 10 meters walking distance when comparing controls to prosthesis users (non-MPK vs. controls $p < 0.001$; MPK vs. control $p < 0.025$).	--
		Significant differences could be shown for the covered walking distance in the 6-minutes walking test when comparing controls to prosthesis users (non-MPK vs. controls $p < 0.001$; MPK vs. control $p < 0.002$).	--
		No differences regarding walking time on the 10-meters, steps taken on the 10 meters and covered walking distance on the 6 minutes walk test could be detected when comparing the MPK group and the non-MPK group.	--
		No differences regarding walking time on the 10-meters, steps taken on the 10 meters and covered walking distance on the 6 minutes walk test could be detected when comparing the MPK group and the non-MPK group.	0

Author's Conclusion

"Persons walking with a non-MPK have an increase in cortical brain activity during ambulation when compared to healthy controls. Furthermore, individuals using an MPK joint require less prefrontal cortical brain activity to ambulate when compared to those using a non-MPK joint. Findings suggest that prescription of microprocessor-controlled prosthetic knees in relatively young and healthy users of TF and KD prostheses may help to reduce reliance on executive functions during gait." (Möller et al, 2019)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.