

## Reference

Engel K, Brüggemann G-P, Heinrich K, Potthast W, Liebau C.

Institute of Biomechanics and Orthopedics, German Sport University Cologne, Am Sportpark Müngersdorf 6, 50933 Cologne, Germany.

# Do counteracting external frontal plane moments alter the intraarticular contact force distribution in the loaded human tibiofemoral joint?

Munich: Elsevier 2015, The Knee 22: 68–72


## Products

### Genu Arthro

## Major Findings

With Genu Arthro:

- **A higher knee flexion, leads to an increase of the medial compartment loading (varus and valgus alignment)**
- **Counteracting abduction moments have a higher positive influence on the force in the medial knee joint compartment than adduction moments for the lateral compartment.**
- **The same level of counteracting moment in valgus alignment results in a much higher shift of the COP than in varus alignment.**




The applied counteracting abduction moments led to higher redistributions than the applied adduction moments.

## Population

Ex vivo: 6 cadaver legs (4 female, 2 male) (59-101 yrs)  
SP1-3: varus joint  
SP4-5: neutral  
SP6: valgus joint

## Study Design

Interventional, comparative:



One cylinder on top of the simulator was used for the up and down movements of the sliding bar. The steel cylinder with the femur was attached to a ball joint fixed to this sliding bar. A predefined range of motion from 10° to 40° flexion of the knee joint was measured on the lateral side of the legs using a goniometer. A reference trial without any applied frontal plane moments was captured of all joints. With Genu Arthro counteracting abduction moments of 2.5 and 5 Nm were applied to the varus joints, adduction moments to the valgus joint, respectively. The neutral aligned joints were tested under all conditions. Each condition was repeated five times. The negative sign indicates applied abduction moments.

## Results

Functions and Activities						Participation
Biomechanics – Static measures	Biomechanics – Gait analysis	X-Ray	EMG	Functional tests	Clinical effects	Satisfaction

Category	Outcomes	Results for Genu Arthro	Sig.*
Biomechanics – Static measures	Knee flexion angle	increasing knee flexion angles, the forces in the medial compartment increase, even for the knee joint aligned in varus	n.a.
		<b>Medio-lateral force ratio (&gt;1→med. knee joint comp.)</b>	
		<b>Flexion</b>	
		<i>SP1</i>	<i>SP2</i>
		<i>SP3</i>	<i>SP4</i>
		<i>SP5</i>	<i>SP6</i>
		10°	1.8 5.1 1.5 0.8 1.3 0.2
		20°	2.1 5.4 1.7 1.8 1.5 0.3
		30°	2.8 5.9 2.2 2.7 2.2 0.5
		40°	3.2 6.3 2.5 3.9 3.0 0.6
		The joints aligned in varus (SP1–3) and neutral (SP4-5) show higher loading in the medial compartment, whereas more load is transferred through the lateral compartment in the knee joint aligned in valgus (SP6).	
	Medio-lateral force-ratio	The applied counteracting abduction moments lead to higher redistributions than the applied adduction moments.	n.a.
	Shift of the COP (centre of pressure)	The largest COP shift effect of 4.3 mm to the medial side is found in the joint in valgus alignment when a 5 Nm adduction moment is applied. By comparison, the highest COP shift for the lateral side is 1.8 mm.	n.a.

\* no difference (0), positive trend (+), negative trend (–), significant (++)/(--), not applicable (n.a.)

---

## Author's Conclusion

“The current study demonstrated that direct intraarticular measurements could be used to study the relationship between knee joint compartment loading and mechanically induced external moments. The presented findings have clinical relevance and offer novel insights to guide development and optimization of mechanical aids for the treatment of the osteoarthritic knees. Future of load redistributing knee orthotics should consider adaptive approaches to apply frontal plane moments dependent on the movement phase since the intraarticular load is related to the knee joint flexion angle.” (Engel et al. 2015)

© 2014, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.