

Reference

Bouwsema H, van der Sluis C, Bongers R

Center of Human Movement Sciences, University of Groningen, Groningen

Learning to Control Opening and Closing a Myoelectric Hand

American Congress of Rehabilitation Medicine 2010; 91:1442-6


Products

Virtual hand – PAULA; Myoelectric simulator; Table-top hand (acts like Sensor Hand Speed)

Major Findings

- **Prosthetic users differ in learning capacity which determines time needed to learn how to use myoelectric prosthesis.**
- **Acquired control of a myoelectric hand is irrespective of the type of device used for training (PAULA/ simulator/ table-top hand)**
- **PAULA software is as effective as tabletop hand and prosthetic simulator.**

Peak velocity after the training period




Graph shows peak velocities of opening and closing the hand reached in the post-test (after the training period) for the high capacity learners (HCL) and low capacity learners (LCL) plotted for each of the velocity conditions – slow, comfortable and fast. High-capacity learners could make a good distinction between the 3 different velocity conditions, whereas low-capacity learners could not make this distinction.

Population

Subjects: 34 able-bodied participants
Previous: none
Amputation causes: none
Mean age: 21 years
Mean time since amputation: none

Study Design

A randomized study:



After entering into the study, the subjects were randomized into three groups based on type of the training they will receive. On the first day a pretest was conducted. Afterwards, the subject's control of the hand was trained on 3 consecutive days either by using virtual hand, tabletop hand or prosthetic simulator. After the last training session on the 3rd day, a posttest was administered to determine the level of

skill after the training. The pretest and the posttest test were the same and consisted of 2 parts: the participant was asked to first provide a maximum myoelectric signal for at least 2 seconds (this was repeated 5 times) and, second, to open and close the hand to the maximal aperture on 3 different velocities at command. Participants were asked to control hand opening and closing at the slowest speed possible, at a comfortable speed, and at the highest speed possible. All velocities were executed 3 times in a random order. When the hand was not fully opened or closed, the participants were corrected and instructed again.

Results

Body Function		Activity			Participation	Others	
Mechanics	Pain	Grip patterns / force	Manual dexterity	Activities of daily living (ADL)	Satisfaction and Quality of life (QoL)	Training	Technical aspect

Category	Outcomes	Results for training with PAULA vs simulator vs table-top hand:	Sig.*
Training	Peak and mean velocity	Both peak velocity and mean velocity showed the same main effects.	0
	Number of peaks	A large effect of the velocity conditions showed that in the slow condition the most peaks occurred, whereas in the fast condition the fewest number of peaks were shown.	0

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

"In conclusion, learned control of a myoelectric hand does not depend on the type of training (with a virtual hand, an isolated hand, or a prosthetic simulator). Prosthetic users may differ in learning capacity, and this should be taken into account when choosing the appropriate type of control for each patient." (Bouwsema et al. 2010)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.