

Reference

Andreas Kannenberg¹, MD (GER), PhD, Russell Lundstrom¹, MS, Karl D. Hibler², MA, Shawn Swanson Johnson³, OTR/L

¹Department of Clinical Research & Services, Otto Bock Healthcare LP, Austin, Texas

²Bradenton, Florida

³SSJ Rehab Services LLP, Houston, Texas

Differences in Two Multiarticulating Myoelectric Hands for Facilitating Activities of Daily Living in Individuals with Transradial Amputation: A Cross-Sectional Study

Journal of Prosthetics and Orthotics (2022); 00:00-00


Products

Bebionic vs. i-Limb

Major Findings

With bebionic compared to i-Limb:

- **No significant differences in ADL ease and usefulness between bebionic, i-Limb and historical data for Michelangelo**
- **Higher ease and usefulness scores than previously reported for conventional myoelectric hands**



Population

Subjects:

bebionic-group:

10 transradial amputees (n = 5 male, n = 5 female)

i-Limb group:

10 transradial amputees (n = 9 male, n = 1 female)

Previous prosthesis:

bebionic group:

i-Limb (n = 2), Greifer (n = 1), ETD-powered hook (n = 2), body-powered (n = 3), passive hand (n = 1), none (n = 1)

i-Limb group:

Amputation causes: i-Limb (n = 3), Sensor Hand (n = 1), Greifer (n = 1), ETD-powered hook (n = 1), Body-powered (n = 2), none (n = 1), unknown (n = 1)

bebionic group:
Congenital deformity (n = 3), Trauma (n = 6), other (n = 1)

i-Limb group:
Congenital deformity (n = 3), Trauma (n = 4), Cancer (n = 2), Infection/Sepsis (n = 1)

Mean age (± SD) [years]: bebionic group:
37.4 ± 14.2

i-Limb group:
50.4 ± 17.6


Mean time since Amputation (± SD) [years]: bebionic group:
16.1 ± 19.6

i-Limb group:
16.1 ± 19.6

MFCL: n.a.

Study Design

Observational study design:



Bebionic and i-Limb users were asked about demographics e.g.: age, sex, years of prosthetic use, amputation side and etiology of amputation. Following this, patients were asked to answer a hybrid outcome measure that combined the modified Orthotics and Prosthetics User Survey–Upper Extremity Functional Status (OPUS-UEFS) and the Prosthetic Upper Extremity Functional Index (PUFI).

Results

Body Function		Activity			Participation	Others	
Mechanics	Pain	Grip patterns / force	Manual dexterity	Activities of daily living (ADL)	Satisfaction and Quality of life (QoL)	Training	Technical aspect

Category	Outcomes	Results						Sig.*
		Bebionic			i-Limb			
	OPUS-UEFS ease score	M ¹	IQR ¹	Mean ± SD ¹	M ¹	IQR ¹	Mean ± SD ¹	
Activity, Mobility, Activities of Daily Living (ADLs)	<i>All activities (23):</i>	36	19.5-43.2	32.5 ± 13.5	30.5	23.75-44.5	34.9 ± 14.9	0
	<i>Bimanual activities (14):</i>	23	12.25-29.0	21.1 ± 9.5	21.5	18-28	23.6 ± 8.2	0
	<i>Monomanual activities (9):</i>	12	7-17.25	11.4 ± 6.0	10	5-16.5	11.3 ± 7.4	0

Category	Outcomes	Results		Sig.*
		Bebionic	i-Limb	
Activity, Mobility, Activities of Daily Living (ADLs)	Mean number of ADLs per usefulness of the prosthesis rating category per prosthetic hand	Mean ± SD	Mean ± SD	
	<i>Not useful:</i>	9.7 ± 4.0	8.9 ± 5.0	0
	<i>Useful:</i>	3.8 ± 3.1	6.1 ± 2.6	0
	<i>Very useful:</i>	9.2 ± 3.7	7.2 ± 4.4	0
	Mean number of ADLs per way-of-prosthesis-use rating category per prosthetic hand	Mean ± SD	Mean ± SD	
	<i>Both hands, prosthesis used actively to grasp</i>	10.7 ± 2.9	9.8 ± 3.0	0
	<i>Bot hands, prosthesis used passively to stabilize</i>	2.1 ± 2.5	2.2 ± 1.8	0
<i>With assistance of residual limb</i>	1.8 ± 2.3	1.3 ± 1.4	0	
<i>Second hand alone</i>	6.8 ± 2.7	8.2 ± 2.4	0	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

¹ M: Median, IQR: Interquartile-range, Mean: mean value, SD: Standard deviation

Author's Conclusion

“The differences in overall ADL [activity of daily living] ease and usefulness of the prosthesis between the i-Limb and bebionic hands were clinically negligible. Ease and usefulness scores were higher than previously reported for conventional myoelectric hands. Interestingly, the availability of more grip types in bebionic and i-Limb did not result in greater ease or usefulness than previously reported for the Michelangelo hand with fewer grip types. However, the multiarticulating hands showed different activity ease profiles that they facilitate. Thus, clinicians should have access to all advanced prosthetic hands to be able to match their patients’ functional needs with the differential functional ease profiles of these hands. Future research that compares all available multiarticulating hands using performance-based and patient-reported outcomes is warranted to further guide clinicians’ and payers’ decision making.” (Kannenberget al. 2022)

© 2022, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended

to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.