

## Reference

Rohlmann, A., Zander, T., Graichen, F., Bergmann, G.

# Effect of an orthosis on the loads acting on vertebral body replacement

Clinical Biomechanics 28, 2013: 490-494.

<http://dx.doi.org/10.1016/j.clinbiomech.2013.03.010>

## Products

**Lumbo TriStep (LTS), Hyperextension orthosis (HEO, medi 3C)**


## Major Findings

With Lumbo TriStep:

→ **The average resultant force on the vertebral body replacement (VBR) for 26 activities was reduced**

- by 9% with Lumbo TriStep (LTS)
- by 19% with hyperextension orthosis (HEO)
- The force reduction is usually more pronounced for activities performed during sitting

## Load changes due to orthosis use



Changes of max. resultant force on vertebral body replacement VBR. The values are related to the situation without an orthosis which was regarded as 100%. The results of the 5 patients are compared.

## Population

Subjects: 5 patients with a severe fracture of L1 or L3 vertebral body (4 male, 1 female)

Age: 62 to 71 years

Measurement: Telemeterized vertebral body replacement (VBR) were implanted. The implant allows the measurement of 6 load components acting on it.

Implantation date: 09/2006 - 07/2008


Time between implantation and measurement: 150 to 774 days

Intervention: For several activities during standing, sitting and walking, implant loads were measured with and without an orthosis.

**Study Design**

Descriptive, in-vivo:

**Measurements were performed in one session**



Orthoses:


LTS: Lumbo TriStep: stabilizing orthosis with mobilizing function

HEO: Hyperextension brace: immobilizing orthosis (thoracic and/or lumbar)

**Results**

Functions and Activities						Participation
Biomechanics – Static measures	Biomechanics – Gait analysis	X-Ray	EMG	Functional tests	Clinical effects	Satisfaction


Category	Outcomes	Results for Lumbo TriStep and hyperextension orthosis
Functional tests	Max. resultant force on VBR	The average resultant force on the vertebral body replacement (VBR) for all 26 activities measured was reduced <ul style="list-style-type: none"> <li>- by 9% with Lumbo TriStep (LTS)</li> <li>- by 19% with Hyperextension brace (HEO)</li> </ul>
Functional tests	Max. resultant force on VBR while walking	Changes of maximum resultant force on a vertebral body replacement (VBR) due to an orthosis during walking. A negative value indicates an unloading of the VBR. Considerable inter- and intra-individual variations were observed.



Changes of maximum resultant force on the vertebral body replacement due to an orthosis for walking. The values are relative to the situation without an orthosis which was regarded as 100%. The results for 5 patients (WP1-WP5) are compared.


Biomechanics – Static measures	Biomechanics – Gait analysis	X-Ray	EMG	Functional tests	Clinical effects	Satisfaction
--------------------------------	------------------------------	-------	-----	------------------	------------------	--------------

Category	Outcomes	Results for Lumbo TriStep and hyperextension orthosis
Functional tests	Max. resultant force on VBR while standing	<p>Changes of maximum resultant force on a vertebral body replacement (VBR) due to an orthosis while standing (15 exercises).</p> <p>A negative value indicates an unloading of the VBR.</p> <p>Considerable inter- and intra-individual variations were observed.</p>



Changes of max. resultant force on VBR due to an orthosis for 15 different activities while standing. The median values and the ranges are shown. For LTS n=5, for HEO n=4.

Functional tests	Max. resultant force on VBR while sitting	<p>Changes of maximum resultant force on a vertebral body replacement (VBR) due to an orthosis while sitting (10 exercises).</p> <p>A negative value indicates an unloading of the VBR.</p> <p>Considerable inter- and intra-individual variations were observed.</p>
------------------	---	---



Changes of max. resultant force on VBR due to an orthosis for 10 different activities while sitting. The median values and the ranges are shown. For LTS n=5, for HEO n=4.

---

## Author's Conclusion

“The forces on a VBR and thus on the anterior column of the spine are on average slightly reduced when wearing a LTS brace and more pronounced due to a hyperextension orthosis. However, large inter- and intra-individual variations exist. Therefore, from the biomechanical point of view, no clear recommendation to wear an orthosis can yet be given since the clinically relevant reduction of the implant force is unknown.” (Rohlmann et al. 2013)

© 2014, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.