
Reference

Kutzner I, Kuther S, Heinlein B, Dymke J, Bender A, Halder AM, Bergmann G.

Julius Wolff Institute, Charité – Universitaetsmedizin Berlin, Augustenburger Platz 1, 13353 Berlin, Germany.

The effect of valgus braces on medial compartment load of the knee joint– in vivo load measurements in three subjects

J Biomech 2011. 44: 1354–1360.

Products

Genu Arthro (Otto Bock)

MOS Genu (Bauerfeind)

Major Findings

With Genu Arthro (GA) compared to MOS Genu (MOS) and wearing no orthosis:

→ **Stiffness (measured by spring constant) of MOS is 145% higher than of GA**

→ **Wearing GA or MOS has no significant influence on walking speed and stride length during level walking.**

→ **Overall: forces (Fmed & Fz) are reduced by both orthoses in 88% of 52 measurements. MOS leads for almost all conditions to lower forces than GA.**

→ **Level walking:**

GA (8° valgus): Fmed: 7% lower
Fz: 2-3% lower

MOS (0°, 4°, 8° valgus): Fmed: 9-30% lower
Fz 2-9% lower

→ **Stair ascent:**

GA (0°, 8° valgus): Fmed: 2-9% lower

MOS (0°, 8° valgus): Fmed: -2% higher - 26% lower


→ **Stair descent:**

GA (0°, 8° valgus): Fmed: 5-7% lower
Fz: 3-7% lower

MOS (0°, 8° valgus): Fmed: 2-24% lower
Fz 6-16% lower

Caution: The author herself calls this study a case report with a limited number of 3 subjects, which does not allow general conclusions. Differences in 21% of 52 force measurements were significant. Moreover: the methods do not explain why an additional valgus angle of 4° was used only for MOS during level walking.

Reduction of the medial force (Fmed) during stair ascent (orthosis compared to wearing no orthosis)




Population

Subjects: 3 male
 Mean age: 65 ± 4.6 years
 Mean body mass: 98.3 ± 3.3 kg
 Inclusion criteria: medial osteoarthritis
 Mechanical axis angle: 1.7 ± 1.25° varus

Study Design

Interventional, comparative:



The braces were first fitted to the leg in a neutral position to examine whether the brace itself already has an influence on joint loading. After performing the activities with the brace in neutral position, additional valgus angles of 4° (for MOS and level walking only) and 8° (MOS and GA) were adjusted. Three activities of daily living were investigated: walking at a self-selected speed on level ground, ascending stairs, and descending stairs. An implanted, instrumented tibial tray was developed to measure the 6 components (3 forces and 3 moments) of the knee contact loading in vivo. Additionally, the stiffness of the braces in the frontal plane was determined by a testing machine with a maximum test load of 100N.

Results

Functions and Activities						Participation
Biomechanics – Static measures	Biomechanics – Gait analysis	X-Ray	EMG	Functional tests	Clinical effects	Satisfaction

Category	Outcomes	Results for Genu Arthro (GA) & MOS Genu (MOS)	Sig.*
Biomechanics – Static measures	Stiffness	The spring constant of the MOS is 145% higher than that of GA .	n.a.

Category	Outcomes	Results for Genu Arthro (GA) & MOS Genu (MOS)					Sig.*
Biomechanics – Gait analysis	Fmed 1 st peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (4°) vs without	MOS (8°) vs without	
		1% higher	7% lower	10% lower	18% lower	23% lower	
		-	+	++	++	++	
	Fmed 2 nd peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (4°) vs without	MOS (8°) vs without	
		4% higher	7% lower	9% lower	24% lower	30% lower	
		-	++	+	++	++	
	Fz	With GA (0°) Fz increases by 2-3%.					-
		During usage of GA (8°) , Fz decreases by 2-3%					+
		Fz is reduced by 2-9% with the MOS (0°, 4°, 8°) .					+
	Walking speed	No significant differences.					0
Stride length	No significant differences between walking with or without orthosis.					0	
Functional tests	Stair ascending – Fmed 1 st peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (8°) vs without		
		2% lower	6% lower	2% higher	11% lower		
		+	+	-	+		
	Stair ascending – Fmed 2 nd peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (8°) vs without		
		6% lower	9% lower	5% lower	26% lower		
		+	+	+	++		
	Stair ascending - Fz	No significant differences for all conditions.					0
	Stair descending – Fmed 1 st peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (8°) vs without		
		5% lower	7% lower	12% lower	24% lower		
		+	+	+	+		
Stair descending – Fmed 2 nd peak	GA (0°) vs without	GA (8°) vs without	MOS (0°) vs without	MOS (8°) vs without			
	7% lower	6% lower	2% lower	17% lower			
	+	+	+	++			
Stair descending - Fz	With GA (0°, 8°) Fz is reduced by 3-7%					+	
	Fz is decreased by 6-16% with MOS (0°, 8°) 2nd peak of MOS 8° (-16%) is significant.					+	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

“Since the patient's acceptance of the brace and the wearing comfort is of major importance, the amount of applicable external valgus moment is limited. Whereas no major discomfort was reported by the subjects when wearing the GA brace, discomfort was reported when walking with the MOS brace in 8° valgus. Since the chosen valgus settings of 8° with the MOS brace would probably not have been tolerated for a long duration by the subjects, medial load reductions of more than 25% cannot be expected permanently. ...

... Due to the variability between the subjects, the authors suggest that valgus braces should only be used if a patient reports pain relief. The unloading effect of braces must furthermore be compared to other more comfortable conservative methods such as forearm crutches, laterally wedged shoes or weight reduction." (Kutzner et al. 2011)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.