

---

## Reference

Cheesborough J, Smith L, Kuiken T, Dumanian G

Neural Division of Plastic and Reconstructive Surgery, Northwestern University, Chicago, Illinois

# Targeted Muscle Reinnervation and Advanced Prosthetic Arms

Semin Plast Surg 2015;29:62–72.

---

## Products

**Targeted Muscle Reinnervation and Coapt Pattern Recognition control myoelectric prosthesis**


---

## Major Findings

The effect of Targeted Muscle Reinnervation (TMR) and Coapt pattern recognition control on myoelectric prosthetic use:

→ **Similar functional performance of the shoulder disarticulation and the transhumeral side of a bilateral amputee, despite poorer expectation for the higher-level amputation side**

### Clothespin Relocation Test



During functional tests such as the clothespin relocation test the patient demonstrated similar performance between his left - shoulder disarticulation ( $60.6 \pm 11.5s$ ) and right - transhumeral sides ( $59.7 \pm 10.6s$ ), despite the difference in amputation level. The results suggest the intuitiveness of control with TMR and pattern recognition control, as higher-level amputation would otherwise be expected to provide poorer performance in functional tasks.


---

## Population

Subjects:	one male bilateral amputee (left shoulder disarticulation and right transhumeral amputation)
Amputation etiology:	trauma
Age at amputation:	43 years
Age at TMR	45 years
Previous prosthesis:	hybrid prosthesis on his transhumeral side, which included a passive (non-moving) elbow and a myoelectric hook; myoelectric prosthesis after TMR surgery on his shoulder disarticulation side
Intervention prosthesis:	bilateral fitting with Coapt pattern-recognition myoelectric prostheses

## Study Design

### Case report:



A 43 year old male sustained a severe electrical burn injury and required a left side amputation at the shoulder disarticulation level and a right side amputation at the transhumeral level. Two years after the injury, the TMR surgery was first performed on the patient's left side and four months later on his right side. The patient was initially fitted with a hybrid prosthesis on his transhumeral side and with a myoelectric prosthesis on his shoulder disarticulation side. After one year of prosthetic use, the patient was bilaterally fitted with Coapt pattern-recognition myoelectric prostheses as intervention devices.

## Results

Body Function		Activity			Participation	Others	
Mechanics	Pain	Grip patterns / force	Manual dexterity	Activities of daily living (ADL)	Satisfaction and Quality of life (QoL)	Training	Technical aspect

Category	Outcomes	Results for pattern recognition myoelectric prosthetic use after TMR:	Sig.*
Pain	Self-reported	Five months after the TMR surgery, the patient reported complete resolution of his neuroma pain bilaterally. Some occasional phantom limb pain was reported.	n.a.
Manual dexterity	Box and Blocks	Similar performance of the left ( $11.0 \pm 1.5$ blocks) and right ( $14.3 \pm 0.3$ blocks) side with pattern recognition prosthesis, despite the difference in amputation level.	n.a.
	Clothespin Relocation Task	Similar performance of the left ( $60.6 \pm 11.5$ s) and right ( $59.7 \pm 10.6$ s) side with pattern recognition prosthesis, despite the difference in amputation level.	n.a.
Activities of daily living (ADL)	Self-reported	Many tasks were easier to perform with the pattern recognition controlled myoelectric prosthesis after TMR: eating with a fork, drinking from a water bottle, carrying a laundry basket, yard work, as well as placing, retrieving and replacing items from a refrigerator.	n.a.

\* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

## Author's Conclusion

"Targeted muscle reinnervation combined with existing and emerging prosthetic technology allows for intuitive control of myoelectric prostheses for amputees at multiple levels. For complex amputees, such as the patient presented in the case example, a strategic and orderly approach to care is essential, understanding that each patient will present unique challenges." (Cheesborough et al., 2015)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary