

Reference

Fantini Pagani CH, Hinrichs M, Brüggemann G-P.

Institute of Biomechanics and Orthopaedics, German Sport University Cologne, Germany.

Kinetic and Kinematic Changes with the Use of Valgus Knee Brace and Lateral Wedge Insoles in Patients with Medial Knee Osteoarthritis

Orthop Res 2012; 30:1125-1132.

Products

Genu Arthro vs Lateral Wedge Insoles

Major Findings

With Genu Arthro compared to lateral wedge insoles and wearing no orthosis:

→ **External knee adduction moments with Genu Arthro at the second peak of the gait cycle were**

18-21% lower than without orthosis

11-15% lower than with lateral wedge insoles


→ **Knee lever arm with Genu Arthro at 20-30% and 70-80% of the stance phase was**

9-12% shorter than without orthosis (20-30% of the stance phase)

19-21% shorter than without orthosis (70-80% of the stance phase)

The effect was twice to three times bigger than with the lateral wedge insoles.

External knee adduction moments at second peak of the gait cycle




Population

Subjects:	10 patients (8 female, 2 male)
Mean age:	57.5 ± 7.1 yrs
Mean body mass:	78.8 ± 12.2 kg
Inclusion criteria:	age over 50 yrs, diagnosis of medial osteoarthritis grade II or III

Study Design

Observational, comparative:



The patients were not informed about the different adjustments of the orthosis.

The moments calculated using the kinematic and GRF data collected during these trials were defined as external knee adduction moments. For the conditions with orthosis (neutral, 4° and 8°), net moments were calculated by subtracting the orthosis moments from the external knee adduction moments.

Results

Functions and Activities						Participation
Biomechanics – Static measures	Biomechanics – Gait analysis	X-Ray	EMG	Functional tests	Clinical effects	Satisfaction

Category	Outcomes	Results for Genu Arthro	Sig.*															
Biomechanics – Gait analysis	External knee adduction moments	The external peak knee adduction moment (first peak) did not differ significantly.																
		The external peak knee adduction moment (second peak) was significantly lower with Genu Arthro with both adjustments:																
		<table border="1"> <thead> <tr> <th>4° valgus vs no*</th> <th>8° valgus vs no</th> <th>Insoles vs no</th> <th>4° valgus vs insoles</th> <th>8° valgus vs insoles</th> </tr> </thead> <tbody> <tr> <td>18% lower</td> <td>21% lower</td> <td>7% lower</td> <td>11% lower</td> <td>15% lower</td> </tr> <tr> <td>++</td> <td>++</td> <td>++</td> <td>++</td> <td>++</td> </tr> </tbody> </table>	4° valgus vs no*	8° valgus vs no	Insoles vs no	4° valgus vs insoles	8° valgus vs insoles	18% lower	21% lower	7% lower	11% lower	15% lower	++	++	++	++	++	
		4° valgus vs no*	8° valgus vs no	Insoles vs no	4° valgus vs insoles	8° valgus vs insoles												
		18% lower	21% lower	7% lower	11% lower	15% lower												
++	++	++	++	++														
The knee adduction angular impulse was also significantly lower under almost all conditions:																		
<table border="1"> <thead> <tr> <th>4° valgus vs no</th> <th>8° valgus vs no</th> <th>Insoles vs no</th> <th>4° valgus vs insoles</th> <th>8° valgus vs insoles</th> </tr> </thead> <tbody> <tr> <td>14% lower</td> <td>18% lower</td> <td>7% lower</td> <td>8% lower</td> <td>12% lower</td> </tr> <tr> <td>++</td> <td>++</td> <td>++</td> <td>+</td> <td>++</td> </tr> </tbody> </table>	4° valgus vs no	8° valgus vs no	Insoles vs no	4° valgus vs insoles	8° valgus vs insoles	14% lower	18% lower	7% lower	8% lower	12% lower	++	++	++	+	++			
4° valgus vs no	8° valgus vs no	Insoles vs no	4° valgus vs insoles	8° valgus vs insoles														
14% lower	18% lower	7% lower	8% lower	12% lower														
++	++	++	+	++														
Orthosis moment	8° vs 4° valgus: 32% higher (first peak), 30% higher (second peak)	++																
Net knee adduction moment	8° vs 4° valgus: 14.7% lower (first peak), 14.8% lower (second peak)	++																
Net knee adduction angular impulse	8° vs 4° valgus: 22.2% lower	++																

Category	Outcomes	Results for Genu Arthro			Sig.*
	Knee lever arm	The knee lever arm at 20-30% of the stance phase was lower with both adjustments:			
		4° valgus vs no	8° valgus vs no	Wedges vs no	
		9% lower	12.6% lower	5.4% lower	
		+	++	+	
		The knee lever arm at 70-80% of the stance phase was significantly lower with both adjustments:			
		4° valgus vs no	8° valgus vs no	Wedges vs no	
		19.3% lower	21.4% lower	7.1% lower	
		++	++	++	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)
no = without orthosis

Author's Conclusion

"In conclusion, the knee orthosis tested in this study and 4° laterally wedged insoles were effective in reducing knee lever arm in the frontal plane, knee adduction moment, and possibly joint load. However, a small effect size was observed with the insoles. The knee orthosis was more effective than 4° wedged insoles in reducing the external knee adduction moment." (Fantini Pagani et al. 2012)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.