ottobock.

GENIUM



Clinical Study Summaries

This document summarizes clinical studies conducted with the Genium. The included studies were identified by a literature search made on PubMed and within the journals Der Orthopäde, JPO Journal of Prosthetics and Orthotics, Orthopädie-Technik and Technology & Innovation.

Table of content:

Overview table
Summaries of categories
Level walkingp 5
Stairsp7
Ramps, Hillsp9
Uneven ground, Obstacles
Safety p 13
Activity, Mobility, ADL
Preference, Satisfaction, QoL
Summaries of individual studies
Copyright

Genium Study Summaries

The summaries are organized in three levels depending on the detail of information. The overview table (Level 1) lists all the relevant publications dealing with a particular product (topic) as well as researched categories (e.g. level walking, safety, activities, etc). By clicking on individual categories, a summary of all the literature dealing with that category will open (Level 2).

For those interested to learn more about individual studies, a summary of the study can be obtained by clicking on the relevant reference (Level 3).

Reference			Category								
		Functions and Activities									Environment
Author	Year	<u>Level</u> walking	<u>Stairs</u>	<u>Ramps,</u> <u>Hills</u>	<u>Uneven</u> ground Obstacles	Cognitive demand	Metabolic Energy Consumption	<u>Safety</u>	Activity Mobility ADL	Preference Satisfaction QoL	Health economics
<u>Lura</u>	2017		х								
<u>Hahn</u>	2016	х	х	x	x	x	x	х	x		
<u>Highsmith</u>	2016								x		
<u>Highsmith</u>	2016		х					х	x	x	х
<u>Highsmith</u>	2016			x							
Huppert*	2016	x	х	x				х	x	x	
<u>Bell</u>	2016			x							
<u>Schalk**</u>	2015								x	x	
<u>Lura</u>	2015	x		x							
<u>Aldridge</u> Whitehead	2014		x								
<u>Highsmith</u>	2014									x	
Highsmith	2014	x			x			х			

<u>Highsmith</u>	2014		х	х	х						
<u>Kannenberg</u>	2013	x	х	x	х			x	x		
<u>Bellmann</u>	2012	х	х	х				x			
<u>Bellmann</u>	2012		x								
Blumentritt	2012	x	x	x				x			
<u>Kampas</u>	2011	x	x	x	x			x			
Total numbe	r: 18	8	11	10	5	1	1	8	6	4	1
*Review											
** Case study											

2 Summaries of categories

On the following pages you find summaries of categories researched in several studies (e.g. level walking, stairs, etc.). At the end of each summary you will find a list of reference studies contributing to the content of the particular summary.

Level Walking

Major Findings	With Genium compared to C-Leg:						
	 → Increased toe clearance with more consistent maximum knee angle at 64° across different walking velocities → Step length symmetry is improved by up to 60% → Increased toe clearance when walking with small steps Maximum knee flexion angle in swing phase is increased by 11% 						
	→ Reliable swing phase release even when walking with small steps Swing phase is reliably released in 95% of all small steps						
	Difficulty of walking backwards is improved by 26% Safe loading when stepping backwards was in all tests possible						
	Reduced asymmetry of step length with Genium compared to C-Leg						
	Slow walking speed Mid walking speed						

Bellmann	et al.	(2012)
----------	--------	--------

Clinical Relevance The main aim of a prosthesis is the restoration of function. For lower extremities the most important function is ambulation. It has influence on the mobility grade of the subject, the participation of life and, therefore, general quality of life. Furthermore, a natural gait pattern is pursued, since it prevents the sound side from higher or inappropriate loads due to compensatory movements. Gait asymmetries can contribute to secondary diseases such as osteopenia or arthritis. **Summary** Walking velocities were investigated by Bellmann et al. (2012) as well as by Highsmith et al. (2014) at varying speeds and varying distances. No differences between Genium and C-Leg could be observed. However, subjects had the tendency to rate the perceived exertion lower with Genium compared to C-Leg (Highsmith et al. 2014). Regarding biomechanical gait analysis, several improvements with Genium were documented. The maximum knee angle in swing phase is not influenced by gait speed and stays constant at 64° (Bellmann et al. 2012 & Blumentritt et al. 2012). Therefore, adequate toe clearance is ensured even when walking slowly. Over the range from very slow to normal walking velocities, peak knee flexion angle in swing phase is increase by up to 7° with Genium compared to C-Leg. In stance phase, peak knee flexion angle is increased by 2° during slow and normal speed with Genium compared to C-Leg. Knee angles are due to accommodation, training and use

	of Genium closer to the intact limb and therefore a more physiological gait pattern is achieved (Lura et al. 2015). Furthermore, it was observed that, with Genium, the prosthetic knee is in a pre-flexed position of 4° at heel strike. The pre-flex leads to reduced breaking forces which was proved by a decreased ground reaction force on the prosthetic side (Bellmann et al. 2012 & Blumentritt et al. 2012). In conclusion the improved gait characteristics of Genium resulted in a more normalized, anatom- ic movement pattern.
	Walking with small steps is improved with Genium; not only is the maximum knee angle in swing phase increased compared to C-Leg, but also the swing phase is reliably released with Genium in 95% of all steps. With C-Leg, swing phase is re- leased in 75% of all the steps (Bellmann et al. 2012 and Blumentritt et al. 2012).
	When walking backwards, safe loading of the prosthesis is possible with Genium (Blumentritt et al. 2012). Moreover, subjects rated walking backwards as less difficult to perform with Genium than with C-Leg (Kannenberg et al. 2013).
	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. Variation of gait speed exhibit the highest number of sensitive confounders in the functional benefits category. However, the investigated variables failed to exhibit classifying power (e.g. variation of gait speed impacted by mobility grade presented a high significance ($p < 3x10^{-26}$) but a very low r^2 (0.13)).
References of summarized studies	Bellmann, M., Schmalz, T., Ludwigs, E., & Blumentritt, S. (2012). Immediate effects of a new microprocessor-controlled prosthetic knee joint: a comparative biomechan- ical evaluation. Archives of physical medicine and rehabilitation, 93(3), 541–549. doi:10.1016/j.apmr.2011.10.017
	Blumentritt, S., Bellmann, M., Ludwigs, E., & Schmalz, T. (2012). Zur Biomechanik des mikroprozessorgesteuerten Prothesenkniegelenks Genium. Orthopädie- Technik, 01(12), 24–35.
	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M. J., Kahle, J. T., Lura, D. J., Lewandowski, A. L., Quillen, W. S., & Kim, S. H. (2014). Stair ascent and ramp gait training with the Genium knee. Technology & Innovation, 15(4), 349–358. doi:10.3727/194982413X13844488879267
	Kampas, P., & Seyr, M. (2011). Technologie und Funktionsweise des Genium- Prothesenkniegelenks. Orthopädie-Technik, 12(11), 898–903.
	Kannenberg, A., Zacharias, B., Mileusnic, M., & Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f
	Lura, D. J., Wernke, M. M., Carey, S. L., Kahle, J. T., Miro, R. M., & Highsmith, M. J. (2015). Differences in knee flexion between the Genium and C-Leg microprocessor knees while walking on level ground and ramps. Clinical Biomechanics, 30(2), 175–181. doi:10.1016/j.clinbiomech.2014.12.003

↑ Back to overview table

Stairs

Major Findings

With Genium compared to C-Leg:

- → Difficulty of ascending and descending stairs is decreased by 34% and 10%
- → 70 80% of subjects are able to ascend stairs with a reciprocal step-overstep strategy
- \rightarrow Loading of sound side is reduced by 10% during ascent
- → Movement of sound side is within range of healthy subjects during ascent Movement of hip and knee joints on sound side are reduced by 34% and 33%

Percent of subjects by stair ascent strategy with Genium



Stair ascent strategy was assessed observing 14 subjects (Aldridge Whitehead et al. 2014).

Clinical Relevance Stair ambulation is an activity that is important for amputees with an activity level ranging from K2 to K4. Being able to ascend and descend stairs is a requirement to participate in daily life. Evaluation of stair ascent includes stair ascent strategy, use of handrail and/or use of an assistive device. Biomechanical assessment is conducted to determine load on the joints and joint angles and to compare them to values measured in healthy subjects.

Summary

The assessment of stair ascent strategy was conducted by multiple groups. Bellmann et al. (2012) reported that 80% of subjects were able to ascend stairs reciprocally after only one day of using Genium. These findings were confirmed by Highsmith et al. (2014), reporting 70% of subjects, and by Aldridge Whitehead et al. (2014), reporting 72%. Lura et al. (2017) found that from the subjects who were able to use a reciprocally step-over-step strategy to ascend stairs, 41 % preferred it over step-to-step strategies, compared to C-Leg users in which preference was only 5%. Furthermore, an improvement of the Stair Assessment Index (SAI) from 5 to 11 points (Aldridge Whitehead et al. 2014) respectively from 6 to 11 points (Highsmith et al. 2016) was achieved by switching from C-Leg or an NMPK to Genium. A score of 11 points represent a reciprocal stair ascent strategy with hand rail or assistive device use on a scale where a score of 13 points is the maximum. The reciprocal stair ascent strategy is accomplished with Genium by means of on an additional function, activated by a backward movement of the prosthesis after lifting the foot (Kampas et al. 2011).

	Although the duration of a stride during stair ascent is longer with Genium than with C-Leg, a clear approximation to the movement pattern of healthy subjects was observed when using Genium. The movements of the knee and hip joints on the contralateral side were decreased and the loading of the contralateral knee joint was reduced with Genium compared to C-Leg (Bellmann et al. 2012 & Blumentritt et al. 2012). Moreover, hip and knee flexion during swing phase were increased and therefore toe clearance during stair ascent increased (Aldridge Whitehead et al. 2014). The peak flexion angle and swing period increased significantly with Genium (Lura et al. 2017). All these changes led overall to a decreased between-limb difference with Genium compared to C-Leg.
	Besides improvement in stair ascent strategy and gait characteristics, subjects re- ported that ascending and descending stairs was perceived as less difficult with Genium than with C-Leg (Kannenberg et al. 2013).
	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. Descent and ascent from stairs presented a very clear responsiveness in 38.3% and 63.10% of subject's perception category, respectively. Reciprocal stair ascent exhibit the highest number of sensitive confounders (25) for the advanced manoeuvres category. However, none of the factors qualified as predictor for performance.
References of summarized studies	Aldridge Whitehead, J. M., Wolf, E. J., Scoville, C. R., & Wilken, J. M. (2014). Does a Microprocessor-controlled Prosthetic Knee Affect Stair Ascent Strategies in Persons With Transfemoral Amputation? Clinical Orthopaedics and Related Re- search [®] . doi:10.1007/s11999-014-3484-2
	Bellmann, M., Schmalz, T., Ludwigs, E., & Blumentritt, S. (2012). Immediate effects of a new microprocessor-controlled prosthetic knee joint: a comparative biomechan- ical evaluation. Archives of physical medicine and rehabilitation, 93(3), 541–549. doi:10.1016/j.apmr.2011.10.017
	Bellmann, M., Schmalz, T., Ludwigs, E., & Blumentritt, S. (2012). Stair ascent with an innovative microprocessor-controlled exoprosthetic knee joint. Biomedizinische Technik. Biomedical engineering, 57(6), 435–444. doi:10.1515/bmt-2011-0029
	Blumentritt, S., Bellmann, M., Ludwigs, E., & Schmalz, T. (2012). Zur Biomechanik des mikroprozessorgesteuerten Prothesenkniegelenks Genium. Orthopädie- Technik, 01(12), 24–35.
	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M. J., Kahle, J. T., Lura, D. J., Lewandowski, A. L., Quillen, W. S., & Kim, S. H. (2014). Stair ascent and ramp gait training with the Genium knee. Tech- nology & Innovation, 15(4), 349–358. doi:10.3727/194982413X13844488879267
	Highsmith, M.J., Klenow, T.D., Kahle, J.T., Wernke, M.M., Carey, S.L., Miro, R.M., Lura, D.J. (2016). Effects of the Genium microprocessor knee system on knee mo- ment symmetry during hill walking. Technology & Innovation, 18: 151-157.
	Kampas, P., & Seyr, M. (2011). Technologie und Funktionsweise des Genium- Prothesenkniegelenks. Orthopädie-Technik, 12(11), 898–903.
	Kannenberg, A., Zacharias, B., Mileusnic, M., & Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f
	Lura, D. J., Wernke, M. M., Carey S. L., Kahle, J. T. Miro, R. M. & Highsmith, M. J. (2017). Crossover study of amputee stair ascent and descent biomechanics using

Genium and C-Leg prostheses with comparison to non-amputee control. Gait & Posture, 58, 103-107.

A Back to overview table

Ramps / Hills

Major Findings	With Genium compared to C-Leg:
	ightarrow Ease of ascending and descending ramps is improved by 24% and 17%
	→ Improved toe clearance during ramp ascent and descent
	Knee angle in swing phase increased by 14° during ramp ascent
	Knee angle in swing phase increased by 13° during ramp descent
	→ Improved symmetry of knee moment while ascending and descending
	ramps
	Ascent: significantly lower at slow and fast walking speed
	Descent: significantly lower at slow and self-selected walking speed



Ramp ambulation is perceived as less difficult with Genium

Kannenberg et al. (2013)

Clinical Relevance Similar to stairs, ramps and hills need to be navigated by amputees with a wide range of activity levels to be able to participate in daily life. Biomechanical assessment is conducted to determine joint angles. The maximum knee flexion angle in swing phase was of special interest. To minimize the risk of stumbling, it is required to have an increased maximum flexion angle relative to level walking to ensure an adequate foot clearance.

Summary

Maximum knee angle was increased by 7° when ascending and by 8° when descending a ramp with Genium compared to C-Leg (Bellmann et al. 2012 and Blumentritt et al. 2012). Increased maximum knee angle in swing phase leads to an increased foot clearance which further decreases the risk of stumbling.

Besides increased safety, patients reported that ascending and descending ramps is less difficult to perform with Genium than with C-Leg (Kannenberg et al. 2013) and also that it's more comfortable. When ascending ramps, less focal pressure near the anterior aspect of the hip was experienced (Highsmith et al. 2014).

Highsmith et al. (2014) observed that when subjects with Genium descended a ramp with only a slight decline such as 5°, a walking pattern similar to level walking characterized by two knee flexion peaks was facilitated. This feature of Genium is enabled through an adapted resistance in stance phase; it allows for a flexion angle which is higher than the maximum knee angle when level walking (17°). In comparison to C-Leg, peak knee flexion angle during descending a slope of 5° was with Genium in swing phase increased by up to 8° and during stance phase increased by up to 4° over a variety of gait velocities from slow to fast. Peak knee flexion an-

	gles of the prosthetic side are therefore closer to peak knee flexion angles of the intact leg and therefore a more normalized, anatomic movement pattern is achieved with Genium (Lura et al. 2015). Similar results were reported by Bell et al. (2016) where additionally more knee flexion at initial heel strike and swing phase were observed as well as faster walking speed by 0.1 m/s compared to other MPKs. Also the knee moment while ascending and descending ramps was more symmetrically with Genium (Highsmith et al. 2016). The values were significant at slow and fast walking speed when ascending a ramp and at slow and self-selected walking speed when descending.
	Furthermore, when walking on ramps swing phase release occurs even when the prosthesis is in a flexed and loaded position (Kampas et al. 2011). Overall, a more natural gait pattern on ramps is achieved with Genium.
	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. Descent from ramps presented a very clear responsiveness in 59 % of subject's perception. Ascent (57.4%) and standing on ramps (76.4%) presented also a very clear responsiveness as performance indicators. However, none of the factors qualified as predictor for performance.
References of summarized studies	Bell, E.M., Pruziner, A.L., Wilken, J.M., Wolf, E.J. (2016). Performance of conv tional and X2 [®] prosthetic knees during slope descent. Clin Biomech, 33: 26–31.
	Bellmann, M., Schmalz, T., Ludwigs, E., & Blumentritt, S. (2012). Immediate effects of a new microprocessor-controlled prosthetic knee joint: a comparative biomechanical evaluation. Archives of physical medicine and rehabilitation, 93(3), 541–549. doi:10.1016/j.apmr.2011.10.017
	Blumentritt, S., Bellmann, M., Ludwigs, E., & Schmalz, T. (2012). Zur Biomechanik des mikroprozessorgesteuerten Prothesenkniegelenks Genium. Orthopädie- Technik, 01(12), 24–35.
	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M.J., Klenow, T.D., Kahle, J.T., Wernke, M.M., Carey, S.L., Miro, R.M., Lura, D.J. (2016). Effects of the Genium microprocessor knee system on knee mo- ment symmetry during hill walking. Technology & Innovation, 18: 151-157.
	Highsmith, M. J., Kahle, J. T., Lura, D. J., Lewandowski, A. L., Quillen, W. S., & Kim, S. H. (2014). Stair ascent and ramp gait training with the Genium knee. Technology & Innovation, 15(4), 349–358. doi:10.3727/194982413X13844488879267
	Kampas, P., & Seyr, M. (2011). Technologie und Funktionsweise des Genium- Prothesenkniegelenks. Orthopädie-Technik, 12(11), 898–903.
	Kannenberg, A., Zacharias, B., Mileusnic, M., & Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f
	Lura, D. J., Wernke, M. M., Carey, S. L., Kahle, J. T., Miro, R. M., & Highsmith, M. J. (2015). Differences in knee flexion between the Genium and C-Leg microprocessor knees while walking on level ground and ramps. Clinical Biomechanics, 30(2), 175–181. doi:10.1016/j.clinbiomech.2014.12.003

Uneven Ground / Obstacle Course

Major Findings

With Genium compared to C-Leg:

- → Clinically relevant increase in perceived ease for walking on unknown or uneven terrain.
- → Clinically relevant increase in perceived ease for stepping over or stepping on minor obstacles.
- → Stair climbing training improves the ability to cross obstacles since same unique function is applied

Differences in perceived ease and safety on uneven ground for Genium compared to C-Leg



Perceived difficulty and safety was rated by experienced C-Leg users after 3 months of Genium use (Kannenberg et al. 2013).

Clinical Relevance

Summary

Besides ambulating on stairs and ramps, overcoming obstacles and walking on different types of terrain are requirements to be able to participate in life. These activities are together with level walking a solid base for an amputee to ambulate independently. The ability to walk on uneven terrain is usually measured by the time required to complete a course and the perceived exertion.

Highsmith et al. (2014) studied the impact of different prostheses when amputees completed a course on sloping terrain over trimmed grass, sand, rocks and small roots. They could not observe any difference in gait speed between Genium and C-Leg; however, amputees showed a trend to rate perceived exertion lower when using Genium. When asked about the difficulty in an additional study, amputees rated walking on uneven and unknown terrain as clinically relevant less difficult and safer with Genium than with C-Leg. Furthermore, clinically relevant increase in perceived ease for stepping over and stepping on minor obstacles was reported (Kannenberg et al. 2013).

In a study which investigated stair ascent training using Genium, subjects believed that stair climbing practice improved their ability to cross obstacles (Highsmith et al. 2014). It is confirmed in a technical report that obstacles can be crossed with the same function as used for ascending stairs, activated by a backward movement of the prosthesis after lifting the foot. Using this function, the obstacle is crossed with the prosthetic limb first. Conventionally, crossing an obstacle is only possible through either stepping with the sound limb first or through a circumduction movement (Kampas et al. 2011).

	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. Crossing obstacles presented a very clear responsiveness in 65.9 % of subject's perception. However, none of the factors qualified as predictor for performance.
References of summarized studies	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M. Jason; Kahle, Jason T.; Lura, Derek J.; Dubey, Rajiv V.; Carey, Stephanie L.; Quillen, William S.; Mengelkoch, Larry J. (2014): Short and Mid- Distance Walking and Posturography With a Novel Microprocessor Knee. In: Tech- nology & Innovation 15 (4), S. 359–368. DOI: 10.3727/194982413X13844488879302.
	Highsmith, M. J., Kahle, J. T., Lura, D. J., Lewandowski, A. L., Quillen, W. S., & Kim, S. H. (2014). Stair ascent and ramp gait training with the Genium knee. Tech- nology & Innovation, 15(4), 349–358. doi:10.3727/194982413X13844488879267
	Kampas, P., & Seyr, M. (2011). Technologie und Funktionsweise des Genium- Prothesenkniegelenks. Orthopädie-Technik, 12(11), 898–903.
	Kannenberg, A., Zacharias, B., Mileusnic, M., & Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f

↑ Back to overview table

Safety

Major Findings

With Genium compared to C-Leg:

- → 60% of the examined activities of daily living (ADLs) are rated as safer above a threshold considered to be clinically relevant
- → Movement control in all three backward directions improved by up to 10%





Safety was assessed for 45 activities (Kannenberg et al. 2013).

Clinical Relevance	Safety aspects of the prosthesis are highly relevant for the patients. Since the fear of falling can have a negative impact on activities of daily living as well as on participation, perceived safety is regarded as an important factor for quality of life of an amputee. Information about perceived safety when performing different activities is gathered through a questionnaire. Balance tests are conducted to obtain objective information about the patients' ability to react in potential falling situations.
Summary	The results of an activity of daily living (ADL) questionnaire show, that 60% of the examined ADLs were rated as safer with Genium compared to C-Leg with a comparative rating above the threshold considered to be a clinically relevant difference (25% of the maximum possible difference). The other 40% of the examined ADLs also showed improved safety, but below this threshold. Especially the category 'Family and Social Life' with 83% of ADLs rated as safer, and 'Mobility and Transportation' with 63% of ADLs rated as safer were clearly in favour for Genium (Kannenberg et al. 2013).
	Highsmith et al. examined 2016 the safety with the four square step test, a test of dynamic balance that clinically assesses the person's ability to step over objects forward, sideways, and backwards. The test could be finished by 9% faster with Genium than with C-Leg.
	Other studies tested the influence of Genium when standing on a decline. It was observed that with Genium not only was the body weight evenly distributed on both legs, but also the body posture was more relaxed than with C-Leg. The latter can most likely be explained by the reduced activity required to keep up the posture and the decreased hip moments which were found in patients using Genium. Furthermore, observed that the postural sway on the prosthetic side is reduced when standing on a decline using Genium compared to C-Leg (Bellmann et al. 2012 and Blumentritt et al. 2012).

	Highsmith et al (2014) investigated postural stability measured by a balance system. They could show that posterolateral directional stability over the sound side is im- proved with Genium compared to C-Leg.
	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. The category Safety within Functional Benefits presented a very clear responsiveness of 49.94%. However, none of the factors qualified as predictor for performance.
References of summarized studies	Bellmann, M., Schmalz, T., Ludwigs, E., & Blumentritt, S. (2012). Immediate effects of a new microprocessor-controlled prosthetic knee joint: a comparative biomechan- ical evaluation. Archives of physical medicine and rehabilitation, 93(3), 541–549. doi:10.1016/j.apmr.2011.10.017
	Blumentritt, S., Bellmann, M., Ludwigs, E., & Schmalz, T. (2012). Zur Biomechanik des mikroprozessorgesteuerten Prothesenkniegelenks Genium. Orthopädie- Technik, 01(12), 24–35.
	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M. Jason; Kahle, Jason T.; Lura, Derek J.; Dubey, Rajiv V.; Carey, Stephanie L.; Quillen, William S.; Mengelkoch, Larry J. (2014): Short and Mid- Distance Walking and Posturography With a Novel Microprocessor Knee. In: Tech- nology & Innovation 15 (4), S. 359–368. DOI: 10.3727/194982413X13844488879302.
	Highsmith, M.J., Klenow, T.D., Kahle, J.T., Wernke, M.M., Carey, S.L., Miro, R.M., Lura, D.J. (2016). Effects of the Genium microprocessor knee system on knee mo- ment symmetry during hill walking. Technology & Innovation, 18: 151-157.
	Kampas, P., & Seyr, M. (2011). Technologie und Funktionsweise des Genium- Prothesenkniegelenks. Orthopädie-Technik, 12(11), 898–903.
	Kannenberg, A., Zacharias, B., Mileusnic, M., & Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f

↑ Back to overview table

Activity / Mobility / Activity of daily living (ADL)

Major Findings	With Genium compared to C-Leg:
	→ 53% of the examined activities of daily living (ADLs) show a clinically relevant decrease in difficulty → 75% of ADLs in category 'Eamily and Social Life' show a clinically relevant
	decrease in difficulty
	→ 58% of ADLs in category 'Mobility and Transportation' show a clinically relevant decrease in difficulty
	 → Mobility increased by 5% → Physical functional performance scores significantly higher with Genium
	Clinically relevant decrease in difficulty to perform acitivities of daily living with Genium
	47% 53%
	Less difficult with Genium, below clinically relevant threshold

Perceived difficulty was assessed for 45 activities (Kannenberg et al. 2013).

Clinical Relevance	With an Activity of daily living (ADL) questionnaire, information is gathered about self-care activities such as functional mobility, dressing, eating and personal hygiene, as well as activities for living independently in a community, such as shopping, housework and transportation. The questionnaire is used as a tool to measure the general independence of patients.				
Summary	The results of an activity of daily living (ADL) questionnaire demonstrate that 53% of the examined ADLs were rated to show a clinically relevant decrease in perceived difficulty with Genium compared to C-Leg. The other 47% of examined ADLs showed a trend to be rated as less difficult. Especially in the categories, 'Family and Social Life' and 'Mobility and Transportation', the results were clearly in favour of Genium: 75% and 58% of the examined ADLs, respectively, were rated as less difficult to perform (Kannenberg et al. 2013). This improvement in the activities of daily living was confirmed by a study by Highsmith et al. (2016). The Perceived function and safety in three of five ADL domains improved significantly. The other two domains showed no difference. This study group also showed that the overall mobility increased by 5%, as measured with the Amputee Mobility Predictor.				
	In a case study, a bilateral transfemoral amputee was fitted with Genium with osseo- integrated prosthesis fixation (OPF). Compared to the situation before OPF, the patient was able to perform four more activity categories of daily living listed in lower extremity functional scale (LEFS) (any of your usual work, housework, etc; usual hobbies, recreational or sporting; getting into or out the bath; walking between rooms) (Schalk et al. 2015).				

	Functional performance differences were assessed between Genium, C-Leg and intact knees by Highsmith et al. 2016 (b) using the Continuous-Scale Physical Functional Performance-10 (CS-PFP10). Subjects with Genium presented a significant (p<0.05) score improvement for upper-body flexibility, balance and endurance, compared to C-Leg users. Compared to non-amputees, Genium users presented a significantly lower score only for endurance domain.
	In a retrospective, cross-sectional cohort analysis from Hahn et al. 2016, clinically important factors on performance using Genium were analysed based on 899 trial fittings. While none of the factors qualified as predictor for performance, toileting was found to be the most responsive indicator for subjects' perception. The complexity of the task may have been insufficiently considered so far.
References of summarized studies	Hahn, A., Lang, M., Stuckart, C. (2016). Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings. Medicine, 95 (45):e5386.
	Highsmith, M. J. Kahle, J. T., Miro, R. M. Cress, M. E., Lura, D. J., Quillen, W. S. Carey, S.L., Dubei, R. V., Mengelkoch, L. J. (2016)(b). Functional performance differences between the Genium and C-Leg prosthetic knees and intact knees. J Rehabil Res Dev. 2016;53(6):753–66.
	Highsmith, M.J., Klenow, T.D., Kahle, J.T., Wernke, M.M., Carey, S.L., Miro, R.M., Lura, D.J. (2016). Effects of the Genium microprocessor knee system on knee mo- ment symmetry during hill walking. Technology & Innovation, 18: 151-157.
	Kannenberg, A., Zacharias, B., Mileusnic, M., Seyr, M. (2013). Activities of Daily Living: Genium Bionic Prosthetic Knee Compared With C-Leg. JPO Journal of Prosthetics and Orthotics, 25(3), 110–117. doi:10.1097/JPO.0b013e31829c221f
	Schalk, S.A.F., Jonkergouw, N., van der Mer, F., Swaan, W.M., Aschoff, H.H., van der Wurff, P. (2015). The Evaluation of Daily Life Activities after Application of an Osseointegrated Prosthesis Fixation in a Bilateral Transfemoral Amputee: A Case Study. Medicine (Baltimore), 94(36): e1416.

A Back to overview table

Preference / Satisfaction / Quality of Life (QoL)

Major Findings	With Genium compared to C-Leg:
	→ Improved prosthesis-related quality of life
	→ 45% of prosthetic-related scales are improved
	→ 70% of activities, relevant to the physical performance aspects of Genium, are improved
	→ 80% preferred Genium

Percent of scales showing improvement in prosthetic function and quality of life with Genium



In total 41 items divided in 9 scales were assessed (Highsmith et al. 2014).

Definition/Clinical Relevance	Satisfaction and quality of life can be measured to determine the general well-being of a person. They are all very meaningful parameters to investigate, since they have the most direct impact on the amputee's well-being. They are influenced by other categories and can therefore be seen as a summary of possible activities, independence and perceived safety. A common outcome measure in prosthetic research is the Prosthesis Evaluation Questionnaire (PEQ), a questionnaire with a total of 84 items. Several selected items are further used in 9 subscales.
Summary/Discussion	Comparing all items of the Prosthesis Evaluation Questionnaire (PEQ) in aggregate, the Genium was rated to be improved compared to C-Leg. 4 out of 9 scales were rated as improved, namely perceived response, social burden, utility, and well-being. In addition, the scales appearance and sounds, showed a trend to be improved with Genium. Furthermore, 10 items of special interested were looked at separately. 7 out of these 10 were improved with Genium compared to C-Leg; comfort standing, ability to walk in close spaces, go down stairs, walking up and down steep hills, walk on slippery surfaces and satisfaction with walking (Highsmith et al. 2014). In one publication by Highsmith et al. (2016) the preference for Genium was 80%.
	In a case study, a bilateral transfemoral amputee was fitted with Genium with osseo- integrated prosthesis fixation (OPF). Life habits were recorded with LIFE-H ques- tionnaire for daily activities and social roles and showed increases in five and eight respectively out of 12 categories in comparison to the situation before OPF (Schalk et al. 2015).

References of summarized studies	Highsmith, M.J., Klenow, T.D., Kahle, J.T., Wernke, M.M., Carey, S.L., Miro, R.M., Lura, D.J. (2016). Effects of the Genium microprocessor knee system on knee mo- ment symmetry during hill walking. Technology & Innovation, 18: 151-157.
	Highsmith, M. J., Kahle, J. T., Miro, R. M., Lura, D. J., Dubey, R. V., Carey, S. L., . Mengelkoch, L. J. (2014). Perceived differences between the Genium and the C- Leg microprocessor prosthetic knees in prosthetic-related function and quality of life. Technology & Innovation, 15(4), 369–375. doi:10.3727/194982413X13844489091297
	Schalk, S.A.F., Jonkergouw, N., van der Mer, F., Swaan, W.M., Aschoff, H.H., van der Wurff, P. (2015). The Evaluation of Daily Life Activities after Application of an Osseointegrated Prosthesis Fixation in a Bilateral Transfemoral Amputee: A Case

Study. Medicine (Baltimore), 94(36): e1416.

↑ Back to overview table

3 Summaries of individual studies

On the following pages you find summaries of studies that researched Genium. You find detailed information about the study design, methods applied, results and major findings of the study. At the end of each summary you also can read the original study authors' conclusions.

Reference	Lura DJ, Wernke MM, Care	ey SL, Kahle JT, Miro RM, H	ighsmith MJ.			
	School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.					
	Crossover study of amputee stair ascent and descent biomechanics using Genium and C-Leg prostheses with comparison to non-amputee					
	control					
	Gait & Posture 2017; 58: 1	03-107.				
Products	Genium vs C-Leg					
Major Findings	With Genium compared to	C-Leg:				
	angle and swing period → Ability and preference → With Genium, 41% of s pattern prefer it to step prefer step-over-step (vious gait pattern.	d increased significantly). to use a step-over-step ga ubjects who were able to p-to-step (ST); with C-Leg, SOS), likely due to comfor	it increased significantly. use a step-over-step (SOS) only 5% subjects would t and stability of their pre-			
	Stair ascent	step-over-step (SOS)* gai	t, mean values			
	Genium (GP) and contralateral (GC, n=17); C-Leg (CP) and contralateral (CC, n=4);					
	(n =	subjects who were able to use SO	S gait).			
	140 120 Lagrage Flexion (Degrees) 100 100 100 100 00 100 100 100	2,2 1,8 1,4 1,4 1,4 0,6 0,2 -0,2 -0,2 -0,2	1,2 1,2 1,2 0,8 0,0 0,6 0,4 0,0 0,4 0,2 0,2 0,2 0,2 0,2 0,2 0,2 0,2			

*"[...] since participants were not asked to perform the ST gait if they preferred the SOS gait, and not all participants were able to perform the SOS gait, each grouping did not have the same number of trials. For the descent trials, since only one participant preferred the ST gait for descent with the Genium there was not a sufficient sample to perform statistical analysis with for study outcomes for ST descent" (Lura et al., 2017)

Population	Subjects:	20 unilateral, transfemoral amputees, 5 non- amputees
	Previous prosthesis:	C-Leg
	Amputation causes:	70% trauma, 20% malignancy, 10% vascular disease
	Mean age:	Patients: 46.5 yrs (±14.2 yrs), Controls 57.2 (±15.7 yrs)
	Mean time since amputation:	17.7 yrs (±15.6 yrs)
	MFCL:	КЗ

Study Design

Interventional, randomized crossover design:



Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results (compared to non-amputees)	
		C-Leg	Genium
Stairs	Descent		
	Peak flexion angle		
	prosthetic to contralateral side		
	contralateral side to non-amp.		
	Peak extension moment		
	contralateral to prosthetic side contra-	++	++
	lateral side to non-amp.	++	++
	Prosthetic side to non-amp.	0	0
	Swing duration		
	prosthetic to contralateral side	+	+
	prosthetic side to non-amp.	+	+
	contralateral side to non-amp.	+	+
	C-Leg to Genium	+	
	Ascent (step-over-step pattern)		
	Peak flexion angle		
	prosthetic to contralateral side		++
	prosthetic side to non-amp.		0
	contralateral side to non-amp.	0	0
	Genium to C-Leg		++

Category	Outcomes	Results (compared to non-ampute	
		C-Leg	Genium
	Peak extension moment		
	prosthetic to contralateral side		
	prosthetic side to non-amp.		
	Genium to C-Leg prosthetic		+
	Genium to C-Leg contralat.		+
	Swing duration		
	prosthetic to contralateral side	++	++
	prosthetic side to non-amp.	++	++
	contralateral side to non-amp.	0	0
	Genium to C-Leg		++

Author's Conclusion

"Use of the Genium knee enabled the majority of the participants to use a reciprocal SOS gait pattern for stair ascent, increased knee flexion during swing phase of stair ascent, and generally contributed to a more symmetric gait. However, there was not a significant change in gait parameters for participants while descending stairs, and the swing duration while using the Genium for stair ascent was marginally longer than while using the C-Leg." (Lura et al., 2017)

A Back to overview table

Reference	Hahn A, Lang M, Stuckart C.
	Analysis of clinically important factors on the performance of advanced hydraulic, microprocessor-controlled exo-prosthetic knee joints based on 899 trial fittings
Products	Genium®, C-Leg®, mechanical knee joints
Major Findings	 Overall responsiveness exceeded 90% after 1 week trial fitting. Genium responders span a wide range of demographic and epidemiologic characteristics. Investigated variables failed to exhibit classifying power. Eligibility decisions for the patient based on such variables only (e.g. mobility grade, comorbidities, BMI and over 50 more) might thus be unjustified. Single and multiple regression models detected influences of mobility grade, residual limb conditions, socket type, use of liners, vacuum technology, dynamic foot and others on likelihood to improve performance indicators. Trial fittings may pose a very appropriate method to investigate individuals' potential to benefit from Genium. Bias in the sample is estimated to be below 5%. Trial outcomes vary in sensitivity. Variation of gait speed, perception of toileting and reciprocal stair ascent exhibit the highest number of sensitive confounders in their respective categories (functional benefits (FB), subject perception (SP) and advanced maneuvers (AM), respectively). C-Leg walkers profit from an upgrade to Genium in advanced maneuvers and subjects' perception categories.



Responsiveness to Functional Benefits (FB) variable subgroup

Population

899

Subjects:

Amputation level:

Mean age:

MFCL:

Previous prosthetic knee: C Leg (689), mechanical hydraulic knees (38), pneumatic knees(22), other polycentric (19), 4-bar knees (15), brake knees (9), locked knees (2). Amputation causes: 68.9% trauma, 15.4% tumor, 6% vascular disease. 80.1% TF, 18.9% KD. 49.0 <u>+</u> 12.9 y. First prosthesis since: 21.2 <u>+</u> 15.6 y. 12.5% K2, 64.1% K3, 22.8% K4.

Retrospective, cross-sectional cohort analysis.



Data from routine trial fittings was retrieved from customer support service, from 272 prosthetic clinics between 2011 and 2015. Influence of clinical variables (i.e. mobility grade) on performance indicators (functional befits (FB), subjects perception (SP), advanced manoeuvring (AM)) were analysed using single and multiple linear and logistic regression modelling. FB by prosthetist and SP by patient were assessed on 5 point Likert scales. The upper two ratings of the Likert scale were defined as being respondent. Advanced manoeuvres were categorized by the prosthetist non-ordinally and dichotomized for analysis. The number of statistically relevant confounders were identified in each performance-indicating category. Clinical significance was derived by regression estimate (e).

Results

Functions an	d Activities					Participation			Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health Economics
Category		Outco	omes		Results f	for Genium	®		Sig.*
Category Level walking		Harm (FB)	onization of	gait patterr	FB total r very cle clear: 4 # of de	% nders: 10	++		
C (F R (F V A	Capa (FB, S	bility to vary SP)	r gait speed	FB total r very cle clear: 3 # of det	++				
				SP total r very cle clear: 3 # of de	++				
				Impacted e: 0.36	++				
	Redu (FB)	ction of ove	rall effort	FB total r very cle clear: 4 # of de	% nders: 13	++			
	Walki AM)	ng backwai	ds (SP,	SP total r very cle clear: 2 # of de	% nders: 5	++			
				AM responsiveness safety: 92.28% # of detected sensitive confounders: 16				++	
		"Door	r Test" (AM)		AM respo	onsiveness s	afety: 88.3	5 %	++

Category	Outcomes	Results for Genium®	Sig.*
	Abrupt change of walking direction	# of detected sensitive confounders: 13	
Stairs	Descent (SP, AM)	SP total responsiveness: 66.82 % Very clear: 38.30 % Clear: 28.52% # of detected sensitive confounders: 13	++
		AM responsiveness safety: 70.52 % # of detected sensitive confounders: 20	++
	Ascent (reciprocal) (SP, AM)	SP total responsiveness: 88.64% very clear: 63.10% clear: 25.54 % # of detected sensitive confounders: 12	++
		AM responsiveness safety: 32.81 % # of detected sensitive confounders: 25	++
Ramps, Hills	Descent (SP, AM)	SP total responsiveness: 85.78% very clear: 59% clear: 26.78 % # of detected sensitive confounders: 6	++
		AM responsiveness safety: 70.86 % # of detected sensitive confounders: 21	++
	Ascent (SP, AM)	SP total responsiveness: 86.03% very clear: 57.40% clear: 28.63 % # of detected sensitive confounders: 9	++
		AM responsiveness safety: 74.53 % # of detected sensitive confounders: 14	++
	Standing (SP)	SP total responsiveness: 95.69% very clear: 76.60% clear: 19.09 % # of detected sensitive confounders: 11	++
Uneven ground, Obstacles	Crossing Obstacles (SP, AM)	With prosthetic side first SP total responsiveness: 92.36% very clear: 65.90% clear: 26.46% # of detected sensitive confounders: 8	++
		With contralateral side first SP total responsiveness: 85.32% very clear: 54% clear: 31.32% # of detected sensitive confounders: 5	++
	Stepping on Obstacles (SP)	# detected sensitive confounders: 10	++
Cognitive Demand	Divided Attention (FB)	FB total responsiveness: 98.05% very clear: 56.57% clear: 41.48% # of detected sensitive confounders: 18	++
	Dual Tasking (SP)	# detected sensitive confounders: 14	++
Metabolic Energy Consumption	Carrying objects with visual obstruction (SP, AM)	SP total responsiveness: 82.83% very clear: 51.30% clear: 31.53%	++

Category	Outcomes	Results for Genium®	Sig.*
		# detected sensitive confounders: 9	
		AM responsiveness safety: 87.74% # detected sensitive confounders: 11	++
	Carrying heavy loads (SP, AM)	SP total responsiveness: 82.40% very clear: 51% clear: 31.40% # detected sensitive confounders: 6	++
		AM responsiveness safety: 72.73% # detected sensitive confounders: 14	++
Safety	Safety (FB)	FB total responsiveness: 97.04% very clear: 49.94% clear: 47.10% # detected sensitive confounders: 14	++
Activity, Mobility, Activities of Daily Living (ADLs)	Change of mobility grade	FB total responsiveness: 96% very clear: 57.10% clear: 38.90% # detected sensitive confounders: 14	++
	Toileting (SP)	# detected sensitive confounders: 18	++

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

"Responders to Genium trial fittings span a wide range throughout the entire investigated variables. None of those variables nor their combination seem to **qualify as predictor for an individual response to a performance indicator**. This is confirmed for mobility grade and further includes age, etiology, residual limb conditions, and comorbidities. BMI fails to exhibit statistical significance. Decision making processes that rely on those variables without appropriately considering the subjects' individual potential and capabilities do not seem to be supported by these findings. As no data is available supporting such approaches, the denial of access to advanced technology based on such variables may indeed be questionable. A threshold value for walking capacity cannot be excluded and may pose a **component of a possible predictor**. Future research may consider a minimum walking capacity as a component of a predictive instrument.

Toileting was identified as the most responsive indicator in the subject's perception. Difficulties associated with this specifically demanding task may be insufficiently considered and may play a more important role when deciding upon the appropriate prosthetic components. Future protocols for trial fittings may consider limiting the number of performance indicators to those with high differentiating power. Subjects having previously been fitted with C-Leg show benefits when fitted with Genium. Most of these benefits can be found in perception and advanced manoeuvres among which is stairs ascent. Liners and the use of a higher dynamic response foot further contribute to a better utilization of functional benefits." (Hahn et al., 2016).

↑ Back to overview table

Reference	Highsmith MJ, Kahle JT, Miro RM, Cress ME, Lura DJ, Quillen WS.						
	School of Physical Therapy and Rehabilitation Sciences, University of South Florida, Tampa, FL.						
	Functional performance differences between the						
	Genium and C-Leg prosthetic knees and intact						
	knees						
	Journal of Rehabilitation Research and Development 2016; 53(6): 753-766.						
Products	Genium vs C-Leg (vs Nonamputees)						
Major Findings	 With Genium compared to C-Leg: → Genium scores significantly higher than C-Leg in the upper-body flexibility (UBF), balance(BAL) and endurance (END) domains → Genium scores higher than C-Leg upper-body and lower-body strength scores 						
	Genium scores higher than C-Leg in Physical Functional Performance 90						
	80 [++] [++] [++]						
	50 50 ■C-Leg						
	O 30 Image: Solution of the second seco						

UBS = upper-body strength; LBS = lower-body strength; UBF = upper-body flexibility; BAL = balance; END = endurance; ++ p < 0.05 (significant)

UBF

BAL

END

LBS

Population	Subjects:	20 unilateral transfemoral amputees (AMP), 5 nonamputee controls (NAMP)
	Previous prosthesis:	C-Leg
	Amputation causes:	Trauma (75%), Malignancy (20%), PVD (5%)
	Mean age:	46.5 ± 14.2 yrs (AMP); 57.2 ± 15.7 yrs (NAMP)
	Mean time since amputation:	17.7 ± 15.6 yrs
	MFCL:	MFCL 3

UBS

10 0

Total score

Interventional, randomized crossover design:



Results

Functions and Activities Level Stairs Ramps, Uneven Cognitive Energ					Participation			Environment	
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Category	Outcomes	Results for Genium vs C-Leg	Sig.*
Activity, Mobility, Activities of Daily Living (ADLs)	CS-PFP10 (Continuous-Scale Physical Functional Perfor- mance 10)	The total score showed a trend of improvement by 7.4%.	+
		The upper-body flexibility score was significantly improved by 7%.	++
		The balance score was significantly improved by 7.6%.	++
		The endurance score was significantly improved by 8.4%.	++
		The upper-body and lower-body strength scores showed a trend of improvement (+5.4% and +8.1%).	+
Category	Outcomes	Results for Genium vs nonamputees	Sig.*
Activity, Mobility, Activities of Daily Living	CS-PFP10 (Continuous-Scale Physical Functional Perfor-	Nonamputees scored higher in all five domains but only significantly higher in	+
(ADLS)	mance 10)	the endurance domain (+22.4%).	++
Category	Outcomes	Results for C-Leg vs nonamputees	Sig.*
Activity, Mobility, Activities of Daily Living (ADLs)	CS-PFP10 (Continuous-Scale Physical Functional Perfor- mance 10)	Nonamputees scored significantly higher in the total score (-24.4%) and lower-body strength (-27.6%), upper- body flexibility (-13.4%), balance (-27.1%) and endurance (-28.9%) do- mains but not significantly higher in the domain upper-body strength.	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "There were no significant differences in functional UBS between nonamputees and persons with TFA regardless of knee condition. Compared with the C-Leg, Genium use improved the UBF, BAL, and END domains of functional performance, likely because of improved confidence, willingness to lift and carry greater mass, and

ability to move faster during activity. These benefits may be technologically due to the incorporation of a faster processing speed and axial load data assisting in regulating knee resistance and offering new functions such as stance locking and backward stepping. In the LBS, UBF, BAL, and END domains, C-Leg use resulted in significantly lower scores compared with nonamputees. Genium use significantly reduced the magnitude of impairment. The only domain in which persons with TFA performed significantly lower than nonamputees regardless of knee condition was the END domain. In terms of total CS-PFP10 performance, C-Leg use resulted in significantly lower function compared with nonamputees, whereas Genium use was not significantly different from nonamputees. Nonetheless, regardless of knee condition, persons with TFA did not equal or surpass nonamputees in any functional performance domain, suggesting room for improvements in TFA integrated functional performance. Further, the CS-PFP10 test was able to detect statistically significant differences of small effect size between prosthetic knee conditions, which should be interpreted with caution because the test has not been formally assessed in persons with TFA." (Highsmith et al., 2016)

<u>A Back to overview table</u>

Reference	Highsmith MJ, Kahle JT, Wernke MM, Stephanie LC, Miro RM, Lura DJ, Sutton BS.						
	School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.						
	Effects of the Genium knee system on functional						
	level, stair ambulation, perceptive and economic						
	outcomes in transfemoral amputees						
	Technology and Innovation 2016; 18: 139-150						
Products	Genium vs C-Leg						
Major Findings	With Genium compared to C-Leg:						
	 → The quality of stair ascent and descent improved significantly → Mobility and functional level improved significantly → Perceived function and safety in ADLs was as good as with C-Leg or significantly better → 80% preferred Genium 						
	Preference for prosthetic knee joint						
	80%						
	چ 70%						
	tig 60%						
	រម្មវ័ជ្ ភក្តិ 50%						
	б де 40%						
	W 30%						

Population	Subjects:	20 unilateral, transfemoral amputees
	Previous prosthesis:	C-Leg
	Amputation causes:	70% trauma, 20% malignancy, 10% vascular disease
	Mean age:	$46.5 \pm 14.2 \text{ yrs}$
	Mean time since amputation:	17.7 yrs ± 15.6 yrs
	MFCL:	K3

■Genium ■C-Leg

60% 50% 40% 30% 20% 10% 0%

Interventional, randomized crossover design:



Results

Functions and Activities Level Stairs Ramps, Uneven Cognitive					Participation			Environment	
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Outcomes	Results for Genium	Sig.*
SAI (Stair Assessment Index)	SAI stair <i>ascent</i> score improved signifi- cantly from 6 to 11 points (median).	++
	SAI stair <i>descent</i> median score was 11 for both knee joints. The mean score signifi- cantly improved by 9%.	++
	Stair ascent and descent completion time did not differ significantly.	0
Four Square Step Test	Time to complete the test significantly de- creased by 9%.	
AMP (Amputee mobility predictor)	Mobility increased significantly by 5%.	++
Step activity derived functional level	The functional level significantly increased by 6%.	++
ADL survey	Perceived function and safety in three of five ADL domains improved significantly. The other two domains showed no difference.	++
Survey	80% of the participants preferred Genium.	++
ICER (Incremental cost- effectiveness ratio)	The ICER for reimbursing Genium ranges from \$6,000 to \$6,522 per unit of functional increase assuming a \$30,000 intervention cost difference.	n.a.
	Outcomes SAI (Stair Assessment Index) Four Square Step Test AMP (Amputee mobility predictor) Step activity derived functional level ADL survey Survey ICER (Incremental cost- effectiveness ratio)	OutcomesResults for GeniumSAI (Stair Assessment Index)SAI stair ascent score improved signifi- cantly from 6 to 11 points (median).SAI stair descent median score was 11 for both knee joints. The mean score signifi- cantly improved by 9%.Stair ascent and descent completion time did not differ significantly.Four Square Step TestTime to complete the test significantly de- creased by 9%.AMP (Amputee mobility predictor)Mobility increased significantly by 5%.Step activity derived functional levelThe functional level significantly increased by 6%.ADL surveyPerceived function and safety in three of five ADL domains improved significantly. The other two domains showed no difference.SurveyICER (Incremental cost- effectiveness ratio)The ICER for reimbursing Genium ranges from \$6,000 to \$6,522 per unit of functional increase assuming a \$30,000 intervention cost differ- ence.

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion	"Accommodation and use of the Genium knee system comp proved stair walking performance, multi-directional stepping perceived function. Genium was also preferred compared to high functioning community ambulators with unilateral trans Finally, Genium is a more costly microprocessor knee syste patients, is worth funding due to significant differences in fu	pared with C-Leg im- g, functional level, and o C-Leg in this group of femoral amputation. m but, in this group of unctional performance
	with activities of daily living." (Highsmith et al., 2016).	Back to overview table

Reference	Highsmith MJ, Klenow TD, Kahle JT, Wernke MM, Carey SL, Miro RM, Lura DJ.						
	School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA. Effects of the Genium microprocessor knee system on knee moment symmetry during hill walking.						
	Technology and Innovation 2016; 18: 151-157.						
Products	Genium vs C-Leg						
Major Findings	With Genium compared to C-Leg:						
	→ The degree of asymmetry for knee moment during hill walking is lower with Genium						
	→ The degree of asymmetry for knee moment d with Genium						
	The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment during the selection of th	g speed ected walking speed ng hill walking at					
	The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durir self-slected walking s	g speed ected walking speed ng hill walking at speed					
	The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0.07	g speed ected walking speed ng hill walking at speed					
	→ The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durir self-slected walking s 0,08 0,07 E 0,06	g speed ected walking speed ng hill walking at speed					
	→ The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0,07 0,06 0,05	g speed ected walking speed ng hill walking at speed					
	 → The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0,07 0,06 0,05 0,04 	g speed ected walking speed ng hill walking at speed					
	→ The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durir self-slected walking s 0,08 0,08 0,07 0,06 0,06 0,05 0,04 0,03	g speed ected walking speed ng hill walking at speed					
	→ The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0,07 0,07 0,06 0,05 0,03 0,02	g speed ected walking speed ng hill walking at speed Genium					
	 The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0,07 0,06 0,05 0,04 0,03 0,02 0,01 	g speed ected walking speed ng hill walking at speed Genium					
	The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durir self-slected walking s 0,08 0,07 0,06 0,05 0,04 0,02 0,01 0	g speed ected walking speed ng hill walking at speed Genium C-Leg					
	The degree of asymmetry for knee moment d with Genium Uphill: significantly lower at slow and fast walking Downhill: significantly lower at slow and self-sele Asymmetry of knee moment durin self-slected walking s 0,08 0,07 0,06 0,05 0,04 0,05 0,03 0,02 0,01	g speed ected walking speed ng hill walking at speed Genium C-Leg					

A positive value indicates a greater knee moment on the sound side when ascending ramps, a negative value a greater knee moment on the prosthetic side and a value of zero perfect symmetry.

Population	Subjects: Previous prosthesis: Amputation causes: Mean age: Mean time since amputation: MFCL:	20 unilateral, transfemoral amputees C-Leg 70% trauma, 20% malignancy, 10% vascular disease 46.5 ± 14.2 yrs 17.7 ± 15.6 yrs K3
	WI CE.	NO

Interventional, randomized crossover design:



Results

Functions and Activities				Participation			Environment		
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Category	Outcomes	Results for Genium	Sig.*
Ramps, Hills	Motion analysis Ramp ascent	The degree of asymmetry for knee moment was significantly lower at slow and fast walking speed and not significantly different at self-selected walking speed.	
	Motion analysis Ramp descent	The degree of asymmetry for knee moment was significantly lower at slow and self- selected walking speed and not significantly different at fast walking speed.	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

"Accommodation and use of the Genium knee system, compared with C-Leg, improved knee moment symmetry in slow speed walking up and down a five degree ramp. Additionally, the Genium improved knee moment symmetry when walking downhill at comfortable speed. At fast walking speed, variance in knee moment symmetry was lower when using Genium. These results were found in a sample of high functioning persons with unilateral transfemoral amputation; however, the results likely have application in other patients who could benefit from more consistent knee function, such as older patients and others who have slower walking speeds." (Highsmith et al., 2016)

<u>A Back to overview table</u>

Reference	Huppert L, Mileusnic M, Hahn A.
	Otto Bock Healthcare Products GmbH; Brehmstraße 16, 1110 Vienna.
	Das Genium-Prothesenkniegelenk – ein
	Überblick über die wissenschaftliche Evidenz
	(Genium prosthetic knee joint–Overview of scientific evidence)
	Orthopädie Technik 2016. 4: 44-49.
Products	Genium vs C-Leg
Major Findings	With Genium compared to C-Leg:
	→ Level walking:
	→ More physiological walking due to increased knee flexion angle during
	standing and swing phase while walking on level ground and ramps.
	→ Reduced impact forces through 4° "Prefiex" at initial neel contact
	\rightarrow Reduction of step length gait asymmetry while level walking by 40-60%.
	→ Correct swing initiation of 95% of the subjects when walking with small steps with Genium instead of 75% with C-Leg.
	→ Safe detection of walking backwards and therefore reliable blocking of the swing phase release.
	→ Improved stair ambulation:
	→ 70-80% of the patients used step-over-step strategy for stair ascent.
	→ Range of motion (ROM) of the hip and knee joint of the sound side was
	reduced by one third and is nearly equivalent to able bodied persons.
	→ More balanced and safer standing on ramps.
	→ Activities of daily living (ADLs) showed a clinically relevant decrease in perceived difficulty (53% of ADLS) and gain in safety (60% of ADLS).
	→ Quality of life (QoL) is significantly improved including 4 out of 9 scales of Prosthetic Evaluation Questionnaire (PEQ)
	Genium makes activities of daily living safer
	Sports
	Moving around in a confined place
	Stepping in a sidewalk curb
	Stepping over minor obstacles
	Stepping on minor obstacles like
	Walking up stairs
	Walking down stairs
	Walking up ramps 0 = no difference
	Walking down ramps 1 = safer 2 = much safer
	Walking with different speeds
	Chasing a child
	Walking in a crowd
	Walking on uneven terrain
	Walking a dog
	0 0,5 1 1,5 2
	Safer with Genium
	Kannenberg et al., 2013.

Study Design

Systematic review:

Nine publications were identified comparing Genium to C-Leg with each including on average between 10 and 20 transfemoral participants. The following table lists topics that were reviewed in this overview including the number of supporting studies:

Results	Number of studies
Level walking	4
Stairs	3
Ramps, Hills	3
Safety	4
ADLs	1
Quality of Life (QoL)	1
Health	4

Results

Functions and Activities				Participation			Environment		
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Category	Outcomes	Results for Genium vs C-Leg	References
Level walking	Forces at initial heel contact	The 4° "Preflex" at initial heel contact reduces impact forces, thus protecting the body.	[1,2]
	Knee flexion during standing phase	Increased knee flexion angle during stance phase up to 2° with Genium while walking very slow, slow or fast.	[5]
	Adaptive swing phase control	Maximum knee flexion angle was 64°, which ensures toe clearance at different gait velocities.	[1,2]
	With Genium, the knee flexion increased signifi- cantly at very slow, slow and fast walking speed compared to C-Leg. These angles are nearly equivalent to those of able bodied persons.		[5]
		Adding more weight on the prosthetic foot (like heavy shoes), led to higher knee flexion angles.	[3]
		At 95% of small steps, swing was initiated cor- rectly through adaptive swing phase control of Genium. With C-Leg the percentage was only 75%.	[1,2]
	Asymmetry of step length	Asymmetry of step length was reduced by 40-60%, depending on gait velocity.	[1,2]
Stairs	Stair ascent strategy	70-80% of the patients could use step-over-step strategy to ascent stairs with Genium.	[4,6,7]
	Range of motion (ROM)	Compensations in terms of ROM of the hip and knee joint on the sound side were reduced by about one third, which is nearly equivalent to an able bodied person.	[6]
Ramps, Hills	Maximum knee flexion during stance phase	During ramp descent at slow and fast walking speed the knee flexion angle increased signifi- cantly with Genium.	[5]

Category	Outcomes	Results for Genium vs C-Leg	References
	Maximum knee during swing phase	7° to 8° higher knee flexion angle during ramp ascent and descent with Genium compared to C-Leg.	[1,2,5]
	Standing on a 10 degree ramp for 3 minutes	 Higher loading of the affected side up to 86%. Sagittal knee flexion moment on the prosthetic side increased by 92%. Reduction of postural sway of the prosthetic side. 	[1,2]
Safety	Stumbles and falls	 The risk for stumbling or falling can be reduced through: Better toe clearance through higher knee flexion. [1,2,5] Initiation of the swing phase, while making small steps. [1,2] Walking backwards detection, thus blocking the swing phase release. [2] 	[1,2,5]
	ADL questionnaire	60% of ADLs showed a clinically relevant gain in safety.	[8]
Activity, Mobility, Activities of Daily Living (ADLs)	ADL questionnaire	53% of ADLs showed a clinically relevant decrease in perceived difficulty. Especially ascending and descending stairs and ramps as well as walking backwards improved significantly.	[8]
Preference, Satisfaction, Quality of Life (QoL)	Prosthetic Evaluation Questionnaire (PEQ)	 4 out of 9 scales were rated significantly higher: Perceived Response Social Burden Utility Well-being 	[9]
		"Appearance" and "Sounds" had the tendency to be rated higher, but not significantly.	
		 3 out of 9 scales were unchanged: Ambulation Frustration Residual Limb Health 	
Author's Conclusion	"Erste wissenschaf um biomechanisch cherheit und Einfac sowie den Einfluss mit dem Genium se werden können. Vo	ttliche Studien, welche die Leistungsversprechung , funktionell und hinsichtlich des subjektiven Zug chheit in der Durchführung von Aktivitäten des täg auf die Lebensqualität überprüften, liefern Hinweis elbst im Vergleich zum C-Leg weitere Gebrauchsvo on Bedeutung sind diese Ergebnisse insbesonde	gen des Geni- jewinns an Si- jlichen Lebens se darauf dass orteile realisiert re im Hinblick

<u>A</u> Back to overview table</u>

apparates." (Huppert, 2016)

auf die Sicherheit des Anwenders sowie deren Schutz ihres gesamten Bewegungs-

Bell EM, Pruziner AL, Wilken JM, Wolf EJ.

Department of Rehabilitation, Walter Reed National Military Medical Center, Bethesda, MD, 20889, USA.

Performance of conventional and X2[®] prosthetic knees during slope descent.

Clin Biomech, 2016; 33: 26-31.

Products	Genium X2®					
Major Findings	 With Genium X2[®] (X2) compared to Mechanical knee (MECH) (Total Knee & Mauch, Össur) and Standard microprocessor controlled knee (MPK) (C-Leg, Otto Bock & Rheo Knee, Össur): → Walking speed is faster by 9.6% (+0.1 m/s) compared to MPK. 					
	→ Step length is longer with X2 by 1.6-16.9% and leads to a more harmonized movement in slope descent.					
	→ Through the stance flexion resistance, walking with the X2 was more com- parable to able-bodied individuals.					
	Initial knee flexion (0% GC): 326.7% more flexion Max knee swing flexion (50-100% GC): 21.2% more flexion					
	→ A heightened use of the intact limb for support in descent could be indi- cated by the significant increase of the max. support moment flexion with X2 compared to MPK:					
	<u>0% GC</u> : 26.3% higher <u>35-75% GC</u> : 21.2% higher					
	→ With X2, the prosthetic limb was utilized and loaded more normative. Ther- fore, the first vertical impact maximum (0-30%) increased up to 13.2%.					

<u>Please note</u>: The percentage differences were calculated between the published Median values.



Population	Subjects:	amputees	
	Previous prosthetic knee:	Mechanical knee (n=8)	→ Total Knee (Össur) → Mauch (Össur)
		Standard MPK (n=13)	→ C-Leg (Otto Bock) → Rheo Knee (Össur)
	Amputation causes:	Trauma	
	Mean age:	32.7 yrs (± 5.3 yrs)	
	Time since amputation:	≥ 2 years	
	MFCL:	K4	

Study Design

Interventional, pre- to post-test design:



Use of handrails influenced what data were available for biomechanical analysis. As such, data from participants who self-selected to use handrails for support were not included in temporal-spatial, kinematic, or kinetic analyses.

	Temporal-spatial	<u>Kinematic</u>	<u>Kinetic</u>
MECH	n=4	n=4	n=3
MPK	n=12	n=12	n=11

Results

Functions and Activities					Participation			Environment	
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Category	Outcomes	Results for Genium X2				
Ramps, Hills	Temporal-spatial	Walking speed is faster:				
		X2 vs. MPK	X2 vs. MECH			
		9.6% faster	8.2% faster			
		++	+			
		Step length is longer with X2:				
		X2 vs. MPK	X2 vs. MECH			
		16.9% longer	1.6% longer			
		++	+			
		No significant differences were found for stance time.				
	Kinematic	The initial knee flexion (0% GC) increased with X2 compared to MPK (significantly) and MECH:				
		X2 vs. MPK	X2 vs. MECH			
		326.7% higher	63.3% higher			
		++	+			

Outcomes	Results for Genium X2				
	The max. knee swing flex creased:	ion (50-100% GC) in-			
	X2 vs. MPK 21.2% higher	X2 vs. MECH 8.83% higher			
	 No significant differences were found for: Ankle excursion (0-100% GC) Hip excursion (0-100% GC) 				
Kinetic	The max. support momen with X2 compared to the	t (0-30% GC) increased previous knee:			
	X2 vs. MPK 26.3% higher ++	X2 vs. MECH 240.6% higher +			
	The max. support moment (35-75% GC) increased:				
	X2 vs. MPK 21.2% higher ++	X2 vs. MECH 147.5% higher +			
	Due to more prosthetic limb utilization and norm tive loading, the first vertical impact maximum (0-30%) increased:				
	X2 vs. MPK 13.2% higher	X2 vs. MECH 12.3% higher			
	++	+			
	 No significant differences were found for: Max. braking force (0-30%) Max. propulsive force (35-75%) 				
	Kinetic	OutcomesResults for Genium X2The max. knee swing flex creased:X2 vs. MPK 21.2% higher ++No significant differences w • Ankle excursion (0- • Hip excursion (0-10)KineticThe max. support momen with X2 compared to the X2 vs. MPK 26.3% higher ++X2 vs. MPK 26.3% higher ++The max. support momen X2 vs. MPK 21.2% higher ++Due to more prosthetic lin tive loading, the first verti (0-30%) increased:X2 vs. MPK 13.2% higher ++No significant differences w • Max. braking force • Max. propulsive for • Consert duration line			

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "The aim of the current study was to determine if use of the X2[®] improves overall slope descent mechanics by assessing self-selected technique of decent, and subsequent changes in temporal-spatial outcomes and joint mechanics. Although this analysis finds normalization of some temporal-spatial outcomes and joint mechanics were likely achieved due to the use of stance flexion resistance with the X2® device allowing for improved control lowering the body when both leading and trailing, some values continued to deviate from those of able-bodied individuals. Nevertheless, decreased reliance on handrail use as MECH users descended in the X2® suggest improved function and perhaps greater confidence in the device possibly reducing the risk of falling. Furthermore, overall reductions in intact limb loading and more symmetric loading at impact could indicate more normative loading patterns and a possible reduction of intact limb overuse during downslope walking. Reducing compensatory gait strategies during slope descent, perhaps through use of the X2®, could thus play a role in mitigating longer-term overuse injuries commonly associated with TFA." (Bell et al., 2016)

A Back to overview table

Reference Schalk SAF, Jonkergouw N, van der Mer F, Swaan WM, Aschoff HH, van der Wurff P.

Military Rehabilitation Centre Aardenburg, Korte Molenweg 3, 3941 PW, Doorn, The Netherlands.

The Evaluation of Daily Life Activities after Application of an Osseointegrated Prosthesis Fixation in a Bilateral Transfemoral Amputee: A Case Study

Medicine (Baltimore) 2015. 94(36): e1416.

Products Genium with osseointegrated prosthesis fixation (OPF) Major Findings With Genium with OPF compared to previous situation without a prosthesis: → The patient was able to perform four more activities (in total 11 out of 20) of Lower extremity functional scale (LEFS)

→ LIFE-H (life habits questionnaire) increased in daily activities and social roles in five and eight respectively out of 12 categories.



More activities possible with Genium with OPF

Population	Subjects:	One bilateral, transfemoral amputee (male)
	Previous:	No prosthesis, because several attempts to control
		the problems related to rotation and fixation of the
		socket failed.
	Amputation causes:	Trauma (injured by an improvised explosive device
		(IED) during a military mission)
	Stump length:	13 cm (right), 19 cm (left)
	Age:	21 years
	Time since amputation:	4 years

Study Design

Case study:



Results

Functions and Activities				Participation			Environment		
Level walking	Stairs	Ramps, Hils	Uneven ground, Obstacles	Cognitive demand	Energy	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health Economics

Category	Outcomes	Results for Genium with OPF				
Activity, Mobility, Activities of Daily Living (ADLs)	Lower extremity functional scale (LEFS)	The LEFS is a questionnaire containing 20 activities whic subject ranks according to difficulty (0 - extreme difficultie to 4 - no difficulties).				
		 The subject was capable of performing four more activities with Genium with OPF compared to no prosthesis. Any of your usual work, housework, Usual hobbies, recreational or sporting Getting into or out the bath Walking between rooms 	+			
		 Three activities were rated to be easier without prosthesis (LEFS-score = 4) than with Genium (LEFS-score = 3) Lifting an object (e.g. bag) from the floor Performing light activities around your home Performing heavy activities around your home As statement about this trend, the patient admitted that he may have overestimated the expected capabilities of OPF. 				
		No differences between Genium with OPF and no prosthesis for four activities that are capable: • Getting into or out of a car • Walking two blocks • Sitting for one hour • Rolling over in bed	0			
Preference, Satisfaction, Quality of Life (QoL)	LIFE-H (life habits ques- tionnaire)	This LIFE-H questionnaire is divided in two domains activities" and "social role") and 12 categories (where education was not applicable for this patient).	("daily eby			
		Daily activities: LIFE-H increased with Genium with OPF in five categories (recreation, community life, mobility, housing and fitness).	+			

Outcomes	Results for Genium with OPF	Sig.*
	LIFE-H decreased with Genium with OPF in three categories (work, relationships, personal care). This could be a result of the social change in 2013, when the patient moved from his parent's house to his own place.	-
	Social roles: LIFE-H increased with Genium with OPF in eight categories (recreation, community life, mobility, housing, communication, personal care, fitness and nutrition).	+
	The patient rated the situation without prosthesis as being better than with Genium in two categories (work, relationships). This could again be a result of the social change in 2013, when the patient moved from his parent's house to his own place.	-
	Outcomes	OutcomesResults for Genium with OPFLIFE-H decreased with Genium with OPF in three categories (work, relationships, personal care). This could be a result of the social change in 2013, when the patient moved from his parent's house to his own place.Social roles: LIFE-H increased with Genium with OPF in eight categories (recreation, community life, mobility, housing, communication, personal care, fitness and nutrition).The patient rated the situation without prosthesis as being better than with Genium in two categories (work, relationships). This could again be a result of the social change in 2013, when the patient moved from his parent's house to his own place.

Author's Conclusion

"In conclusion, although several studies clearly demonstrate the benefits of using an OPF, studies describing the long-term effects are lacking. In this specific case we conclude that the quality of life improved through the use of an OPF. However, OPF might not be the appropriate device for every individual with TFA, due to varying bone compositions and co-morbidities." (Schalk, 2015)

↑ Back to overview table

Reference	Lura DJ, Wernke MM, Carey SL, Kahle JT, Miro RM, Highsmith MJ.							
	School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.							
	Differences in knee flexion between the Genium							
	and C-Leg microprocessor knees while walking							
	on level ground and ramps							
	Clinical Biomechanics 2015; 30(2): 175-181.							
Products	Genium vs C-Leg							
Major Findings	With Genium compared to C-Leg:							
	 → More physiological movement pattern → Increased knee flexion angle in both stance and swing phase due to accommodation, training and use of Genium during level walking and walking on ramps by up to 7° during swing phase when level walking by up to 2° during stance phase when level walking by up to 8° during swing phase when descending a ramp by up to 4° during stance phase when descending a ramp by up to 9° during swing phase when ascending a ramp 							
	Increased swing phase knee flexion angle during 5° ramp							
	ascent							
	75							
	65							
	5 55							
	- Jest							
	genium							
	35							
	25							
	slow normal fast							

Peak knee flexion angle was measured at slow, normal and fast walking speed.

Population	Subjects:	20 unilateral, transfemoral amputees
	Previous prosthesis:	C-Leg
	Amputation causes:	70% trauma, 20% malignancy, 10% vascular
		disease
	Mean age:	46.5 yrs (±14.2 yrs)
	Mean time since amputation:	17.7 yrs (±15.6 yrs)
	MFCL:	КЗ
	Mean age: Mean time since amputation: MFCL:	disease 46.5 yrs (±14.2 yrs) 17.7 yrs (±15.6 yrs) K3

Interventional, randomized crossover design:



Results

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results for Genium compared to C-Leg			
		Very slow	Slow	Normal	Fast
Level Walking	Peak knee flexion angle in swing phase	++ 7 °	++ 6°	++ 6 °	+
	Peak knee flexion angle in stance phase	+	++ 2 °	++ 2 °	+
	0.5 kg ankle weight at- tached to each leg Peak knee flexion angle in swing phase	++ 6°	++ 5°	++ 4°	++ 3°
	0.5 kg ankle weight at- tached to each leg Peak knee flexion angle in stance phase	++ 1 °	+	+	++ 2°
Ramps, Hills	5° slope Descent Peak knee flexion angle in swing phase		++ 8°	++ 7°	++ 3°
	5° slope Ascent Peak knee flexion angle in swing phase		+	++ 9°	-
	5° slope Descent Peak knee flexion angle in stance phase		++ 3°	+	++ 4°
	5° slope Ascent Peak knee flexion angle in stance phase		+	+	+

Category	Outcomes	Results for Genium compared to C-Leg				
		Very slow	Slow	Normal	Fast	
-	10° slope		+	+	+	
	Descent					
	Peak knee flexion angle					
	in swing phase					
	10° slope		+	++ 8°	0	
	Ascent					
	Peak knee flexion angle					
	in swing phase					
	10° slope		+	+	+	
	Descent					
	Peak knee flexion angle					
	in stance phase					
	10° slope		+	+	+	
	Ascent					
	Peak knee flexion angle					
	in stance phase					

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "Accommodation, training and use of the Genium were found to produce increased knee flexion compared with the C-Leg in both the stance and swing phases of gait. This increased knee flexion is clinically significant as it better recreates a normalized, anatomic movement pattern. The knee flexion angle of the non-amputated leg was not significantly affected by use of the Genium relative to the C-leg. Control subjects typically had the greatest knee flexion, followed by the amputees' sound side, and then prosthetic side of the subjects with the Genium and C-Leg respectively. This shows that Genium use increases stance and swing knee flexion angles compared with the C-Leg, but improvements are still possible, especially in certain walking conditions such as when walking uphill." (Lura et al., 2014)

↑ Back to overview table

Reference	Aldridge Whitehead JM, Wolf EJ, Scoville CR, Wilken JM.						
	Extremity Trauma and Amputation Center of Excellence, Department of Orthopae- dics and Rehabilitation, Brooke Army Medical Center, Ft Sam Houston, TX, USA.						
	Does a Microprocessor-controlled Prosthetic Knee Affect Stair Ascent Strategies in Persons with Transfemoral Amputation?						
	Clinical Orthopaedics and Related Research 2014; 472(10):3093-3101						
Products	X2 vs C-Leg and Total Knee						
Major Findings	With X2 compared to C-Leg:						
	 strategy → Greater range of motion (ROM) on prosthetic side throughout gait cycle of stair ascent: 62° more knee flexion during swing phase 30° more knee flexion at initial contact 22° more hip flexion during swing phase 						
	Improved stair ascent strategy with X2						
	10 8 8 E-Leg and						
	G Total Knee See 4 See See See See See See See See Se						
	0 Prosthesis						

Population

Subjects:14 transfemoral amputeesPrevious prosthesis:12 C-Leg (MPCK), 2 Total Knee (NMPCK)Amputation causes:not reportedMean age:31.1 yrsMean time since amputation:≥ 6 monthsMFCL:K3 – K4 (independent ambulators)

Interventional, pre- to post-test design:



Results	lts
---------	-----

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results for X2 compared to C-Leg and Total Knee				
Stairs	Self-selected stair ascent		Step-to	Skip-step	Step-over	
	strategy	X2	14.25%	14.25%	71.5%	n.a.
		C-Leg and Total Knee	71.5%	21.5%	7%	n.a.
	Stair Assessment Index (SAI)	SAI score improved from 5 to 11 points (medi- an).				n.a.
		While using ascent 60% 30% used the handrai	g X2 during 6 of subjec one handra Is.	step-over-sto ts used two h ail and 10% c	ep stairs andrails, lid not use	n.a.
		for X2 (stair function) compared to C-Leg (no stair function)				Sig.*
	Motion Analysis	62° more prosthetic knee flexion during swing phase.				++
		22° more p swing pha	prosthetic se.	limb hip flex	tion during	++
		30° more p ground co	prosthetic ntact.	knee flexion	at initial	++
		Prosthetic during pus	limb peal sh-up was	k hip power increased b	generated y 143%.	

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "Stair ascent can be difficult for individuals with transfemoral amputation because of the loss of knee function. Most individuals with transfemoral amputation use either a step-to- step or skip-step strategy because it allows the leading intact limb to do the majority of work. A new microprocessor-controlled knee (X2[®]) uses flexion/extension resistance to allow step-over-step stair ascent. We compared stair ascent strategies and joint mechanics as individuals with transfemoral amputation ascended stairs using their conventional prosthetic knee and the novel X2[®] device. Most participants self-selected a step-to-step stair ascent strategy while using their conventional device and a step-over-step strategy while using the X2[®] device. Participants were

more symmetrical while using the X2[®] than the conventional device to include more similar peak knee and hip flexion during swing and peak hip power generation during push-up when comparing between limbs. Although the X2[®] resulted in greater prosthetic knee ROM and fewer between-limb differences than the conventional knee, stair ascent gait deviations still persisted compared to individuals without amputation. Peak knee flexion during swing while using the X2[®] device was the only prosthetic limb measure that was not different from individuals without amputation. Correlational analysis revealed that greater X2[®] knee flexion during initial contact and swing was associated with greater prosthetic limb hip power during pull-up and push-up/early swing, respectfully [sic]." (Aldridge Whitehead et al., 2014)

<u>A Back to overview table</u>

Reference	Highsmith MJ, Kahle JT, Miro RM, Lura DJ, Dubey RV, Carey SL, Quillen WS, Mengelkoch LJ. School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.						
	Perceived differer	nces between th	e Genium und				
	the C-Leg microp	rocessor prosthe	etic knees in				
	prosthetic-related	l function and q	uality of life				
	Technology and Innovation 20	014; 15(4):269-375.					
Products	Genium vs C-Leg						
Major Findings	Genium compared to C-Leg:						
	 → Prosthesis-related quality of life is improved → 45% of prosthetic-related scales are improved → 70% of activities, relevant to the physical performance aspects of Genium, are improved 						
	Percent of scales showing i and quality of life with Genin	mprovement in prosthetic um	function				
	33%	45%	■ higher with Genium				
			tend to be higher with Genium				
	229	%	■ unchanged				
Population	Subjects: Previous prosthesis:	20 unilateral, transfemora	al amputees				
	Amputation causes:	70% trauma, 20% malig	nancy, 10% vascular				
	Mean age:	disease 46.5 yrs (±14.2 yrs)					
	Mean time since amputation:	17.7 yrs (±15.6 yrs)					
	MFCL:	K3 - K4 (ambulate withou community)	It assistive device within				

Interventional, randomized crossover design:



Results

Activities							Participation	Environment
Level walking	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results for Genium compared to C-Leg	Sig.*
Preference,	Prosthesis Evaluation	Comparing all 41 items:	
Satisfaction,	Questionnaire (PEQ)	Genium was rated to be improved	++
Quality of Life (QoL)			
		4 out of 9 scales were rated higher:	
		Perceived Response	++
		Social Burden	++
		Utility	++
		Well-Being	++
		2 out of 9 scales had the tendency to be rated	
		higher:	
		Appearance	+
		Sounds	+
		3 out of 9 scales were unchanged:	
		Ambulation	0
		Frustration	0
		Residual Limb Health	0
		7 out of 10 chosen items, relevant to the physi-	
		cal performance aspects of Genium, were	
		improved:	
		Comfort standing	++
		Ability to walk in close spaces	++
		Go down stairs	++
		Walk up steep hill	++
		Walk down steep hill	++
		Walk on slippery surfaces	++
		Satisfaction with walking	++

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

"In this sample, several PEQ scales were rated as being significantly improved following Genium use. These included Perceived Response, Social Burden, Utility, and Well-Being. Movement analysis data from other studies may explain improvements in the Perceived Response scale. Further, the Genium was rated as improving specific mobility activities including walking on stairs, slopes, and slippery surfaces as well as increasing comfort during standing. The entire PEQ, which represents patient-perceived prosthetic-related function and quality of life, was rated as a significant improvement when participants utilized the Genium MPK. This study presents additional data to make future comparisons in persons with lower limb amputation who are assessed utilizing the PEQ in its ordinally scaled format. Further research is needed to more emphatically confirm the perceived improvements with objective functional and movement analyses. Finally, the PEQ does not assess perceived differences in nonambulatory functional tasks, preference, or safety." (Highsmith et al., 2014)

A Back to overview table

Reference Highsmith MJ, Kahle JT, Lura DJ, Dubey RV, Carey SL, Quillen WS, Mengelkoch LJ.

School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.

Short and mid-distance walking and posturography with a novel microprocessor knee

Technology and Innovation 2014; 15(4):259-368.

Products	Genium vs C-Leg
Major Findings	With Genium compared to C-Leg:
	 → Movement control in all three backward directions improved by up to 10% → Walking velocities during short- and mid-distance are maintained whereas levels of perceived exertion tend to decrease

Improved limits of stability with Genium



Limits of stability (LOS) was measured by Biodex Balance SD system.

Population	Subjects:	20 unilateral, transfemoral amputees
	Previous prosthesis:	C-Leg
	Amputation causes:	70% trauma, 20% malignancy, 10% vascular
		disease
	Mean age:	46.5 yrs (± 14.2 yrs)
	Mean time since amputation:	17.7 yrs (± 15.6 yrs)
	MFCL:	K3 - K4 (ambulate without assistive device within community)

Interventional, randomized crossover design:



Results

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results for Genium compared to C-Leg	Sig.*
Level Walking	75 m self-selected walk- ing speed (SSWS)	Walking velocity increased by 2% (1.17 m/s vs 1.15 m/s).	+
		Perceived exertion was rated lower.	+
	75 m fastest possible walking speed (FPWS)	No difference in walking velocity. Perceived exertion was rated lower.	0 +
	6 m fastest possible walking speed (FPWS)	No difference in walking velocity. Perceived exertion was rated lower.	0 +
Uneven Ground,	38 m fastest possible	No difference in walking velocity.	0
Obstacle Course	walking speed (FPWS), sloping terrain over trimmed grass, sand, rocks, and small roots	Perceived exertion was rated lower.	+
Safety	Postural stability and limits of stability (LOS) both measured by Biodex	No differences in postural stability.	0
	Balance SD system	Improved by 10% in backward sound side direction	++
		Improved by 10% in backward direction	+
		Improved by 9% in backward amputated side direction	+
		Decreased by 6% in forward sound side direction	-
		Improved by 8% in forward direction	+
		Decreased by 12% in forward amputated side direction	
		Time to complete LOS test tended to be decreased by 2%.	+

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion

'During short to mid-distances, the Genium knee sustains the walking speed improvements realized by the C-Leg with a trend toward decreased levels of perceived exertion. In terms of directional control, TFAs demonstrate multidirectional impairment compared with nonamputees. However, C-Leg use results in improved anterolateral directional control compared with Genium, possibly due to the toe load requirement needed to initiate swing phase knee flexion. Conversely, Genium use results in a trend of movement control improvements in all three rearward directions compared with the C-Leg., (Highsmith et al. 2014)

<u>A Back to overview table</u>

Highsmith MJ, Kahle JT, Lura DJ, Lewandowski AJ, Quillen WS, Kim HS.

School of Physical Therapy & Rehabilitation Sciences, University of South Florida, Tampa, FL, USA.

Stair ascent and ramp gait training with the Genium knee

Technology and Innovation 2014; 15(4):349-258.

Products	Genium
Major Findings	With Genium:
	→ 70% of subjects are able to climb stairs reciprocally → Ability to cross obstacles improves when subject are able to climb stairs reciprocally
	→ Focal pressure near the anterior aspect of the hip during ramp ascent is experienced as decreased compared to other prosthesis

Stair climbing technique with Genium



Population	Subjects: Previous prosthesis: Amputation causes: Mean age: Mean time since amputation:	20 transfemoral amputees not reported trauma not reported not reported
	MFCL:	K3 – K4 (unlimited community ambulators)

Study Design

Technical report about stair ascent and ramp gait training with Genium

Results									
Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics
Category	1	O	utcomes		Results fo	or Genium	1		Sig.*
Stairs		Observations			70% of subjects demonstrated ability to climb stairs reciprocally.				n.a.
Ramps, Hills		Observations			Ascent: less focal pressure was experienced by subjects near the anterior aspect of the hip.				n.a.
	Descent: On 5° declines, more active walkers utilized a stepping strategy similar to level walk- ing with two sagittal knee flexion peaks.				n.a. k-				
Uneven Ground, Subject's report Obstacle Course			Subjects believe that stair climbing practice improved the ability to cross obstacles				n.a.		
* no differ	rence (0), p	ositive trend	l (+), negativ	e trend (–)	, significant	(++/), n	ot applicable	e (n.a.)	

Author's Conclusion "Technological developments in assistive devices continue to outpace rehabilitation strategies to maximize their utilization and implementation. Clinical rehabilitation techniques remain limited. This technical note presents strategies for training the patient with transfemoral amputation in how to utilize the reciprocal stair ascent and ramp gait functions of the Genium knee. Additional training suggestions for further advanced training with these skills are also discussed. Functional training strategies introduced here were specifically used with the Genium knee in high-functioning patients. Therefore, they may not be appropriate for all patients with transfemoral amputation based on component or functional level. Thus, clinical judgment and patient goals are vital in the decision of whether or not to include such training in the course of an amputee's therapy. We maintain that ramp and stair training in a broader context may be functionally important even if a patient indicates these obstacles are not often encountered in their usual routines. This training is important because it is difficult to determine when daily activities require out-of-the-ordinary settings or present unanticipated challenges. Supervised practice and familiarity may improve safety by decreasing fall risk, should the situation arise." (Highsmith et al., 2014)

↑ Back to overview table

Otto Bock HealthCare GmbH, Duderstadt, Germany.

Activities of Daily Living: Genium Bionic Prosthetic Knee Compared with C-Leg

Journal of Prosthetics & Orthotics 2013; 25(3):110-117.

Products	Genium vs C-Leg							
Major Findings	With Genium compared to C-Leg:							
	→ 60% of Activities of Daily Living (ADLs) show a cl safety → 53% of ADLs show a clinically relevant decrease	inically relevant gain in in difficulty.						
	→ Ease of ascending and descending stairs is improved	oved by 34% and 10%						
	\rightarrow Ease of ascending and descending ramps is imp	roved by 24% and 17%						
	\rightarrow Ease of walking backwards is improved by 26%							
	Clinically relevant improvement in safety of activities living with Genium	of daily						
	40%	Safer with Genium						
		above clinically relevant						
		threshold						
	60%	threshold Safer with Genium, below clinically relevant threshold						
	60%	threshold Safer with Genium, below clinically relevant threshold						
	60% Safety was assessed for 45 activities.	threshold Safer with Genium, below clinically relevant threshold						
Population	60% Safety was assessed for 45 activities. Subjects: 10 unilateral, transfemora Previous prosthesis: C-Leg Amputation causes: trauma	threshold Safer with Genium, below clinically relevant threshold 						
Population	60% Safety was assessed for 45 activities. Subjects: 10 unilateral, transfemora Previous prosthesis: C-Leg Amputation causes: trauma Mean age: 36.7 yrs (±10.2 yrs)	threshold Safer with Genium, below clinically relevant threshold						
Population	60% Safety was assessed for 45 activities. Subjects: 10 unilateral, transfemoral Previous prosthesis: C-Leg Amputation causes: trauma Mean age: 36.7 yrs (±10.2 yrs) Mean time since amputation: 12.5 yrs (±9.6 yrs) MFCL: 40% K3, 60% K4	threshold Safer with Genium, below clinically relevant threshold al amputees						
Population Study Design	60% Safety was assessed for 45 activities. Subjects: 10 unilateral, transfemoral Previous prosthesis: C-Leg Amputation causes: trauma Mean age: 36.7 yrs (±10.2 yrs) Mean time since amputation: 12.5 yrs (±9.6 yrs) MFCL: 40% K3, 60% K4 Interventional, pre- to post-test design:	threshold Safer with Genium, below clinically relevant threshold al amputees						
Population Study Design	60% Safety was assessed for 45 activities. Subjects: 10 unilateral, transfemora Previous prosthesis: C-Leg Amputation causes: trauma Mean age: 36.7 yrs (±10.2 yrs) Mean time since amputation: 12.5 yrs (±9.6 yrs) MFCL: 40% K3, 60% K4 Interventional, pre- to post-test design: C-Leg Genium	threshold Safer with Genium, below clinically relevant threshold al amputees						

Results										
Activities									Participation	Environment
Level walking	Stairs	Ramps Hills	Uneven ground, Obstacle	Cog den	nitive nand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics
Category			Outcomes			Results fo	r Genium	compared	o C-Leg	Sig.*
Level Walking			ADL Questionnaire, separate for each MPCK			Difficulty of creased by Walking at less difficul	of walking y 26%. varying spo It.	backwards	s was de- o be rated a	++ s +
Stairs			ADL Questi separate for	onnaire [·] each N	, IPCK	Difficulty of stairs was	of ascendi decrease	ng and des d by 34% a	cending nd 10%.	++
Ramps, Hills ADL Questionnaire, separate for each MPCK					, IPCK	Difficulty of slopes wa	of ascendi s decreas	ng and des ed by 24% a	cending and 17%.	++
Uneven Ground, ADL Questionnaire, Obstacle Course separate for each MPCK				, IPCK	Walking on tends to be	uneven ar rated as le	nd unfamiliar ess difficult.	ground	+	
Safety			ADL Questi comparative	onnaire	1	60% of AD gain in sat 40% of AD ly relevant of In the cate of ADLs sh safety. In the cate tion' 63% of vant gain i	PLs showe fety. DLs showed gain in safe egory 'Fam howed a c egory 'Mot of ADLs sl in safety.	d a clinical a trend tow ty. hily and Soc linically relation bility and Tr nowed a cli	y relevant ard a clinica cial Life' 83 evant gain i ansporta- nically rele	n.a. ו- א in
Activity Mo Activities o (ADLs)	bility, f daily living	1	ADL Questi comparative	onnaire	1	53% of AD decrease i 47% of AD cally releva In the cate of ADLs sh crease in p In the cate tion' 58% of vant decrease	PLs showed in perceive PLs showed the showed egory 'Fam howed a c perceived egory 'Mot of ADLs sh	d a clinical ed difficulty l a trend tow e in perceive nily and Soo linically rela difficulty. pility and Tr nowed a cli	ly relevant vards a clini- rards a clini- rards a clini- transporta- nically rele- riculty	n.a. %
* no differe	ence (0), po	sitive tre	end (+), nega	ative tre	nd (–),	significant ((++/), no	t applicable	(n.a.)	
Author's C	Conclusion		"This study	showed	l that th	e new techr	ological fu	nctions of th	ne Genium F	lionic Pros.

Author's Conclusion "This study showed that the new technological functions of the Genium Bionic Prosthetic Knee for transfemoral amputees of MFCL 3 and 4 lead not only to clear biomechanical benefits compared with the C-Leg but also to a further improvement in the subjective perception of safety and perceived difficulty of many ADLs." (Kannenberg et al., 2013)

A Back to overview table

Bellmann M, Schmalz T, Ludwigs E, Blumentritt S.

Department of Research, Otto Bock HealthCare, Duderstadt, Germany

Immediate effects of a new microprocessorcontrolled prosthetic knee joint: a comparative biomechanical evaluation

Archives of physical medicine and rehabilitation 2012; 93(3):541-549.

Products	Genium vs C-Leg
Major Findings	With Genium compared to C-Leg:
	 → Step length symmetry is improved by up to 60% → Knee flexion angle in swing phase is nearly constant at 63.8° across walking velocities → Loading of sound side is reduced by 10% when climbing stairs → Knee flexion is increased by 14% during ramp ascent and descent, leading to an increased toe clearance
	More symmetrical step length with Genium



Asymmetry of step length was determined by measuring difference in step length between prosthetic limb and sound limb when subjects walked at slow, medium and fast velocities.

Popula	tion
--------	------

Subjects:11 unilateral, transfemoral amputeesPrevious prosthesis:C-LegAmputation causes:91% trauma, 9% tumourMean age:37 yrs (± 10.2 yrs)Mean time since amputation:12.5 yrs (range from 3 – 34 yrs)MFCL:45% K3, 55% K4

Interventional, pre- to post-test design:



Results

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Category	Outcomes	Results for Genium	Slow v	Mid v	Fast v
Level Walking	Motion analysis	Mean walking velocity	0	0	0
		Results for Genium Slow v Mid v F Mean walking velocity 0 0 Asymmetry of step length decr ++ d	decr +		
		External knee flexion mo- ment at weight ac- ceptance	incr +	incr +	incr ++
		External hip flexion mo- ment at weight ac- ceptance	incr –	0	decr +
		Vertical ground reaction forces on prosthetic side at weight acceptance			
		 With stance flexion 	incr – –	decr +	decr +
		 Without stance flexion movement 	0	decr +	decr ++
		Horizontal ground reac- tion forces on prosthetic side at weight ac-			
		With stance flexion movement	decr +	decr ++	decr ++
		 Without stance flexion movement 	decr ++	decr +	decr ++
		Results for Genium			Sig*
		Maximum knee flexion ang constant at 63.8° with Gen ties. With C-Leg the mean angle increased by 14.6° v increased by 1 m/s.	n.a.		
	Motion analysis when walking with small	Mean maximum knee fle phase increased by 11%	++		
	steps	Genium failed to switch int 4.9% of steps, whereas C swing phase mode in 24.7	to swing pha -Leg failed t % of steps.	n.a.	

Category	Outcomes	Results for Genium	Slow v	Mid v	Fast v
Stairs	Motion analysis descending	Mean sagittal external l the prosthetic side at m 15°.	++		
	Motion analysis ascending	The mean duration of a er with the Genium (usin method) than with C-Le method). A clear approximation to t healthy subjects was obse	 n.a.		
		Range of motion (ROM) and prosthetic knee and with Genium (step-over ROM of the contralatera ventional method).	ralateral vas lower d) than C-Leg (con-	++	
		Loading of the contralar reduced by 10%.	int was	++	
		Mean loading of the contr slightly increased.	-		
Ramps, Hills	Motion analysis de- scending and ascend- ing a ramp (10° de-	Increased mean knee fl ascending (14% higher) higher) a ramp.	++		
	cline)	Sagittal external knee f effected side increased ing.	lexion mome by 11% whe	ents on the n descend-	++
		Decreased contralateral v forces at weight acceptan	+		
Safety	Motion analysis and sway measurements	The affected side was a loads.	86% higher	++	
	when standing on a ramp (10° decline)	Sagittal external knee f prosthetic side increase			
		Reduced moments of the 69%.	hip centre of	rotation by	+
		Reduced postural sway	of the prost	hetic side.	++

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "For the rehabilitation of lost lower limbs and the simulation of natural gait, increasingly complex technologies are required. The functions of the Genium knee, such as stance function, stairs function, and improved functionality for walking on level ground and ramps, provide further approximation to a more natural gait and physiologic loading pattern immediately after transition to this knee joint. This offers the user increased functionality and clinically relevant advantages in everyday situations like standing with an automatically lockable knee joint, walking on level ground at different velocities and step lengths, and walking up and down ramps and stairs. This increased functionality is represented in this study by direct comparison with the C-Leg, whose benefit has been demonstrated in numerous scientific studies. Additional investigations are indicated to evaluate if improvements in biomechanical parameters identified by this study can be confirmed after a longer period of adaptation to the Genium or if further advancements can be detected. It would also be of interest to evaluate if above-knee amputees with a lower activity level would benefit from the Genium functions." (Bellmann et al., 2012) <u>A Back to overview table</u>

Bellmann M, Schmalz T, Ludwigs E, Blumentritt S.

Research Biomechanics, Otto Bock HealthCare GmbH, Göttingen, Germany.

Stair ascent with an innovative microprocessorcontrolled exoprosthetic knee joint

Biomedizinische Technik 2012; 57(6):435-444.

Products	Genium vs C-Leg						
Major Findings	With Genium compared to C-Leg:						
	 → 80% of subjects were abl using Genium for one da → Movements of hip and kr and are within a normal r 	e to climb stairs with step-over-step y nee on the sound side are reduced by range	pattern after 34 and 33%				
	Decreased movement on s	ound side with Genium during stair					
	ascent						
	90						
	80						
	70						
	ି 60 ତୁ 60		■ C-Leg (step- by-step)				
	ociji al		■ Genium (step				
	40		over-step)				
	30						
	20 ——— hip	knee	_				
Population	Subjects:	10 unilateral, transfemoral amputees					
	Previous prosthesis:	C-Leg					
	Amputation causes:	90% trauma, 10% tumour					
	Mean time since amputation:	$\frac{30 \text{ yrs}}{10.4 \text{ yrs}} = \frac{1000 \text{ z}}{10.4 \text{ yrs}}$					
	MFCL:	40% K3, 60% K4					

Study Design

Interventional, pre-to post-test design:



Activities								Participation	Environmen
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics
Category		O	utcomes		Results fo	or Genium C-Leg (ste	ı (step-over p-by-step)	-step) com-	Sig.*
Stairs		M	otion analysis	S	80% of subjects were able to climb stairs with step-over-step pattern after one day using Ge- nium. With C-Leg subjects ascended stairs using the step-by-step strategy.				
						Duration of step cycle was increased by 31%.			
					On the sound side, movement of the hip was decreased by 34% and movement of the knee was decreased by 33%.				
					Decreased flexion angles of hip and knee on the prosthetic side when contact is made with the step.				
					Increased extension in hip on the prosthet- ic side at the end of the extension phase.				
					Joint powe by 14%.	er of the co	ontralateral k	nee decrease	d +

Author's Conclusion

"The study shows that a definite approximation of a natural gait and reduction of stress to the residual musculoskeletal system is achieved with the Genium prosthetic knee joint. Movements and motion sequences are nearly normal on the prosthetic side as well as the contralateral side, and the stress acting on the contralateral knee joint also tends to be reduced. The stump can help support the motion sequence and can operate in its original function as hip joint extensor in this gait situation. These improvements are apparent after just a short adaptation phase to this function." (Bellmann et al., 2012)

↑ Back to overview table

Blumentritt S, Bellmann M, Ludwigs E, Schmalz T.

Otto Bock HealthCare, Research Department, Duderstadt, Germany.

Zur Biomechanik des mikroprozessorgesteuerten Prothesenkniegelenks Genium

(Biomechanics of the microprocessor controlled knee prosthesis Genium)

Orthopädie-Technik 2012; 01:24-35.

Products	Genium vs C-Leg							
Major Findings	 With Genium compared to C-Leg: Improved toe clearance based on a constant knee angle of 64° prevents stumbling and therefore the risk of falling is reduced The safety potential when walking backwards as well as when walking with small steps is improved Improved toe clearance during ramp ascent and descent Knee angle in swing phase increased by 14% during ramp ascent 							
	75							
	70							
	E 65							
	55 Geniu							
	50							
	Gait velocity [m/s]							
Population	Subjects:11 unilateral, transfemoral amputeesPrevious prosthesis:C-LegAmputation causes:not reportedMean age:36.7 yrs (± 10.2 yrs)Mean time since amputation:3 – 34 yrsMECL:K2 – K4							

Interventional, pre- to post-test design:



Results

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics

Outcomes	Results for Genium compared to C-Leg		
Motion analysis when walking at slow, normal, and fast self-selected walking velocity	Maximum knee angle in swing phase is not influenced by walking velocity and stays con- stant at 64°, as physiological. Knee flexion at heel strike is constantly at 4°. More symmetrical step length. Decreased ground reaction force at the begin- ning of stance phase.	n.a.	
Motion analysis when walking with small steps	When walking with small steps, in 95% of all steps swing phase was released. With C-Leg swing phase was released in 75% of steps.	n.a.	
Walking backwards	Prosthesis never switched to swing phase resistance and therefore safe loading was al- ways possible.	n.a.	
Motion analysis ascend- ing	73% of subjects were able to ascend stairs reciprocally. With C-Leg subjects ascended stairs using the step-by-step strategy.	n.a.	
Motion analysis descend- ing	All subjects descended stairs reciprocally with Genium and C-Leg. The motion cycle is similar.	n.a.	
Motion analysis ascend- ing (10° decline)	Knee angle was increased by 14%.	n.a.	
Motion analysis descend- ing (10° decline)	Knee angle was increased by 13%.	n.a.	
Motion analysis and sway measurements when standing on a horizontal	Standing on horizontal ground: No difference.	n.a.	
ground or ramp (10° decline)	Standing on decline: Body weight was distributed evenly on both legs. Body posture was more relaxed. Reduced activity to keep up posture. Decreased hip moments by 76%.		
	OutcomesMotion analysis when walking at slow, normal, and fast self-selected walking velocityMotion analysis when walking with small stepsWalking backwardsWalking backwardsMotion analysis ascend- ingMotion analysis descend- ing (10° decline)Motion analysis and sway measurements when standing on a horizontal ground or ramp (10° decline)	OutcomesResults for Genium compared to C-LegMotion analysis when walking at slow, normal, and fast self-selected walking velocityMaximum knee angle in swing phase is not influenced by walking velocity and stays con- stant at 64°, as physiological. Knee flexion at heel strike is constantly at 4°. More symmetrical step length. Decreased ground reaction force at the begin- ning of stance phase.Motion analysis when walking with small stepsWhen walking with small steps, in 95% of all steps swing phase was released. With C-Leg swing phase was released in 75% of steps.Walking backwardsProsthesis never switched to swing phase resistance and therefore safe loading was al- ways possible.Motion analysis ascend- ing73% of subjects were able to ascend stairs reciprocally. With C-Leg subjects ascended stairs using the step-by-step strategy.Motion analysis descend- ing (10° decline)All subjects descended stairs reciprocally with Genium and C-Leg. The motion cycle is similar.Motion analysis descend- ing (10° decline)Knee angle was increased by 14%.Motion analysis ascend- ing (10° decline)Standing on horizontal ground: No difference.Motion analysis ad sway measurements when standing on a horizontal ground or ramp (10° decline)Standing on decline: Body weight was distributed evenly on both legs. Body posture was more relaxed. Reduced activity to keep up posture. Decreased action to keep up posture.	

"With the Genium knee prosthesis a new generation of microprocessor controlled knees is available to the supply team. This joint offers the technical requirements for the expansion of mobility for amputees in daily life and this is exceeding the dimensions known until now. Affected aspects are function as well as safety. The designed advantages of the specific functions – sitting, standing, walking – and the intuitive switch between these function states could be proved with measured data and were confirmed by patient reports. A reliable and precise switch between stance phase resistance and swing phase resistance plays the key role in safe functioning of the prosthesis. This switch in functionally is necessary and at the same time safety-critical. The high reliability of C-Leg could be proved for Genium too. The safety potential of Genium is even improved when amputees walk backwards or when they walk with small steps." (Blumentritt et al.. 2012)

<u>A Back to overview table</u>

Otto Bock Healthcare Products, Vienna, Austria.

Technologie und Funktionsweise des Genium-Prothesenkniegelenks

(Technology and functionality of the Genium knee prosthesis)

Orthopädie-Technik 2011; 12:898-903.

Products	Genium						
Major Findings	With Genium:						
	 → Allows for a symmetry → Increased safety re → Stair function enable 	etrical and phys sults from assu les reciprocally	iological gait patte red toe clearing at stair climbing stra	ern any walking velocity itegy			
			3				

Stair and Obstacle Function: (1) Release, (2) Knee flexion, (3) Unloaded extension movement, (4) Loaded extension movement. (Kampas *et al* 2011)

Populatio	n	No	t applicable						
Study Des	sign	Technical Report							
Results									
Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consump- tion	Safety	Activity, Mobility, ADLs	Preference, Satisfac- tion, QoL	Health economics
Category									
Level Wall	king	Pre	e-Flex		 Knee Facilit Full gr earlier Peak f ward r tion ef 	angle at h ated yielc round cor r stage. force of h movemen ffect.	neel strike in a ling. ntact of the fo eel-strike is c t which resul	a 4° flexion p ot is reached converted inte ts in a shock	osition. d in an o a for- absorp-
		Ad Co	aptive Yieldi ntrol	ng-	 Maxim 17°. When angle 	num knee walking o higher th	flexion angle on a decline, an 17° possil	in stance ph maximum kn ble.	ase is ee flexion

Category		
	Swing phase release by dynamic stability control (DSC)	 Swing phase is released depending on the position of the point of applied force on the foot and not on the amount of load. → Swing phase is released even when the subject is carrying additional load or if he supports himself on a shopping cart. Forward rotation of the prosthesis. → Release of swing phase when stepping backwards is avoided. Straight knee joint. Lower leg is in a forward tilted position. Minimum load.
	Adaptive swing phase control	 Pre-setting of the maximum knee angle is 65° (adjustable). Knee angle is constant even when walking at varying speeds. If the maximum knee angle is not reached when walking slowly, the extension movement is slowed down and foot clearance is assured.
Stairs	Release of function	Through a backward movement of the prosthesis after lifting the foot.
	Knee flexion	A low flexion resistance of the prosthesis enables that knee flexion is initiated through a hip flexion movement.
	Extension before ground contact (unloaded)	Extension movement in swing phase is performed faster and therefore the foot can be put on the next step in a controlled way. Flexion resistance is blocked.
	Extension after ground contact (loaded)	Extension resistance is increased further, so that the knee reaches a straight position in a controlled manner. Flexion resistance stays blocked.
Ramps, Hills	Adapted resistance	Allows for an increased flexion angle in stance phase compared to level walking, where it is lim- ited to 17° knee flexion.
	Release of swing phase	 Occurs even when the prosthesis is in a flexed and loaded position. Flexion is opened only continuously so that the prosthesis allows only for smooth motion clear- ance. Additional foot clearance is provided
Uneven Ground, Obstacle Course	Crossing Obstacles	Obstacles can be crossed using the same function as for stair ascent.
Safety	Dynamic stability control (DSC)	Allows for a safe switch between stand and swing phase resistance.
Activity, Mobility, Activities of daily living (ADLs)	Functions	 Activated intuitively through a motion pattern. The joint is releasing, blocking or changing the resistance continuously between conditions. After finishing conducting the activity the joint switches back automatically to base mode. Following functions can be activated: Stair and obstacle

Category

	Standing Sitting
MyModes	 Activated through repeated bouncing or through remote control Five MyModes for special activities can be pro- grammed.

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

Author's Conclusion "The manufacturer is convinced that Genium prosthetic system has a significant impact on rehabilitation and disability compensation based on optimized physiological walking including dynamic stability control, implementation of functions and five MyModes, as well as other in the article mentioned improvements. The prior target is to minimize the difference between the natural body function of the sound locomotor system and the artificial replacement by the prosthesis during walking but also during other activities of daily living." (Kampas & Seyr 2011)

↑ Back to overview table

Copyright:

© 2018, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This document contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this document.