

## Reference

Major M, Stine R, Heckathorne C, Fatone S and Gard S

Northwestern University Prosthetics-Orthotics Center, Northwestern University  
Feinberg School of Medicine, Chicago, USA

# Comparison of range-of-motion and variability in upper body movements between transradial prosthesis users and able-bodied controls when executing goal-oriented tasks

Journal of NeuroEngineering and Rehabilitation 2014, 11:132.

## Products


### Myoelectric prosthesis

## Major Findings

With myoelectric prosthesis users compared to able-bodied controls:

- **Shoulder and trunk movements are common compensatory motions in prosthesis users.**
- **Increased variability in movement suggests that prosthesis users do not stick to a defined motor strategy.**
- **Kinematic repeatability may increase with prosthesis experience.**

### Average range of motion for carton pouring task



Upper body range of motion (RoM) was analysed on able-bodied controls and myoelectric transradial prosthesis users during execution of carton pouring task (lifting a carton, located at midline of the body, and emptying the liquid contents into a jar on the contralateral side with minimal spilling). Results indicate that prosthesis users demonstrate a significant increase in shoulder abduction, trunk transverse rotation, trunk lateral flexion and trunk forward flexion than able-bodied subjects.

## Population

Subjects:	6 able-bodied controls 7 myoelectric transradial prosthesis users
Prosthesis:	System Electric Hand, MyoHand VariPlus Speed Hand, Transcarpal Hand, Motion Control Hand, i-Limb Ultra Revolution, i-Limb Ultra and i-Limb Hand
Amputation causes:	4 traumatic, 3 congenital
Mean age:	able-bodied individuals - $35 \pm 11$ year prosthetic users - $49 \pm 18$ years
Mean time since amputation:	$9.5 \pm 11.0$ years

## Study Design

Observational study:

Participants were requested to execute five goal-oriented tasks while seated (carton pouring, page turning, food cutting, lifting and transferring weighted object, lifting and transferring tray). Able-bodied controls and prosthesis users performed these tasks using their non-dominant and prosthetic limb, respectively.

## Results

Body Function		Activity			Participation	Others	
Mechanics	Pain	Grip patterns / force	Manual dexterity	Activities of daily living (ADL)	Satisfaction and Quality of life (QoL)	Training	Technical aspect

Category	Outcomes	Results for myoelectric prosthesis users compared to able-bodied controls:	Sig.*
Mechanics	Goal orientated tasks: <ul style="list-style-type: none"> <li>• carton pouring</li> <li>• page turning</li> <li>• food cutting</li> <li>• lifting and transferring weighted object</li> <li>• lifting and transferring tray</li> </ul>	The majority of prosthesis users were unable to routinely execute food cutting and page turning tasks.	n.a.
		<b>Prosthesis users demonstrated a significant increase in shoulder abduction, trunk transverse rotation, trunk lateral flexion, and trunk forward flexion RoM when executing carton pouring, lifting and transferring tasks.</b>	--
		No difference in shoulder and elbow flexion/extension RoM was observed.	0
		<b>Kinematic variability was high for prosthetic users.</b>	--
		<b>Kinematic repeatability was low for prosthetic users.</b>	--

\* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

## Author's Conclusion

"Transradial prosthesis users utilize shoulder abduction and trunk movement as compensatory motions to execute goal-oriented tasks, and the majority of these motions are accompanied by increased kinematic variability when compared to able-bodied controls. The average repeatability of upper body kinematics was positively associated with prosthesis experience. As these dynamics may be necessary to compensate for the absence of active distal DoFs (degrees of freedom) in the prosthetic arm, transradial prosthesis users may benefit from dedicated training that: 1) encourages optimization of these dynamics to facilitate execution of ADLs, and 2) fosters adaptable but reliable motor strategies." (Major et al. 2014)

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.