

Reference

Ernst M, Altenburg B, Bellmann M, Schmalz T.

Research Biomechanics, CR&S, Otto Bock HealthCare GmbH, Göttingen, Germany.

Standing on slopes – how current microprocessor-controlled prosthetic feet support transtibial and transfemoral amputees in an everyday task

Journal of NeuroEngineering and Rehabilitation (2017) 14:117.

Products

(Meridium, Elan, Proprio, TSA, Raize) vs conventional prosthetic feet

Major Findings

→ Only Meridium

- Joint angles and joint torques are closest to non-amputees for
 - Standing on an upward slope of 10°
 - Standing on a downward slope of 10°
- Autoadaptive dorsiflexion stop and sufficient range of motion improve symmetric loading
 - Clear superiority for Meridium compared to other microprocessor-controlled feet (MPFs)

→ With microprocessor-controlled prosthetic feet (MPFs) compared to conventional prosthetic feet:

- Full adjustment of the ankle joint improves symmetry of vertical ground reaction forces
- Compensatory posture necessary for transtibial and transfemoral amputees, when prosthetic foot has no automatic ankle angle adaptation





Figure 1: Differences in ankle angles when standing on a downward slope (10°) are illustrated for Meridium, other MPFs, conventional prosthetic feet and non-amputees.

Population	Subjects:	4 unilateral transtibial amputees (TT) 4 unilateral transfemoral amputees (TF) 20 non-Amputees (control group)
	Previous prosthesis foot:	Conventional prosthetic feet (Non-MPF)
	Amputation causes:	not reported
	Mean age:	4 TT: 56.2 yrs ± 12 yrs; 4 TF: 44.5 yrs ± 3 yrs 20 non-Amputees: 22.5 yrs ± 3 yrs
	Mean time since amputation:	> 3 yrs
	MFCL:	K3 and K4

Study Design Interventional, crossover design:

Transtibial (N = 4) and transfemoral* (N = 4) amputees



The order of wearing the MPFs was randomized for each subject. The graph shows an example, where Meridium is selected as the third MPF.

* Transfemoral amputees were not equipped with the Raize foot, which reduced the number of data collection session from 6 to 5.

Results

Activities								Participation	Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consumption	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health economics
Category	Outcomes		Results for prosthetic TT and TF vs. non-Amp		Sig.*	Results for sound TT and TF vs. non-Amp		Sig.*	
Level walking	Ankle torque		<i>Positive values: Dorsiflexion; Negative values: Plantarflexion</i>		0	TT:		++	0
			No sig. differences for all feet.			Elan: +0.34 ± 0.08			
	Knee torque		<i>Positive values: Knee extension; Negative values: Knee flexion</i>			TT:			
		No sig. differences for all feet.		..		No sig. differences for all feet.		0	
		Elan: -0.01 ± 0.06							
Level walking	Hip torque		<i>Positive values: Hip flexion; Negative values: Hip extension</i>		..	TT:		0	0
			No sig. differences for all feet.			Proprio: +0.07 ± 0.06			
			No sig. differences for all feet.						

Category	Outcomes	Results for prosthetic TT and TF vs. non-Amp	Sig.*	Results for sound TT and TF vs. non-Amp	Sig.*
Ramps, Hills	Ankle torque	<i>Positive values: Dorsiflexion; Negative values: Plantarflexion</i>			
	Down (10°)	TT:		TF:	
		Elan: -0.10 ± 0.08	--	Everyday Feet: +0.43 ± 0.13	++
		Proprio: -0.10 ± 0.15	--		
		TSA: +0.04 ± 0.02	++	Elan: +0.39 ± 0.08	++
		Raize: +0.06 ± 0.03	++	Proprio: +0.39 ± 0.03	++
		TF:		TF:	
		Elan: -0.05 ± 0.07	--	Everyday Feet: +0.42 ± 0.14	++
		TSA: +0.04 ± 0.04	++		
		Up (10°)	TT:		
		Everyday foot: +0.62 ± 0.15	++		
		Elan: +0.23 ± 0.01	++		
		Proprio: +0.46 ± 0.02	++		
		Raize: +0.52 ± 0.15	++		
		TF:		No sig. differences for all feet.	0
		Everyday foot: +0.66 ± 0.07	++		
		Elan: +0.48 ± 0.09	++		
		Proprio: +0.52 ± 0.04	++		
		TSA: +0.61 ± 0.08	++		
	Knee torque	<i>Positive values: Knee extension; Negative values: Knee flexion</i>			
	Down (10°)	TT:		TF:	
		Everyday feet: -0.16 ± 0.04	--	Elan: -0.07 ± 0.11	--
		Elan: -0.17 ± 0.04	--		
		Proprio: -0.16 ± 0.06	--		
		Raize: -0.03 ± 0.07	--		
		TF:			
		Everyday feet: -0.21 ± 0.28	--		
		Elan: -0.21 ± 0.05	--		
		Proprio: -0.24 ± 0.03	--		
		TSA: -0.09 ± 0.09	--		
	Up (10°)	TT:			
		Elan: +0.26 ± 0.04	++		
		Proprio: +0.38 ± 0.06	++		
		TF:		No sig. differences for all feet.	0
		Everyday feet: +0.29 ± 0.07	++		
		Elan: +0.30 ± 0.06	++		
		Proprio: +0.31 ± 0.09	++		
		TSA: +0.31 ± 0.04	++		

Category	Outcomes	Results for prosthetic TT and TF vs. non-Amp	Sig.*	Results for sound TT and TF vs. non-Amp	Sig.*
	Hip torques	<i>Positive values: Hip flexion; Negative values: Hip extension</i>			
	Down (10°)	TT:			
		Proprio: +0.09 ± 0.04	++		
		TSA: +0.14 ± 0.10	++		
		Raize: +0.12 ± 0.05	++		
		TF:		No sig. differences for all feet.	0
		Everyday feet: +0.26 ± 0.12	++		
		Elan: +0.1 ± 0.12	++		
		Proprio: +0.24 ± 0.1	++		
		TSA: +0.13 ± 0.04			
	Up (10°)	TT:		TF:	
		Raize: +0.15 ± 0.14	++	Everyday feet: +0.05 ± 0.05	++
* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)					

Author's Conclusion

"A prosthetic foot that combines both key features – an auto-adaptive dorsiflexion stop and sufficient ROM to completely adapt to inclinations - enables lower limb amputees to stand on slopes in an almost natural manner. The biomechanical parameters indicate that this concept is superior to conventional passive feet or feet which provide only one key design feature such as a sufficient ROM. Finally, the results indicate that both, TT and TF amputees, benefit from such a foot.." (Ernst et al, 2017)

© 2018, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.