

## Reference

Grabowski AM, D'Andrea S.

Integrative Physiology Department, University of Colorado Boulder, Boulder, CO.

# Effects of a powered ankle-foot prosthesis on kinetic loading of the unaffected leg during level-ground walking

J Neuroeng Rehabil. 2013 Jun 7;10:49.

## Products

### BiOM (Bionic powered ankle-foot prosthesis)

## Major Findings


With BiOM compared to conventional passive-elastic feet (Passive) and non-amputees (Control):

#### → Reduction of the resultant ground reaction forces (GRF) acting on the unaffected

Decrease by 2.1% to 10.7% (for walking speeds from 0.75 to 1.5 m/s) with BiOM compared to Passive

#### → BiOM decreased the external adduction moment (EAM) on the unaffected knee

Reduction by 20.6% and 12.2% (for walking speeds of 1.5 and 1.75 m/s, respectively) with BiOM compared to Passive




The average resultant first peak GRF on the unaffected leg for all walking speeds (0.75 – 1.5) was calculated; BiOM presented a reduced resultant GRF compared to passive feet.

## Population

Subjects:	Seven unilateral, transtibial amputees (Seven male) Seven non-amputees (Control)
Previous prosthetic feet:	Flex-Foot, Ossur (3); Axion, Otto Bock (1); Venture, College Park (1); Renegade, Freedom Innovations (1); Limb Logic, Ohio Willow Wood (1)
Amputation causes:	Trauma
Mean age:	Amputees: 45 ± 6 yrs Control: 48 ± 7 yrs
Mean time since amputation:	21.6 ± 11.6 yrs
MFCL:	K3

## Study Design

Interventional, pre- to post design:



Subjects with an amputation completed two experimental walking sessions; one using their own passive-elastic foot and one using the powered ankle-foot prosthesis (acclimation session of at least 2 hours). Non-amputee subjects completed one experimental session.

Each subject walked at 0.75, 1.00, 1.25, 1.50, and 1.75 m/s, while the stiffness and power delivery of the powered prosthesis was adjusted so that prosthetic ankle angle at toe-off and net positive mechanical work matched average biological ankle data.

## Results

Functions and Activities					Participation			Environment	
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic energy consumption	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health Economics
Category	Outcomes		Results for BiOM vs Passive vs Control					Sig.*	
Level Walking	Unaffected leg: 1 <sup>st</sup> peak ground reaction force (GRF) [N/kg]		<b>For all walking speeds from 0.75 to 1.5 m/s, amputees with BiOM showed a significant decrease of the 1<sup>st</sup> peak of GRF by 2.1% to 10.7% when compared to Passive.</b>					++	
			There were no significant differences between BiOM and Control for the 1 <sup>st</sup> peak GRF while walking with different walking speeds.					0	
	Unaffected leg: GRF rate [N/kg/s]		There were no significant differences between BiOM, Passive and Control while walking, except for one specific walking speed: <b>For 1.25 m/s, Control showed significant lower GRF rates than amputees with BiOM (+30.3%) or Passive (+49.1%).</b>					--	
	Unaffected leg: 1 <sup>st</sup> peak external knee adduction moment (EAM) [Nm/kg]		<b>For the two fastest walking speeds (1.5 and 1.75 m/s), the BiOM reduced the 1<sup>st</sup> peak of the EAM significantly by 20.6% and 12.2%, respectively, when compared to Passive.</b>					++	
			Apart of that, there were no significant differences between BiOM, Passive and Control while walking with different walking speeds.					0	
	Unaffected leg: knee EAM rate [Nm/kg/s]		There were no significant differences between BiOM, Passive and Control in knee EAM rate for all walking speeds.					0	

\* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

---

## Author's Conclusion

"A passive-elastic prosthesis cannot emulate normative biological function during the stance phase of walking; thus people with a lower-extremity amputation employ compensatory mechanics and have a higher incidence of musculoskeletal injury, specifically knee osteoarthritis in their unaffected leg. A biomimetic prosthesis could mitigate the risk of knee osteoarthritis by decreasing unaffected leg forces and knee moments. In this investigation, we found that when people with a unilateral transtibial amputation due to trauma and K3 level of amputation used a powered ankle-foot prosthesis during level-ground walking over a range of speeds, they reduced the peak resultant force and knee adduction moment on their unaffected leg compared to when they used their own passive-elastic prosthesis. At the walking speed closest to preferred, subjects with an amputation using a powered ankle-foot prosthesis reduced their unaffected peak knee EAM by over 20%. A significant reduction in peak knee EAM has the potential to decrease the risk of knee osteoarthritis. Based on these results, we conclude that a biomimetic powered ankle-foot prosthesis could potentially limit musculoskeletal stress to the contralateral leg during walking, thus decreasing the risk of secondary injury in people with a lower-extremity amputation." (Grabowski and D'Andrea, 2013)

© 2018, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.