

The biomechanical response of persons with transfemoral amputation to variations in prosthetic knee alignment during level walking

Journal of Rehabilitation Research and Development (2016); 53(6), 1089-1106

DOI: 10.1682/JRRD.2014.12.0311. [Open Access](#).

With 3R95 in different knee alignment conditions (bench, 2cm anterior, and 2cm posterior):

→ **Anterior alignment creates instability and triggers multiple compensatory strategies to coordinate stance-phase knee flexion compared to the bench alignment at the time of opposite toe off**

- decreased loading rate of the prosthesis
- decreased step length (~6% smaller, $p=0.023$)
- increased trunk flexion (factor ~1.6x higher, $p=0.007$)
- affected limb in a more vertical posture (~73% more upright, $p=0.002$)
- significantly decreased extension moment at the knee joint (~38% smaller, $p=0.02$) → corresponds to a decrease in the margin of stability
- significantly increased early-stance hip extension moment (~46% higher, $p=0.019$) (interpretation: instability caused by anterior alignment must be compensated using hip extensors to maintain knee joint stability)

→ **Posterior alignment has fewer impact on patients regarding stance-phase, knee joint control compared to the bench alignment**

- significantly longer terminal double support phase (~7% longer, $p = 0.007$)
- extension moment at the knee joint significantly increased at toe off (factor ~1.3x higher, $p<0.001$) → corresponds to increase in the margin of stability

→ **hip torque strategy alone might be insufficient to counteract the effect of the alignment perturbation**





Figure 1. Mean affected-side joint kinematics at opposite toe off during level walking. * thigh angle > 0: distal end of thigh is anterior to proximal end. ** value < 0: extension. (--) significant difference to bench alignment ($p<0.025$).

Population	Subjects:	11 participants (9 = male, 2 = female)
	Amputation Level:	unilateral transfemoral (n = 10), unilateral knee disarticulation (n = 1)
	Previous prosthesis:	n.a.
	Amputation causes:	trauma (n = 8), cancer (n = 3)
	Mean age:	46 ± 16 years
	Mean time since amputation:	n.a.
	MFCL:	K3 or higher

Study Design Single-blinded, pseudo-randomized crossover design:



*BASE: knee bench-aligned, ANT: 2 cm anterior knee translation, POST: 2 cm posterior knee translation

First, the subjects underwent a physical examination where residual limb length, hip range of motion and muscle strength were examined. The participants could practice for about 15 minutes on the treadmill, or until comfortable and select their personal walking speed.

All subjects started with the BASE condition. To accommodate to the BASE condition, the subjects could walk overground for <5 minutes. Next, the participants walked at a self-selected speed on a level treadmill surface for 2.5 minutes. After 2 minutes, quantitative gait data was recorded for 30 seconds. Subsequently, the protocol was repeated with the same constant speed for the 2 other alignment conditions (ANT, POST) in a randomized order. Prior to each condition, they were again offered an accommodation time of <5 minutes.

For the anterior knee alignment, the knee was translated 2 cm in anterior direction (ANT). Similar in posterior direction (POST). All subjects completed all three alignment conditions.

Results

Functions and Activities						Participation			Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic Energy Consumption	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health Economics

Category	Outcomes	Results for 3R95	Sig. ^a (vs. BASE)
----------	----------	------------------	---------------------------------

Category	Outcomes	Results for 3R95			Sig. ^a (vs. BASE)	
Level Walking ^b	Spatiotemporal characteristics	<u>BASE and POST:</u> slightly longer step length on the affected side compared to the intact side				
		<u>ANT:</u> significantly decreased step length on the affected side compared to BASE				
			BASE	ANT	POST	
		Intact Limb Step length [cm]	52 ± 7	54 ± 5	53 ± 6	0
Affected Limb Step length [cm]	54 ± 8	51 ± 7	55 ± 8	--		
Stance phase [ms]	874 ± 115	859 ± 134	884 ± 123	0		
<u>Initial double support phase:</u>						
- significantly longer for ANT compared to BASE						
- shortest overall for POST						
<u>Terminal double support phase:</u>						
- significantly shorter for ANT						
- significantly longer for POST						
		BASE	ANT	POST		
Initial double support [ms]	219 ± 61	238 ± 75	207 ± 59	--		
Single support [ms]	442 ± 41	426 ± 50	450 ± 48	0		
Terminal double support [ms]	214 ± 33	196 ± 31	229 ± 33	--		
Affected-side joint kinematics [°] ^c	Opposite Toe Off	<u>ANT compared to BASE:</u>				
		- persistent significant increase in trunk flexion & anterior tilt				
		- thigh: subjects maintained their affected limb in a significantly more vertical posture				
			BASE	ANT	POST	
Ankle	-3.9 ± 2.4	-3.2 ± 2.5	-3.9 ± 2.5	0		
Knee	-4.8 ± 3.0	-5.4 ± 2.5	-5.0 ± 2.7	0		
Hip	19.6 ± 5.7	21.1 ± 5.6	19.0 ± 5.7	0		
Thigh	3.0 ± 2.3	0.8 ± 3.1	2.8 ± 2.3	--		
Pelvis	16.1 ± 5.9	19.6 ± 6.9	15.6 ± 6.0	--		
Trunk	1.2 ± 2.9	3.8 ± 3.6	0.6 ± 3.1	--		
Knee Break (late-stance initiation of knee flexion)						
<u>ANT compared to BASE:</u>						
- persistent significant increase in trunk flexion & anterior tilt						
	BASE	ANT	POST			
Ankle	7.1 ± 4.0	7.7 ± 4.2	7.1 ± 3.8	0		

Category	Outcomes	Results for 3R95			Sig. ^a (vs. BASE)
	Knee	0.6 ± 2.5	-0.2 ± 2.1	0.9 ± 2.5	0
	Hip	-5.6 ± 6.9	-4.5 ± 6.3	-6.0 ± 6.5	0
	Thigh	-20.7 ± 3.2	-22.2 ± 3.3	-21.1 ± 3.2	0
	Pelvis	17.0 ± 6.1	19.1 ± 6.6	17.0 ± 6.2	--
	Trunk	1.2 ± 2.8	3.1 ± 3.2	1.1 ± 3.1	--

Affected-side joint kinetics and GRFs^d

Opposite Toe Off

ANT compared to BASE:

- extension moment at the knee joint significantly decreased
- significantly increased internal hip extension moment
- fore/aft braking force significantly decreased

POST compared to BASE:

- extension moment at the knee joint significantly increased
- no sig. difference in the hip extension or flexion moment

	BASE	ANT	POST	
Knee Moment [Nm/kg]	-0.13±0.08	-0.08 ± 0.07	-0.30±0.08	--
Hip Moment [Nm/kg]	0.26 ± 0.16	0.38 ± 0.19	0.23 ± 0.14	--
Fore/aft GRF [N/kg]	-0.05±0.04	-0.02 ± 0.04	-0.06 ± 0.04	--
Vertical GRF [N/kg]	0.84 ± 0.08	0.84 ± 0.08	0.87 ± 0.09	0

Knee Break (late-stance initiation of knee flexion)

ANT compared to BASE:

- propulsive force significantly increased

POST compared to BASE:

- propulsive force significantly decreased

	BASE	ANT	POST	
Knee Moment [Nm/kg]	0.009 ± 0.03	0.008 ± 0.03	-0.004 ± 0.04	0
Hip Moment [Nm/kg]	-0.69 ± 0.43	-0.60 ± 0.36	-0.75 ± 0.49	0
Fore/aft GRF [N/kg]	-0.0004 ± 0.01	0.011 ± 0.01	-0.01 ± 0.01	--
Vertical GRF [N/kg]	0.38 ± 0.21	0.39 ± 0.22	0.35 ± 0.21	0

^a no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)

^b Results shown as Mean ± standard deviation (SD);

self-selected walking speed: 0.8 ± 0.2 m/s

^c Positive angle indicates the distal end of the thigh is anterior to the proximal end.

^d Positive values indicate: external knee flexion moment, internal hip extension moment, anterior-directed fore/aft GRF, and a superior-directed vertical GRF.

Author's Conclusion

“In summary, subjects responded to a destabilizing alignment perturbation by increasing their hip extension moment during early stance phase, which confirmed our hypothesis that persons with transfemoral amputation rely on a hip extensor strategy to maintain knee joint stability. However, subjects also decreased the rate at which they loaded their prosthesis, decreased their step length, increased their trunk flexion, and maintained their limb in a more vertical posture at the time of opposite toe off, suggesting that a hip torque strategy alone was perhaps insufficient to counteract the effect of the alignment perturbation. Although we hypothesized that subjects would respond to a stabilizing alignment by increasing their hip flexion moment to initiate knee break during late stance, our results did not support this hypothesis. Instead, knee break was delayed for the posterior alignment condition. From a clinical perspective, these results suggest that a bias toward posterior alignment may have fewer implications for the patient in terms of knee joint control. However, it is possible that an overly stable knee joint may also have negative consequences for swing phase. To our knowledge, this is the first systematic study of prosthetic knee joint alignment in which lower-extremity and trunk kinematics were investigated in conjunction with sagittal-plane hip and knee kinetics to understand the mechanisms of prosthetic knee joint control. Accordingly, we believe the results from this study may provide new insights into the biomechanical mechanisms contributing to the sagittal stability and control of a passive prosthetic knee joint.” (Koehler-McNicholas et al., 2016)

Author's Affiliation

¹Rehabilitation Engineering Research Program, Minneapolis VA Health Care System, Minneapolis, MN

²Northwestern University Prosthetics-Orthotics Center for Education and Research, Department of Physical Medicine and Rehabilitation, Feinberg School of Medicine, Northwestern University, Chicago, IL

³Department of Biomedical Engineering, McCormick School of Engineering and Applied Science, Northwestern University, Evanston, IL

⁴Jesse Brown VA Medical Center, Chicago, IL

⁵Rehabilitation Institute of Chicago, Chicago, IL

© 2025, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.