

Reference

Atsushi Makimoto^{1,2}, Yoko Sano^{1,2}, Saturo Hashizume¹, Akihiko Murai¹, Yoshiyuki Kobayashi¹, Hiroshi Takemura^{1,2}, and Hiroaki Hobara¹

Ground Reaction Forces During Sprinting in Unilateral Transfemoral Amputees

Journal of Applied Biomechanics (2017); 33:406–409.

DOI: 10.1123/jab.2017-0008

Products

3S80, 1E90

Major Findings

With 3S80 knee with 1E90 runner foot comparing prosthetic and intact limb:

→ **Asymmetrical loading of the limbs due to significantly smaller forces and impulses in the prosthetic limb than the intact limb during sprinting.**

- Bilateral differences in peak ground reaction forces for the vertical (26%), braking (54%), propulsive (27%) and medial components (50%) ($p < 0.01$, large effect size).

→ **Altered braking propulsion strategy in the prosthetic limb compared to the intact side.**

- The prosthetic limb exhibits significantly reduced braking impulses (-52% , $p < 0.01$) but comparable forward propulsive impulse relative to the intact limbs.
- 25% earlier transition during stance from braking to propulsion in the prosthetic limb (zero fore-aft shear^d, $p < 0.01$) – indicating a shortened braking phase relative to the intact limb.


→ **The study supports the development of limb-specific rehabilitation and training strategies.**

Population

Subjects:	9 sprinters (5 male, 4 female)
Amputation level:	unilateral transfemoral
Prosthetic components:	Knee: 3S80 (9), Foot: 1E90 (9)
Amputation causes:	Trauma (5), Cancer (2), Infection (1), Sarcoma (1)
Mean time since amputation:	12.8 ± 8.3 years (range 3.5–29 years)
Mean age:	32.9 ± 11.7 years
Height:	1.61 ± 0.09 m
Mean body mass:	55.9 ± 9.2 kg without prosthesis
Personal sprinting record:	17.3 ± 2.3 s for 100 m
Running experience:	range: 0.9–3.2 years
MFCL:	not reported

Study Design

Prospective study design:



The participants all used the same knee units and running-specific prosthesis. They were asked to perform maximal sprinting on a 40m runway while kinetic and kinematic data were acquired using force plates and an optical motion capture system. Between trials, the participants rested for an adequate period. Four successful steps of each limb were analysed for each participant.

Results

Activities		Participation				Body function			Other
Sprinting, running, jumping	Other sports	Leisure / recreational sports	Competitive sports	Paralympic sports	Preference, satisfaction, QoL	Biomechanics (kinematics / kinetics)	Clinical (metabolic / performance)	Medical (pain, injuries)	Technical aspects / alignment

Category	Outcomes	Results for 3S80, 1E90 ^a	Sig. ^b																								
Biomechanics (kinematics / kinetics)	Stance Time (s)	Significantly longer stance duration on the prosthetic limb.																									
		<table><tr><th>Prosthetic limb</th><th>Intact limb</th><th>Difference</th></tr><tr><td>0.140 (±0.013)</td><td>0.127 (±0.015)</td><td>-12.6% ^c</td></tr></table>	Prosthetic limb	Intact limb	Difference	0.140 (±0.013)	0.127 (±0.015)	-12.6% ^c	++																		
Prosthetic limb	Intact limb	Difference																									
0.140 (±0.013)	0.127 (±0.015)	-12.6% ^c																									
	Ground Reaction Force (GRF) (N/BW)	Significantly lower peak forces in the prosthetic leg during sprinting in all measured directions compared to the intact limb (p<0.01, large effect size), except for the peak lateral force.																									
		<table><tr><th>Force component</th><th>Prosthetic limb</th><th>Intact limb</th><th>Difference</th></tr><tr><td>Peak vertical</td><td>2.86 (±0.67)</td><td>3.86 (±0.58)</td><td>25.9%</td></tr><tr><td>Peak braking</td><td>-0.44 (±0.14)</td><td>-0.95 (±0.20)</td><td>53.9%</td></tr><tr><td>Peak propulsive</td><td>0.40 (±0.10)</td><td>0.55 (±0.07)</td><td>27.2%</td></tr><tr><td>Peak medial</td><td>0.23 (±0.07)</td><td>0.46 (±0.12)</td><td>50.0%</td></tr><tr><td>Peak lateral</td><td>-0.08 (±0.05)</td><td>-0.15 (±0.17)</td><td>48.3%</td></tr></table>	Force component	Prosthetic limb	Intact limb	Difference	Peak vertical	2.86 (±0.67)	3.86 (±0.58)	25.9%	Peak braking	-0.44 (±0.14)	-0.95 (±0.20)	53.9%	Peak propulsive	0.40 (±0.10)	0.55 (±0.07)	27.2%	Peak medial	0.23 (±0.07)	0.46 (±0.12)	50.0%	Peak lateral	-0.08 (±0.05)	-0.15 (±0.17)	48.3%	-- -- -- -- 0
Force component	Prosthetic limb	Intact limb	Difference																								
Peak vertical	2.86 (±0.67)	3.86 (±0.58)	25.9%																								
Peak braking	-0.44 (±0.14)	-0.95 (±0.20)	53.9%																								
Peak propulsive	0.40 (±0.10)	0.55 (±0.07)	27.2%																								
Peak medial	0.23 (±0.07)	0.46 (±0.12)	50.0%																								
Peak lateral	-0.08 (±0.05)	-0.15 (±0.17)	48.3%																								
	Zero fore-aft shear ^d (%)	The transition from braking to propulsion during stance happens 25.4 % earlier (p<0.01) in the prosthetic limb compared to the intact limb.																									
		<table><tr><th>Prosthetic limb</th><th>Intact limb</th><th>Difference</th></tr><tr><td>37.10 (±2.09)</td><td>49.73 (±2.88)</td><td>25.4%</td></tr></table>	Prosthetic limb	Intact limb	Difference	37.10 (±2.09)	49.73 (±2.88)	25.4%	--																		
Prosthetic limb	Intact limb	Difference																									
37.10 (±2.09)	49.73 (±2.88)	25.4%																									
	Impulses (Ns/BW)	Significantly lower braking impulse in the prosthetic limbs than in intact limbs.																									
		<table><tr><th>Component</th><th>Prosthetic limb</th><th>Intact limb</th><th>Difference</th></tr><tr><td>Braking impulse</td><td>-0.01 (±0.00)</td><td>-0.02 (±0.01)</td><td>52.2%</td></tr></table>	Component	Prosthetic limb	Intact limb	Difference	Braking impulse	-0.01 (±0.00)	-0.02 (±0.01)	52.2%	--																
Component	Prosthetic limb	Intact limb	Difference																								
Braking impulse	-0.01 (±0.00)	-0.02 (±0.01)	52.2%																								
		No significant differences in vertical, propulsive, medial or lateral impulses between limbs.	0																								

^a results shown as Mean (± Standard Deviation); relative difference expressed in percentage of intact limb (%INT)

^b no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)
significance set at p<0.05; trends set at 0.1>p>0.05

^c stance time difference: value stated in paper, calculation error possible

Category	Outcomes	Results for 3S80, 1E90 ^a	Sig. ^b
----------	----------	-------------------------------------	-------------------

^d Zero fore-aft shear: The time point during stance at which the anterior–posterior ground reaction force crosses zero, indicating the transition from braking to propulsion.

Author’s Conclusion	“Peak forces in vertical, braking, propulsive, and medial directions were significantly greater in intact limbs than those in prosthetic limbs, whereas there were no significant differences in peak lateral force between limbs. Further, significantly less braking impulses were observed in prosthetic limbs than in intact limbs; however, the other measured impulses were not different between limbs. Therefore, the results of the present study suggest that limb-specific rehabilitation and training strategies should be developed for transfemoral amputees wearing running-specific prostheses.” (Makimoto et al., 2017)		
Author’s Affiliation	¹ National Institute of Advanced Industrial Science and Technology ² Tokyo University of Science		

©2026, Otto Bock HealthCare Products GmbH (“Otto Bock”), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a ‘fair use’ of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond ‘fair use’, you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.