

Reference

Francois Genêt^a, Axel Ruetz^b, Rania Belmahfoud^c, Isabelle Loiret^d, Caroline Navarre^e, Isabelle Noizette^f, Laurent Thefenne^g, Léo Borrini^h, Frédéric Charlateⁱ, Claire Delbrouck^j, Frédéric De Lucas Vasquez^k, Stéphane Vigier^l, Guillaume Bokobza^m, Vincent Moiziardⁿ, Virgile Pinelli^o, Brice Lavrard^p, Frank Braatz^q

Impact of a microprocessor-controlled knee-ankle-foot orthosis in community ambulators with quadriceps insufficiency fitted with an SCO: a randomized crossover trial

Annals of physical and rehabilitation medicine, vol. 69,1 102057. 22 Nov. 2025, doi:10.1016/j.rehab.2025.102057

Products

C-Brace

Major Findings

With C-Brace compared to previous Stance-Control Orthosis (SCO):

→ Significantly improved mobility and endurance

- Significantly improved PLUS-M T-score by 21.3% (+9.9 points) with C-Brace ($p < 0.001$)
- Significantly improved distance covered during the 6MWT by 19.5% (+65.9 meters) with C-Brace ($p < 0.001$)

→ Lower risk of falling

- Significantly improved Activity specific balance confidence (ABC) scale by 52.5% (+28.8 points) with C-Brace ($p < 0.001$)
- Reduced risk of falling with C-Brace, as the ABC score was $> 80\%$




Figure 1: Improvements in the PLUS-M Mobility score, 6-minute walking test (6MWT) and the Activity specific balance confidence scale for C-Brace compared to SCO in the Per Protocol (PP) group. Legend: * $p < 0.001$.

→ Significantly improved participation for important functional activities


- Significantly improved Patient-Specific Functional Scale (PSFS) value by 4.1 points with C-Brace ($p < 0.001$)
- Improvements in ADLs were considered clinically relevant as they were difficult to manage with the SCO from the participants perspective

→ Significantly improved quality of life

- Significantly improved EQ-5D-5L health utility score by 27.2% (+0.19) with C-Brace ($p < 0.001$)
- C-Brace achieved values for the EQ-5D-5L utility score (0.88) reaching close to values reported for French cohort (0.905) and German cohort (0.88) in the literature compared to SCO (0.692)

- **Significant reduction in the use of walking aids with C-Brace ($p = 0.005$) in outdoor conditions**
- **87% of the participants preferred the C-Brace over their SCO**

Population	Subjects:	Intention to Treat (ITT): 38 adults (12 female) Per protocol (PP): 30 adults (10 female)
	Etiology:	SCO users with quadriceps insufficiency ITT: Polio (20), Trauma (4), SCI (2), Other (12) PP: Polio (16), Trauma (3), SCI (1), Other (10)
	Previous SCO	ITT: E-Mag (30), SPL-Basko (7), NEURO TRONIC (1) PP: E-Mag (25), SPL-Basko (4), NEURO TRONIC (1)
	Time since SCO use:	At least 3 months
	Mean age:	ITT: 52.3 ± 12.8 years PP: 50.8 ± 11.3 years
	MFCL:	ITT: 56.7% K2 ($n = 22$); 43.7% K3 ($n = 16$) PP: 56.7% K2 ($n = 17$); 43.7% K3 ($n = 13$)
	Contralateral side (PP):	86.7% no deficiency, 6.7% orthopaedic shoe, 3.3% orthopaedic insoles, 3.3% orthopaedic pad
	*SCI=Spinal Cord Injury, Polio=Poliomyelitis	



*Several rehabilitation sessions are provided after C-Brace fitting. The protocol imposed a minimum of 4h of rehabilitation.

For either SCO or C-Brace, the follow-up period after fitting lasted 2 to 3 months. At the end of each follow-up period, data was collected. After the C-Brace follow-up period and re-fitting with their SCO, a 1h rehabilitation session was required.

A total of 38 participants were enrolled (ITT = 38) and randomized in two arms. 37 participants used the C-Brace at least once. 30 participants completed all assessments without any major deviation from the protocol (PP = 30, C-Brace/SCO = 13 and SCO/C-Brace = 17)

Results

Functions and Activities						Participation	Environment
Biomechanics – Static Measurement	Biomechanics – Gait analysis	X-Rays	EMG	Functional tests	Clinical effects	Satisfaction	Health Economics

Category	Outcomes	Results for C-Brace vs. SCO (PP group)	Sig. ¹				
Functional tests	6-min walk test (6MWT) [m]	Significantly longer distance covered during the 6MWT by 19.5% (+65.9 meters) with C-Brace (p < 0.001):	++				
		<table><tr><th>C-Brace (PP)</th><th>SCO (PP)</th></tr><tr><td>404.3 ± 126.4 m</td><td>338.4 ± 108.5 m</td></tr></table>	C-Brace (PP)	SCO (PP)	404.3 ± 126.4 m	338.4 ± 108.5 m	
C-Brace (PP)	SCO (PP)						
404.3 ± 126.4 m	338.4 ± 108.5 m						

Category	Outcomes	Results for C-Brace vs. SCO (PP group)	Sig. ¹
Clinical effects - Mobility	Simplified Activities-specific Balance Confidence scale (ABC scale)	Significantly higher ABC scale by 52.5% (+28.8 points) with C-Brace (p < 0.001):	++
		C-Brace (PP)	SCO (PP)
		83.6 ± 17.6	54.8 ± 21.5
	Mobility Questionnaire PLUS-M (primary outcome measure)	Significantly higher PLUS-M mobility T-score by 21.5% (+9.9 points) with C-Brace compared to SCO (p < 0.001):	++
		C-Brace (PP)	SCO (PP)
		55.9 ± 8.5 points	46 ± 6.4 points
	Use of Walking Aids	Similar walking aids habits for C-Brace and SCO in indoor conditions for 80% of the participants.	0
		Significantly reduced use of walking aids in outdoor conditions with C-Brace (p = 0.005):	++
		• 30% (9/30) of the participants did not require aids anymore with C-Brace	
	Psychological Impact of Assistive Device scale (PIADS)	Improvement in the global score on the PIADS scale by 1.8 points.	n.a.
Satisfaction	Use of orthosis	70% of participants (21/30) reported daily use of C-Brace and 60% (18/30) use it more than 8h per day. The SCO was used daily by 63% (19/30) of the participants and over 8h per day by 53% (16/30).	n.a.
	Satisfaction questionnaire - QUEST 2.0	Significantly improved global satisfaction by 12.5% with C-Brace (p < 0.001):	++
		C-Brace (PP)	SCO (PP)
		4.5 ± 0.4	4.0 ± 0.6
		Higher satisfaction for sub score device (+0.6 points) and the services provided (+0.4 points) with C-Brace than with SCO.	n.a.
		Most important items reported were safety, effectiveness and weight after C-Brace use. After SCO use the main items were safety, effectiveness and comfort.	n.a.
		Preference	87% of participants (26/30) preferred the C-Brace over the SCO. Only 13% preferred the SCO.
Satisfaction – Quality of Life	Quality of life (EQ-5D-5L utility score)	Significantly higher EQ-5D-5L utility score by 27.2% with C-Brace (p < 0.001):	++
		C-Brace (PP)	SCO (PP)
		0.880 ± 0.106	0.692 ± 0.296
	EQ-VAS health	Significantly higher EQ-VAS health score by 21.6% with C-Brace (p = 0.002):	++
	C-Brace (PP)	SCO (PP)	
	76.3 ± 17.4	63 ± 21.8	

Category	Outcomes	Results for C-Brace vs. SCO (PP group)	Sig. ¹				
Satisfaction - Participation	Patient-Specific Functional Scale (PSFS)	Significantly improved PSFS score by 4.1 points with C-Brace (p < 0.001):	++				
		<table><tr><th>C-Brace (PP)</th><th>SCO (PP)</th></tr><tr><td>7 ± 2.6</td><td>2.9 ± 1.8</td></tr></table>	C-Brace (PP)	SCO (PP)	7 ± 2.6	2.9 ± 1.8	
		C-Brace (PP)	SCO (PP)				
		7 ± 2.6	2.9 ± 1.8				
Activities focused on mobility challenges encountered in daily life, such as:	n.a.						
<ul style="list-style-type: none">• Playing with children in the garden• Walking without walking aids• Walking downstairs and upstairs• Walking on ramps etc. <p>The improvements with C-Brace were considered clinically relevant, as the activities are difficult to manage with SCO from the participants own perspective.</p>							

¹ no difference (0), positive trend (+), negative trend (–), significant (++/–), not applicable (n.a.); significance set at $p < 0.05$; trends set at $0.1 > p > 0.05$.

Author's Conclusion

"This multicenter and international randomized crossover clinical trial asked 38 community ambulators with a quadriceps insufficiency to test and compare 2 knee-ankle-foot orthoses: their SCO and the C-Brace. Our results show that the C-Brace significantly improved mobility, endurance, confidence, participation, satisfaction, psychological adjustment, and quality of life in this population. Moreover, the C-Brace led to a decrease in the use of walking aids when walking outdoors, even though safety has been reported as the most important satisfaction criterion for participants. In all, community ambulators requiring the use of KAFO for walking could greatly benefit from the use of the C-Brace orthosis to improve their outdoor mobility and facilitate completion of daily activities. Further studies including people with bilateral quadriceps insufficiency are an interesting prospect to assess the possible advantages of the C-Brace for this population." (Genêt et al., 2026)

Author's Affiliation

^a Department of Physical Medicine and Rehabilitation, Neuro-Orthopedics Unit (UPOH), Raymond Poincaré Hospital, Garches, U1179 END-ICAP, INSERM, UFR Simone Veil-Santé, Versailles Saint-Quentin-en-Yvelines University, Paris Saclay University, Montigny-Le Bretonneux, IPS Foudation, Paris, France

^b Katholisches Klinikum Koblenz-Montabaur, Rudolf-Virchow-Straße 7, Koblenz 56073, Germany

^c CRRF La Chataigneraie, 48 rue de la Convention, Paris 75015, France

^d Centre Louis Pierquin, Institut Régional de Médecine Physique et de Réadaptation de Nancy, UGECAM Nord-Est, 75 boulevard Lobau 54042, Université de Lorraine, DevAH, Nancy F-54000, France

^e Centre de Rééducation Fonctionnelle de Salies du Béarn, 3 boulevard Saint-Guily, Salies-de-Béarn 64270, France

^f CSSR LADAPT, 14 allée de la terrasse BP 70011, Thionville 57101, France

^g Hôpital d'Instruction des Armées Laveran, 34 boulevard Alphonse Laveran, Marseille 13384, France

^h Hôpital National D'Instruction des Armées Percy, 101 avenue Henri Barbusse BP 406, Clamart 92141, France

ⁱ Centre Jacques Calve - Fondation Hopale, 72 esplanade Parmentier, Berck 62600, France

^j Pôle de Readaptation de Cornouaille, 61 Rue de Tregunc, Concarneau 29900, France

^k CRRF Hôpital Leopold Bellan, 16 rue de l'Aqueduc, Paris 75010, France

^l Clinique SMR Le Mont Veyrier, 475 route des Menthonnex, Argonay 74370, France

^m Centre de readaptation La Tourmaline, 31 boulevard Salvador Allende, Saint-Herblain 44800, France

ⁿ Centre Bouffard Vercelli, 334 rue Diego Velasquez, Perpignan 66000, France

^o Service de Chirurgie Vasculaire, CHU Hôpital de Rangueil, 1 avenue du Professeur Jean Poulhes, Toulouse 31400, France

^p Institut Robert Merle d'Aubigne, 2 Rue Emilion Michaut et Lucien Rabeux, Valenton 94460, France

^q Universitätsmedizin Göttingen, Robert-Koch Strasse 40, Göttingen 37075, German

©2026, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.