

Harmony vs other socket systems

Limb Volume Fluctuation

Major Findings

With vacuum assisted socket system (VASS) compared to other socket systems:

→ **Prevents volume loss of the residual limb during activity**


Residual limb volume increased 3.7% with VASS

Residual limb volume decreased 6.5% with suction socket system (SSS)

→ **Limb volume change is more positive than with pin suspension system (PSS) during activity**

→ **Net gain of limb volume is achieved with under-sized, neutral and over-sized socket design after activity**

Increased stump volume with VASS after walking



Volume of residual limb before and after 30 minutes of walk was measured for vacuum-assisted socket system (VASS) and suction socket system (SSS). (Board et al., 2001)

Clinical Relevance

One of the most important benefits of the VASS is that limb volume is maintained throughout the day. With other suspension systems, the residual limb loses volume as fluid gets pushed out of the limb during stance phase. As a result, the socket fit is less consistent and users need to remove prosthesis multiple times/day while attempting to manage volume changes with additional sock layers.

Summary

VASS shifts the fluid balance in the residual limb to one of maintenance or a slight gain by driving less fluid out of the limb during stance phase and drawing more fluid in during swing phase. An increase in the residual limb volume by 3.7% after 30 min of walking was observed with VASS, whereas the residual limb volume decreased by 6.5% with the suction socket system (Board et al., 2001). A case report demonstrated that 2 out of 3 subjects showed an increase in limb fluid volume during activity by 1.2%, respectively 0.4%. Moreover, 3 out of 3 subjects showed a more positive limb volume change during activity compared to pin suspension system (PSS) (Sanders et al., 2011).

A study investigating the effect of different socket sizes measured a change of the fluid balance of the residual limb towards a net gain for all socket sizes during walking. The improved fluid balance is independent of socket size. Nevertheless a tight

socket is crucial to make the VASS work as well as for safety reasons. (Goswami et al., 2003).

Beil et al. (2002) proposed that drawing more fluid in based on the 27% increase of negative pressure in swing phase is probably most responsible for volume maintenance. The drop in pressure is hypothesized by Street (2006) to be the result of the anchored liner. Thus tissues elongate, tissue pressure drops and therefore fluid is drawn into the limb.

References of summarized studies

Board, W. J., Street, G. M., & Caspers, C. (2001). A comparison of trans-tibial amputee suction and vacuum socket conditions. *Prosthetics and Orthotics International*, 25(3), 202–209. doi:10.1080/03093640108726603

Goswami, J., Lynn, R., Street, G., & Harlander, M. (2003). Walking in a vacuum-assisted socket shifts the stump fluid balance. *Prosthetics and Orthotics International*, 27(2), 107–113.

Kahle, J. T., Orriola J.J., Johnston W., & Highsmith, M. J. (2014). The effects of vacuum-assisted suspension on residual limb physiology, wound healing, and function: A systematic review. *Technology & Innovation*, 15(4), 333–341.

Klute, G. K., Berge, J. S., Biggs, W., Pongnumkul, S., Popovic, Z., & Curless, B. (2011). Vacuum-Assisted Socket Suspension Compared With Pin Suspension for Lower Extremity Amputees: Effect on Fit, Activity, and Limb Volume. *Archives of physical medicine and rehabilitation*, 92(10),

Sanders, J. E., Harrison, D. S., Myers, T. R., & Allyn, K. J. (2011). Effects of elevated vacuum on in-socket residual limb fluid volume: case study results using bioimpedance analysis. *Journal of rehabilitation research and development*, 48(10), 1231–1248.

Street, G. M. (2007). Vacuum Suspension and its Effects on the Limb. *Orthopädie-Technik*, 04.

© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.