

Benchmarking the performance of Reduced Representation Methylation Sequencing (RRMS)

To benchmark **RRMS** performance, we performed human **RRMS** on five replicates of a metastatic melanoma cell line and its normal pair for a male individual (COLO829/COLO829_BL) and a triple-negative breast cancer cell-line pair (HCC1395/HCC1935_BL), and compared the results with those obtained using Reduced Representation Bisulfite Sequencing (**RRBS**). **RRBS** is a method used to obtain genome-wide methylation analysis without the need to sequence the whole genome: CGIs are enriched through restriction digest and size selection followed by bisulfite conversion of methylated bases and next-generation sequencing. The bisulfite workflow is complex, time consuming and is also imprecise, as it does not specifically target any promoter region.

Library preparation

Genomic DNA was extracted from human cells grown in culture using the QIAGEN Genra Puregene Cell Kit (for gDNA extraction from blood or tissues, please see the [Nanopore Community](#) for more details) and quantified using the Qubit fluorimeter. For **RRBS**, samples were sent to Diagenode for processing and sequencing. For **RRMS**, sample preparation and sequencing was performed in-house. Firstly, DNA was sheared using a Covaris g-TUBE (Fig. 1).




Figure 1. Fragment-length distribution of **RRMS** libraries after g-TUBE fragmentation.

Fragmented DNA was then prepared for sequencing using the Ligation Sequencing Kit, loaded on MinION flow cells and run on a GridION MK1. Off-target reads are ejected from the pore after approximately 1 second of sequencing and consequently have a read length of ~500 bp. On-target reads are accepted and sequence as normal, generating read lengths ~6 kb (as a result of the g-TUBE shearing, Fig. 2a). To optimise output, flow cells were washed twice, after 30 and 64 hours, and fresh library was re-loaded each time (Fig. 2b).




Figure 2. a) Read-length distribution and b) output from a typical **RRMS** run

Basecalling and analysis

Basecalling and calling of 5mC modifications was performed using *guppy 6.1.5* (dna_r9.4.1_450bps_modbases_5mc_cg_sup) and output as unmapped BAM files in batches of 4,000 reads. BAM files were concatenated using *samtools* and mapped to the human reference (hg38) using *minimap2* maintaining per read modification information.

Modbam2bed was used to compute 5mC frequencies for all genomic CpG positions. Methylation frequencies from forward and reverse strands were aggregated enabling *modbam2bed --aggregate* option. The aggregated output was used for downstream analysis. All CpGs with < 10 overlapping reads discarded using a custom script. Finally, *DSS* was used to compute differentially methylated regions between tumour and normal pairs. CNV profiles were computed from ejected reads using *QDNASeq*.

Please be aware that a new chemistry together with a new basecaller

(<https://github.com/nanoporetech/dorado>) has been made available starting from early 2023. As a result, the recommended workflow for basecalling has been updated.

Modbam2bed has been also replaced by *modkit* (<https://github.com/nanoporetech/modkit>).

For more detailed information, please refer to the [RRMS protocol](#) accessible within the community resources. **This upgrade is not expected to influence results showed in this document.**

Results

RRMS generated ~35x average coverage across target regions and results in high confidence methylation calls for 7.3 - 8.5 million CpGs per sample, covering >90% of targeted promoters, CGIs, shores and shelves (NB if >50% of CpGs within a given feature were called then the feature is classified as being covered), and coverage across target regions is uniform and reproducible. **RRBS**, in comparison, yields 1.7 - 2.5 high-confidence calls per sample (Fig. 4).




Figure 4. Typical performance of **RRMS** run. Panel a) shows the average -fold coverage obtained for on-target ROI is >30x. Panel b) shows the number of CpG sites where the methylation status could be determined with high confidence (i.e. >10 spanning reads).

Methylation frequencies called by **RRMS** and **RRBS** are highly similar for CpGs covered by both technologies ($R > 0.967$) (Fig. 5).




Figure 5. Comparison of CpG calls between **RRMS** and **RRBS** shows that where a CpG is covered by both technologies the methylation calls are highly correlated.

With **RRMS** we were able to detect ~62 Mb of differentially methylated regions (DMRs) between tumour and normal pairs, of which a high proportion overlapped with cancer census genes. In comparison **RRBS** yielded ~20 Mb of DMRs. When investigating the DMRs further we identified that the de-novo methyltransferase DNMT3A promoter is clearly methylated in both tumours but not in their normal counterparts. These methyltransferases are essential for establishing and maintaining normal levels of methylation, and so their dysregulation can contribute to cancer development. Therefore, **RRMS** can provide key information for tumour characterisation as well as methylation status (Fig. 6).




Figure 6. Differences in methylation between samples can be determined using **RRMS**. The methylation status of a 6.5 kb region of chromosome 2 is compared for two tumour/normal cell line pairs and clear differences in the methylation patterns can be observed between the normal (BL) and tumour samples.

Finally, as an additional benefit of **RRMS**, it is possible to call copy-numbers variants (CNV) across the whole genome without any additional sample preparation (i.e. you get CNV calling “for free”), by using the reads which are ejected by AS during the course of an **RRMS** run. We observe that the calls made using the ejected reads from **RRMS** are highly correlated to the calls made using short-read WGS and ONT WGS (Fig. 7).




Figure 7. CNV calling using **RRMS**. a) The reads which are ejected during the course of an **RRMS** run can be used to call copy-number variation across the whole genome. b) The CNV calls made using **RRMS** are highly correlated to the calls made using short-read WGS and ONT WGS.

Conclusion

RRMS is a method that combines direct methylation detection and adaptive sampling to generate >7 million high-confidence CpG calls, covering all annotated CpG islands, shores, shelves and annotated promoters within the human genome, with minimal upfront library preparation. It enables the detection of differentially methylated regions between samples as well as whole-genome CNV calling. In addition, it can be extended to include custom targets like marker SNPs or SVs. **RRMS** therefore provides key information with potential application to tumour characterisation and methylation status. Combined with its ease of use and ability to scale to a high number of samples, **RRMS** is well suited to investigating methylation differences in large cohorts, as well as providing deeper insights into the mechanisms behind diseases such as cancer and monitoring tumour progression.