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
Figs. 1&2 A 1.2mL unit dose syringe is twisted firmly onto the IndiSpense syringe. While supporting plunger of unit dose syringe with the non-dominant hand, depress plunger of IndiSpense and fill unit dose syringe to desired level.

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Description:

Ultra-Etch is a 35% phosphoric etchant solution with an optimum viscosity. The flowability allows precise placement, including occlusal grooves, yet is viscous enough to prevent migration. Ultra-Etch is self-limiting in its depth of etch (average depth of 1.9μ with 15 second etch). This etchant contains no glycerin and is designed for rapid and complete removal upon rinsing.

Indications for Use:

Used for etching enamel and dentin prior to applying bonding adhesives and placing composite restorations or dental sealants.

Directions for Use:

- For 1.2 syringe
 - Remove Luer cap.
 - Securely attach working tip of choice (Blue Micro™ or Inspiral™ Brush tip).
 - Verify etchant flow prior to applying intraorally.
 - Verify the color and consistency of the product flow before applying. Initially a small amount of clear liquid may express from the syringe; this is normal.
- For IndiSpense™ Syringe (Figure 1 & 2)
 - Remove Luer cap from IndiSpense syringe.
 - Attach 1.2 syringe to the male threads of the IndiSpense syringe.
 - Depress IndiSpense plunger while guiding 1.2mL syringe plunger to desired fill.
 - Separate syringes and re-cap IndiSpense syringe with Luer cap.
 - Securely attach working tip (Blue Micro™ or Inspiral Brush tip) to 1.2mL syringe.
- Direct placement instructions

Precautions and Warnings:

- Isolate tooth (teeth).
 - Use pumice, disk or diamond bur on uncut enamel.
 - Rinse and dry prepared area.
 - Verify etchant flow prior to applying intraorally.
 - Apply etchant to enamel and dentin (15 seconds)
 - Rinse thoroughly, dry and proceed per adhesive manufacturer's instructions.
 - Restorative phosphoric acid cleansing
 - Apply etch to inside of porcelain veneer or crown for 5-10 seconds.
 - Rinse and dry
- Precautions and Warnings:
- Carefully read and understand all instructions before using.
 - Keep out of reach of children.
 - If product comes in contact with any soft tissue, IMMEDIATELY rinse area with copious amounts of water.
 - Near pulp exposures should be treated with a protective base prior to placing etchant.
 - Re-cap, disinfect and wipe syringe with an intermediate level disinfectant between uses. If syringe cover is used, remove tip, re-cap, and discard syringe cover
 - Dispose of tip and empty syringe properly.
 - Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips. To avoid cross contamination, do not re-use tips.
 - Store at room temperature.
 - Do not use after expiration date noted on container.
- NOTE: For SDS and additional information about using Ultra-Etch, please go to www.ultradent.com

Ultra-Etch™
35% Phosphoric Acid

Manufactured by:

Ultradent Products, Inc.
505 West Ultradent Drive (10200 South)
South Jordan, Utah 84095 USA

Manufactured in the USA from globally sourced materials

 **ULTRADENT**
PRODUCTS, INC.

800.552.5512 |  ultradent.com

© 2020 Ultradent Products, Inc. All Rights Reserved. 1009279AR01 111020