IFU of UltraSeal XT Hydro

Description:

UltraSeal XT Hydro is a 53% filled, light-cure, radiopaque, methacrylate-based, thixotropic resin sealant. It is hydrophilic and has a self-adhesive quality. UltraSeal XT Hydro chemistry provides the option of a visual verification for marginal retention with the use of a UV light, upon placement and at recall visits. The UV light is not included in the UltraSeal XT Hydro kit.

Accessories list

No.

Indications for Use:

Use UltraSeal XT Hydro for prophylactic sealing of pits and fissures.

Contraindications

Confirm that the patient will not be allergic to this product.

For patients or users with allergy concerns, refer to common product allergen document available at www.ultradent.com. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.

Warnings and Precautions

- Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin, if contact occurs, thoroughly rinse skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, have patient consult their physician.
- For daily use, store at room temperature. Refrigerate bulk storage. Bring refrigerated products to room temperature before using.
- Redirect overhead light, or use UV protective light shield, to prevent premature polymerization.
- Re-cap immediately following use to avoid polymerization in syringe.
- Tips are single use. To avoid cross-contamination, do not re-use tips.

Directions for Use:

UltraSeal XT Hydro sealant –Pit and fissure sealant

Preparation:

- Attach Inspiral™ Brush tip or Blue Micro™ tip to Ultra-Etch 35% etchant.
 - Attach Inspiral Brush tip to UltraSeal XT Hydro sealant syringe. If desired, the Inspiral Brush tip bristles length can be adjusted by grasping and moving inward or outward.
 - Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.
 - To prepare tooth prior to sealants use one of the following methods:
 - an air abrasion unit with appropriate particles
 - a thin diamond fissure/prep bur
 - air polishing prophylaxis system
 - Consepsis™ Scrub slurry
 - Prophy paste
 - Flour of pumice
- NOTE: Use high magnification when cleaning grooves.
- Apply etch to fissures for 20 seconds.
 - Thoroughly rinse.
 - Lightly air dry the tooth to remove visible moisture. Do not desiccate the tooth.

Placement:

- Express a small drop of UltraSeal XT Hydro sealant into the grooves. Work the resin down into the grooves with the brush tip.
 - Light cure with VALO™ LED curing light on Standard Power for 10 seconds or 3 seconds on Xtra Power mode.
 - Use an explorer to ensure pits and fissures are sealed and adjust occlusion if necessary.
- NOTE: As an option, a black light may be used for visual verification of marginal retention. Either the VALO curing light with the Black Light Lens specifically designed to filter for UV light, or a small black light, is available through Ultradent. For best results, redirect overhead dental light when using black light.

Processing

To avoid cross-contamination, wipe syringe with an intermediate level disinfectant between uses and use a disposable barrier sleeve. If a barrier sleeve is not used, syringe should be treated as single use.

Storage condition

Store in room temperature.

UltraSeal XT™ Hydro

Manufactured by:
Ultradent Products, Inc.
505 West Ultradent Drive (10200 South)
South Jordan, Utah 84095 USA
Manufactured in the USA from globally sourced materials

ULTRADENT
PRODUCTS, INC.

8 0 0 . 5 5 2 . 5 5 1 2 |  [ultradent.com](http://www.ultradent.com)

© 2020 Ultradent Products, Inc. All Rights Reserved. 1009467AR01 121020

Manufacturing date

Refer to product label

Shelf life

18 months

Explanation for the symbol on the label

Catalog number

Lot number

Do not re-use to avoid cross contamination

See instructions

Use by date

For professional use only

Manufacturer

Revised date:

2020/7/16

Registrant/Manufacturer: Ultradent Products, Inc.
Address: 505 West 10200 South, South Jordan, UT 84095, USA
Contact: +1-888-863-5883

Agent/altersales: Ultradent Products, China
Address: Room 102, No.34 building, No.67 Gongye Dadao Bei Road, Haizhu District, Guangzhou 510250, P.R.China
Contact:020-84156506

Registration number: 国械注进20162632765
Product Standard number: YZB/USA 0684-2016
Product name: Pit and Fissure Sealant and Conditioner Model/.

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