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Porcelain Etch & Silane Manual

Description:

Ultradent Porcelain Etch etchant is a buffered, gelled hydrofluoric acid to be judiciously used for intraoral or extraoral etching of porcelain. Ultradent Silane is a mono-component solution for application to porcelain following etching with Porcelain Etch etchant.

Indications:

Porcelain Etch etchant intends for etching the porcelain material intraorally and extraorally. It is also used for pretreating the porcelain before bonding.

Contraindications:

Confirm that patient has no known allergies to treatment materials.

Important Note: Hydrofluoric acid gives off a vapor that irritates respiratory passageways and on a chronic basis can damage them. Ultradent Porcelain Etch etchant is buffered, gelled, and dispensed in a way that minimizes these problems. Instructions must be followed to ensure removal of any harmful material. Silane is for bonding after etching in preparation for bonding resin. If proper precautions are followed, porcelain may be safely etched intraorally for bonding as follows:

Procedure:

1. Freshen the porcelain surface to be repaired with a high speed diamond bur. If a previous repair is present, remove all old composite to fresh porcelain surface.
 2. Remove luer lock caps from syringes. For optimum control, hold syringes in the palm of the hand.
 3. Isolate/protect adjacent tissues and restorations by applying a bead or layer of EtchArrest™ acid barrier/neutralizer medium with the Black Micro™ tip (Fig. 1). An alternative isolation is rubber dam (DermaDam™) and rubber dam sealant (OraSeal™).
 4. Twist an Inspirail Brush tip securely onto syringe. The Inspirail™ Brush tip ensures controlled flow from the tip.
 5. Before applying, dispense a small sample of etchant from tip onto a mixing pad to ensure evenness and to become familiarized with flow rate. If there is excess resistance to flow, DO NOT proceed. Attach a new tip and test again before using on patient. Dispose of removed tip immediately (according to government regulations) to avoid contact.
 6. Express an even coating of Porcelain Etch etchant onto prepared porcelain (Fig.2) using the Inspirail Brush tip. Etch porcelain 90 seconds. (note Fig.3 and 4 showing before and after SEMs)
 7. Assistant should hold vacuum tip adjacent to the surface while the dentist dispenses the Porcelain Etch etchant. This ensures removal of any irritating vapors while acid reacts on porcelain (90 seconds). To neutralize the Porcelain Etch etchant, blend EtchArrest acid barrier/neutralizer into the etchant until the amber color of the etch is no longer identifiable.
 8. Remove neutralized acid by placing the vacuum tip next to the surface and then rinsing with water spray. Rinse thoroughly and dry. Porcelain should have a frosted appearance. Repeat etching if not.
 9. Apply Ultra-Etch™ etchant for 5 seconds to remove porcelain salts and debris formed by hydrofluoric etchant. Rinse and dry.
 10. Apply Silane with a Mini Brush tip. Let evaporate for 60 seconds. If not completely dry after 60 seconds, air dry with oil-free air.
 11. Place bonding resin according to product instructions.
- Important Note:** Testing shows highest bonding values when etched porcelain is silanated with Ultradent Silane followed by Peak™ Universal Bond bonding agent in the following manner: Apply Silane to porcelain using Mini Brush tip. Gently brush the Silane into the porcelain for 60 seconds using small strokes. Air dry and proceed with application of Peak Universal Bond bonding agent.


Fig. 1 Protective bead of EtchArrest is placed on soft tissues adjacent to the tooth surfaces.


Fig. 2 Carefully express an even coating of Ultradent Porcelain Etch to prepared porcelain.


Fig. 3 SEM of diamond-ground porcelain.


Fig. 4 Note increased surface area and undercuts afforded by etch with Porcelain Etch for 90 seconds.

Precautions:

1. Caution: Chemical may set in tip after a few hours. To avoid clogging, remove tip and re-cap with Luer Lock Caps between each use. Always verify material flow prior to applying intraorally. If resistance is met, replace tip and re-check. This does NOT apply to dual-barrel syringes and mixing tips.
2. Operator, assistant, and patient must wear protective eyewear.
3. Replace luer lock caps on syringes when not in use. Tips are disposable.
4. We recommend Ultradent Syringe Covers for a protective barrier to prevent cross-contamination.
5. Clinician and patient must wear protective eyewear with side shields when handling or applying this product.

General Precautions:

1. Review instructions, precautions, and SDS before beginning treatment. Use only as directed..
2. For Professional use only.
3. Keep products out of heat/sunlight.
4. Avoid skin exposure to resins.
5. Isolate strong chemicals to area of treatment.
6. Test flow of materials from syringe and tip before using intraorally.
7. Never force syringe plungers.

Please visit www.ultradent.com for product technical instruction.

Ultradent syringes have an expiration date stamped on the side of the syringe consisting of one letter and three numbers. The letter is a lot number used for manufacturing purposes and the three numbers are the expiration date. The first two numbers are the month, and the third number is the last number of the year.

Transport and Storage condition:

Store in room temperature.

Manufacturing date / expired date:

Find it in the product label

Explanation for the symbols in the product label:

See instructions

Do not re-use to avoid cross contamination.

Catalog number

Use by date

Lot number

Manufacturer

For professional use only

Recycle

Revision date:

June 29, 2021

Registrant/Manufacturer: Ultradent Products, Inc.

Address: 505 West Ultradent Drive (10200 South), South Jordan, Utah 84095, USA
Contact: +1-888-863-5883

Agent/aftersales: Ultradent Products, China

Address: Room 102, No.34 building, Fenghuang Chuangyi, No.67 Gongye Dadao Bei Road, Haizhu District, Guangzhou 510250, P.R.China

Contact: 020-84156506

Registration number/ Product technical requirement Number: 国械注进20162172247

Product Name: Porcelain Etch and Silane

Model: 405, 406, 410, 406-1, 410-1 (details see below table)

Model	Description	Configuration
406	Porcelain Etch Syringe 2pk	2x1.2mL
405	Porcelain Etch kit	2x1.2mL (Porcelain Etch)
		2x1.2mL (Saline)
		20 x Black Mini Brush tips
		20 x Inspirail Brush tips
410	Silane Syringe 2pk	2x1.2mL
406-1	Porcelain Etch Single Kit	1x1.2mL+2xInspirail Brush tips
410-1	Saline Single Kit	1x1.2mL+2xBlack Mini Brush tips

Porcelain Etch & Silane

Manufactured by:
Ultradent Products, Inc.
505 West Ultradent Drive (10200 South)
South Jordan, Utah 84095 USA
Manufactured in the USA from globally sourced materials

ULTRADENT
PRODUCTS, INC.

800.552.5512 | ultradent.com
© 2022 Ultradent Products, Inc. All Rights Reserved. 1009280AR03 052322