Jiffy Diamond Strips and Proximal Saw

Instructions For Use

The user instructions for your device are provided in both electronic and printed format. If you wish to have a full printed copy of instructions, please call 800.552.5512.

Scope of this information:

The following safety and reprocessing information applies to the finishing strips offered by Ultradent.

Intended use:

Jiffy Finishing Strips and Saw are for multiple uses in dental treatments. They are used for interproximal polishing, reduction, and cleaning as well as finishing of dental enamel and dental materials such as composites, metal alloys, ceramics and dental resins.

They are intended for manual use by trained technicians and dentists.

Transport, storage, and treatment of new instruments prior to first use:

There are no particular conditions necessary for transport.

Please store all original packaged instruments in a clean and dry place at room temperature.

New instruments must be cleaned, disinfected, and sterilized according to preparation instructions before use.

Disposal:

All our instruments are free of hazardous substances. Only the requirements for protection against infection and contamination need to be considered for disposal.

Instructions for use:

1. Users and assistants are advised to wear mask, goggles, and gloves. The treatment unit and environment are also to be appropriately disinfected.
2. Strips must be cleaned, disinfected, and sterilized before each use.
3. Processing extraneous materials in the mouth, such as filling materials, can release nanoparticles of these. Depending on the starting material, these may be bioactive. Thus suction and possibly other protective measures such as a dental rubber dam, etc. are recommended.
4. Please work with as little pressure as possible and do not twist the instrument.
5. If the instrument jams remove the instrument carefully without twisting it and check for damage before using it again.
6. Dull and damaged instruments must no longer be used. Please check the instruments prior to each use. Signs of damage are kinks or tears in the strip.

Note regarding the obligation to report:

According to the relevant statutory regulations, serious incidents which occur in conjunction with a product from us must be reported to Ultradent Products, Inc. at 505 West Ultradent Drive (10200 South), South Jordan, Utah 84095 USA and to the responsible national authorities of the country in which the user works.

In terms of the law, incidents are subject to mandatory reporting if the use of the instrument directly or indirectly caused or could have caused the death or serious damage to health of the patient, user or another person (source: the German Safety Plan for Medical Devices [MPSV], Section 2).

Reprocessing Instructions

Basic legal and regulatory information:

These preparation instructions were produced in accordance with the EU medical device regulation (MDR/2017), DIN EN ISO 17664: 2004 and the recommendations of the commission for hospital hygiene and infection prevention at the Robert Koch Institute (RKI): "Infection prevention in dentistry – requirements for hygiene" (2006) and "Hygiene Requirements for the Reprocessing of Medical Devices" (2012).

Risk assessment and classification of the product groups according to the guideline cited above:

Due to their area of application and abrasive surface finishing strips are to be classified as semi-critical B for polishing, with particular requirements for preparation. They are technically suitable for preparation in most types of ultrasonic baths, cleaning and disinfection equipment, and autoclaves. Cleaning and disinfection equipment can be insufficient for sole use in cleaning under some circumstances due to the aforementioned special considerations. (Please refer to critical process steps as well.) Our instruments are not suited for chemical or hot air sterilization.

Critical process steps:

The post-cleaning check is to be considered a critical process step (critical control point). The rough surfaces from diamond coating or serration of instruments result in a risk that tissue remains and other impurities are not removed completely. This can reduce the effectiveness of disinfection and sterilization. Therefore, drying of impurities is to be avoided, and the post-cleaning check should be performed with particular care using aids such as a magnifier or microscope.

Expertise and protection of the person performing preparation:

The strips must only be prepared by persons with the necessary expertise. The requirements for this are defined in the relevant regulations and statutes.

Gloves, mask, safety glasses, and protective clothing must be worn to avoid contact with potentially infectious material.

Suitable media for cleaning and disinfection:


Only certified media produced according to the relevant standards may be used. The manufacturers of these cleaning and disinfection materials check whether a product is fundamentally usable for rotary instruments. Thus all media for which the manufacturer indicates a corresponding approval can be used. They should contain corrosion protection agents. Cleaning and disinfection media, which contain strong acids and bases or hydrogen peroxide and sodium hypochlorite, are not suited for use.

Preparation of used instruments:

All product groups must be prepared prior to each use according to the following steps.

Cleaning and Sterilization (Reprocessing) of Used Devices

<p>Scope</p>	<p>These instructions are applicable to all Ultradent Jiffy Diamond Strips and Interproximal Saws and are classified as semi-critical B. After every use, the devices need:</p> <ol style="list-style-type: none"> 1. Cleaning 2. Disinfection 3. Sterilization 4. Release 5. Storage <p>NOTE: Personal protection equipment including safety glasses and gloves should be used when reprocessing these devices.</p>
<p>1. Cleaning</p>	<p>Manually rinse devices under running water. Use a brush to remove any tissue and rinse again.</p>
<p>2. Chemical Disinfection</p>	<p>All Jiffy Diamond Strips and Interproximal Saws can be disinfected using a non-corrosive pH-neutral chemical disinfection bath. Follow manufacturer's directions for specific contact times and concentrations. Strips and saws can also be cleaned by using a disinfection solution in an ultrasonic bath (max. 45°C, 113F). To avoid damage, use holding devices to prevent the devices from touching the surrounding materials. After disinfection, rinse the devices with distilled water and dry thoroughly.</p> <p>NOTE: Cleaning and disinfection machines or chemo-thermal disinfection machines are NOT recommended.</p> <p>Control: Prior to sterilization, all devices must be checked, preferably under magnification, for cleanliness and proper functioning. Discard worn, bent, or damaged instruments. If necessary, repeat the cleaning and disinfection cycle.</p>
<p>3. Sterilization</p>	<p>All devices which have come into contact with mucous membrane, blood, or other tissues in the mouth, must be sterilized after disinfection. Put devices into suitable trays or posts then autoclave using the cycle pre-programmed to sterilize the devices at 134°C, 273 F /2.1 bar pressure and for an exposure time of 5 minutes.</p>
<p>4. Release</p>	<p>Clinics working with a QM (Quality Management) system must document the preparation and the devices must be released for use or storage.</p>
<p>5. Storage</p>	<p>The devices should be stored at room temperature in a clean, dry environment. Used devices must be stored in a manner that precludes contamination of man or instrument. All devices need to be processed according to the above instructions prior to the next use.</p>

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