Thermo Clone™

VPS Impression Material

English	3
Français	9
Español	15
Português	21

- Thermo Clone™ Fast Set Putty**
- Thermo Clone™ Fast Set Heavy Body**
- Thermo Clone™ Heavy Body**
- Thermo Clone™ Fast Set Medium Body**
- Thermo Clone™ Medium Body**
- Thermo Clone™ Fast Set Bite Registration**
- Thermo Clone™ Fast Set Light Body**
- Thermo Clone™ Light Body**
- Thermo Clone™ Fast Set Superlight Body**
- Thermo Clone™ Fast Set Clear Bite Registration**
- Thermo Clone™ Alginate Alternative**

DESCRIPTION/INDICATIONS:

Thermo Clone VPS Impression Material is a super hydrophilic vinyl polysiloxane impression material and is an exciting new advancement in prosthodontics and restorative dentistry.

Thermo Clone VPS Impression Materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation. These properties promote immediate inflow to the primary crowns, cast support or preparation limits even in the moist oral environment. Thermo Clone VPS impression materials are available in 1:1 50 ml and 5:1 380 ml cartridges.

1. Thermo Clone VPS Fast Set Putty is to be used as preliminary materials for:
 - Two-step putty-wash impression technique
 - One-step putty-wash impression technique
 - Two-step putty-wash impression technique using a foil (plastic putty spacer)
 - One-step putty impression technique for forming functional periphery
2. Thermo Clone VPS Impression Material Heavy Body in both Standard and Fast Set are to be used as heavy body materials for:
 - One-step impression technique (simultaneous technique) using single or dual viscosities
 - Two-step impression technique using dual viscosities
 - Functional impressions
3. Thermo Clone VPS Impression Material Medium Body in both Standard and Fast Set are to be used for:
 - Making impressions over fixed/removable restorations and implants (i.e.: transferring impression posts and bridge components).

4. Thermo Clone VPS Fast Set Bite Registration is used for:
 - Making intraoral or extraoral occlusal bite registrations for fixed or removable restorations and implants.
5. Thermo Clone VPS Impression Material Light Body in both Standard and Fast Set are used for:
 - Making impressions over fixed/removable restorations and implants (i.e.: transferring impression posts and bridge components)
6. Thermo Clone VPS Impression Material Superlight Body Fast Set is used for:
 - Making impressions over fixed/removable restorations and implants (i.e.: transferring impression posts and bridge components)
7. Thermo Clone VPS Fast Set Clear Bite Registration Material is a syringeable addition-curing elastomeric bite registration for:
 - Making accurate occlusal records. It can be used for making accurate occlusion registrations.
 - The transparency of the bite registration material enables direct visual inspection of the correct bite.
8. Thermo Clone VPS Alginate Alternative is addition-curing, elastomeric impression materials intended to be used as an alternative to traditional alginate materials (e.g., for anatomical impressions). Available in 1:1 50 ml cartridges suitable for the Ultradent dispensing impression gun and in 5:1 380 ml cartridges for use in most commercially available other automatic dispensing and mixing machines.

Indications/Techniques

Thermo Clone VPS Alginate Alternative is useful as alginate replacements for making impressions for:

- Opposing dentitions
- Fabricating temporary crowns and bridges
- Anatomical models
- Fabricating simple removable prosthetic restorations
- Producing models for the construction of trays for whitening solutions, bites and surgical guides.
- Orthodontic work
- Case study models
- Fabricating removable retainers and splints

Model Preparation

Impressions with Thermo Clone VPS Alginate Alternative may poured immediately after disinfection of the impression up to several weeks with standard dental stone (class III – V). Impressions can also be poured several times. A surfactant solution is not required.

Disposal

The empty cartridges / mixing tips / applicator tips and set material can be disposed in the household garbage.

PROCEDURE:

CARTRIDGES (1:1 50 ml)

1. Load the cartridge in the dispensing gun.
2. Immediately prior to taking the impression remove the cartridge cap and discard.
3. Gently squeeze the dispenser to bleed enough material from the cartridge to ensure even flow of the base and catalyst. (About ¼ inch of material).
4. Attach the mixing tip.
5. Squeeze the dispenser in order to extrude the product.
6. After use, store cartridges on end with the used mixing tips attached to protect the remaining materials from contamination.
7. Before using the cartridge again, remove the filled mixing tip and check the cartridge outlet openings to ensure that no polymerized materials have sealed the openings. Clear any clogged openings. Attach a new mixing tip and dispense a small amount of impression material until equal amounts are extruded at the same rate.

CARTRIDGES (5:1 380 ml)

1. Prior to first use unscrew the blue plug counter clock wise and pull it out, check the orifices for blockage.
2. Only use the dynamic mixing tip that comes with the new 380 ml cartridge.
3. The dynamic mixing tip must be tightened into the proper position before any impression material is dispensed. Make sure that you see a homogeneous mix extruding from the mixing tip.
4. When reusing a previously used cartridge make sure that the machine plungers are retracted and then advanced when inserting and beginning the dispensing process.
5. Follow the instructions for your 5:1 mixing/dispensing machine.

- At least 3 cm of the mixed material must be discarded each time the plunger is restarted or when using the cartridge in different mixing units.
- Prior to the next use remove the previously utilized dynamic mixer, turning it to the left. Ensure that the dispensing orifices are not blocked and then attach a new dynamic mixing tip.

TIPS:

Use Ultradent Impression Mixing Tips attached to the impression cartridge with intraoral tip to mix material and inject around preparation.

PLEASE NOTE:

- To ensure optimal impressions, Thermo Clone VPS Impression Material should be stored at 23°C (74°F). Working and setting time will be affected by deviation from this temperature.
- Ideally it should be kept at room temperature (18-25°C) or (64-77°F).
- Should gloves be worn, test them for compatibility prior to mixing with a sample of Thermo Clone VPS Impression Material to be mixed. Certain types of gloves can prevent the polymerization (e.g. latex gloves) or impair the adhesive bond (e.g. nitrile gloves). Using vinyl gloves is recommended.
- Avoid contaminating Thermo Clone VPS Impression Material with chemicals that could inhibit the set (e.g. acrylic and methacrylate residues, latex rubber and sulfur compounds).

Cleaning and Disinfection of Ultradent Impression Gun:

- Remove the cartridge.
- Remove any uncured paste using a cloth soaked in alcohol.
- The Ultradent Impression Gun is reusable if properly disinfected.
- Ultradent recommends exchanging the gun if there are visible traces of use.
- Disinfection can be carried out by immersing impressions in a commercially available cold sterilization solution.
- Ultradent recommends the use of glutaraldehyde-based solutions.

DISINFECTION:

The impressions should be rinsed under cold running water after removal from the mouth. Soaking for 30 minutes in liquid disinfectants such as 2% glutaraldehyde

(acidic), 0.5% - 1% sodium hypochlorite, 0.1% - 1.5% povidone iodine, halogenated phenol, or 0.13% glutaraldehyde (neutral) will not affect the surface or dimensions.

CASTING:

The impression may be poured immediately after removal from the mouth following disinfection. Surface tension will be reduced if the impression is briefly washed out with a surfactant. Rinse thoroughly afterwards with water.

Ultradent™ Impression Materials are compatible with all dental stones, epoxy resins, and polyurethanes.

CONTRA-INDICATIONS:

- There are no known contra-indications provided when the product is used as directed. There are no known harmful reactions or secondary effects on patients and/or dental personnel using vinyl polysiloxane impression materials.
- Do not use Thermo Clone VPS Fast Set Clear Bite Registration and Thermo Clone VPS Fast Set Bite Registration for conventional impressions. The high final hardness values achieved severely hamper withdrawal of the impression and may lead to inadvertent extraction of teeth.

IMPORTANT NOTES AND PRECAUTIONS:

The temperature of Thermo Clone VPS Impression Material affects working time. Ideally it should be kept at room temperature between (18°-23°C) or (64°-74°F). Even within those boundaries, changes in temperature will affect the working time: You may refrigerate to increase working and setting time. Clear any clogged openings and dispense a small amount of impression material until equal amounts are extruded at the same rate. This will ensure even dispensing of fresh base and catalyst pastes. For even more precise application, Ultradent™ Intraoral Tips are available.

- Certain chemical agents like superoxol or hydrogen peroxide may interfere with the set of Thermo Clone VPS Impression Materials, care should be taken to remove these residual materials for impression site. Latex gloves (or anything that contains sulfur) and hand lotions can inhibit or contaminate setting of all Thermo Clone VPS Impression Materials.
- Should gloves be worn, test them for compatibility prior to mixing with a sample of Thermo Clone VPS Impression Material to be mixed. Certain types of gloves

can prevent the polymerization (e.g. latex gloves) or impair the adhesive bond (e.g. nitrile gloves). Using vinyl gloves is recommended.

3. Proper care should always be taken to ensure good preparation, cleanliness, and moisture control including hemorrhage and saliva in order to achieve an adequate impression.
4. There is no bubbling or gasing with Thermo Clone VPS Impression Materials.
5. Disinfection can be accomplished with any appropriate disinfectant following the original manufacturer's procedures.
6. Multiple casts can be poured from an impression and Thermo Clone Impression Material is compatible with all gypsum products, electroplating and refractory casts for up to 15 days.
7. Temperature affects all set times. All testing for working and setting times was done at 74°F/ 23°C and/or 95°F/ 35°C.
8. Shelf life is at least 1.5 years.
9. Putty:

Do not wear latex gloves or hand lotions while mixing the putty. Wash hands and mix equal amounts of catalyst and base. Use the scoops to dispense equal amounts of both base and catalyst. Do not interchange the colored base or catalyst lids or scoops for the hand mixed putty. In order to ensure optimal adhesion of the two materials during the two step putty wash technique both impression steps must take place one immediately after the other. Close containers after use and ensure that the lids are not interchanged. Contamination of the base with catalyst in the container will make the material unusable. Knead the base and catalyst for 45 seconds until the material is homogeneous in color. Prior to making the wash impression the putty impression must be cleaned with water, dried, trimmed as needed.

CAUTIONS:

- Do not use Thermo Clone VPS Impression Materials and Thermo Clone VPS Bite Registration Materials after the expiration date.
- Do not leave residual bite registration material in the sulcus or oral cavity.
- Do not swallow Thermo Clone VPS Impression Materials and Thermo Clone VPS Bite Registration Materials. If swallowed: consult a medical doctor if problems arise or persist.
- Avoid contact with the eyes. If there is accidental contact with the eyes: rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

- The products do not normally cause allergic reactions; however, for sensitive persons, an allergic reaction to the product cannot be ruled out.
- Use of products containing active sulfur, aluminum chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide bite registration materials etc.) in conjunction with the Thermo Clone VPS Impression Material products will interfere with the setting reaction of the vinyl polysiloxane material. Use of these materials requires the area to be rinsed thoroughly to remove all residues.
- Do not allow the material to enter the sewerage system or waterways to avoid an environmental hazard.
- Avoid contact with clothing, as the material cannot be removed by dry cleaning.

COMPATIBILITY:

Thermo Clone VPS Impression Material is addition-reaction type silicone impression material. It may not be combined with standard or condensation type silicone materials or polyether materials.

TECHNICAL DATA:

GENERAL PRECAUTIONS:

1. For dental professional use only.
2. Review instructions, precautions, and SDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight.
4. Avoid skin exposure to resins.
5. Isolate strong chemicals to area of treatment.
6. Confirm that patient has no known allergies to treatment materials.
7. Test flow of materials from syringe and tip before using intraorally.
8. Never force syringe plungers.
9. Use Ultradent™ Syringe Covers and/or clean and disinfect syringes between patients.

Product Name	Mixing Ratio	Working Time at 23°C/74°F<	Working Time at 35°C/95°F<	Intraoral setting Time at 35°C/95°F>
Thermo Clone™ Fast Set Putty	1:1	2 min.	N/A	2 min.
Thermo Clone™ Fast Set Heavy Body	5:1	1 min.	N/A	2 min. 15 sec.
Thermo Clone™ Fast Set Heavy Body	1:1	1 min.	N/A	2 min. 15 sec.
Thermo Clone™ Heavy Body	5:1	2 min. 15 sec.	N/A	2 min. 30 sec.
Thermo Clone™ Heavy Body	1:1	2 min. 15 sec.	30 sec.	2 min. 30 sec.
Thermo Clone™ Fast Set Medium Body	5:1	1 min.	30 sec.	2 min. 15 sec.
Thermo Clone™ Fast Set Medium Body	1:1	1 min.	30 sec.	2 min. 15 sec.
Thermo Clone™ Medium Body	5:1	2 min. 15 sec.	1 min.	2 min. 30 sec.
Thermo Clone™ Medium Body	1:1	2 min. 15 sec.	1 min.	2 min 30 sec.
Thermo Clone™ Fast Set Bite Registration	1:1	15 sec.	N/A	55 sec.
Thermo Clone™ Light Body	1:1	2 min. 15 sec.	1 min.	2 min. 15 sec.
Thermo Clone™ Fast Set Light Body	1:1	1 min.	30 sec.	1 min. 15 sec.
Thermo Clone™ Fast Set Superlight Body	1:1	1 min.	30 sec.	1 min. 15 sec.
Thermo Clone™ Fast Set Clear Bite Registration	1:1	15 sec.*	N/A	45 sec.
Thermo Clone™ Alginate Alternative	1:1	1 min 15 sec.	N/A	1 min. 15 sec.
Thermo Clone™ Alginate Alternative	5:1	1 min 30 sec	N/A	1 min. 30 sec.

Note: The total working time and intraoral setting time should be strictly adhered to in order to ensure optimal impressions and bite registrations. Should the temperature of the material deviate noticeably from 23°C (74°F), the total working time and setting times will be affected.

Product Name	Total Setting Time>	Hardness (approx.) Shore	Linear Dimensional Change (approx.) %	Elastic-Recovery Test (approx.) %	Strain in Compression (approx.) %
Thermo Clone™ Fast Set Putty	4 min.	A66	-0.20	99	2.7
Thermo Clone™ Fast Set Heavy Body	3 min. 15 sec.	A62	-0.20	99.7	2.5
Thermo Clone™ Fast Set Heavy Body	3 min. 15 sec.	A62	-0.20	99.5	3.0
Thermo Clone™ Heavy Body	4 min. 45 sec.	A62	-0.20	99.5	3.0
Thermo Clone™ Heavy Body	4 min. 45 sec.	A62	-0.20	99.5	3.0
Thermo Clone™ Fast Set Medium Body	3 min. 15 sec.	A58	-0.20	99.7	3.5
Thermo Clone™ Fast Set Medium Body	3 min. 15 sec.	A58	-0.20	99.7	3.0
Thermo Clone™ Medium Body	4 min. 45 sec.	A58	-0.20	99.7	3.0
Thermo Clone™ Medium Body	4 min. 45 sec.	A58	-0.20	99.7	3.0
Thermo Clone™ Fast Set Bite Registration	1 min. 10 sec.	A85	-0.20	N/A	N/A
Thermo Clone™ Light Body	4 min. 30 sec.	A46	-0.20	99.3	3.5
Thermo Clone™ Fast Set Light Body	2 min. 15 sec.	A46	-0.20	99.3	3.5
Thermo Clone™ Fast Set Superlight Body	2 min. 15 sec.	A46	-0.20	99.3	3.5
Thermo Clone™ Fast Set Clear Bite Registration	1 min.	A62	-0.25	N/A	N/A
Thermo Clone™ Alginate Alternative	2 min. 30 sec.	A43	-0.20	99.5	7.2
Thermo Clone™ Alginate Alternative	3 min.	A45	-0.20	99.5	5

- Thermo Clone™** Putty prise rapide
- Thermo Clone™** Heavy prise rapide
- Thermo Clone™** Heavy
- Thermo Clone™** Medium prise rapide
- Thermo Clone™** Medium
- Thermo Clone™** Enregistrement occlusal à prise rapide
- Thermo Clone™** Light prise rapide
- Thermo Clone™** Light
- Thermo Clone™** Superlight prise rapide
- Thermo Clone™** Enregistrement occlusal transparent à prise rapide
- Thermo Clone™** Alternative à l'alginat

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Ultradent Products, Inc