## Opalescence™ Office Tooth Whitening System

**Description:** Opalescence™ Office is a dentist-controlled, dentist-applied, cosmetic in-office 6% hydrogen peroxide bleaching gel. The 2-barrel jet mixing syringe ensures maximum strength for each application and precise dosing of the activator. One barrel contains sodium fluoride and potassium nitrate, along with a unique chemical activator, which is combined with a second barrel containing hydrogen peroxide. The final hydrogen peroxide concentration is 6%.

Carefully read and understand SDS and all instructions before using.

**Ingrédients:**
Aufhellungs-Gel
Glycerin, Water (Aqua), Hydrogen Peroxide, Silica, Carbomer, PVP, Sodium Fluoride, Sodium Hydroxide, EDTA
Aktivator
Glycerin, Potassium Nitrate, Sodium Hydroxide, Carbomer, Iron Oxide
Combined
Glycerin, Water (Aqua), Hydrogen Peroxide, Silica, Carbomer, PVP, Sodium Hydroxide, Potassium Nitrate, Sodium Fluoride, Iron Oxide, EDTA

**Intended Use:**
Opalescence Office is a professionally applied teeth whitener used to whiten and brighten a patient’s smile. It is applied by the dental professional for bleaching one or more teeth.

**Preliminary Procedure:**
1. Clean and polish teeth approximately two weeks prior to bleaching appointment.
2. Determine and record baseline shade. Use photographs when indicated for baseline documentation.
3. Prior to mixing, bring Opalescence Office to room temperature.
4. Cover any sensitive areas of exposed dentin.

**Direction for Use:**
**Mixing Opalescence in-office whitening:**
1. Check syringes are securely attached.
2. Press the small clear plunger into the middle small clear syringe to rupture internal membrane and combine bleach and activator.
3. Press the plunger of the red syringe in, pushing all contents into the clear syringe.
4. Press the clear plunger completely back into the red syringe.
5. To thoroughly mix activator with bleaching gel, push stems back and forth continually with thumbs and mix a minimum of 50 times rapidly (25 times each side)
6. Press all mixed gel into the RED syringe.
7. Press to separate the two syringes and securely attach the recommended tip (Black Mini™) tip onto the red syringe.

**Isolation - OpalDam™ resin barrier (recommended isolation):**
1. Place Ultradent IsoBlock™ bite block and self-supporting lip/cheek retractors.
2. Remove Luer lock cap from OpalDam resin barrier and securely attach a Black Mini™ tip.
3. Verify flow prior to applying intraorally.
4. Rinse and dry teeth and gums thoroughly.
5. Apply a continuous bead of OpalDam resin barrier along the gingival margin, overlapping approximately 2-3 mm onto the enamel. Begin and finish the bead one tooth beyond the most distal tooth that is being bleached.
7. Express resin through any open embrasures onto the lingual tissue to completely seal and cover exposed papilla.

8. Using the mouth mirror, look up the long axis of the teeth for gaps or minute openings exposing the underlying tissue and re-apply resin.
9. Light cure 20 seconds per arch using a scanning motion. Carefully check resin cure with an instrument.

**Isolation - Rubber Dam:**
1. Apply rubber dam placing ligated clamps bi-laterally and one tooth distally to the teeth to be whitened.

2. Before attaching rubber dam to the frame, apply a bead of rubber dam sealer at the gingival margin beneath the dam.

3. Attach dam to the frame and insert rubber dam at the gingival border. Check for optimum seal.

**Applying Opalescence in-office whitening:**
1. On gage or mixing pad, verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.

NOTE: If resistance continues, DO NOT use that syringe and contact Ultradent Customer Service for replacement.
1. Apply a 0.5-1.0mm thick layer of gel to the clean labial surface of the tooth.
2. Apply gel to remain on the tooth 20 minutes, periodically checking and re-applying areas that have thinned or need replenishing.

3. Suction gel from teeth using the Luer Vac and SST™ (Surgical Suction Tip) or a surgical suction tip. To avoid gel splatter, do not use water while withdrawing gel. When no gel is visible, lightly rinse and air dry. Use caution not to dislodge the isolation barrier or rubber dam seal.
5. Evaluate shade change after each application and monitor patient for sensitivity and irritation. If persistent tooth sensitivity occurs, discontinue use.
6. If using a rubber dam, follow “Clean-Up” Step 1, remove clamps, dam, and frame.
7. Remove rubber dam sealer with spatula, soft toothbrush, floss, suction and firm air/water spray.

8. Repeat steps 1-5 up to two additional times if desired results have not been achieved. Do not exceed 3 applications per visit.

8. If additional in-office whitening is desired and no significant sensitivity is noted, re-schedule and retreat patient as desired. If sensitivity was experienced, repeat treatment when the sensitivity has subsided.

**Clean-up:**
1. Suction gel from teeth using the Luer Vac and SST™ (Surgical Suction Tip) or a surgical suction tip, then thoroughly rinse teeth with an air/water spray and the high volume suction.
2. Using a dental instrument, gently slide the tip beneath the OpalDam resin barrier and lift off. Check and remove any interproximal remnants.

- If using a rubber dam, follow “Clean-Up” Step 1, remove clamps, dam, and frame.
- Remove rubber dam sealer with spatula, soft toothbrush, floss, suction and firm air/water spray.
- NOTES:**
  - Whitening effect may continue for 24-48 hours following bleaching, hence a slight increase in whitening may be realized during this period.
  - A single, or a few teeth can be bleached using this same instruction and applying the bleach only to those specific areas of discoloration.
  - For post-up sensitivity, use potassium nitrate products such as UltraZr™ desensitizing gel immediately after treatment and for 15 minutes to 1 hour per day until sensitivity subsides.
  - Etching enamel or prophly paste polishing before any bleaching procedure is not recommended and does not potentiate the bleaching process, but can cause sensitivity.

**Precautions and Warnings:**
1. Pregnant or nursing women should not use teeth whitening products.
2. For patients or users with allergy concerns, refer to product allergen document available at www.ultradent.com. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.
3. Patients with serious health concerns should consult their primary care provider prior to treatment.
4. To prevent additional irritation, gingival and general oral health should be evaluated prior to initiating treatment.
5. Avoid contact with soft tissue or salivary flow.
6. After eye contact: Rinse opened eye for several minutes under running water.
7. Wear protective glasses when handling this product.
8. Clinician, assistant, and patient must wear protective eyewear with side shields when mixing and applying Opalescence in-office whitening.
9. Restorative materials will not whiten.
10. Wait 7-10 days after bleaching procedures before placing adhesive restorations, as whitening materials interfere with bonding agents.
11. White spots on the enamel may become more apparent during the whitening process and often blend with continued whitening.
12. Teeth can re-stain over time, do touch up teeth whitening treatments as necessary.
13. Monitor progress to prevent whitening beyond the desired shade or translucency.
14. Opalescence in-office whitening should be used within 10 days of mixing. Keep refrigerated between uses.
15. Use only recommended tips.
16. To be only sold to and applied by dental practitioners. This product shall not be used by patients but must be used by a licensed professional only.
17. Not to be used on a person under 18 years of age.
18. Contains sodium fluoride.
19. Please note the expiration date printed on the outer package/syringe.
20. To avoid cross-contamination, wipe syringe with an intermediate level disinfectant between uses and use a disposable barrier sleeve. If a barrier sleeve is not used, syringe should be treated as single-use.
21. Dispose of waste according to local rules, guidelines, and regulations.
22. Report any serious incident to the manufacturer and the competent authority. www.ultradent.com / 1-800-552-5512 / 801-572-4200

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