Description:

UltraEZ™ Desensitizing Gel is a unique, sustained release, 3% potassium nitrate and 0.11% weight by weight fluoride ion, sticky, viscous gel. UltraEZ is designed to be used in a custom-fabricated tray.

Indications:

UltraEZ Desensitizing Gel provides a film-like varnish for sensitive teeth, sealing dentinal tubules of exposed dentin and other areas where post-operative sensitivity is a concern. The product is used either by a dental professional in the office or provided to the patient for home treatment of dentin sensitivity.

Procedure:

1. Read and understand instructions before using or dispensing to patient.
2. Fabricate tray to meet patient's needs. If patient already has a custom bleaching tray, UltraEZ may be applied in the same tray and manner as bleaching gel.

Note: If large embrasures exist from prior periodontal surgery that might present problematic impression removal from the teeth, block out embrasures with OraSeal™ putty.

3. UltraEZ Desensitizing Gel is dispensed in 1.2ml unit dose syringes. A complete arch should use most of the contents of a syringe. If sensitivity is localized, place gel in affected area(s) and a small amount in the molar region to hold tray in place.
4. The length and number of times tray is worn depends on condition, patient, and clinician. Recommended treatment times can range from 15 minutes to 1 hour.

5. Tray should be removed at meal times. Immediately following tray removal, patient should brush and rinse teeth with water to remove all remaining gel. DO NOT swallow solution during rinsing. Wash tray and store in a tray storage case. Fresh gel should be used with each application.
6. Gel may be warmed before placing in tray to avoid thermal sensitivity.

Laboratory Instructions:

1. Pour and trim cast as you would for a standard bleaching tray. Apply 0.5-1.0mm of Ultradent™ LC Block-Out Resin to a clean, dry cast on areas in which patient is experiencing

sensitivity (Figs. 1 & 2). Also apply resin to areas where there are open embrasures, etc. For periodontal patients, clinician may determine that pre-operative impressions are indicated. In such cases, reservoir spacing will be created post-op because of recontoured and apically positioned gingiva.

2. Cure LC Block-Out for ~2 minutes in light curing unit. A hand-held intraoral light Valo™ Curing Light) can be used instead (20-40 seconds per tooth). Wipe off oxygen inhibition layer.

3. With vacuum former (Ultra-Form™ or EconoForm™), heat a Sof-Tray™ 0.035" sheet until it sags ~1 inch. Activate vacuum and adapt softened plastic over mold. Cool and remove model.

4. Cut excess bulk of material away with serrated plastic trimmers.
5. Trim excess tray material, allowing 4-6mm to overlap onto gingival tissue. Contour tray margin to avoid anatomical structures such as frenum attachments, tori, etc.

6. Return tray to model; check tray extensions. Gently flame polish edges, one quadrant at a time (Blazer Micro Torch).

7. While still warm, hold periphery of each segment firmly against model for three seconds with water-moistened finger. If an area is short of the desired length, gently heat and push tray material to the desired location. If this thins the material too much, a new tray should be fabricated.

Precautions:

1. Keep desensitization product away from small children.
2. This product is NOT for nighttime use.
3. Do not use tobacco products or eat while using this product.
4. Not intended for use by pregnant women.


5. If any problems occur, discontinue treatment and notify your dentist.
6. Keep out of heat and/or sunlight at all times.

7. Patients should not swallow large volumes of UltraEZ Desensitizing Gel. Excess material at insertion and residue upon removal should be brushed off and rinsed: DO NOT swallow.


Dentist Instructions & Lab Instructions

EN Construct and dry the stone model, cutting palate and tongue areas out.

Vacuum adapt Sof-Tray™ .035" sheet to model. Cool and separate.

For sensitive roots use Ultradent's LC Block-Out Resin on gingival margin to provide a reservoir space for UltraEZ Desensitizing Gel. It is recommended to lap tray 2-3mm onto gingiva.

General Precautions:

1. For professional use only.
2. Review instructions, precautions, and MSDS before beginning treatment. Use only as directed.
3. Keep products out of heat/sunlight.
4. Confirm that patient has no known allergies to treatment materials.

Patient Instructions


EN Follow your dentist's instructions on how to load gel into your custom bleaching tray. Use no more than 1/3 to 1/2 of syringe.

Brush teeth, then insert tray. Lightly tap tray to adapt tray sides to teeth. Remove tray after 15 minutes to 1 hour or as directed by your dentist.

Remove excess gel with clean finger or soft toothbrush. Rinse twice; do not swallow rinses.


Clean tray with soft brush and cool tap water. Store tray in case provided.

Patient Instructions:

EN Center tray on arch.

Lightly press into place.

Remove outer tray.

After placing tray in mouth, gently press the tray from biting surface of teeth towards the gums to ensure complete coverage of gel.

UltraEZ Desensitizing Gel: Wear for 15 minutes to 1 hour per day or as directed by your dentist.

When finished remove tray and brush teeth.

UltraEZ™ Desensitizing Gel

Flammability

Reactivity

Health

EN- Recycle

EN- Keep away from heat/sunlight

EN- See Instructions

EN- Storage temperature

RxOnly

EN- For professional use only

EN- Use by date

LOT

EN- Lot Number

HAZARD RATING
4 = Severe
3 = Serious
2 = Moderate
1 = Slight
0 = Minimal

For professional use only. Keep out of reach of children.

For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free 1-800-552-5512. Outside U.S. call (801) 572-4200 or visit www.ultradent.com.

For product SDS please see our website:
www.ultradent.com
1-800-552-5512; 801-572-4200

Ultradent syringes have an expiration date stamped on the side of the syringe consisting of one letter and three numbers. The letter is a lot number used for manufacturing purposes and the three numbers are the expiration date. The first two numbers are the month, and the third number is the last of the year.

U.S. federal law restricts the sale of this device by or on the order of a dentist.

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© 2022 Ultradent Products, Inc. All Rights Reserved. 1009739A003 062322