CE

170K200=*, 170K201=*, 170K202=*

DE	Gebrauchsanweisung	3
EN	Instructions for use	39
FR	Instructions d'utilisation	75
IT	Istruzioni per l'uso	111
NL	Gebruiksaanwijzing	147
SV	Bruksanvisning	183
DA	Brugsanvisning	219
NO	Bruksanvisning	255
FI	Käyttöohje	291

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DE

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Table of contents

EN

1	Foreword	40
2	Product description	40
2.1	Components/design	40
2.1.1	Design	40
2.2	Available sizes	40
3	Intended use	43
3.1	Indications for use	43
3.2	Indications	43
3.3	Lifetime	43
3.4	Qualification	43
4	Safety	43
4.1	Explanation of warning symbols	43
4.2	Safety instructions	43
5	Preparing the product for use	44
5.1	Adaptation	44
5.2	Application	72
5.3	Removal	72
6	Cleaning	72
7	Maintenance	72
8	Disposal	73
9	Legal information	73
9.1	Liability	73
9.2	CE conformity	73

1 Foreword

INFORMATION

Date of last update: 2025-02-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the C-Brace interim orthosis 17OK200=*, 17OK201=* and 17OK202=* for the C-Brace joint unit.

The C-Brace interim orthosis 17OK2=* is also referred to below as the CBI, product or orthosis and the C-Brace adapter joint is referred to as the orthotic joint.

2 Product description

2.1 Components/design

2.1.1 Design

The C-Brace interim orthosis 17OK2=* allows the fabrication of a unilateral leg orthosis (KAFO) by assembling the thigh sleeve, lower leg shell, foot component and the screw connection to the orthotic joint. Combined with the C-Brace joint unit, the leg orthosis is complete for the user and functional for temporary use. The position of the anchor screws and the short cap screws determines whether the orthotic joint is for the right or left side. The angular degrees for the orthotic joint are determined with the help of the orthotic joint templates during the fitting of the leg orthosis.

2.2 Available sizes

Thigh sleeve

1

The thigh sleeve consists of the following components:

1. Thigh sleeve
2. Compensation pad
3. Hook-and-loop patches

Thigh sleeve size chart	Reference number	Side	Size	Circumference in cm
Circumference 15 cm from the centre of the patella	17OK202=L/R-S	left/right	S	30 – 45
	17OK202=L/R-M	left/right	M	40 – 55
	17OK202=L/R-L	left/right	L	50 – 65

Ankle-foot orthosis

2

The ankle-foot orthosis consists of the following components:


1. Lower leg shell
2. Foot component
3. Orthotic ankle joint

Ankle-foot orthosis size chart	Reference number	Side	Size	Ankle measurement in mm	Shoe size
Ankle – floor measurement in mm or shoe size	17OK200=L/R-S	left/right	S	70	36 – 41
	17OK200=L/R-L	left/right	L	85	39 – 46

CBI Orthotic joint

The orthotic joint 17OK201=* of the C-Brace interim orthosis has special screw points for the thigh sleeve and the lower leg shell. It is available in 0°, 3°, 6° and 9° angles.

3


17OK201=* scope of delivery		
Item	Quantity [Piece(s)]	Designation
1	1	Upper joint section
2	3	Anchor screw
3	3	Cap screw, long
4	1	Axial washer, medial
5	1	Lower joint section
6	1	Bushing
7	1	Axial washer, lateral
8	1	Joint screw
9	3	Cap screw, short
10	1	Adapter plate
11	3	Oval head screw, lower joint section
12	7	Drill sleeve
13	4	Washer
14	4	Oval head screw, upper joint section

Not included in scope of delivery			
Item	Quantity [Piece(s)]	Designation	Reference number
Not illustrated	1	C-Brace joint unit	17KO1
Not illustrated	1	Service set (5x lateral axial washers, 1x medial axial washer, 1x bushing)	17KF103=16

Tools

4

The tool set 17OK203=N-1 is required to fabricate the CBI.

1. Clamping lever
2. Retaining screw
3. Saw guide
4. Clamp adapter
5. Template set (see fig. 5)

5

The template set 17OK203=N-2 with the angle positions 0°, 3°, 6° and 9° is used to select the correct orthotic joint.

Required tools			
Quantity [Piece(s)]	Designation	Quantity [Piece(s)]	Designation
1	Torx 25, 30, 40	1	Phillips screwdriver
1	Torque wrench	1	Metal saw
1	Measuring devices	1	Knee pivot gauge 743A8
1	Drill bit 2.5 mm, 7 mm	1	Countersink
1	Textile scissors, micro-serrated	1	Hole punch 717L4
1	Grease pen 645C2	1	Marking pen
1	Riveting hammer	1	Double-sided adhesive tape 633D5

3 Intended use

3.1 Indications for use

The C-Brace interim orthosis is intended **exclusively** for orthotic treatment of the lower limbs of users weighing up to 110 kg.

3.2 Indications

- Unilateral or bilateral leg paresis or flaccid paralysis, e.g. due to polio or incomplete paraplegia
 - The physical prerequisites such as muscle status, mobility grade or axis deviations are crucial
- Indications must be determined by the physician.

3.3 Lifetime

The product is designed for a maximum lifetime of **6 months**.

3.4 Qualification

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

 Warning regarding possible risks of accident or injury.

 Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Reaching into the joint mechanism

Risk of pinching due to movement of the joint

- ▶ Do not reach into the joint mechanism.
- ▶ Use caution during assembly and adjustments.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- Keep the product away from open flames, embers and other sources of heat.

⚠ CAUTION

Product contains magnets

Influence on or damage to devices or items due to a strong magnetic field

- Keep the product away from devices and items that are sensitive to magnetic fields. (e.g. pacemakers, electronic devices, swipe cards).
- Pay attention to any manufacturer's information (e.g. for a medical implant).

⚠ CAUTION

Damaged electronic components

Electric shock due to contact with live components

- Check the product for damage (e.g. cable, housing) prior to each use while in a de-energised state.

⚠ CAUTION

Incorrect adjustment or positioning of the product

Risk of injury due to falling

- Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- Use the product on only one patient.
- Observe the maintenance recommendations.

⚠ CAUTION

Exposure of the product to unsuitable environmental conditions

Patient injury, damage, brittleness or destruction due to improper handling

- Do not expose the product to condensing ambient humidity or liquids.
- Do not expose the product to abrasive substances (e.g. sand, dust).
- Do not expose the product to temperatures below **-10 °C** (14 °F) or above **+40 °C** (104 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

Handling and use of the product must be learned before using it for the first time.

Applying and removing, sitting down and standing up as well as walking must be practised.

5.1 Adaptation

⚠ CAUTION

Incorrect or excessively tight application of the product to the body

Risk of pressure points and constriction of blood vessels and nerves due to improper application

- Ensure that the product is applied properly and fits correctly.
- Do not continue to use the product if it fits incorrectly.

NOTICE**Damage due to improper handling (heat treatment, drilling holes)**

Loss of orthosis functionality

- ▶ Do not carry out any heat treatment on the orthosis. The orthosis is not thermoformable.
- ▶ Do not drill any holes in the bottom of the sole or the connecting elements as holes break the fibres and weaken the component.

NOTICE**Use without suitable footwear**

Restriction of functionality due to insufficient stabilisation

- ▶ Only use the orthosis with closed footwear, taking the allowable effective heel height into account.

Selecting footwear

- ▶ Choose a sturdy shoe with laces or hook-and-loop fasteners and a solid toe cap to ensure the optimal effectiveness of the orthosis. The effective heel height should be between 0.5 cm and 1.5 cm.

Adjusting the ankle-foot orthosis

6

Take the insole out of the shoe.

7

INFORMATION: The foot component is selected in combination with the lower leg shell according to the table on page 5.

Make sure that the foot is positioned on the insole or customised insole when taking measurements.

- 1) **Optionally** fabricate an insole and temporarily attach it to the foot component of the orthosis with double-sided adhesive tape.
- 2) The foot component is selected based on the ankle-floor measurement and/or shoe size.

8

INFORMATION: The lateral section of the foot component can be adapted to a limited extent. (see fig. 12)

Optional: Optimise the fit by using padding material.

9

- 1) Place the foot on the insole.
- 2) Mark the position of the lateral ankle on the insole with a pen or adhesive tape.

10

INFORMATION: For selection of the ankle-foot orthosis see the table on page 5.

- 1) Position the insole on the foot component.
- 2) Align the mark on the insole with the orthotic ankle joint.

11

- 1) Position the foot with the insole on the foot component.
- 2) Mark the outline of the insole on the foot component.

12

NOTICE: Do not sand the areas marked in red.

13

- 1) Sand the foot component according to the markings.
- 2) Remove grooves and burrs.

14

NOTICE: If an insole was also used to determine the ankle height, it must still be placed on the foot component of the orthosis.

Position the foot component in the shoe and check the fit.

15

INFORMATION: Position the lower leg shell as close as possible to the patella.


Mark the front height of the lower leg shell.

- 1) Check the fit of the lower leg shell.
- 2) Optionally optimise the fit using additional pads.
- 3) Temporarily secure the lower leg shell to the leg with the strap.
- 4) Check the fit in the flexion cut-out while the strap is closed.

- 1) Position the foot component relative to the lower leg shell.
- 2) Make sure that the insole used is positioned on the foot component.
- 3) Mark the length of the lower leg bar to match the insertion area of the orthotic ankle joint.

INFORMATION: Two guide grooves are provided for the saw to cut the lower leg bar to length.

Prerequisite: The laser marking with the arrow symbol points distally. **NOTICE:** Incorrect alignment of the saw guide will damage the drill holes in the lower leg bar.

Example of shortening the lower leg bar:

The possible positions of the drill holes in the lower leg bar for the insertion area of the orthotic ankle joint are shown in green in this example.

INFORMATION: The mark for cutting the bar section to length is at the level of the groove in this example.

Position the lower leg bar in the saw guide.

20


Secure the lower leg bar in the saw guide with the screw.

21

INFORMATION: The shape of the saw guide enables secure clamping in the vice.

22

Shorten the bar section with the saw.

23

Round off the edges by grinding and fit it into the insertion area of the orthotic ankle joint.

24

- 1) Check the fit of the bar section in the orthotic joint.
- 2) Remove grooves and burrs by grinding.
- 3) Attach the bar section in the orthotic joint using screws.

25

Prerequisite: The orthosis has been placed with the insole in the shoe and applied to the patient.

- 1) Check the fit of the orthosis.
- 2) Mark the position of the strap so it is centred on the lower leg shell and shorten it if necessary.

TIP: The central positioning of the strap on the lower leg shell makes it easier for the user to open and close the hook-and-loop fastener.

26

- 1) Remove the Y-hook-and-loop and cut the strap to length.
- 2) Then reattach the Y-hook-and-loop to the strap with hook-and-loop fastener.
- 3) Secure the orthosis to the leg with the strap.

Marking the knee pivot point

27


Indicate the markings to determine the compromise pivot point according to Nietert.

Mark the knee pivot point.

TIP: Determine the a-p measurement with the 60:40 knee pivot gauge.

Adapting the thigh sleeve

Determine the position for the circumference measurement on the thigh.

INFORMATION: The thigh sleeve is selected according to the table on page 4.

Measure the circumference to select the thigh sleeve.

NOTICE: Only process the carbon up to the visible and palpable mark at most.

32

NOTICE: Only cut the flexible contact surfaces/flaps to the visible and palpable mark at most.

33

Prerequisite: The pad is attached to the thigh sleeve with a hook-and-loop fastener.
Mark the course of the edge with the thigh sleeve applied.

34

Adapt the flexible course of the edge of the thigh sleeve.

Tip: Use micro-serrated scissors to cut the flexible course of the edge.

35

- 1) Mark the length of the strap.
- 2) Remove the Y-hook-and-loop and cut the strap to length.
- 3) Then reattach the Y-hook-and-loop to the strap with hook-and-loop fastener.

- 1) Mark the course of the edge of the thigh sleeve on the carbon.
- 2) Grind the carbon according to the markings.
- 3) Remove grooves and burrs.

Materials required: degreasing cleaner

- 1) Clean the adhesive surface on the thigh sleeve using a degreasing cleaner.
- 2) Apply the self-adhesive hook-and-loop patches (see fig. 1, item 3) according to the course of the pads.

TIP: To ensure good adhesion of the hook-and-loop patches to the flexible sections of the thigh sleeve, allow the cleaner to evaporate for 2 minutes.

Adapt the pad to the course of the edge of the thigh sleeve.

TIP: For better wearing comfort, allow the proximal course of the edge of the pad to protrude by 2 cm and fold over the edge (see Fig. 39).

Fold over the pad and secure it with the hook-and-loop fastener.

40

Mark the areas of the flap to be cut out with a pen.

TIP: Parts of the flap on the thigh sleeve can be cut out to optimise the fit.

41

Cut out the marked areas with a pair of scissors.

TIP: Use micro-serrated scissors.

Selecting the orthotic joint

42

INFORMATION: When putting on the thigh sleeve, make sure there is enough space between the lower edge of the carbon shell and the compromise pivot point.

43

INFORMATION: Use the template set to select the orthotic joint (see fig. 5)

Before positioning the template, align the screw connection surfaces of the thigh sleeve and the lower leg section parallel to each other.

Align the mark on the template at the level of the knee pivot point.


Place the template on the outside or inside of the screw connection surface of the lower leg section.

Prerequisite: The leg of the user is extended and the orthosis shells are correctly aligned. The template lies flat on the screw connection surfaces of the thigh sleeve and the lower leg section.

Determine the required medial-lateral orthotic joint angle by positioning the template.

Positioning the orthotic joint

INFORMATION: The clamp adapter aligns and secures the lower leg bar of the orthotic joint on the screw connection surface in two planes. The lower leg bar is positioned up/down and forwards/backwards in the clamp adapter.

48

INFORMATION: The orthotic joint must be secured in the extended position for the following steps.

Install the locking screw in the orthotic joint.

49

- 1) Loosen the 4 screws of the clamp adapter.
- 2) Slide the clamp adapter onto the screw connection surface to the upper end.

50

Prerequisite: The thigh sleeve and ankle-foot orthosis have been applied to the extended leg of the user. The locking screw (see fig. 3, item 5) is installed in the orthotic joint.

- 1) Position and align the lower leg bar between the clamp adapter and the screw connection surface.
- 2) Secure the position of the lower leg bar by tightening the screws.

51

NOTICE: 3 screw connections are required for final assembly of the lower leg bar with the screw connection surface.

Make sure that the lower leg bar does not come into contact with the carbon due to excess length.

Optional: Shorten the lower leg bar of the orthotic knee joint.

52

INFORMATION: As with the lower leg bar, the saw guide can also be used to shorten the bar of the orthotic joint.

Insert the bar section of the orthotic joint into the saw guide.

53

Secure the bar section in the saw guide with the screw.

54

- 1) Shorten the bar section with the saw.
- 2) Remove grooves and burrs by grinding.

55


- 1) Check the pivot point with a suitable tool (e.g. an Allen key).
- 2) Position the lower leg bar until the desired result is obtained.

INFORMATION: After aligning the orthotic joint pivot point, the position of the orthotic joint lower section is secured in the clamp adapter.


Tighten the screws of the clamp adapter.

Connect the upper joint section to the thigh sleeve.


INFORMATION: The mechanism is released only when the clamping lever is in the vertical position.

The upper joint section (item 2) of the orthotic joint is secured via the clamping lever (item 5), the clamp washer (item 4) and the insert nut (item 3) on the screw connection surface of the thigh sleeve (item 1).

Application example: vertically shift the upper joint section.

- 1) Release the clamping lever.
- 2) Align the slotted hole (item 1) of the clamp washer vertically.
- 3) Shift the upper joint section.
- 4) Adjust the clamping force by screwing the clamping lever in/out.

Application example: horizontally shift the upper joint section.

- 1) Release the clamping lever.
- 2) Align the slotted hole (item 1) of the clamp washer horizontally.
- 3) Shift the upper joint section.
- 4) Adjust the clamping force by screwing the clamping lever in/out.

- 1) Bring the clamp washer and the clamping lever together.
- 2) Then screw the clamping lever into the thread in the thigh sleeve.

Determine the desired position by moving the upper joint section towards the thigh sleeve.

- 1) Adjust the clamping force by screwing in the clamping lever in the released position.
- 2) Secure the specified position of the upper joint section on the thigh sleeve by locking the clamping lever.

INFORMATION: The orthotic joint must be able to move freely to check the pivot point. To do this, the locking screw is loosened to allow the leg to be flexed freely in the orthosis.

Loosen the locking screw.

TIP: Use a suitable auxiliary tool (e.g. Allen key) to project the orthotic joint pivot point onto the body.

- 1) Check the pivot point position of the orthotic knee joint with the leg flexed and extended.
- 2) Optionally adjust the position of the upper joint section until the desired result is achieved.

Mark the position of the upper joint section on the screw connection surface.

Mark the position of the lower leg bar on the screw connection surface.

68

Open the straps and take off the orthosis.

Connecting the orthotic joint to the orthosis

69

Loosen the screws from the orthotic ankle joint and remove the lower leg bar from the foot component.

70

Loosen the joint screw.

71

- 1) Disassemble the orthotic joint.
- 2) Remove the medial and lateral axial washer.

72

Required tools: Drill sleeves (see fig. 3, item 12)

NOTICE: 3 screws are required to secure the lower leg bar of the orthotic joint on/under the screw connection surface of the lower leg shell. The drill sleeves prevent damage to the internal threads during the drilling process in the lower leg bar.

Position the drill sleeves in the internal threads of the lower leg bar.

73

TIP: Secure the lower leg section with the saw guide in the vice for drilling.

74

Required tools: 2.5 mm drill bit

- 1) Using the drill bit, drill all three holes into the carbon of the screw connection surface.
- 2) Remove the drill sleeves.

75

- 1) Loosen the screws of the clamp adapter.
- 2) Remove the lower leg bar and clamp adapter from the screw connection surface.

Tools required: 7 mm drill bit, counterbore

- 1) Ream the drill holes using a 7 mm drill bit.
- 2) Deburr the drill holes with a counterbore.

Materials required: degreasing cleaner, Loctite 241

- 1) Clean the threads of the screw connections using a degreasing cleaner.
- 2) Apply Loctite 241 to all screws.

Tools required: torque wrench, Torx bit 30

- 1) Position the lower leg bar on the screw connection surface.
- 2) Secure all screws with a tightening torque of **8 Nm**.

Remove the remains of the Loctite with a cloth.

NOTICE: 3 or 4 screws are required to attach the upper joint section to the screw connection surface of the thigh sleeve.

Prerequisite: The upper joint section is positioned in the centre of the screw connection surface.

Example 1, application case with 3 screws:

The positions of the drill holes are shown in green in this example.

Prerequisite: The upper joint section is not positioned in the centre of the screw connection surface.

Example 2, application case with 4 screws:

The positions of the drill holes are shown in green in this example.

Required tools: Drill sleeves (see fig. 3, item 12)

NOTICE: The drill sleeves prevent damage to the internal threads during the drilling process in the upper joint section.

Position the drill sleeves in the internal threads of the upper joint section.

Required tools: 2.5 mm drill bit

- 1) Use the drill bit to drill 3 or 4 holes in the carbon of the screw connection surface.
- 2) Remove the drill sleeves.

84

- 1) Remove the clamping lever and the clamp washer.
- 2) Remove the upper joint section from the thigh sleeve.

85

Tools required: 7 mm drill bit, counterbore

- 1) Ream the drill holes using a 7 mm drill bit.
- 2) Deburr the drill holes with a counterbore.

86

INFORMATION: Remove the insert nut for the clamping lever as it is no longer required.

- 1) Screw a screw into the insert nut.
- 2) Remove the insert nut from the thigh sleeve.

87

Required tools and materials: torque wrench, Torx bit 30, degreasing cleaner, Loctite 241

- 1) Clean the threads of the screw connections using a degreasing cleaner.
- 2) Apply Loctite 241 to all screws.
- 3) Position the upper joint section on the screw connection surface and secure it with screws (see fig. 3, item 15) and washers (see fig. 3, item 14).
- 4) Secure all screws with a tightening torque of **8 Nm**.
- 5) Remove the remains of the Loctite with a cloth.

88

INFORMATION: For a better fit of the thigh sleeve, the cavity under the screw connection surface is filled with a compensation pad (see fig. 1, item 2).

Fit the compensation pad into the thigh sleeve by grinding.

89

Use a grease pen to mark the screws over a large area.

90

INFORMATION: When the compensation pad is pressed on, the markings are transferred to the pad via the screws.

Press the compensation pad into the recess.

91

Apply double-sided adhesive tape to the inner surfaces of the compensation pad.

92

INFORMATION: Punching out the compensating material at the marked point creates a free space for the screws.

Use a corresponding hole punch to punch out the marked areas.

93

Glue the compensation pad into the thigh sleeve.

Mounting the orthotic joint

94

Fasten the pads in the thigh sleeve and the lower leg section with a hook-and-loop fastener.

95

Glue the insole to the foot component.

96

Required tools and materials: torque wrench, Torx bit 30, degreasing cleaner, Loctite 241

- 1) Clean the threads of the screw connections using a degreasing cleaner.
- 2) Position the lower leg bar in the insertion area of the orthotic ankle joint.
- 3) Secure all screws with Loctite® 241 and a tightening torque of **8 Nm**.

97

Remove the remains of the Loctite® with a cloth.

98

Tools required: torque wrench, Torx bit 30

NOTICE: Note the order of the axial washers during assembly (see fig. 3).

- 1) Place the medial axial washer (item 13) on the upper joint section.
- 2) Position the lower joint section with the lateral axial washer (item 7) on the upper joint section.
- 3) Assemble the adapter joint with the joint screw (item 6).
- 4) Tighten the joint screw with a tightening torque of **8 Nm**.

Mounting the anchor screws


The positioning of the 3 anchor screws in the upper joint section and the 3 short cap screws in the lower joint section determines the side of the orthotic joint. The anchor bolts are marked with rings. All required screws and anchor screws are included in the scope of delivery.

Detailed information on the use of the anchor screws is part of the C-Brace basic seminar.

Example: left upper joint section, 0°

The following illustration shows the final positions of the anchor screws for the left orthotic joint at a joint angle of 0°.

99

Required tools and materials: torque wrench, Torx bit 40, degreasing cleaner, Loctite 241

- 1) Screw the anchor screw with 4 rings into the threaded hole marked with **0000 L** in the upper joint section (see fig. 99).
- 2) Screw the anchor screw with 3 rings into the threaded hole marked with **000 L** in the upper joint section.
- 3) Screw the anchor screw with 4 rings into the threaded hole marked with **0000 L** in the upper joint section.

INFORMATION: 2 anchor screws are always screwed into the anterior threaded holes in the upper joint section.

100


- 1) Clean the threads of the screw connections using a degreasing cleaner.
- 2) Secure the anchor bolts with Loctite 241 and a tightening torque of **8 Nm**.
- 3) Remove the remains of the Loctite® with a cloth.

Example: left upper joint section, 3°, 6°, 9°

The following illustration shows the final positions of the anchor screws for the left orthotic joint at joint angles of 3°, 6° and 9°.

101

Assembling the C-Brace joint unit

102

Required tools and materials: torque wrench, Torx bit 25, degreasing cleaner, Loctite 241

NOTICE: Use the supplied cap screws (see fig. 3, item 3, item 9) of the orthotic joint to assemble the C-Brace joint unit.

- 1) Clean the threads of the screw connections using a degreasing cleaner.
- 2) Apply Loctite 241 to all screws.

103

- 1) Position the C-Brace joint unit on the orthotic upper joint section.
- 2) Secure all screws with Loctite® 241 and a tightening torque of **7 Nm**.

104

Position the C-Brace lower joint section with the adapter plate on the orthotic lower joint section.

105

- 1) Secure all screws with Loctite® 241 and a tightening torque of **7 Nm**.
- 2) Remove the remains of the Loctite® with a cloth.

Checking the fit

106

Place the orthosis on the user.

107

Check the fit of the orthosis on the user.

108

NOTICE: To adjust the C-Brace joint unit with the user, observe the information in the 647G1337=* instructions for use.

Carry out the static and dynamic fitting with the user.

5.2 Application

NOTICE

Use of a worn or damaged product

Limited effectiveness

- Before each use, check the product for functional reliability and for possible wear or damage.
- Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Open the straps of the thigh and lower leg shells.
- 2) Flex the orthosis and apply it to the thigh and lower leg.
- 3) Insert the orthotic foot component into the shoe.
- 4) Fasten all straps on the orthosis.

5.3 Removal

- 1) Open all straps on the orthosis.
- 2) Remove the shoe from the orthotic foot component.
- 3) Take off the orthosis.

6 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads on a user-specific basis.

- Shorten the maintenance intervals according to the expected loads.

Spare parts are listed in the section "Components/design".

Replacing the axial washers and bushing

Required materials: service set (17KF103=16), T30 Torx key.

> **Prerequisite:** The C-Brace joint unit has been disassembled.

- 1) Loosen the joint screw ⑥ (see fig. 3).
- 2) Disassemble the orthotic joint; remove the bushing ④ and the axial washers ⑬, ⑦.
- 3) Place the bushing with the rounded side in the lower joint section, avoiding canting or damaging the bushing.
- 4) Press the bushing in to the stop.
- 5) Replace the medial axial washer ⑬ and position it on the upper joint section.
- 6) Replace the lateral axial washer ⑦ and position it on the lower joint section.
- 7) Assemble the adapter joint with the joint screw.
- 8) Optimise the movement of the orthotic joint by gradually replacing the lateral axial washers.

- 9) Tighten the joint screw with a tightening torque of **8 Nm**.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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