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Highsmith MJ, Kahle JT, Wernke MM, Carey, SL, Miro RM, Lura DJ, Sutton BS.
School of Physical Therapy & Rehabilitation Sciences, University of South Florida,
Tampa, FL, USA.

Effects of the Genium knee system on functional level, stair ambulation, perceptive and economic outcomes in transfemoral amputees

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
With Genium compared to C-Leg:

- **The quality of stair ascent and descent improved significantly**
- **Mobility and functional level improved significantly**
- **Perceived function and safety in ADLs was as good as with C-Leg or significantly better**
- **80% preferred Genium**

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
Interventional, randomized crossover design:

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“Accommodation and use of the Genium knee system compared with C-Leg improved stair walking performance, multi-directional stepping, functional level, and perceived function. Genium was also preferred compared to C-Leg in this group of highfunctioning community ambulators with unilateral transfemoral amputation. Finally, Genium is a more costly microprocessor knee system but, in this group of patients, is worth funding due to significant differences in functional performance with activities of daily living.” (Highsmith et al., 2016)

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