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Samitier BC, Guirao L, Costea M, Camós JM, Pleguezuelos E.

Hospital de Mataró, Barcelona, Spain.

The benefits of using a vacuum-assisted socket system to improve balance and gait in elderly transtibial amputees

Journal of Prosthetics and Orthotics. Epub.

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Vacuum-assisted socket system* (VASS) vs other socket system

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With VASS compared to previous socket:

MFCL K2 subjects

→ Improvement of prosthesis use by 7.4%

→ Balance increased by 15%

→ Trend towards decreased risk of falling

Trend towards improvements in timed up and go test (TUG) by 15%

MFCL K3 subjects

→ Decreased risk of falling

Timed up and go test (TUG) improved by 22%

→ Walking velocity increased by 15% and mobility by 19%

→ Increased balance by 19%

→ Subjects tend to more satisfied

→ Time of prosthesis use tends to be increased

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
Decreased time to complete TUG with VASS

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
The timed up and go test (TUG) includes standing up from a chair, walking 3 meters, turning around, walking 3 meters, sitting down. The marked part indicated a greater risk of falling assessed from a healthy geriatric population.

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Subjects: 16 unilateral, transtibial amputees
 Previous socket system: not reported
 Amputation causes: 100% peripheral vascular disease
 Mean age: 65 ± 10 yrs
 Mean time since amputation: 5.2 ± 2.2 yrs
 MFCL: 37.5% K2, 62.5% K3

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“In conclusion, the Harmony® P2 & HD is a useful device in dysvascular transtibial amputees over 50 years of age. In our study, use of the VASS improved balance, gait and transfers in patients with MFCL-3 mobility grade and balance and prosthesis use in patients with MFCL-2 activity level. In patients with a lower activity level, the use of an additional distal valve in the socket should be considered.” (Samitier et al. 2014)

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