

## 3R20, 3R36

<b>DE</b> Gebrauchsanweisung .....	8	<b>PL</b> Instrukcja użytkowania.....	60
<b>EN</b> Instructions for use .....	14	<b>HU</b> Használati utasítás .....	66
<b>FR</b> Instructions d'utilisation.....	19	<b>CS</b> Návod k použití.....	72
<b>IT</b> Istruzioni per l'uso .....	25	<b>RO</b> Instrucțiuni de utilizare.....	77
<b>ES</b> Instrucciones de uso .....	31	<b>TR</b> Kullanma talimatı .....	83
<b>PT</b> Manual de utilização.....	37	<b>EL</b> Οδηγίες χρήσης .....	89
<b>NL</b> Gebruiksaanwijzing.....	43	<b>RU</b> Руководство по применению .....	95
<b>SV</b> Bruksanvisning .....	49	<b>ZH</b> 使用说明书.....	101
<b>DA</b> Brugsanvisning .....	55		

1 **i** 5 Lieferumfang/  
Scope of delivery


2 **i** 6 Technische Daten/  
Technical data

3

## i 7.2 Grundaufbau der Prothese/

### Bench alignment of the prosthesis

4


## 7.2 Grundaufbau der Prothese/ Bench alignment of the prostheses


Bench alignment of the prostheses

①


①

②


③


②


①

②

③

5

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### INFORMATION

Last update: 2018-10-24

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

### 1.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

## 2 Product description

The product (3R20, 3R36) has the following key features:

- Polycentric prosthetic knee joint
- Product components for stance phase stability
  - Polycentric kinematics (adjustable extension stop)

- Product components for swing phase control:
  - Extension assist (adjustable spring force)
  - Adjustable friction brake as movement resistance for extension and flexion
- Material: Steel (3R20), Titanium (3R36)

## 3 Intended use

### 3.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

### 3.2 Area of application

#### Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

#### Allowable field of application (3R20, 3R36)

Recommended mobility grade: <b>1 + 2</b>	Everyday prosthesis
Allowable body weight: <b>≤ 100 kg</b>	

#### MG 1: indoor walker

The patient has the ability or the potential to use the prosthesis for the purpose of moving or proceeding slowly on level floors. The amount of time and the distance he/she can walk are seriously limited due to his/her condition.

#### MG 2: Restricted outdoor walker

The patient has the ability or the potential to walk slowly with the prosthesis and to negotiate low environmental obstacles like curbs, single steps or uneven surfaces. The amount of time and the distance he/she can walk are seriously limited due to his/her condition.

### 3.3 Combination possibilities

#### **⚠ CAUTION**

#### **Failure to observe the manufacturer requirements regarding combination possibilities**

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

#### **INFORMATION**

When fabricating a prosthesis for the lower limbs, all of the prosthetic components used have to meet the requirements in regard to the patient's body weight and activity level.

### 3.4 Environmental conditions

#### **⚠ CAUTION**

#### **Use under unallowable environmental conditions**

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

#### **Allowable environmental conditions**

Temperature range: -10 °C – +60 °C

Use - Relative humidity: 0 % – 90 % (Non-condensing)

Storage/Transportation - Relative humidity: 100 %

Tyre pressure: 250 – 1100 mbar

#### **Allowable environmental conditions**

Commercially available, solvent-free cleaning agents

#### **Unallowable environmental conditions**

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

### 3.5 Reuse and service life

#### **⚠ CAUTION**

#### **Reuse on another patient**

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

#### **⚠ CAUTION**

#### **Exceeding the service life**

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the service life defined in this section is not exceeded.

Maximum service life for this product is 3 to 5 years, depending on the patient's level of activity.

### 4 General safety instructions

#### **⚠ CAUTION**

#### **Reaching into the area of the joint mechanism**

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

**⚠ CAUTION****Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

**Signs of changes in or loss of functionality during use**

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

**5 Scope of delivery**

For the scope of delivery, see Figure **1** on page 2.

**6 Technical data**

For the technical data, see Figure **2** on page 2.

**2**: ① Proximal system height, ② Distal system height

**7 Preparation for use****7.1 Information on fabrication of a prosthesis****⚠ CAUTION****Incorrect alignment, assembly or adjustment**

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

**⚠ CAUTION****Initial use of the prosthesis by the patient**

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

**7.2 Bench alignment of the prosthesis****⚠ CAUTION****Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

**INFORMATION**

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

- ▶ To enable the patient to stand safely, align the prosthesis using the alignment images (Figure **3** on page 3) as well as the instructions for use of all prosthesis components used.
- ▶ Additionally, observe Figure **4** on page 4 and the following points for the bench alignment of the prosthesis:
- ▶ ① Turn the adjustable extension stop to the right to the maximum extension position for bench alignment.

- ▶ Use the rotation to the left shown in ② for dynamic fitting if the patient has difficulty initiating the flexion. Note that this reduces the stability of the prosthetic joint.

### 7.3 Optimising the static alignment

- ▶ Check and optimise static alignment on standing patient (see Fig. 5 on page 4).
- ▶ If a measuring device for representing load lines is available, check the distances indicated in Figure 5 on page 4 and correct the prosthesis alignment as necessary.

### 7.4 Optimising during dynamic trial fitting

#### ⚠ CAUTION

#### Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
  - ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.
- ▶ If the patient has difficulty initiating the flexion, turn the adjustable extension stop to the left as shown in Figure 4 on page 4 in ②.

#### Adjusting the extension assist

- ▶ For adjusting the extension assist, note Figure 6 on page 5.
- 6: ① Increasing the spring tension, ② Reducing the spring tension
- ▶ **CAUTION!** The lower leg of the prosthesis must not swing through too far in the dorsal direction at the beginning of the swing phase.  
Adjust the extension assist so that the lower leg of the prosthesis swings through in the dorsal direction according to the patient's requirements, but reaches full extension in a timely manner for the next heel strike.

#### Adjusting the swing phase friction

- ▶ For adjusting the swing phase friction, note Figure 7 on page 5.
- 7: ① Increasing the axle friction, ② Reducing the axle friction

- ▶ **CAUTION!** The prosthetic knee joint has to reach full extension, even at slower walking speeds. Observe the information under the heading "Fine-tuning the settings".
- ▶ Adjust the swing phase friction so that the prosthetic knee joint does not swing too hard against the extension stop and so that the prosthetic foot swings through sufficiently according to the needs of the patient.

### Fine-tuning the settings

- ▶ **CAUTION!** The various settings of the prosthetic knee joint cannot be adjusted entirely independently of each other. When the settings cannot be fully adjusted to the comfort needs of the patient, the settings should be established primarily based on safety aspects.
- ▶ Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.
- ▶ If the extension position is not reached with every step and at any walking speed, reduce the braking effect and/or increase the initial tension of the extension assist.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and the annual safety inspection.  
Advise the patient to check the prosthesis if changes in function occur.

### 7.5 Attaching the cosmetic cover

#### ⚠ CAUTION

#### Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.
- ▶ After affixing the cosmetic elements, verify proper functioning of the prosthesis.
- ▶ To reduce friction and to eliminate noise, apply silicone spray directly onto the contact surfaces of the cosmetic foam cover.

## 7.6 Finishing the prosthesis

- ▶ **CAUTION!** To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long, and by tightening all screw connections. In doing so, note the instructions for use for all prosthesis components regarding torque values and thread lock.

## 8 Use

### 8.1 Information for use

#### NOTICE

#### Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

### 8.2 Cleaning

#### CAUTION

#### Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

#### INFORMATION

- ▶ Clean the product if it gets dirty.
- ▶ For disinfection, only use disinfectants that do not affect the materials of the product. For further information, please consult the manufacturer.

- ▶ Observe the cleaning instructions for all prosthetic components.

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

## 9 Maintenance

#### CAUTION

#### Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

#### Lubrication

Figure **8** on page 6 shows which lubricants to use. The lubricants can be ordered from the manufacturer using the reference numbers (address at the end of this document).

- ▶ Use lubricants on the product sparingly.

**8**: ① Lubrication

## Replacing components

- ▶ If necessary, replace the components shown in Figures  (page 6) and  (page 7). The components can be ordered from the manufacturer using the reference numbers (address at the end of this document).

## 10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

## 11 Legal information

### 11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 11.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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