---

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Bellmann M\*, Köhler T M \*, Schmalz T\*

\*Ottobock SE & Co. KGaA, Clinical Research and Services, Research Biomechanics, Göttingen, Germany.

# Comparative biomechanical evaluation of two technologically different microprocessor-controlled prosthetic knee joints in safety-relevant daily-life situations

Biomedical Engineering / Biomedizinische Technik, 20180026, ISSN (Online) 1862-278X, ISSN (Print) 0013-5585, DOI: <https://doi.org/10.1515/bmt-2018-0026>.

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**C-Leg 4, Rheo Knee XC**

---

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→ Rheo Knee XC reliably switches into the swing phase during small steps (default swing) in every trial observed; C-Leg 4 in approximately 90% of all trials (default stance)

→ Walking backwards is safer with C-Leg 4 due to a reliable and stable flexion resistance

- The joint remains in stance phase mode throughout the gait cycle
- Due to the "default swing" concept used in the Rheo Knee XC, walking backwards may lead to an uncontrolled knee flexion and thus to a fall

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→ C-Leg 4 offers a reliable and stable load bearing while descending stairs and slopes, due to:

- engaging stance phase flexion resistance before ground contact; increasing resistance progressively up to the end of the yielding phase
- significantly higher knee flexion moment ( $p < 0.05$ ) and thus higher load bearing capacity on the prosthetic side compared to the Rheo Knee XC; also leading to less compensation on the contralateral side
- knee joint extension position remains constant at the end of every swing phase, leading to a more precise positioning of the foot

→ Higher knee flexion during swing phase with C-Leg 4 when ascending slopes, which may lead to greater ground clearance

→ Swing phase perturbations are better compensated by C-Leg 4

- Knee joint extension after perturbation executed up to a flexion angle of approximately  $46^\circ$  at ground contact ( $11^\circ$  more than Rheo Knee XC)
- Higher flexion resistance leads to slower knee flexion and enhanced load bearing capacity, enabling more time to prevent a fall

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

Subjects: 6 active, unilateral transfemoral amputees  
 Previous prosthesis: Genium (N = 4), C-Leg 3 (N = 2)  
 Amputation causes: Trauma (all subjects)  
 Mean age: 42.5 ± 8.7 years  
 Mean time since amputation: 25.8 ± 6.0 years  
 MFCL: K3+ and K4

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"Safety against uncontrollable knee flexion and an overall reliable functionality of prosthetic knee joints is the basis for a successful clinical rehabilitation of transfemoral amputees. Safety and performance are the clinically relevant parameters that have been in the scope of this study. The objective biomechanical data measured in this context imply functional and safety-related advantages and disadvantages which can be attributed to the unequal technological concepts. The "default swing" principle used in the Rheo Knee XC offers slight advantages in the reproducibility of the swing phase release when walking with small steps, but at the disadvantage for walking backwards safely. This may lead to an uncontrolled flexion of the knee joint and thus to a fall. The C-Leg provides a reliably stable loadbearing prosthesis. Furthermore, the Rheo Knee XC shows lower safety reserves and requires increased compensation of the remaining locomotor system when walking down stairs, walking on slopes or while recovering from a stumble. These findings suggest that the technological concept used in C-Leg 4 provides an enhanced functional quality and advantages in daily-life situations compared to the technology used in Rheo Knee XC, especially concerning safety-relevant aspects." (Bellman et al., 2018)

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© 2018, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.