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Hahn A, Sreckovic I, Reiter S, Mileusnic M.

Ottobock Healthcare Products, Vienna, Austria.

First results concerning the safety, walking, and satisfaction with an innovative, microprocessor-controlled four-axes prosthetic foot.

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With Meridium compared to previous prosthetic foot:

- Improved walking on level ground (54% of subjects) and uneven ground (82%) as well as ramp ascent (97%) and descent (91%)
- Increased safety and stability while standing and walking for more than 45% of the users.
- The use of Genium seems to significantly increase the preference of Meridium for transfemoral amputees.

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

Subjects:	70 subjects (at 44 international prosthetic clinics) - 77% male, 23% female - 64% Transtibial amputation, 36% Transfemoral amputation or kneedisarticulation
Previous prosthetic foot:	85% Carbon fiber feet, 7% Microprocessor controlled feet, 4% SACH, 2% SAFE, 2% Single axis feet
Amputation causes:	66% Trauma, 13% Cancer, 13% Vascular Disease, 5% Infection, 3% Congenital
Mean age:	45.6 ± 13.7 years
Mean time since amputation:	13.9 ± 13.6 years
MFCL:	63% K3, 37% K4
対象:	70名(世界中、44カ所の義肢クリニックにおいて) - 77% 男性、23% 女性 - 64% 下腿切断、36% 大腿切断または膝離断
以前使用の足部:	85% カーボンフット、7% コンピューター制御足部、4% SACH、2% SAFE、2% 単軸足部
切断原因:	66% 外傷、13% 悪性腫瘍、13% 血管障害、5% 感染、3% 先天性
平均年齢:	45.6歳±13.7歳
切断後の平均年数:	13.9歳±13.6歳
活動レベル:	63% K3、37% K4

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Observational, prospective study:

“For 70 participants (81%), at least one patient questionnaire was completed. Therefore, data on 70 participants were used for analysis. Average Meridium use was 106.3 ± 57.04 days. To investigate effects related to trial duration, stratification for time of use differentiated between those having worn Meridium for less than 70 (46 ± 15 days; 51% participants) and for more than 70 (136 ± 35 days; 49% participants).

Fear of falling and number of falls and stumbles during the observation period were recorded. Comparative evaluations were recorded on 5-point Likert scales indicating “much better with Meridium,” “better with Meridium,” “no difference with Meridium,” “worse with Meridium,” and “much worse with Meridium.” Questions regarding pain, concentration, and exertion used the terms “more” and “less.”

Responders are classified as subjects who experienced “better” or “much better” performance with Meridium.” (Hahn et al, 2018)

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“First fittings with Meridium showed that users perceived benefits and prosthetist recommended the device. Alignment and the use of the software were rated appropriate for routine use. Meridium seems to be preferred by amputees with a preference for natural walking and the requirement of safely and comfortably negotiating uneven terrain and slopes. Subject’s preferences do not correlate with amputation level, age, and mobility grade. In transfemoral amputees, the use of Genium seems to significantly increase the perception of walking-related benefits offered by Meridium. Individual assessment and trial fittings seem to be essential to identify responders to the new foot.” (Hahn et al,2018)

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