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Kannenberg A, Zacharias B, Mileusnic M, Seyr M.

Otto Bock HealthCare GmbH, Duderstadt, Germany.

Activities of Daily Living: Genium Bionic Prosthetic Knee Compared with C-Leg

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Genium vs C-Leg

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With Genium compared to C-Leg:

- **60% of Activities of Daily Living (ADLs) show a clinically relevant gain in safety**
- **53% of ADLs show a clinically relevant decrease in difficulty.**
- **Ease of ascending and descending stairs is improved by 34% and 10%**
- **Ease of ascending and descending ramps is improved by 24% and 17%**
- **Ease of walking backwards is improved by 26%**

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
Clinically relevant improvement in safety of activities of daily living with Genium

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Subjects: 10 unilateral, transfemoral amputees
 Previous prosthesis: C-Leg
 Amputation causes: trauma
 Mean age: 36.7 yrs (±10.2 yrs)
 Mean time since amputation: 12.5 yrs (±9.6 yrs)
 MFCL: 40% K3, 60% K4

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Interventional, pre- to post-test design:

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"This study showed that the new technological functions of the Genium Bionic Prosthetic Knee for transfemoral amputees of MFCL 3 and 4 lead not only to clear bio-mechanical benefits compared with the C-Leg but also to a further improvement in the subjective perception of safety and perceived difficulty of many ADLs." (Kannenberg et al., 2013)

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