Ease of Activities of Daily Living with Conventional and Multigrip Myoelectric Hands

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Michelangelo Hand vs previous myoelectric prostheses 以前使用の筋電義手

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With Michelangelo Hand compared to previous myoelectric prostheses:


- **Perceived ease of use to perform ADLs increased by 35%**
- **Bimanual activities were easier to perform by 31%**
- **Participants used the prosthesis to actively grasp an object in more bimanual activities**

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Perceived Ease of Performing Activities of Daily Living

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Perceived ease to perform 23 activities of daily living was measured with OPUS-UEFS questionnaire. Total OPUS-UEFS score was 35% higher with Michelangelo prostheses, meaning that tasks were easier to conduct with Michelangelo hand.

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
Subjects: 16 subjects
 Previous: 10 Sensor hand speed; 3 Myohand VariPlus Speed; 1 Motion Control Hand, 1 DMC plus Myohand, Elektrogreifer
 Amputation causes: 8 traumas, 6 congenital deformities, 1 cancer and 1 sepsis
 Mean age: 41 ± 14 years
 Mean time since amputation: 12.8 ± 16.1 years

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“These results suggest that the Michelangelo hand offering more grip types and functional hand positions as well as a flexible wrist may be used more actively and for more ADLs. These effects seem to be promoted by a reduced perceived difficulty to perform many ADLs as Michelangelo shortens the functional gap between prosthetic and sound human hands.”

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© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.