

1C56 Taleo Adjust. Ready for everyday life.

The *Taleo Adjust* was designed for active users who navigate a wide range of indoor and outdoor environments, and who place a high value on flexibility in the choice of footwear. At the touch of a push button, users can set the appropriate heel height from 0 up to 7 cm.

Designed for optimal outcomes

- Heel height adjustable carbon foot
- The easy, stepless heel height adjustment from 0 to 7 cm allows users to wear different shoes. The individually optimal prosthetic alignment can be reproduced.
- The seamless rollover from heel strike to toe-off allows easy and comfortable walking without dead spots – depending on the footwear.
- Thanks to three different heel wedges, both the impact stress at heel strike and the gait dynamics can be customised to individual needs.
- Openings in the sole of the footshell prevent water from collecting in the prosthesis.
- Protected against fresh, salt and chlorinated water.

More than a foot. A foundation.



Mobility grade	2, 3
Max. body weight	115 kg
Side	left (L), right (R)
Sizes	22 – 28 cm
Weight without footshell*	680 g
Footshell shape	normal, normal low cut, slim low cut
Footshell colour	beige (4), light brown (15)
Weight with normal footshell*	930 g
Build height with normal footshell*	118 mm
Heel height*	0 – 7 cm

* Reference size is 26 cm.



Easy adjustment of the heel height

Wear a variety of shoes

Proven Taleo Low Profile base spring

1C56. Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Selection of the spring stiffness relative to

1. body weight and activity as well as

Body weight (kg)	Low activity level	Normal activity level
Up to 51	1	2
52 to 58	2	3
59 to 67	3	4
68 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 115	7	-

2. foot size

Stiffness	Foot size						
	1	2	3	4	5	6	7
22						-	-
23							-
24							-
25							
26							
27	-	-					
28	-	-					

no restrictions

Example for ordering with 2C8* footshell:

1C56=R26-4-P/4

Article no.	Side	Size	Stiffness	Colour
1C56	=		-	P /

Side	Size (cm)	Stiffness	Colour	Side
R Right	22, 23, ..., 28	1, 2, ..., 7	4 Beige	R Right
L Left			15 Light brown	L Left

Example for ordering with 2C15* low cut footshell:

1C56=R26-4-P/4N-L

Article no.	Side	Size	Stiffness	Colour	Shape
1C56	=		-	P /	

Side	Size (cm)	Stiffness	Colour	Shape
R Right	22, 23, ..., 26	1, 2, ..., 7	4 Beige	N-L Normal low cut (24-26 cm)
L Left			15 Light brown	S-L Slim low cut (22-23 cm)



The scope of delivery includes the 1C56 foot module, the 2C8 footshell including connection cover or the 2C15 low cut footshell, 2F50 heel wedge set (comprising three different degrees of hardness) and a black Spectra sock.

Recommended combinations*

TT	TF
 4R140 One-way valve + 453A30/40 ProFlex Plus sealing sleeve + 6Y95 Caleo 3D or 4R180 Harmony P4 + 453A30/40 ProFlex Plus sealing sleeve + 6Y512=*-P Uneo Fresh	 7T451=1 SiOCX TF Pro or 6Y110 Skeo Sealing + 21Y21 Click Valve
 4R11=* Quickchange	 3C60* Kenevo, 3R85* Dynion, 3C98* C-Leg, 3B1* Genium or 3B5* Genium X3  4R57=WR Waterproof Rotation Adapter 4R11=* Quickchange
	 3F1=* + 99B120=* Functional cosmesis

* All components are sold separately and are available Ottobock products that are compatible with the 1C56 *Taleo Adjust* foot, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.