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Luchetti M., Cutti AG, Verni G, Sacchetti R, Rossi N.

Department of Psychology, University of Bologna, Italy.

Impact of Michelangelo prosthetic hand: Findings from a crossover longitudinal study

Journal of Rehabilitation Research & Development 2015; 52(2):605-618.

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With Michelangelo compared to conventional myoelectric prosthesis:

→ A higher functionality can be achieved

Score of Box and Block Test increased by 20.8%

Score of SHAP increased by 11.4%

Time needed to perform Minnesota Manual Dexterity Test decreased by 14.8%

→ ADL execution is easier

Reported by 84% of patients

→ Michelangelo hand is more actively used at home

Lateral grip preferred over opposition grip

→ Gesture and posture is more natural

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
Box and Blocks is a manual dexterity test where number of blocks transported from one box to another in 60s is assessed. The users were able to transport 5 blocks more on average with Michelangelo hand. Four out of six participants (67%) had a sore above minimal clinically relevant detectable change (more than 6.5 blocks were transported with Michelangelo hand in 60s than with conventional myoelectric prosthesis).

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Previous prosthesis: tridigital myoelectric prosthesis

Amputation causes: trauma

Mean age: 47 yrs (range: 35-65 yrs)
Mean time since amputation: 15 yrs (range: 4.5-48 yrs)

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Interventional pre- to post-test design:

The subjects were provided with 5 days of occupational therapy after they have been fitted with Michelangelo.

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"Amputation-and prosthetic-related factors, along with psychological factors (e.g., patient coping strategies, attitude, expectations) and social factors (i.e., support of family and friends, reactions of others), need to be screened in the prosthesis fitting process. The present study shows that the M is effective in improving the functional ability and in easing the social interaction of previous active users of a myoelectric prosthesis." (Luchetti et al. 2015)

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© 2014, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.