
3R15, 3R49

DE Gebrauchsanweisung	10	PL Instrukcja użytkowania.....	67
EN Instructions for use	16	HU Használati utasítás	74
FR Instructions d'utilisation.....	22	CS Návod k použití.....	80
IT Istruzioni per l'uso	29	RO Instrucțiuni de utilizare.....	86
ES Instrucciones de uso	35	HR Upute za uporabu	92
PT Manual de utilização.....	42	TR Kullanma talimatı	98
NL Gebruiksaanwijzing.....	48	EL Οδηγίες χρήσης	104
SV Bruksanvisning	55	JA 取扱説明書.....	111
DA Brugsanvisning	61	ZH 使用说明书.....	116

1 **i** 5 Lieferumfang/
Scope of delivery

3R15/3R49 506G1=M8x25

21Y67


21Y113=2

21Y70=N



513D28 (2x)


647H84

21Y66

2 **i** 6 Technische Daten/
Technical data


1 - 2


≤100 kg
≤220 lbs

3R15: 490 g
3R49: 360 g


150°

75 mm

h_{s-prox} : 8 mm

h_{s-dist} : 1 mm

①

Proximale Systemhöhe
Proximal system height




②

Distale Systemhöhe
Distal system height

3

7.2 Grundaufbau der Prothese/ Bench alignment of the prosthesis




Bench alignment of the prosthesis

4


7.2 Grundaufbau der Prothese/ Bench alignment of the prosthesis

Bench alignment of the prosthesis

5

7.2 Grundaufbau der Prothese/ Bench alignment of the prosthesis

6

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INFORMATION

Last update: 2018-12-18

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

This document and the product are intended exclusively for the fabrication of a prosthesis by orthopaedic technicians with technical knowledge of lower limb prosthetic fittings.

1.1 Explanation of warning symbols

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warning regarding possible technical damage.

1.2 Meanings of pictograms in the illustration

	Numbering for the illustrations		Reference to the number of the corresponding section
	Numbering for a defined sequence		Numbering for the parts of an illustration
	Right		Wrong

2 Product description

The product (3R15, 3R49) has the following key features:

- Monocentric prosthetic knee joint
- Product components for stance phase stability:
 - Adjustable, load-dependent braking mechanism to secure the extension position
- Product components for swing phase control:
 - Extensor (spring force adjustable)
- Material: Stainless steel (3R15), Titanium (3R49)

3 Intended use

3.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

3.2 Area of application

CAUTION

Excessive strain on the product

Fall due to breakage of load-bearing components

- ▶ Only use the product according to its allowable field of application.

Allowable field of application (3R15)

Recommended mobility grade: 1 + 2	Everyday prosthesis
Allowable body weight: ≤ 100 kg	

MG 1: indoor walker

The patient has the ability or the potential to use the prosthesis for the purpose of moving or proceeding slowly on level floors. The amount of time and the distance he/she can walk are seriously limited due to his/her condition.

MG 2: Restricted outdoor walker

The patient has the ability or the potential to walk slowly with the prosthesis and to negotiate low environmental obstacles like curbs, single steps or uneven surfaces. The amount of time and the distance he/she can walk are seriously limited due to his/her condition.

Unallowable amputation types

Hip disarticulation, Hemipelvectomy

3.3 Combination possibilities

CAUTION

Failure to observe the manufacturer requirements regarding combination possibilities

Injuries, malfunctions or product damage due to unallowable combination of prosthesis components

- ▶ Based on the instructions for use of all prosthetic components used, verify that they may be combined with each other and are approved for the patient's field of application.
- ▶ Contact the manufacturer with questions.

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

Unallowable combinations

Prosthetic hip joints	All reference numbers
-----------------------	-----------------------

3.4 Environmental conditions

CAUTION

Use under unallowable environmental conditions

Fall due to damaged product

- ▶ Do not expose the product to unallowable environmental conditions (see the table "Unallowable environmental conditions" in this section).
- ▶ If the product was exposed to unallowable environmental conditions, take suitable steps (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Allowable environmental conditions

Temperature range: -10 °C – +60 °C

Use - Relative humidity: 0 % – 90 % (Non-condensing)

Storage/Transportation - Relative humidity: 100 %

Tyre pressure: 250 – 1100 mbar

Commercially available, solvent-free cleaning agents

Unallowable environmental conditions

Storage/Transportation: Mechanical vibrations, Impacts

Hygroscopic particles (e.g. Talcum powder), Dust, Sand, Fresh water, Salt water, Acids, Sweat, Urine

Cleaning agents containing solvents

3.5 Reuse and service life

⚠ CAUTION

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

⚠ CAUTION

Exceeding the service life

Fall due to change in or loss of functionality and damage to the product

- ▶ Ensure that the service life defined in this section is not exceeded.

Maximum service life for this product is 3 to 5 years, depending on the patient's level of activity.

4 General safety instructions

⚠ CAUTION

Reaching into the area of the joint mechanism

Pinching of limbs (e.g. fingers) and the skin due to uncontrolled joint movement

- ▶ Do not reach into the joint mechanism during daily use.
- ▶ Close attention is required during assembly and adjustment tasks.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Changes in functionality can manifest themselves, for example, through a changed gait (swing phase, stance phase), incomplete extension, stiffness and the development of noise.

5 Scope of delivery

The scope of delivery is shown in Figure [1](#) on page 2.

Only product components with reference numbers in the illustration can be reordered separately.

6 Technical data

The technical data are shown in Figure [2] on page 2.

①: Proximal system height (PSH); ②: distal system height (DSH)

7 Preparation for use

7.1 Information on fabrication of a prosthesis

⚠ CAUTION

Incorrect alignment, assembly or adjustment

Injuries due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Initial use of the prosthesis by the patient

Fall due to lack of patient experience, incorrect alignment or incorrect adjustment of the prosthesis

- ▶ For the safety of the patient, use a suitable device (e.g. parallel bars, handrail, wheeled walker) during initial standing and walking.

7.2 Bench alignment of the prosthesis

INFORMATION

The prosthetic components and prosthesis shown in the illustrations are used as examples to illustrate the general process. The instructions for use of the prosthetic components selected for the patient contain detailed information and have to be used for fabrication of the prosthesis.

- ▶ **CAUTION!** To enable the patient to stand safely, align the prosthesis using the alignment images for the sagittal plane (Figure [3] on page 3) and for the frontal plane (Figure [4] on page 3) as well as the instructions for use of all prosthetic components that were used.

- ▶ Additionally, observe the following points for the bench alignment of the prosthesis:

- ▶ [5] on page 4: Mount the prosthetic joint with its extension assist cable as shown in the illustration.

7.3 Optimising the static alignment

- ▶ To optimise the static alignment, position the patient on the measuring device.
- ▶ In the sagittal plane, check the course of the load line relative to the marked alignment reference point and correct this if needed (see Figure [6] on page 4).
- ▶ In the frontal plane, check the course of the load line relative to the marked points (prosthetic foot, prosthetic knee joint, spina iliaca anterior superior) and correct if needed (see Figure [6] on page 4).

7.4 Optimising during dynamic trial fitting

⚠ CAUTION

Adjusting the Settings

Fall due to incorrect or unfamiliar settings

- ▶ Only adapt the settings to the patient gradually.
- ▶ Explain the effects of the adjustments on the use of the prosthesis to the patient.

Restore default settings

- ▶ Prior to the dynamic trial fitting, restore the factory settings as shown in the illustrations if the factory settings have been changed ([7], [8] and [9] on pages 5 and 6).
- ▶ [9] on page 6: ① Adjust gradually, ② Minimal braking action
When resetting, set the axis friction gradually (45°) to minimal resistance to motion ④ and then turn back slightly without creating play ⑤.

Setting the load-related braking action ([10] on page 6)

- ▶ Set the load-related braking action as shown in the illustration.

- ① Increase threshold – Joint brakes at a higher load
- ② Reduce threshold – Joint brakes at a lower load
- ▶ Increase gradually (rotation: approx. 45° – 90°).
- ▶ When the desired effect is reached, continue with "Setting axis friction".

Setting axis friction (11 on page 7)

- ▶ Set the axis friction as shown in the illustration.
 - ① Decrease play – Increase resistance to motion
 - ② Increase play – Decrease resistance to motion
- ▶ Reduce play and increase resistance to motion gradually (rotation: approx. 20°).
- ▶ Set the load-related braking action and axis friction so that the prosthetic knee joint switches from the stance phase to the swing phase and no longer jams.
- ▶ **NOTICE!** To prevent increased wear and tear, do not use axis friction to control the swing phase.

Setting the extension assist (12 on page 7)

- ▶ Set the extension assist as shown in the illustration if necessary.
 - ① Increase spring tension
 - ② Reduce spring tension
- ▶ **CAUTION!** Set the extension assist so that the lower leg of the prosthesis swings through in the dorsal direction according to the patient's requirements, but reaches full extension in time before the next heel strike.

Fine-tuning the settings

- ▶ **CAUTION!** The various settings of the prosthetic knee joint cannot be adjusted entirely independently of each other. When the settings cannot be fully adjusted to the comfort needs of the patient, the settings should be established primarily based on safety aspects.
- ▶ Adapt the settings of the prosthetic knee joint to the patient by means of fine-tuning and exercises.

- ▶ If the extension position is not reached with every step and at any walking speed, reduce the braking effect and/or increase the initial tension of the extension assist.
- ▶ Check the adjustment of the prosthesis settings during normal consultations and the annual safety inspection. Advise the patient to check the prosthesis if changes in function occur.

7.5 Finishing the prosthesis

Put on protective cover

- ▶ When the prosthesis is finished, put on the protective cover as shown in the illustration [13] on page 8.
 - ①: Alternative positioning
 - ②: Recommended positioning
- ▶ **CAUTION!** To avoid product damage and the risk of falling, finish the prosthesis by replacing set screws that are too short or too long, and by tightening all screw connections. In doing so, note the instructions for use for all prosthesis components regarding torque values and thread lock.

7.5.1 Attaching the cosmetic cover

CAUTION

Use of talcum

Risk of injury, damage to the product due to lack of lubrication

- ▶ Do not use talcum on the product or other prosthesis components.

Recommended combinations

Foam cover	3R6, 3S106
------------	------------

- ▶ After affixing the cosmetic elements, verify proper functioning of the prosthesis.
- ▶ To reduce friction and to eliminate noise, apply silicone spray directly onto the contact surfaces of the cosmetic foam cover.

8 Use

8.1 Information for use

NOTICE

Mechanical overload

Impaired functionality due to mechanical damage

- ▶ Check the product for damage prior to each use.
- ▶ Do not use the product if its functionality has been impaired.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

8.2 Cleaning

CAUTION

Use of unsuitable cleaning agents or disinfectants

Impairment of functionality and damage due to incorrect cleaning agents or disinfectants

- ▶ Clean the product only according to the instructions given in this section.
- ▶ For the prosthesis, observe the cleaning instructions for all prosthetic components.

INFORMATION

- ▶ Clean the product if it gets dirty.
- ▶ For disinfection, only use disinfectants that do not affect the materials of the product. For further information, please consult the manufacturer.
- ▶ Observe the cleaning instructions for all prosthetic components.

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

9 Maintenance

CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.

- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ Following an individual period for the patient to get accustomed to the prosthesis, check the settings of the prosthetic joint and adapt them to the patient's requirements again as needed.
- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ As part of the safety inspections, inspect the prosthetic joint for wear and proper functionality. Special attention should be paid to movement resistance, bearings and abnormal noises. Full flexion and extension must be ensured at all times. Readjust as required.

Lubrication

- ▶ **14** on page 8: Do not lubricate or grease the prosthetic joint.

Replacing components


CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.

► Follow the instructions regarding the length of the screws and about how to secure the screws.

- When necessary, remove a worn brake bushing (4B74=20.4X24.5-K) as shown in the illustration  on page 9.
- To reassemble, perform these steps in reverse order.
- In doing so, note the following points:
 -  Ensure the correct assembly of the bushing safety device.
 -  Clean the knee axis with isopropyl alcohol ( 634A58) before assembly.

10 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

11 Legal information

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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