C-Brace joint unit 17K01-2=*

DE	Gebrauchsanweisung (Benutzer)	3
EN	Instructions for use (user)	41
FR	Instructions d'utilisation (utilisateur)	79
IT	Istruzioni per l'uso (utilizzatore)	119
ES	Instrucciones de uso (usuario)	157
PT	Instruções de utilização (utilizador)	195
NL	Gebruiksaanwijzing (gebruiker)	233
TR	Kullanım kılavuzu (kullanıcı)	271
EL	Οδηγίες χρήσης (Χρήστης)	307
JA	取扱説明書(装着者用)	347

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1	Foreword	43
2	Product description	43
2.1	Design	43
2.1.1	Control panel on the joint unit	44
2.2	Function	44
3	Intended use	45
3.1	Indications for use	45
3.2	Conditions of use	45
3.3	Indications	45
3.4	Contraindications	45
3.4.1	Absolute Contraindications	45
3.4.2	Relative Contraindications	45
3.5	Qualification	45
4	Safety	45
4.1	Explanation of warning symbols	45
4.2	Structure of the safety instructions	46
4.3	General safety instructions	46
4.4	Information on the Power Supply/Battery Charging	48
4.5	Battery charger information	48
4.6	Notes on applying the product	49
4.7	Information on Proximity to Certain Areas	49
4.8	Information on Use	50
4.9	Notes on the movement patterns	51
4.10	Notes on the safety modes	53
5	Scope of Delivery and Accessories	54
5.1	Scope of delivery	54
5.2	Accessories	54
6	Charging the battery	55
6.1	Connecting the power supply and battery charger	55
6.2	Connect battery charger to the product	56
6.3	Display of the current charge level	56
6.3.1	Display of battery charge level without additional devices	56
6.3.2	Display of the current charge level during the charging process	57
7	Use	57
7.1	Application	57
7.2	Removal	58
7.3	Movement patterns in basic mode (mode 1)	58
7.3.1	Standing	58
7.3.1.1	Stance function	59
7.3.2	Walking	59
7.3.3	Sitting down	59
7.3.4	Sitting	60

7.3.5	Standing up	60
7.3.6	Walking down stairs	60
7.3.7	Walking up stairs.....	61
7.3.8	Walking up a ramp	61
7.3.9	Walking down a ramp	61
7.3.10	Walking down flat steps	61
7.3.11	Kneeling	62
7.3.12	Intuitive cycling (region-specific)	62
7.4	Switching the product on/off	63
7.5	Bluetooth.....	63
7.5.1	Pairing with a device	63
7.5.2	Switching on Bluetooth	64
7.5.3	Switching off Bluetooth	64
7.6	Recommendations for air travel	64
8	MyModes	64
9	Additional operating states (modes)	64
9.1	Empty battery mode	64
9.2	Mode for charging the product	65
9.3	Safety mode.....	65
9.4	Overheating mode	65
10	Cleaning	65
11	Maintenance	66
12	Legal information	66
12.1	Liability	66
12.2	Local Legal Information	66
12.3	CE conformity	67
12.4	Trademarks	67
13	Technical data	68
14	Appendices	70
14.1	Symbols Used.....	70
14.2	Operating states/error signals.....	70
14.2.1	Status display on the control panel	70
14.2.2	Error while charging the product	74
14.3	Directives and manufacturer's declaration	75
14.3.1	Electromagnetic environment	75

1 Foreword

INFORMATION

Date of last update: 2025-11-17

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Obtain instruction from the qualified personnel in the safe use of the product.
- ▶ Please contact the qualified personnel if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "17KO1-2=* C-Brace joint unit" is referred to as the product/component/orthosis/joint unit below.


These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.


2 Product description

2.1 Design

The fabricated orthosis and the course of the orthosis shell edges can be individually designed. For this reason, the following illustration and the illustration on the cover page of these instructions for use only show one possible version:

2.1.1 Control panel on the joint unit

1. **Long button press:** turns the component on/off (see page 63)
Short button press: queries the status (see page 70)
Query whether the Bluetooth function is switched on (symbol in the button ✱ lights up blue ✱) or switched off (symbol not lit up).
2. Displays the charge level of the integrated battery (see page 56)
3. Maintenance required (see page 70)
4. **Long button press:** turns the Bluetooth function on/off (Switching Bluetooth of the component on/off).
Short button press: if the Bluetooth function is on, the component becomes "visible" for about 2 minutes so it can be detected by another device such as a smartphone.

Under the cover at the end of the control panel: charging receptacle to connect the battery charger (see page 55)

2.2 Function

The product is an orthopaedic technology device that makes everyday activities such as standing and walking possible or easier for the user. If the user's activity level increases in the course of use, adjustments are required. Changes will also have to be made if the user's physical condition worsens and more support from the product is needed.

This product features microprocessor control of the stance and swing phase (SSCO).

The microprocessor uses the measurements of an integrated sensor system as a basis to control a hydraulic unit that influences the extension and flexion resistance of the product.

These sensor data are updated and evaluated 100 times per second. As a result, the behaviour of the product is adapted to the current motion situation (gait phase) dynamically and in real time.

Thanks to the microprocessor-controlled stance and swing phase, the product can be individually adapted to your needs.

The product is adjusted by the O&P professional using the adjustment app.

The "connectgo 560X34-*=" user app (see the "Accessories" section (see page 54), can be used to switch from basic mode to pre-configured MyModes, query information about the component (step counter, charge level, etc.) and adjust the component settings within a certain range.

The product features MyModes for special motion types (e.g. golf, basketball ...). These are pre-set by the O&P professional and can be called up via the user app.

In case of an error in the sensor system, hydraulic control or when the battery is empty, safety mode makes restricted operation possible. Resistance parameters that are predefined by the product are configured for this purpose (see page 65).

The microprocessor-controlled hydraulic unit offers the following advantages

- Approximation of the physiological gait pattern
- Stability while standing and walking
- Adaptation of product characteristics to various surfaces, inclines, gait situations and walking speeds

Essential performance of the product

- Stability in the stance phase

3 Intended use

3.1 Indications for use

The product is intended **solely** for orthotic fittings of the lower limbs.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).

The O&P professional must be consulted about any sports activities.

See the technical data for the permissible environmental conditions (see page 68).

The product is intended **exclusively** for use on **one** user. Use of the product by another person is not approved by the manufacturer.

3.3 Indications

- Unilateral or bilateral lower limb paresis or flaccid paralysis, e.g. due to post-polio syndrome, traumatic paresis, incomplete paraplegia.
- Physical prerequisites such as muscle status, joint mobility and possible axis deviations are crucial and proper control of the orthosis must be guaranteed.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations
- The existing muscle strength of the hip extensors and flexors must permit the controlled swing-through of the limb (compensation using the hip is possible).

3.4 Contraindications

3.4.1 Absolute Contraindications

- Flexion contracture in the knee and/or hip joint in excess of 10°
- Knee varus/valgus malposition in excess of 10°
- Severe spasticity
- Frequent spasticity
- Body weight over 125 kg
- Orthoprosthesis

3.4.2 Relative Contraindications

- Uncontrolled moderate spasticity
- Infrequent spasticity
- Leg length discrepancy in excess of 15 cm

3.5 Qualification

The product may only be fitted by certified qualified personnel.

4 Safety

4.1 Explanation of warning symbols

 WARNING	Warning regarding possible serious risks of accident or injury.
 CAUTION	Warning regarding possible risks of accident or injury.
 NOTICE	Warning regarding possible technical damage.

4.2 Structure of the safety instructions

WARNING

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

WARNING

Use of the product while operating a vehicle

Risk of accident due to unexpected behaviour of the product.

- ▶ All users are required to observe their country's national and state driving laws when operating vehicles with the orthosis. For insurance purposes, drivers should have their driving ability examined and approved by an authorised test centre.
- ▶ The leg on which the orthosis is worn may not be used to control the vehicle or its peripheral components (e.g. clutch pedal, brake pedal, gas pedal...).

WARNING

Use of damaged power supply unit, adapter plug or battery charger

Risk of electric shock due to contact with exposed, live components.

- ▶ Do not open the power supply unit, adapter plug or battery charger.
- ▶ Do not expose the power supply unit, adapter plug or battery charger to extreme loading conditions.
- ▶ Immediately replace damaged power supply units, adapter plugs or battery chargers.

CAUTION

Failure to observe warning/error signals

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (Warnings/error signals) and corresponding change in damping settings must be observed.

CAUTION

Penetration of dirt and moisture into components of the orthosis

- > Falling due to unexpected behaviour of the orthosis as a result of malfunction.
- > Falling due to breakage of load-bearing components.
- ▶ Ensure that neither solid particles, foreign bodies nor liquids penetrate the orthosis components.
- ▶ The orthosis components are protected against splashed water from all directions.
- ▶ However, the orthosis components are not protected against submersion, jets of water and steam.
- ▶ If water has penetrated the components of the orthosis, remove the straps and pads if possible and allow the components to dry.
- ▶ In case of a malfunction after drying, the orthosis must be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.

- ▶ Always close the charging receptacle with the protective cap after disconnecting the charging cable.

CAUTION

Independent manipulations on the joint unit and orthosis components

Falling due to the breakage of load-bearing components or malfunction of the orthosis.

- ▶ Manipulations on the joint unit and orthosis components other than the tasks described in these instructions for use are not permitted.
- ▶ The battery may only be handled by Ottobock authorised, qualified personnel (no replacement by the user).
- ▶ The joint unit and orthosis components may only be opened and repaired by authorised, qualified Ottobock personnel.

CAUTION

Use of the product when battery charge level is too low

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Check the current charge level before use and charge the product if required.
- ▶ Note that the operating time of the product may be reduced at low ambient temperatures or due to ageing of the battery.

CAUTION

Mechanical stress on the product

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

CAUTION

Signs of wear and tear on the product components

Falling due to damage or malfunction of the product.

- ▶ Regular service inspections (maintenance) are mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty.

CAUTION

Use of unapproved accessories

- > Falling due to product malfunction as a result of reduced interference resistance.
- > Interference of other electronic devices due to increased emissions.
- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the sections "Scope of delivery" (see page 54) and "Accessories" (see page 54).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Clean the product with a damp cloth only (fresh water).

4.4 Information on the Power Supply/Battery Charging

CAUTION

Charging the product with damaged power supply unit/charger/charger cable

Falling due to unexpected behaviour of the product caused by insufficient charging.

- ▶ Check the power supply unit, charger and charger cable for damage before use.
- ▶ Replace any damaged power supply unit, charger or charger cable.

CAUTION

Charging while wearing the product

- > Falling due to walking and getting caught on the connected battery charger.
- > Falling due to unexpected product behaviour because of changed damping behaviour.
- ▶ Only charge the product in a sitting position while wearing it.

NOTICE

Use of incorrect power supply unit/battery charger

Damage to product due to incorrect voltage, current or polarity.

- ▶ Use only power supply units/battery chargers approved for this product by Ottobock (see instructions for use and catalogues).

4.5 Battery charger information

NOTICE

Penetration of dirt and humidity into the product

Lack of proper charging functionality due to malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

NOTICE

Mechanical stress on the power supply/battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Do not subject the power supply/battery charger to mechanical vibrations or impacts.
- ▶ Check the power supply/battery charger for visible damage before each use.

NOTICE

Operating the power supply unit/charger outside of the permissible temperature range

Lack of proper charging functionality due to malfunction.

- ▶ Only use the power supply unit/charger for charging within the allowable temperature range. The section "Technical data" contains information on the allowable temperature range (see page 68).

NOTICE

Independent changes or modifications carried out to the battery charger

Lack of proper charging functionality due to malfunction.

- ▶ Have any changes or modifications carried out only by Ottobock authorised, qualified personnel.

4.6 Notes on applying the product

CAUTION

Foreign objects between the leg and orthosis shells

Pressure points on the leg due to foreign objects at the contact points between the leg and orthosis shells.

- ▶ Smooth out wrinkles in the padding material and clothing.
- ▶ Check the leg for pressure points.

CAUTION

Improper fit of the orthosis

Falling/skin irritation due to insufficient support/hold by the orthosis shells.

- ▶ Take the orthosis off immediately and put it back on.
- ▶ Observe the instructions for applying and removing.

CAUTION

Pinching skin in the area of the closures

Injuries and swelling due to disruption of blood circulation to the skin in the area of the closures.

- ▶ Do not over-tighten the closures when applying the product.

CAUTION

Volume fluctuations of the leg or problems with the fit of the orthosis

Injuries, friction and local pressure caused by improper fit of the orthosis shells (too firm/too loose).

- ▶ Sensitivity problems and skin damage require special attention with regard to the fit. Check skin areas daily.
- ▶ Consult a physician or O&P professional in case of any signs of skin damage, even if they are only minor.
- ▶ If the improper fit is the result of weight gain or loss, have new thigh shells and lower leg shells fabricated based on a new plaster cast.
- ▶ Check the leg for pressure points.

4.7 Information on Proximity to Certain Areas

CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Maintaining a minimum distance of **30 cm** from HF communication devices is therefore recommended.

CAUTION

Operating the product in very close proximity to other electronic devices

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Do not operate the product in the immediate vicinity of other electronic devices.
- ▶ Do not stack the product with other electronic devices during operation.

- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for people (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, ...).

If remaining in these areas cannot be avoided, make sure to at least walk/stand in a safe manner (e.g. by using a handrail or holding onto another person).

- ▶ When passing through theft prevention systems, body scanners or metal detectors, watch out for unexpected changes in the damping behaviour of the product.
- ▶ In general, monitor the product for unexpected changes in the damping behaviour when electronic or magnetic devices are in the immediate vicinity.
- ▶ When using the product in an airplane, pay attention to unexpected changes in the damping behaviour of the product

CAUTION

Entering a room or area with strong magnetic fields (e.g. magnetic resonance tomographs, MRT (MRI) equipment...)

- > Falling due to unexpected restriction of the product's range of motion caused by metallic objects adhering to the magnetised components.
- > Irreparable damage to the product due to the effect of a strong magnetic field.
- ▶ Take off the product before entering a room or area with strong magnetic fields and store the product outside this room or area.
- ▶ Damage to the product caused by exposure to a strong magnetic field cannot be repaired.

CAUTION

Remaining in areas outside the allowable temperature range

Falling due to breakage of load-bearing components or malfunction of the product.

- ▶ Avoid remaining in areas with temperatures outside of the permissible range (see page 68).

4.8 Information on Use

CAUTION

Improper mode switching

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Ensure that you stand securely during all switching processes.
- ▶ Verify the changed damping characteristics after switching and observe the feedback from the acoustic signal emitter.
- ▶ Switch back to basic mode once the activities in MyMode have been completed.
- ▶ Take the weight off the product and correct the switching, if required.

CAUTION

Risk of pinching in the joint flexion area

- > Injury due to pinching of body parts.

- > Damage to garments due to pinching in the mechanism of the joint unit or medial support.
- ▶ Ensure that no body parts or garments are in this area when flexing the product.

CAUTION

Overheating of the hydraulic unit due to uninterrupted, increased activity (e.g. extended walking downhill)

- > Falling due to unexpected behaviour of the product because of switching into overheating mode.
- > Burns due to touching overheated components.
- ▶ Be sure to pay attention when pulsating vibration signals start (Warnings/error signals). They indicate the risk of overheating.
- ▶ As soon as these pulsating vibration signals begin, you have to reduce the activity level so the hydraulic unit can cool down.
- ▶ Full activity may be resumed after the pulsating vibration signals stop.
- ▶ If the activity level is not reduced in spite of the pulsating vibration signals, this could lead to the hydraulic element overheating and, in extreme cases, cause damage to the product. In this case, the product has to be inspected by an authorised Ottobock Service Centre. The O&P professional is your contact person.

CAUTION

Overloading due to unusual activities

- > Falling due to unexpected product behaviour as the result of a malfunction.
- > Falling due to breakage of load-bearing components.
- > Skin irritation due to defects on the hydraulic unit with leakage of liquid.
- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, extreme sports (free climbing, parachuting, paragliding, etc.), sports activities that involve jumping, sudden movements or fast sequences of steps (e.g. basketball, badminton, sporty riding).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures your personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately by an O&P professional.

CAUTION

Improper use of the “Intuitive cycling” function

Falling due to unexpected product behaviour caused by changed damping behaviour.

- ▶ The “**Intuitive cycling**” function may only be used by experienced cyclists.
- ▶ After getting off the bicycle, check whether the extension and flexion resistances for walking and standing are set.

4.9 Notes on the movement patterns

CAUTION

Improper use of the stance function

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Make sure that you are standing safely when using the stance function and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Make sure you have been instructed in the proper use of the stance function by an O&P professional or therapist. Information on the stance function see page 59.

CAUTION

Walking up stairs

- > Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.
- > Falling due to breakage of load-bearing components.
- ▶ Always use the handrail when walking up stairs and place most of the sole of your foot on the stair surface.
- ▶ Always step up with the less affected leg, followed by the leg with the product.
- ▶ If there is a connection to a foot component, support by bouncing with the product extended while walking up stairs is not allowable and therefore has to be avoided.
- ▶ Particular caution is required when carrying children up the stairs.

CAUTION

Walking down stairs with moveable orthosis foot component

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and place most of the sole of your foot on the stair surface.
- ▶ Observe the warnings/error signals (see page 70).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warning and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

CAUTION

Walking down stairs with rigid orthosis foot component

Falling due to foot being placed incorrectly on stair as a result of changed damping behaviour.

- ▶ Always use the handrail when walking down stairs and roll over the edge of the step with the middle of the shoe.
- ▶ Observe the warnings/error signals (see page 70).
- ▶ Be aware that resistance in the flexion and extension direction can change in case of warnings and error signals.
- ▶ Particular caution is required when carrying children down the stairs.

CAUTION

Walking down stairs and ramps

Falling due to unexpected, high stance phase damping during the transition from a level surface to stairs or ramps.

- ▶ Pay attention to the changed product behaviour.
- ▶ Verify changed stance phase damping before walking on stairs or ramps.

CAUTION

Lack of support from the product when walking down stairs

Falling due to unexpected product behaviour caused by changed damping behaviour.

- ▶ Check whether corresponding flexion of the joint is possible before walking down stairs. If this is not possible, switch back to basic mode by either using the user app or turning the product off/on.

CAUTION

Improper use of the MyMode “Training mode”

Falling due to unexpected product behaviour caused by changed damping behaviour.

- ▶ Make sure that you are standing safely when using this MyMode and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Note that the knee joint is locked in the flexion direction in this MyMode.
- ▶ Make sure you have been instructed in the correct use of this MyMode by the O&P professional and/or therapist. Switch back to basic mode once the activities in this MyMode have been completed.

CAUTION

Improper use of the MyMode “Freeze position”

Falling due to unexpected product behaviour caused by changed damping behaviour.

- ▶ Make sure that you are standing safely when using this MyMode and check the lock of the knee joint before placing your full weight on the orthosis.
- ▶ Note that the knee joint is locked in the both the flexion and extension direction in this MyMode.
- ▶ Make sure you have been instructed in the correct use of this MyMode by the O&P professional and/or therapist.
- ▶ Switch back to basic mode once the activities in this MyMode have been completed.

CAUTION

Lack of support from the product while sitting down

Falling due to unexpected product behaviour caused by changed damping behaviour.

- ▶ Check whether corresponding flexion of the joint is possible before sitting down. If this is not possible, switch back to basic mode by either using the user app or turning the product off/on.

CAUTION

Quickly pushing the hip forward with the orthosis extended

- > Falling due to unexpected activation of the swing phase.
- ▶ Note that the joint may flex unexpectedly when the hip is pushed forward quickly while the orthosis is extended.
- ▶ Therefore, familiarise yourself with swing phase activation in such situations under secure conditions (e.g. while holding on to parallel bars) and with the instruction of trained, qualified personnel.

CAUTION

Changing the gait pattern influences activation of the swing phase

- > Falling due to unexpected activation of the swing phase.
- ▶ Changing the gait pattern can influence activation of the swing phase. Adjustment by the O&P professional has to be repeated in this case.

4.10 Notes on the safety modes

CAUTION

Safety mode cannot be activated due to malfunction caused by water penetration or mechanical damage

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Do not continue using the defective product.
- ▶ Consult the O&P professional promptly.

CAUTION

Safety mode cannot be deactivated

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ If you cannot deactivate safety mode by recharging the battery, a permanent malfunction has occurred.
- ▶ Do not continue using the defective product.
- ▶ The product must be inspected immediately by an O&P professional.

CAUTION

Safety signal occurs (ongoing vibration)

Falling due to unexpected behaviour of the product because of changed damping behaviour.

- ▶ Observe the warnings/error signals (see page 70).
- ▶ Do not continue using the product after the safety signal has been emitted.
- ▶ The product must be inspected by an authorised Ottobock Service Centre. In this case, the contact person is your prosthetist.

CAUTION

Using the product in safety mode

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ The warnings/error signals (see page 70) have to be observed.
- ▶ Particular caution is necessary when using a bicycle without a freewheel (with a fixed gear).

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 757L16-4 power supply
- 1 pc. 4E50 Battery Charger for C-Leg
- 1 pc. 646C107 Bluetooth PIN card
- 1 pc. Orthotic passport
- 1 pc. 17KO1-2=L C-Brace joint unit left or 17KO1-2=R C-Brace joint unit right
- 1 pc. instructions for use (user)

5.2 Accessories

- USB charging adapter: 757L43
To connect the 757L43 USB charging adapter to the respective charger, observe the information in the USB charging adapter instructions for use.
- **“connectgo” user app: 560X34-*=***
for download from the app stores (Apple App Store, Google Play, etc.). To do so, either enter the search terms “Ottobock”, “connectgo” or scan the QR code.
Further information about the app and how it works can be found either in the link in the description of the app stores or in the installed app.

6 Charging the battery

The following points must be observed when charging the battery:

- Use the 757L16-4 power supply / 757L43 charging adapter and 4E50* battery charger to charge the battery.
- With uninterrupted walking, the capacity of the fully charged battery is sufficient for **at least 18 hours**. It lasts **about 2 days** with average use.
- We recommend charging the product every day when used on a daily basis.
- The rechargeable battery should be charged for at least **3 hours** prior to initial use.
- For the maximum operating time with one battery charge, turning the product off when it is not being used is recommended.
- Switching the MyModes and changing setting parameters using the user app is not possible during the charging process.
- After disconnecting the battery charger, the product is in the same state it was in before connecting the battery charger. For example, if the product was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place (see fig. 1).
- 2) Connect the round, **four-pin** plug of the charging cable to the **OUT** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 3) Connect the round, **three-pin** plug of the power supply to the **12 V** receptacle on the battery charger so that the plug locks into place (see fig. 2).
INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.
- 4) Plug the power supply unit into the wall socket.

- The green LED on the back of the power supply (see fig. 3) and the green LED on the battery charger light up.
- If the green LED on the power supply and the green LED on the battery charger do not light up, there is an error (see page 74).

6.2 Connect battery charger to the product

- 1) Open the charging receptacle cover.
- 2) Connect the charging plug to the charging receptacle of the product.

INFORMATION: Make sure to insert it in the right direction!

- All displays light up alternately for a test.
 - After the displays are tested, a short beep is emitted followed by three short vibration signals to indicate that the charging process is started.
 - If the **!** symbol lights up, this means an error was identified during the self-test (see page 70).
 - To check the charge level while the battery charger is connected, briefly press the **!** button on the control panel.
- 3) Disconnect the product after the charging process is complete.
 - After disconnecting the battery charger, a vibration signal is emitted followed by a short beep, and the current status is displayed for approx. 5 seconds (see page 70).
 - When the symbol in the **!** button is lit up green **!**, the product is switched on and ready for operation.
 - 4) Close the charging receptacle cover.

INFORMATION

The product retains its previous state after disconnecting the battery charger. For example, if the product was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger. A descending sequence of beeps is emitted when the battery charger is disconnected while the product is switched off .

INFORMATION

No display after connecting the battery charger

If no symbols light up on the control panel after connecting the battery charger, the battery may be deep discharged. Leave the battery charger connected for at least 15 minutes and check the charge level while charging by disconnecting/connecting the battery charger.

INFORMATION

Depending on the battery charge level, the battery charger can warm up considerably during the charging process. This is not a malfunction.

6.3 Display of the current charge level

6.3.1 Display of battery charge level without additional devices

The current charge level of the installed battery can be displayed by briefly pressing the **!** button:

Sym- bol	Charge level
	Charge level 67 % to 100 %
	Charge level 34 % to 67 %
	Charge level 10 % to 34 %
	Charge level 5 % to 10 %
	Charge level 0 % to 5 %

6.3.2 Display of the current charge level during the charging process

After connecting the battery charger or pressing the  button while the battery charger is connected, the current charge level is shown by an animated symbol on the control panel in addition to the display on the battery charger (, , ).

Control panel	Battery charger	
	   	Battery is being charged, battery charge level is less than 34%
		Battery is being charged, battery charge level is 34% to 50%
	   	Battery is being charged, battery charge level is 50% to 67%
	The on time of the yellow LED gets shorter as the charge level increases. It only flashes briefly at the end of the charging process.	Battery is being charged, battery charge level is 67% to 99%
		Battery is fully charged

7 Use

INFORMATION

Before each use, check the product for functional reliability and for possible wear or damage.

Handling and use of the product must be learned before using it for the first time.

Applying and removing, sitting down and standing up as well as walking must be practised.

In general, major discomfort should not be experienced when wearing the product. The affected limb should be checked periodically for areas of excessive pressure. If pressure points are observed, stop using the product and consult the O&P professional as soon as possible.

7.1 Application

INFORMATION

The individual approach to put on and take off the product should be developed in cooperation with the orthopaedic technician and/or therapist.

- 1) Open all orthosis shell closures.
- 2) Remove the shoe.

- 3) Sit on the front edge of a chair.
- 4) Flex the orthotic joint.
- 5) Insert the foot into the foot component. While doing so, position the heel and lower leg in the shell.
- 6) Slightly extend the leg and apply the orthosis to the thigh and lower leg.
- 7) Fasten the ankle closure if applicable.
- 8) Fasten the lower leg shell closure.
- 9) Fasten the thigh shell closure.
- 10) Retighten the top closure.
- 11) Put on the shoe.
- 12) Stand up from the chair and retighten all closures.
- 13) Check the correct fit of the orthosis.

⚠ CAUTION

Pinching skin in the area of the closures

Injuries and swelling due to disruption of blood circulation to the skin in the area of the closures.

- Do not over-tighten the closures when applying the product.

7.2 Removal

INFORMATION

Before sitting down, check whether adequate resistance is provided for support while sitting down. If a MyMode with high or very low stance phase flexion resistance is activated, the user app must be used to switch back to basic mode before sitting down. Basic mode can also be activated by switching the product off and back on again.

- 1) Sit on a chair.
- 2) Open the orthosis shell closures.
- 3) Take off the orthosis.
- 4) Fasten the orthosis shell closures.
- 5) Store the orthosis and charge the battery if possible.

7.3 Movement patterns in basic mode (mode 1)

INFORMATION

Movement noises of the joint unit and other joints (medial support, ankle joint)

When using orthotic joints, servo motor, hydraulic, pneumatic or brake load dependent control functions can cause movement noises. This kind of noise is normal and unavoidable. It generally does not indicate any problems. If movement noise increases noticeably during the lifecycle of the orthosis, the orthosis should be inspected by the orthotics and prosthetics professional (O&P professional) immediately.

7.3.1 Standing

A stable stance is assured thanks to the static alignment of the orthosis.

The knee joint of the product has no locking function, so that slow bending is possible under flexion load. To restore the stable stance position, place the leg back under the body and put weight on the heel.

The stance function can be enabled by the O&P professional. Further information on the stance function can be found in the following section.

7.3.1.1 Stance function

INFORMATION

To use this function, it needs to be enabled by the O&P professional. It also has to be activated using the user app.

The stance function is a functional supplement to the basic mode. This function makes it easier for the user to stand on an inclined surface for a longer time. The joint is fixed in the flexion direction at a flexion angle between 5° and 65°.

The stance function must be enabled by the O&P professional. In addition, the type of joint locking (deliberate/intuitive) has to be established by the O&P professional. The locking type cannot be changed via the user app.

Intuitive locking of the joint

The intuitive stance function recognises any situation that puts strain on the orthosis in the flexion direction but where flexion is not permitted. Examples of this include standing on uneven or sloping surfaces. The joint is always locked in the flexion direction when the leg with the orthosis is not fully extended and is kept still for a brief moment.

The joint is not locked when the above conditions are met and a sitting position is assumed.

Intuitive unlocking of the joint

► Upon forward or backward rollover or extension, the flexion resistance is immediately switched from high to stance phase resistance again.

Manual joint lock

► Flex the joint between 5° and 60°.

→ Now the blocked joint can be loaded in the flexion direction.

Releasing the manual joint lock

► The Manual Stance function is automatically deactivated again by extending the knee or by repositioning the leg (e.g. taking a step).

7.3.2 Walking

Initial attempts at walking with the product always require instruction from trained, qualified personnel.

The hydraulics stabilise the knee joint in the stance phase and release the knee joint in the swing phase so that the leg can swing forward freely.

Switching to the swing phase requires a rollover of the entire foot.

An acoustic feedback signal when initiating the swing phase can be activated by the O&P professional.

INFORMATION

Changing the gait pattern influences activation of the swing phase. Adjustment by the O&P professional has to be repeated in this case.

7.3.3 Sitting down

INFORMATION

Before sitting down, check whether adequate resistance is provided for support while sitting down. If a MyMode with high or very low stance phase flexion resistance is activated, the user app must be used to switch back to basic mode before sitting down. Basic mode can also be activated by switching the product off and back on again.

The resistance in the product's knee joint while sitting down ensures even lowering into the sitting position.

- 1) Place both feet side by side at the same level.
- 2) While sitting down, weight should be distributed evenly between both legs and the arm supports used if available.
- 3) Move the buttocks in the direction of the back support and lean the upper body forward.

7.3.4 Sitting

INFORMATION

While sitting, the knee joint also switches to energy saving mode. This energy saving mode is activated regardless of whether the sitting function is activated or not.

When a sitting position is maintained for more than two seconds (i.e. the thigh is close to horizontal and there is no load on the leg), the product switches the resistance in the extension direction to a minimum.

If the sitting function has been enabled by the O&P professional and switched on via the user app, the resistance in the flexion direction is also reduced.

7.3.5 Standing up

INFORMATION

If a MyMode with a very high extension resistance (component remains flexed while standing up) or very low flexion resistance (no support as expected) is activated, it is necessary to switch back to basic mode. Basic mode can also be activated by switching the product off and back on again.

- 1) Place the feet at the same level.
- 2) Lean the upper body forward.
- 3) Place the hands on arm supports, if available.
- 4) Stand up with support from the hands while keeping weight evenly distributed on the feet.

7.3.6 Walking down stairs

INFORMATION

Before walking down stairs, check whether adequate resistance is provided for walking down stairs. If a MyMode with high or very low stance phase flexion resistance is activated, the user app must be used to switch back to basic mode before walking down stairs. Basic mode can also be activated by switching the product off and back on again.

INFORMATION

Activating additional support when going down stairs and ramps


The O&P professional can configure additional support for walking down stairs and ramps using the adjustment app. With this setting, the flexion resistance increases depending on the knee angle. This support can be turned off/on in the user app.

This function must be consciously practised and executed. Only when the sole is properly positioned can the product react correctly and permit controlled flexion.

- 1) Hold the handrail with one hand.
- 2) Position the leg with the product on the step so that the foot projects halfway over the edge of the step.
→ This is the only way to ensure a secure rollover.
- 3) Roll the foot over the edge of the step.
→ This flexes the product slowly and evenly at the knee joint.
- 4) Place the foot of the less affected leg onto the next step.

7.3.7 Walking up stairs

Walking up stairs one step at a time

- 1) Hold the handrail with one hand.
- 2) Place the foot of the less affected leg onto the first step.
- 3) Pull up the leg with the product.

Walking up stairs step-over-step

The product does not have an active drive system to support walking up stairs step-over-step. This is possible given certain physical prerequisites (residual muscular function present in the affected leg) and with corresponding training.

7.3.8 Walking up a ramp

- 1) Hold the handrail with one hand.
- 2) Place the foot of the less affected leg onto the ramp.
- 3) Move the leg with the product forward.

7.3.9 Walking down a ramp

INFORMATION

Activating additional support when going down stairs and ramps

The O&P professional can configure additional support for walking down stairs and ramps using the adjustment app. With this setting, the flexion resistance increases depending on the knee angle. This support can be turned off/on in the user app.

Use the handrail if possible.

Slight slope (< 5–10 %)

The movement pattern is equivalent to walking on the level. A swing phase may be initiated.

Moderate/steep slope (> 5–10 %)

The movement pattern is similar to walking down stairs. Place a load on the orthosis, permit knee flexion against the flexion resistance and roll over the forefoot.

7.3.10 Walking down flat steps

A (possibly surprising) swing phase may be triggered when walking down a flat step (such as a curb). The user can either use the stance phase flexion resistance or initiate a swing phase.

7.3.11 Kneeling

The individual approach to kneeling down and getting up again should be developed in cooperation with the O&P professional and/or therapist. Kneeling down is supported by increased flexion resistance, which permits controlled flexion of the knee joint.

7.3.12 Intuitive cycling (region-specific)

CAUTION

Balance problems when using the “Intuitive cycling” function

- > Falling while getting on and/or attempting to ride away.
- > Falling while getting off or finding a secure standing position.
- ▶ The “**Intuitive cycling**” function may be used only by experienced cyclists.

The “**Intuitive cycling**” MyMode allows the use of a bicycle without exiting basic mode.

Note the prerequisites for switching.

Prerequisites for activating the “Intuitive cycling” MyMode

- The bicycle must have the foot position (pedal axle) behind the knee during pedalling. Switching for recumbent bicycles is not possible.
- The bicycle must be equipped with a freewheel.
- The user must be in the sitting position.
- The sitting position must not be too high, otherwise the knee is extended during the pedaling movement, ending the MyMode.
- The sitting position must not be too low. Note the permissible flexion range of the knee joint.
- The feet have to be positioned on the pedals.
- Pedaling movements must be possible.

Activating MyMode “Intuitive cycling”

- 1) Take a seat on the bicycle.
- 2) Put the feet on the pedals.
- 3) Perform pedalling movements or activate the “**2.Intuitive cycling**” MyMode using the user app.
 - After a few pedalling movements, these are recognised by the knee joint and a short beep signal is produced. If this signal is not produced, the prerequisites for activating this MyMode were not met.
 - The short beep signal is produced periodically at intervals during the pedalling movement until the resistances in the flexion and extension direction have been reduced to the extent that the knee joint moves freely.
 - This MyMode (**2.Intuitive cycling**) is shown in the overview in the user app as the active MyMode.

Deactivating MyMode “Intuitive cycling”

- ▶ From the sitting position, either extend the knee or take the foot off the pedal and put it on the floor. The foot has to be ahead of the knee joint when it is set on the floor.
 - This is recognised by the orthosis and a long beep signal is produced. If this signal is not produced, either repeat the process or switch to the “**1.Basic Mode**” MyMode using the user app. It is also possible to switch to basic mode by switching the product off/on with a long press of the  button on the control unit.

- The basic mode is once again shown as the active MyMode in the overview in the user app.

INFORMATION

You must wear a bicycle helmet to ensure safety while cycling.

In addition, the bicycle must have a freewheel function and it is not permissible to secure the shoes on the pedals (with clips, click mountings or the like).

7.4 Switching the product on/off

In certain cases, e.g. for storage or transportation, the product can be purposely switched off.

⚠ CAUTION

Using the product while switched off

Falling due to unexpected product behaviour because of changed damping behaviour.

- ▶ Prior to using the product, briefly press the  button on the control panel to check whether it is switched on. The button symbol lights up green  when the product is switched on.

Switching off

- 1) Briefly press the  button on the control panel to check whether the product is switched on. The symbol in the button lights up green  and the current charge level is displayed (see page 70).
- 2) With the product switched on, press and hold the  button on the control panel until the symbol in the button turns off and a descending series of beeps (shut-down melody)  is emitted.

Switching on

- ▶ Press the  button on the control panel to switch the product on.
- A long vibration signal followed by a short beep is emitted, and the current status is displayed for about 5 seconds (see page 70).
- When the symbol in the  button is lit up green , the product is switched on and operational.
- Basic mode is activated after switching on.

INFORMATION

No display after switching on

If the  symbol does not light up after a long press of the  button, the battery may be deep discharged. Charging for at least 15 minutes is required in this case.

7.5 Bluetooth

A wireless connection to the component can be established with a mobile device and the user app/adjustment app via the Bluetooth function (see "Accessories" section (see page 54)). With these apps, settings can be configured on the component, current component data can be read out (battery charge level, steps walked, etc.) and error messages can be displayed.

Bluetooth must be switched on at the component in order to establish a connection.

The establishment of a connection is secured with a Bluetooth PIN which is provided on the enclosed Bluetooth PIN card. This Bluetooth PIN should not be shared.

7.5.1 Pairing with a device

- 1) Start the relevant app on the mobile device.
- 2) With the component switched on, briefly press the  button on the control panel to activate recognition (visibility) of the Bluetooth connection for 2 minutes.

- 3) Add the component in the app.

Information: To add the component in the app, follow the app instructions for use and the app instructions on the screen.

7.5.2 Switching on Bluetooth

- ▶ With the Bluetooth function deactivated, press and hold the  button on the control panel until a short vibration signal is emitted and the symbol in the button lights up blue .
- Bluetooth is switched on.
- To check whether the Bluetooth function was activated correctly, query the status by pressing the  button (see page 70).

7.5.3 Switching off Bluetooth



- ▶ With the Bluetooth function activated, press and hold the  button on the control panel until a vibration signal is emitted and the symbol in the button turns off.
- Bluetooth is switched off.
- To check whether the Bluetooth function was deactivated correctly, query the status by pressing the  button (see page 70).

7.6 Recommendations for air travel

Observing the following points is recommended prior to air travel or on board the airplane:

- Take along the 647F558 orthotic passport so you can show it on request or in case of questions.
- Turn off the component's Bluetooth function on board the aircraft if required (Switching Bluetooth of the component on/off).
- Take along the corresponding power supply adapter according to your destination. The power supply is suitable for connecting to alternating voltages from 100 V to 240 V at a network frequency of 50 Hz to 60 Hz.

8 MyModes

MyModes are intended for specific types of movement or posture (e.g. golf, basketball, ...). The O&P professional can enable and configure these in addition to basic mode (mode 1) using the adjustment app. The MyModes can be switched using the user app.

Settings can also be adjusted via the user app.

The parameters of the **"Training mode"** and **"Freeze position"** MyModes are preconfigured and cannot be changed.

9 Additional operating states (modes)

9.1 Empty battery mode

Beeps and vibration signals are emitted if the available battery charge level drops below 5 % (see page 70). Prior to switch-off, the flexion resistance is set to the safety mode values.

INFORMATION

The product retains its previous state after disconnecting the battery charger. For example, if the product was switched off before connecting the battery charger, it will also be switched off after disconnecting the battery charger. A descending sequence of beeps is emitted when the battery charger is disconnected while the product is switched off .

9.2 Mode for charging the product

The product is non-functional during charging.

To switch to basic mode, the battery charger must be disconnected from the product when the battery is charged.

9.3 Safety mode

INFORMATION

Change in flexion/extension resistance in safety mode

In safety mode, walking is possible with restrictions. Please note the possible change in flexion/extension resistance.

The product automatically switches to safety mode if a critical fault occurs (e.g. failure of a sensor signal). Safety mode remains in effect until the error has been rectified.

Flexion resistance preconfigured by the O&P professional is activated in safety mode, corresponding at least to the stance phase resistance. The extension resistance is low and cannot be changed. This makes limited walking possible for the user and allows the user to sit down, even though the sensor system is not active.

The switch to safety mode is indicated by the red  symbol on the control panel and by beeps and vibration signals immediately prior to switching (see page 70).

The current safety mode can be reset by switching the product off/on (see page 63) and by connecting then disconnecting the battery charger. If the product switches into safety mode again or at a later time, this means a permanent error exists. The product must be inspected immediately by an O&P professional.

9.4 Overheating mode

When the hydraulic unit overheats due to uninterrupted, increased activity (e.g. extended walking downhill), the flexion resistance is increased along with the rising temperature in order to counteract the overheating. When the hydraulic unit cools down, the product switches back to the settings that existed before the overheating mode.

Overheating mode is not activated in the MyModes.

Overheating mode is indicated by a long vibration every 5 seconds and additionally by the orange  symbol lighting up.

The following functions are deactivated in overheating mode:

- Sitting function
- Switching to a MyMode
- Changing the product settings

10 Cleaning

- 1) Clean the product with a damp cloth and mild soap when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

Orthosis shell padding

- 1) Remove the padding materials from the orthosis shells.
- 2) Fasten all hook-and-loop closures, if any.

- 3) If the padding materials "623P3 terry cloth padding fabric" or "623F62 SpaceTex padding fabric" are used, hand wash these in warm water at 30 °C/86 °F with a standard mild detergent.
If other padding materials are used, observe the care symbols and cleaning instructions for those materials.
- 4) Rinse out detergent residues thoroughly.
- 5) Allow to air dry. Do not expose to sources of direct heat (e.g. sunlight, stove or radiator).
- 6) After the padding material is completely dry, reattach it to the orthosis shells on the correct sides.

11 Maintenance

Regular maintenance (service inspections) is mandatory in the interest of your own safety and in order to maintain operating reliability and protect the warranty, maintain basic safety and the essential performance characteristics, and ensure safety in regards to EMC.

When maintenance is due, this is indicated by feedback after connecting the battery charger (see the section "Operating states/error signals", see page 70).

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

Components required for maintenance or repair:

The orthosis, battery charger and power supply unit.

17KO1-2=* C-Brace joint unit

Regular maintenance (service inspections) must be carried out at intervals of 24 months or 2 million steps, whichever comes first.

Additional orthosis components used, e.g. joints

Observe the service intervals and maintenance instructions for all orthosis components used.

12 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

12.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

12.2 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off

and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s)..

Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population.

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

12.3 CE conformity

Only applicable for “17KO1-2=L C-Brace joint unit”/“17KO1-2=R C-Brace joint unit”

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The full text of the EU Declaration of Conformity is available at the following Internet address: <http://www.ottobock.com/conformity>

12.4 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

The Bluetooth® word mark and logo are registered trademarks of Bluetooth SIG, Inc. and any use of such trademarks by Otto Bock Healthcare Products GmbH is under license.

13 Technical data

Ambient conditions	
Transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage in the original packaging (≤3 months)	-20 °C/-4 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Long-term storage in the original packaging (>3 months)	-20 °C/-4 °F to +25 °C/+77 °F Max. 93% relative humidity, non-condensing
Transport and storage between applications (without packaging)	-25 °C/-13 °F to +35 °C/95 °F Max. 93% relative humidity, non-condensing +35 °C/95 °F to +70 °C/158 °F Water vapour pressure up to 50 hPa
Operation	-10 °C/+14 °F to +40 °C/+104 °F Relative humidity 15% to 93%, non-condensing, at a water vapour pressure up to 50 hPa Air pressure: 606.3 hPa (up to 4000 m without pressure equalisation)
Maximum temperature that can be reached on the surfaces of the orthosis during operation	+44 °C/+111 °F
Time for warming to the operating temperature after storage between applications, from -25 °C/-13 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Time for cooling to the operating temperature after storage between applications, from +70 °C/+158 °F at an ambient temperature of +20 °C/+68 °F	30 minutes
Charging the battery	+10 °C/+50 °F to +40 °C/+104 °F

General	
Reference number	17KO1-2=L/C-Brace joint unit left/17KO1-2=R C-Brace joint unit right
Joint unit weight [g/oz]	Approx. 1000/35
Maximum user body weight [kg/lbs]	125/276
Product lifetime [years]	6
Information on the product's ruleset and firmware version	Accessible using the user app

Data communication	
Wireless technology	Bluetooth 5.0 (Bluetooth® Low Energy)
Distance range	Approx. 10 m/32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Data rate (over the air)	Up to 2 Mbps
Maximum output power (EIRP):	+4 dBm (~2.5 mW)

Battery charger	
Reference number	4E50*

Battery charger	
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F Max. 93% relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F Max. 93% relative humidity, non-condensing
Input voltage	12 V 
Service life	8 years

Power supply unit	
Reference number	757L16-4
Type	FW8001M/12
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10% to 95% relative humidity, non-condensing
Operation	0 °C/+32 °F to +50 °C/+122 °F Max. 95% relative humidity Air pressure: 70–106 kPa (up to 3,000 m without pressure equalisation)
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V 

Rechargeable battery of the product	
Battery type	Li-ion
Battery specifications	3070 mAh/11.05 Wh
Charging cycles (charging and discharging cycles) after which at least 80 % of the original battery capacity remains available	500
Charge level after 1 hour charging time	30 %
Charge level after 2 hours charging time	50 %
Charge level after 4 hours charging time	80 %
Charge level after 8 hours charging time	Fully charged
Product behaviour during the charging process	The product is non-functional
Operating time of the product with new, fully charged battery at room temperature	At least 18 hours of uninterrupted walking Approx. 2 days with average use

14 Appendices

14.1 Symbols Used

	In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with regulations in your country can be harmful to health and the environment. Please observe the instructions of your national authority pertaining to return and collection procedures.
	Manufacturer
	Type BF applied part
	Complies with the requirements of the Radiocommunications Act (AUS)
IP54	Protected against dust, protected against splashed water
	Declaration of conformity according to the applicable European directives
	Serial number (21)YYYYWWNNNN YYYY – year of manufacture WW – week of manufacture NNN – sequential number

	Medical device
	UDI number (Unique Device Identifier)
	Article number
	Data matrix code
	Global Trade Item Number
	Caution, hot surface
	Please note the instructions for use
	Limits for temperature
	Limits for atmospheric pressure
	Limits for relative humidity

14.2 Operating states/error signals

The product indicates operating states and error messages with symbols on the control panel and through beeps and vibration signals.

14.2.1 Status display on the control panel

The current status of the component is displayed on the control panel for 5 seconds in the following cases:

- The  button on the control panel was pressed briefly.
- The component was switched on by pressing the  button.
- The battery charger was disconnected from the component.
- Battery charger was connected to the component.
- An error was identified during use.

INFORMATION

No status display due to battery deep discharge

If the status of the component is not displayed on the control panel, the battery may be deep discharged. In this case it must be charged for at least 15 minutes before the status can be queried again.

Symbol on the control panel	Beep signal	Vibration signal	Event	Required action
All symbols on the control panel light up alternately	—	—	Test of the indicators (LEDs) after connecting the battery charger	Check whether all symbols (LEDs) light up alternately and in the corresponding colours. If a symbol (LED) does not light up in a colour, the product should be inspected by an O&P professional. If no symbols (LEDs) light up, the battery may be deep discharged. Leave the battery charger connected for at least 15 minutes and then repeat this test by disconnecting/connecting the battery charger.
	—	—	The product is switched on and ready for operation	
	1 x short	1 x long and 1 x short	The product was switched on by pressing the  button on the control panel.	
	—	1 x long at interval of approx. 5 seconds	Hydraulics overheated	Reduce activity.
	—	—	A self-test error was detected when connecting the battery charger.	- Conduct the self-test again by connecting/disconnecting the battery charger. - If the  symbol lights up again, the product must be checked by an O&P professional.

Symbol on the control panel	Beep signal	Vibration signal	Event	Required action
	30 x long	30 x long	Severe error/indication of safety mode activation (see page 65) Switch to safety mode if possible	Walking possible with restrictions. Please note the possible change in flexion/extension resistance. 1. Attempt to reset this error by switching the product off/on (see page 63). 2. If the beep/vibration signal recurs, attempt to reset this error by connecting/disconnecting the battery charger. 3. If the beep/vibration signal recurs, use of the product is prohibited. The product must be checked by an O&P professional.
	—	—	Charge level 10 % to 34 %	
	—	—	Charge level 34 % to 67 %	
	—	—	Charge level 67 % to 100 % Indicates that the battery is fully charged during the charging process.	
	—	—	Battery is being charged, charge level is less than 34 %	
	—	—	Battery is being charged, charge level is 34 % to 67 %	
	—	—	Battery is being charged, charge level is 67 % to 99 %	
	3 x long	3 x long	Charge level between 5 % and 10 %	Charge battery soon. Remaining operating time approx. 1 hour.
	5 x long	5 x long	Charge level between 0 % and 5 %	Charge battery immediately; the product will be switched off after the next warning signal.
	10 x long	10 x long	Charge level 0 % After the beep and vibration signals, the product switches to empty battery mode and then switches off.	Charge the battery.

Symbol on the control panel	Beep signal	Vibration signal	Event	Required action
	4 x short at intervals of approx. 65 sec. (continuously)	—	Charging the battery outside the allowable temperature range	Check whether the specified environmental conditions for charging the battery were met (see page 68).
	—	—	Maintenance required because the maintenance date has already been reached or exceeded.	You should visit the O&P professional soon. They will forward the component to an authorised Ottobock Service Center.
	—	—	Bluetooth activated	
—	1 x short	—	Battery charger is connected	
—	1 x long	—	When getting off the bicycle and reactivating basic mode	After getting off the bicycle, check whether the extension and flexion resistances for walking and standing are set.
—	1 x short at periodic intervals	—	Detection and activation of intuitive cycling mode	The beep signal is emitted until the resistances in the extension and flexion directions are reduced to a minimum.
—	—	3 x short	Charging mode started (3 sec. after connecting battery charger)	
—	1 x short	1 x short	Mode switching or changing adjustment parameters is performed using the user app. If the parameter " Volume " is set to "0" in the user app, only a vibration signal is emitted.	

Symbol on the control panel	Beep signal	Vibration signal	Event	Required action
—		—	The product switches itself off. This occurs automatically in the following cases: - The ⓘ button on the control panel was pressed and held for longer than approx. 3 seconds while the product was switched on. - After the  symbol lights up. - After disconnecting the battery charger, if the product was already switched off before connecting the battery charger.	- Charge the battery. - If desired, switch on the product with the ⓘ button.
—	—	Ongoing	Total failure Electronic control no longer possible. Safety mode active or undetermined valve state. Unknown product behaviour.	- Press and hold the ⓘ button on the control panel until the vibration signal stops (approx. 10 seconds), completely switching off the product. - If the vibration signal recurs after switching on, attempt to reset this error by connecting/disconnecting the battery charger. - If the vibration signal recurs, use of the product is prohibited. The product must be checked by an O&P professional.

14.2.2 Error while charging the product

LED on power supply	LED on battery charger	Error	Resolution
	   ⓘ	Country-specific plug adapter not fully engaged on power supply	Check whether the country-specific plug adapter is fully engaged on the power supply.
		Non-functional socket	Check socket with another electric device.
		Defective power supply	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.

LED on power supply	LED on battery charger	Error	Resolution
	   	No connection between battery charger and power supply	Check whether the charging cable plug is fully engaged on the battery charger.
		Defective battery charger	The battery charger and power supply must be inspected by an authorised Ottobock Service Centre.
	   	Battery is fully charged (or connection with product is interrupted).	Take note of the confirmation signal for differentiation. When the battery charger is connected, a self-test is conducted and confirmed by beeps/vibration signals. When these signals are emitted, the battery is fully charged. If no signals are emitted, the connection to the product is interrupted.
			If the connection to the product is interrupted, an authorised Ottobock Service Centre must inspect the product, battery charger and power supply.

14.3 Directives and manufacturer's declaration

14.3.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notes in the section "Information on proximity to certain areas" (see page 45).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drops	IEC 61000-4-11	0% U_T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U_T ; 1 period and 70% U_T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U_T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, GSM 800/900, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.11 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

Immunity to magnetic fields in close range

Test frequency	Modulation	Interference immunity test level [A/m]
30 kHz	CW	8
134.2 kHz	Pulse modulation 2.1 kHz	65
13.56 MHz	Pulse modulation 50 kHz	7.5

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This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

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This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s).

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A series of horizontal lines for writing, spanning the width of the page.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Legal Manufacturer of Custom Made Device

Legal Manufacturer of C-Brace joint unit 17KO1-2=*
Otto Bock Healthcare Products GmbH
Brehmstraße 16 · 1110 Wien · Austria
T +43-1 523 37 86 · F +43-1 523 22 64
info.austria@ottobock.com · www.ottobock.com