# C-Leg 4

# 1 Bluetooth®

Integrated Bluetooth<sup>®</sup> technology permits straightforward communication with the joint. Bluetooth<sup>®</sup> can be deactivated if necessary.

### 2 Knee angle sensor

The knee angle sensor measures the flexion angle and the angular velocity of the joint.

# 3 Battery and electronics

A lithium-ion battery provides the energy required to control the knee joint. It is located directly in the rotation axis of the C-Leg 4. Charging is recommended overnight when used by the patient on a daily basis. The integrated microprocessor coordinates all measurement and control processes.

## **4** Inertial motion unit (IMU)

The gyroscope and the acceleration sensors allow the position and acceleration of the C-Leg 4 in space to be measured. Prosthetic Control is based on motion analysis and additional force determination.

## 5 Hydraulic unit

The hydraulic unit controls the C-Leg 4. It generates movement resistance for flexion and extension during the stance and swing phase.

#### 6 New carbon frame design

In order to withstand the varied demands of everyday life, the frame is made of carbon – an especially strong, high-grade and lightweight material. The frame houses and protects the electronics, the hydraulic unit and the battery. Thanks to its new design, the C-Leg 4 has a lower system height and therefore a larger fitting spectrum.

## Charging

The receptacle for the charger is located on the back of the joint and protected by a cover.



ottobock.

## Advantages

- Walking on difficult surfaces like sand, gravel and stones
- Choice between intuitive and manual stance function
- User can walk backwards confidently without initiating the swing phase
- Smart control via the Cockpit app for Android smartphones

# Technical data

	212	
MOBIS	2, 3, 4	
Max. body weight	136 kg	
Knee flexion angle	130° without flexion stop <sup>1</sup>	
Weight of the knee joint	1,235 g (pyramid adapter) 1,240 g (threaded connector)	
Frame material	Carbon	
Moisture protection	Weatherproof	
Colours available	Volcano shadow, Desert pearl	
Tube adapter	2R57*, 2R67	

<sup>1</sup> Flexion stop reduces the knee flexion angle by 8°

#### System solution



