


Getting Started: Single Cell ATAC

Introduction



The Chromium Single Cell ATAC Solution provides a comprehensive, scalable approach to determine the regulatory landscape of chromatin in hundreds to thousands of cells in a single sample. Nuclei suspensions are incubated in a Transposition Mix that includes a Transposase, which enters the nuclei and preferentially fragments the DNA in open regions of the chromatin. Using a microfluidic chip, nuclei are partitioned into nanoliter-scale Gel Beads-in-emulsion (GEMs). A pool of ~750,000 10x Barcodes is used to separately and uniquely index the transposed DNA of each individual nucleus. Libraries are generated and sequenced, and 10x Barcodes are used to associate individual reads back to the individual partitions, and thereby, to each individual nucleus. For more information about the assay, click [here](#).

Sample Preparation Guidelines

Preparing clean, viable single nuclei suspension is essential for obtaining high quality single cell ATAC data.

To prepare nuclei suspension from cells (fresh or frozen), consult this [Demonstrated Protocol](#). The protocol includes optional DNase treatment to remove ambient DNA from the samples.

To prepare nuclei suspension from tissue (fresh or frozen), consult this [Demonstrated Protocol](#). The protocol describes nuclei isolation from embryonic mouse brain.

Counting	Nuclei Concentration	Quality
After nuclei isolation, a cell viability stain can be used for staining and counting nuclei. If debris makes nuclei counting challenging, a fluorescent dye that stains nucleic acids may enhance counting accuracy.	<p>155-7,700 nuclei/μl Nuclei concentration should be adjusted to be within the optimal range for the targeted nuclei recovery.</p> <p>Use an automated cell counter or a hemocytometer with a live/dead stain to determine cell/nuclei count.</p> <p>Consult Guidelines for Accurate Target Cell Counts.</p> <p>Measuring nuclei count is critical for accurate target cell recovery, as detailed in this Technical Note.</p>	<p>During nuclei isolation, the outer cell membrane must be lysed without comprising the nuclear membrane to prevent nuclei leakage. A lysis timecourse is recommended to find the optimal lysis time for a sample.</p> <p>Good quality nuclei suspension should have intact nuclear membranes and should be free of clumps and debris.</p>

Challenging Sample Types

[Additional cleanup methods](#), like DNase treatment, filtering, FACS, and/or density centrifugation (listed below), may be used when isolating nuclei from complex and challenging sample types, such as tissues.




Sample Type	DNase Treatment	Filtering	FACS	Density Centrifugation
Cell	Yes	Yes	Yes	Yes
Nuclei	No	Yes	Yes, with caution	Yes
Limited Nuclei Yield	No	Yes, with a small volume filter	No	No
Abundant Nuclei Yield	No	Yes	Yes	Yes

Materials Required and Step Overview

Consult the [user guide](#) for a consumables and equipment list validated by 10x Genomics for executing the Chromium Single Cell ATAC protocol.

To enable seamless experimental planning, a breakdown of the protocol steps, along with execution times and stop & store points are listed below.

10x Genomics how-to [video series](#) provides a visual demonstration of the workflow.

	Steps	Timing	Stop & Store
2 h	Nuclei Isolation		
	Dependent on Cell Type		~1-2 h
	Step 1 – Transposition		
	1.1 Prepare Transposition Mix	10 min	
	1.2 Isothermal Incubation	30 min	
4 h	Step 2 – GEM Generation & Barcoding		
	2.1 Prepare Master Mix	10 min	
	2.2 Load Chromium Next GEM Chip H	10 min	
	2.3 Run the Chromium Controller	18 min	
	2.4 Transfer GEMs	3 min	
	2.5 GEM Incubation	45 min	 15°C ≤ 18 h or -20°C ≤ 1 week
	Step 3 – Post GEM Incubation Cleanup		
	3.1 Post GEM Incubation Cleanup – Dynabeads	35 min	
	3.2 Post GEM Incubation Cleanup – SPRIselect	15 min	 4°C ≤ 72 h or -20°C ≤ 2 weeks
6 h	Step 4 – Library Construction		
	4.1 Sample Index PCR	40 min	
	4.2 Post Sample Index Double Sided Size Selection – SPRIselect	20 min	 4°C ≤ 72 h or -20°C long-term
	4.3 Post Library Construction QC	60 min	

Sequencing

The Chromium Single Cell ATAC protocol produces Illumina®-ready sequencing libraries. The recommended library sequencing depth & run parameters are summarized in the user guide and are also available on the [10x Genomics Support website](#).

	Read 1	i7 Index	i5 Index	Read 2
Purpose	Transposed DNA	Sample Index	10x Barcode	Transposed DNA
Length (bp)	50	8	16	50*

*Can be 1 bp shorter to accommodate NovaSeq 100 cycle kits

Consult the [Sequencing Metrics & Base Composition Technical Note](#) for representative data generated using various sequencers validated by 10x Genomics.

© 2022 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

10xGenomics.com

Contact:
support@10xgenomics.com
 10x Genomics
 6230 Stoneridge Mall Road
 Pleasanton, CA 94588 USA

