

## DEMONSTRATED PROTOCOL

# Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols

with Feature Barcode technology

## Overview

Cell surface proteins can be labeled using a specific protein binding molecule, such as an antibody conjugated to a Feature Barcode oligonucleotide. This protocol provides guidance for antibody-oligonucleotide conjugation and outlines cell surface protein labeling for use with Single Cell RNA sequencing protocols with Feature Barcode technology for cell surface protein. This protocol also provides guidance for enriching labeled cells using Fluorescence Activated Cell Sorting (FACS).

Consult Demonstrated Protocol Cell Multiplexing Oligo Labeling for Single Cell RNA Sequencing Protocols with Feature Barcode technology (CG000391) for guidance on cell surface protein and Cell Multiplexing Oligo labeling,

## Additional Guidance

Consult Demonstrated Protocol Cell Preparation Guide (Document CG00053) for Tips & Best Practices on handling cells and Technical Note Guidelines on Accurate Target Cell Counts (Document CG000091) for determining accurate cell counts.

**Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.**

## Preparation – Buffers

### Buffers

Prepare fresh, maintain at 4°C

PBS + 1% BSA

Alternatively, Cell Staining Buffer from BioLegend can be used for labeling.

PBS + 0.04% BSA

PBS + 2% FBS (Only if performing FACS enrichment of labeled cells)

## Specific Reagents & Consumables

### For Antibody-Oligonucleotide Conjugation

Vendor	Item	Part Number
Abcam	Oligonucleotide Conjugation Kit	ab218260
IDT	Custom DNA Oligos (see Table 1)	-
-	100 µg Purified Azide-free Antibody (1 mg/ml)	-

### For Cell Surface Protein Labeling

Vendor	Item	Part Number
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TotalSeq™ Antibody-Oligonucleotide Conjugates*	-
	Cell Staining Buffer	420201
	Antibodies (Fluorophore)†	-
If using FACS for enriching labeled cells		
Thermo Fisher Scientific	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
	Trypan Blue Stain (0.4%)	T10282
Millipore Sigma	Bovine Serum Albumin In DPBS (10%) (alternative to Thermo Fisher product)	A1595
Miltenyi Biotec	MACS BSA Stock Solution (alternative to Thermo Fisher product)	130-091-376
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV
VWR	Fetal Bovine Serum (FBS)	97068-085
<b>Equipment</b>		
Thermo Fisher Scientific	Countess II FL Automated Cell Counter	AMAQAF1000
	Countess II FL Automated Cell Counting Chamber Slides	C10228

This list may not include some standard laboratory equipment.



\*TotalSeq™-B for: Single Cell 3' v3 and v3.1 protocol with Feature Barcode technology for Cell Surface Protein

\*TotalSeq™-C for: Single Cell 5' v1, v1.1 & v2 protocol with Feature Barcode technology for Cell Surface Protein & Antigen Specificity.

†Choose different clones than antibody-oligonucleotide conjugates

## Protocol Overview


### Option A. Custom Conjugated Antibodies

Follow manufacturer's instructions (Abcam Oligonucleotide Conjugation Kit) for conjugation and purification


### Option B. Preconjugated Antibodies

BioLegend TotalSeq™-B  
OR  
BioLegend TotalSeq™-C

## 1. Label Cells Prepare Antibody Mix & FACS Antibody Pool (if performing FACS enrichment) as described in the Labeling Protocol



## 2. Wash Cells





## Cell Surface Protein Labeling Protocol

This protocol was optimized using TotalSeq-B/C antibody-oligonucleotide conjugates from BioLegend. The labeled cells were enriched by FACS (see Appendix).



Use distinct antibody clones for FACS and cell surface protein labeling. Optimize working concentration of each of the antibodies used.

### 1. Label Cells

This protocol was demonstrated using  $0.2\text{--}1 \times 10^6$  cells. Wash cells according to the appropriate 10x Genomics Demonstrated Protocol for the cell type being prepared.

All steps can be performed in 5-ml centrifuge tubes.

#### Prepare Antibody Mix Supernatant:

- Add appropriate/manufacture's recommended amount of antibody-oligonucleotide conjugates to a 1.5-ml microcentrifuge tube.
- Centrifuge the mix at **14,000 rcf** for **10 min** at **4°C**.
- Transfer the supernatant (containing Antibody Mix) to a new tube and maintain at **4°C**.



Use TotalSeq-B for Single Cell 3' v3 and v3.1 (single & dual index) protocol and TotalSeq-C for Single Cell 5' v1, v1.1 & v2 protocol with Feature Barcode technology for Cell Surface Protein.

#### Prepare FACS Antibody Pool:

- Add appropriate/manufacture's recommended amount of fluorophore antibodies to a 1.5-ml microcentrifuge tube on ice.
- Gently pipette mix and maintain at **4°C**. Avoid light exposure.

- Transfer cells to a new 5-ml tube and add chilled PBS + 0.04 % BSA for a total 1 ml volume.
- Centrifuge cells at **4°C**. Use of swinging-bucket rotor is recommended for higher cell recovery. Centrifugation speed and time depends upon the sample type. Larger or fragile cell types may require slower centrifugation speeds. Use the following table for guidance.

**Table 2.** Sample Type Specific Centrifugation Conditions:

Sample Types	Speed (rcf)	Time (min)
Cell Lines	300	5
PBMCs	400	5
Dissociated Tumor Cells	150	10

- Remove the supernatant without disturbing the pellet.
- Resuspend cell pellet in **50 µl** chilled PBS + 1% BSA.
- Add **5 µl** Human TruStain FcX. Gently pipette mix.
- Incubate for **10 min** at **4°C**.
- Add the prepared Antibody Mix supernatant. If also performing FACS enrichment, add FACS antibody pool.
- Add chilled PBS + 1% BSA to the cells to bring the total volume to **100 µl**. Gently pipette mix 10x (pipette set to 90 µl).
- Incubate for **30 min** at **4°C**. If using FACS antibodies, incubate without light exposure.

Recommended incubation temperature for most sample types is **4°C**. However, incubation temperature is sample type dependent and should be chosen accordingly.

### 2. Wash Cells

Thorough washing of cells post labeling is critical to obtain high-quality data. Optimization of centrifugation speed/time may be needed based on cell type.

- Wash 1:** Wash by adding **3.5 ml** chilled PBS + 0.04% BSA to the cells from step 1i. Gently pipette mix.
- Centrifuge at **4°C**. Centrifugation speed and time depends upon the sample type. Use Table 2 for guidance.
- Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet.



Leaving behind excess supernatant may cause non-specific binding, which may result in increased background reads during sequencing.

- Resuspend the pellet in **100 µl** room temperature PBS and transfer to a new 5-ml tube. Incubate for **5 min** at **room temperature**. Proceed to step e (Wash 2).

Step d can be excluded for the fragile cells (see Cell Washing for Fragile Cell Types in Appendix) to maximize cell health and recovery. If not performing step d, proceed directly to step e (Wash 2) and transfer the cells to a new tube after adding **3.5 ml** PBS + 1% BSA in step e.




- Wash 2:** Using a pipette tip, resuspend the pellet or cells in **3.5 ml** chilled PBS + 1% BSA.
- Centrifuge at **4°C**. Centrifugation speed and time depends upon the sample type. Use Table 2 for guidance.
- Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet.
- Washes 3 & 4: Repeat e - g 2x** for a total of three washes.
- OPTIONAL For enrichment of labeled cells by FACS:** Based on starting concentration and assuming ~50% cell loss, add an appropriate volume chilled PBS + 2% FBS (including a dead cell marker) to obtain a final cell concentration of  $5\text{--}10 \times 10^6$  cells/ml and proceed to FACS (see FACS Guidance).  
After FACS, determine cell concentration and viability using a Countess II Automated Cell Counter or a hemocytometer and proceed **immediately** to relevant Chromium Single Cell RNA Sequencing protocols with Feature Barcode technology (see References).
- If not performing FACS:** Based on starting concentration and assuming ~50% cell loss, add an appropriate volume chilled PBS + 1% BSA to obtain a concentration of  $700\text{--}1,200$  cells/µl. Determine cell concentration and viability using a Countess II Automated Cell Counter or a hemocytometer and proceed **immediately** to relevant Chromium Single Cell RNA Sequencing protocols with Feature Barcode technology (see References).

## Appendix

### Illustrative Overview of Wash Steps

To eliminate non-specific binding with comparable efficiency, wash steps may be performed in 5-ml microcentrifuge tubes using indicated buffer volumes. Non-specific binding contributes to increased background reads during sequencing.

#### Wash Steps in 5-ml Microcentrifuge Tubes



\*Exclude these steps for fragile cells (See Cell Washing for Fragile Cell Types) and proceed directly to Wash 2 and transfer the cells to a new tube after adding 3.5 ml buffer.

### Cell Washing for Fragile Cell Types

PBS resuspension and room temperature incubation (step 2d) during cell wash can be excluded for fragile cell types. The exclusion of this step may yield slightly lower Cell Surface Protein library metrics, but will preserve cell health and recovery. Step 2d can be excluded in the following sample types:

- Cell types with low viability (<70%)
- Sample types that are known to degrade quickly
- Samples with very low cell input (<200,000 cells)

### FACS Guidance

Enrich labeled cells using FACS prior to library generation to enable identification of rare subpopulations.

#### FACS Cell Collection

It is recommended to collect FACS enriched cells in up to 20% FBS to maintain cell viability. Cells should be collected either in **20 μl** volume in the collection tube/plate (96-well plate) or in **200 μl** volume in a 1.5-ml tube.

The sort stream should be adjusted so that the cell-droplet falls into the collection buffer. Sorted cells must be counted and viability measured before proceeding to the 10x Genomics Single Cell protocols. If necessary, the collected cells may be concentrated by centrifugation at **350 rcf** at **4°C** and by removing the supernatant.

Cell loss during FACS is common. Optimize the protocol steps accordingly.



Once sorting is complete, proceed **immediately** to relevant Chromium Single Cell RNA Sequencing protocols with Feature Barcode technology (see References).


## Illustrative Overview of Antibody-Oligonucleotide Conjugate Capture



Antibody-oligonucleotide conjugate capture by protocol specific Gel Bead primers is illustrated below.

### Illustrative Overview of Antibody-Oligonucleotide Conjugate Capture

Single Cell 3'  
v3 & v3.1\* - Cell  
Surface Protein  
(CG000185,  
CG000206 &  
CG000317)




\*single & dual index libraries

TotalSeq-B Antibody-oligonucleotide conjugate



Single Cell 5' v1,  
v1.1 & v2 - Cell  
Surface Protein  
(CG000186,  
CG000208 &  
CG000330)



TotalSeq-C Antibody-oligonucleotide conjugate

## Compatible Primers

To generate Cell Surface Protein libraries with TotalSeq-B/C, use the primers indicated below.

Antibody-Oligonucleotide Conjugate	Compatible Primers	10x PN
TotalSeq-B	Feature cDNA Primers 2	2000097
TotalSeq-C	SC5' Feature cDNA Primer/Feature cDNA Primer 4	2000119/ 2000277



If generating libraries with TotalSeq™-A, an additive primer is required for successful amplification (not provided by 10x Genomics). See [Cite Seq Protocols](#) and [BioLegend](#) for details.

## References

This protocol provides guidance for antibody-oligonucleotide conjugation and outlines cell surface protein labeling for use with:

- Chromium Single Cell V(D)J Reagent Kits with Feature Barcode technology for Cell Surface Protein User Guide (CG000186)
- Chromium Next GEM Single Cell V(D)J Reagent Kits v1.1 with Feature Barcode technology for Cell Surface Protein User Guide (CG000208)
- Chromium Next GEM Single Cell 5' v2 (Dual Index) with Feature Barcode technology for Cell Surface Protein & Antigen Specificity User Guide (CG000330)
- Chromium Single Cell 3' Reagent Kits v3 with Feature Barcode technology for Cell Surface Protein User Guide (CG000185)
- Chromium Next GEM Single Cell 3' Reagent Kits v3.1 with Feature Barcode technology for Cell Surface Protein User Guide (CG000206)
- Chromium Next GEM Single Cell 3' v3.1 (Dual Index) with Feature Barcode technology for Cell Surface Protein User Guide (CG000317)

© 2021 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at <http://www.10xgenomics.com/legal-notices>, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**Contact:**  
[support@10xgenomics.com](mailto:support@10xgenomics.com)  
10x Genomics  
6230 Stoneridge Mall Road  
Pleasanton, CA 94588 USA

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

