

TECHNICAL NOTE

Assay Scheme & Configuration of Chromium Single Cell DNA Libraries

Introduction

The Chromium Single Cell DNA Reagent Kits protocol (Document CG000153) produces Single Cell DNA libraries ready for Illumina sequencing. Using the “Cell Bead-Gel Bead” (CBGB) technology, which involves a two-step microfluidics process to partition individual cells into Cell Beads, followed by partitioning of each Cell Bead with a 10x Gel Bead in a Gel Bead-in-Emulsion (GEM), single cell DNA barcoding is achieved. During library preparation, sequence components essential for Illumina sequencing and downstream data analysis are incorporated into the final library construct.

An overview of the Single Cell DNA workflow is presented in Figures 1-3. A single cell or nuclei suspension is combined with Cell Bead reagents in Chromium Chip C and run on the Chromium Controller or Chromium Single Cell Controller to generate Cell Beads. The Cell Beads are collected and processed to lyse the cells and denature the genomic DNA (gDNA), while retaining the gDNA within the Cell Bead. The gDNA inside the Cell Bead is now accessible for downstream amplification and barcoding (Figure 1).

Gel Beads that each deliver oligonucleotides consisting of a 10x Barcode sequence (Figure 2) are then combined with the Cell Beads and a reaction mix containing enzyme in Chromium Chip D (Figure 3). The chip is run on the Chromium Controller or Chromium Single Cell Controller to generate Cell Bead-Gel Bead (CBGB) GEMs, which are droplets containing one Cell Bead and one Gel Bead.




Figure 1. Single cell suspensions are partitioned into Cell Beads. Cell Beads are collected and processed to lyse the cells and denature the gDNA, while retaining the gDNA within the Cell Bead.




Figure 2. Schematic of Gel Bead oligo sequence. Each Gel Bead comprises of many oligos containing (i) a partial Illumina Read 1 (R1) sequence (22 bp) and (ii) a 16 bp 10x Barcode. The 10x Barcode sequence is unique to each Gel Bead.




Figure 3. Cell Beads are partitioned into GEMs containing the components necessary for barcoding and amplifying DNA from single cells.

Once partitioned, the Gel Bead dissolves, releasing Gel Bead oligonucleotides into the aqueous environment of the GEM. The GEMs are isothermally incubated to amplify gDNA using random hexamers supplied in the Reaction Mix (Figure 4). During the isothermal incubation, the Gel Bead oligonucleotides are ligated to the amplified DNA to produce barcoded fragments ranging from a several hundred basepairs to a few kilobases.

The GEMs are then “broken”, pooling the barcoded DNA molecules from the single cell partitions (Figure 4). During library construction, Read 2 is added during Adaptor Ligation. Illumina P5 and P7 sequences and i7 sample index are added during Sample Index PCR. The final library fragments contain a 10x Barcode and the genomic insert, flanked by sequences used in Illumina amplification (P5, P7) and indexing (i7). Detailed information regarding the sequences added during library preparation is provided in Figure 5.


Conclusion

This Technical Note presents an overview of the Single Cell DNA assay scheme and configuration of the sequencing library after completing the protocol outlined in the Chromium Single Cell DNA Reagent Kits User Guide (Document CG000153). The workflow and sequences presented here provide additional insight into Chromium Single Cell DNA libraries and may serve as a reference to customize library preparation.

References

- Chromium Single Cell DNA Reagent Kits User Guide (Document CG000153)

Inside individual GEMs



Pooled amplified DNA processed in bulk




Figure 4. Assay scheme for Single Cell DNA Library Preparation.





Protocol Step 3.5 – GEM Incubation	
Gel Bead Oligo (PN-2000033)	 <p>5' -CTACACGACGCTCTCCGATCT-NNNNNNNNNNNNNNNN-3'</p>
Hexamer	<p>Hexamer</p> <p>5' -NNNNNN-3'</p>
GEM Incubation Product	 <p>5' -CTACACGACGCTCTCCGATCT-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-3'</p> <p>3' -GATGTGCTGCGAGAAGGCTAGA-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-5'</p>
Protocol Step 5.2 – Adaptor Ligation	
Adapter (Read 2)	<p>Read 2</p> <p>5' -GATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'</p> <p>3' -TCTAGCCTTCTCG-5'</p>
Ligation Product	 <p>5' -CTACACGACGCTCTCCGATCT-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-3'</p> <p>3' -GATGTGCTGCGAGAAGGCTAGA-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-TCTAGCCTTCTCGTGTGCAGACTTGAGGTCAGTG-5'</p>
Protocol Step 5.4 – Sample Index PCR	
Sample Index PCR Primer	<p>Forward Primer: SI-PCR Primer P5-Partial Read 1 (PN-220124)</p> <p>Reverse Primer: Chromium i7 Sample Index i7 – Sample Index – Partial Read 2 (PN-220103)</p> <p>5' -AATGATACGGCGACCACCGA-GATCTACACTCTTCCCTACACGACGCTC-3'</p> <p>5' -CAAGCAGAAGACGGCATACGAGAT-NNNNNNNN-GTGACTGGAGTTCAGACGTGT-3'</p>
Sample Index PCR Product	 <p>5' -AATGATACGGCGACCACCGA-GATCTACACTCTTCCCTACACGACGCTCTCCGATCT-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-AGATCGGAAGAGCACACGTCTGAACTCCAGTCAC-NNNNNNNN-ATCTCGTATGCCGTCTTCTGCTTG-3'</p> <p>3' -TTACTATGCCGCTGTGTGCT-CTAGATGTGAGAAGGAGGTGTGCTGCGAGAAGGCTAGA-NNNNNNNNNNNNNNNN-NNNNNN-Genome_Insert-NNNNNN-TCTAGCCTTCTCGTGTGCAGACTTGAGGTCAGTG-NNNNNNNN-TAGAGCATACGGCAGAGAAGACGAAC-5'</p>

Figure 5. Sequences added during Single Cell DNA library construction. Protocol steps and part numbers refer to Chromium Single Cell DNA Reagent Kits User Guide (Document CG000153).

© 2018 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
 10x Genomics
 7068 Koll Center Parkway
 Suite 401
 Pleasanton, CA 94566 USA

