

## TECHNICAL NOTE


# Guide RNA Specifications Compatible with Feature Barcoding technology for CRISPR Screening

## Introduction

The Single Cell Gene Expression Solution coupled with Feature Barcoding technology enables profiling of gene expression in conjunction with additional cellular features from the same single cell. Specifically, the Single Cell Gene Expression and CRISPR Screening Solution provides a high-throughput and scalable approach to obtain gene expression profiles along with CRISPR-mediated perturbation phenotypes in transduced cells via direct capture of poly-adenylated mRNAs and single-guide RNAs (sgRNAs) (Figure 1). This document provides specifications for engineering Feature Barcoding technology compatible sgRNAs for use with the Chromium Single Cell 3' workflow (Figure 2) for CRISPR screening.

## Chromium Single Cell 3' v3 Gel Beads

In addition to a poly(dT) primer sequence that enables the production of barcoded, full-length cDNA from poly-adenylated mRNA, the Single Cell 3' v3 Gel Beads (Figure 1) also include two additional primer sequences (Capture Sequence 1 and Capture Sequence 2) for direct capture and priming of Feature Barcoding technology compatible sgRNAs present in a cell inside a Gel-bead-in-emulsion (GEM).




**Figure 1.** Schematic of a Chromium Single Cell 3' v3 Gel Bead, including the Capture Sequence 1 and the Capture Sequence 2 for direct capture of sgRNA molecules from a cell inside a GEM.

## Capture Sequence Integration in sgRNA

To enable direct capture, each sgRNA should be engineered to contain either Capture Sequence 1 or Capture Sequence 2, along with a protospacer (Feature Barcode), designed to target gene/s of interest. Two possible locations for integrating the capture sequence in the sgRNA include (1) within the sgRNA hairpin structure, or (2) immediately before the sgRNA termination signal, elongating the 3' end of the sgRNA (Table 1). However, alternate sgRNA integration locations for either of the two capture sequences may be possible depending on the specific application, type of construct used etc.

Performing sgRNA QC by qPCR, NGS or other methods is recommended before proceeding with 10x Genomics Single Cell Solutions. The Chromium Single Cell Gene Expression and CRISPR Screening Solution has been shown to be compatible with sgRNA constructs in several CRISPR applications, including interference (CRISPRi), activation (CRISPRa), and cutting.




**Figure 2.** Transduced cells with integrated sgRNA for the Chromium Single Cell 3' workflow for CRISPR screening.

## Capture Sequence 1


Capture Sequence 1 on Gel Bead: 5'-TTGCTAGGACCGCCTTAAAGC-3'

### Capture Sequence 1 integrated in sgRNA hairpin



5'-NNNNNNNNNNNNNNNNNNNNGTTTAAGAGCTAAGCTGGAAACAGCATAGCAAGTTTAAATAAGGCTAGTCCGTTATCAACTTgcccGCTTAAAGCCGGCTCTAGCAAggcccAAGTGGCACCGAGTCGGTGCTTTTTTT-3'

### Capture Sequence 1 integrated in sgRNA 3'-end




5'-NNNNNNNNNNNNNNNNNNNNGTTTAAGAGCTAAGCTGGAAACAGCATAGCAAGTTTAAATAAGGCTAGTCCGTTATCAACTTgaaaAAGTGGCACCGAGTCGGTGCgCCTTAAAGCCGGCTCTAGCAATTTTTTT-3'

## Capture Sequence 2


Capture Sequence 2 on Gel Bead: 5'-CCTTAGCCGCTAATAGGTGAGC-3'

### Capture Sequence 2 integrated in sgRNA hairpin



5'-NNNNNNNNNNNNNNNNNNNNGTTTAAGAGCTAAGCTGGAAACAGCATAGCAAGTTTAAATAAGGCTAGTCCGTTATCAACTTgcccCCTCACCTATTAGCGGCTAAGggcccAAGTGGCACCGAGTCGGTGCTTTTTTT-3'

### Capture Sequence 2 integrated in sgRNA 3'-end



5'-NNNNNNNNNNNNNNNNNNNNGTTTAAGAGCTAAGCTGGAAACAGCATAGCAAGTTTAAATAAGGCTAGTCCGTTATCAACTTgaaaAAGTGGCACCGAGTCGGTGCcCCTCACCTATTAGCGGCTAAGTTTTTTT-3'

Table 1. Integration of Capture Sequence 1 and Capture Sequence 2 in sgRNA.

© 2018 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at [www.10xgenomics.com/legal-notices](http://www.10xgenomics.com/legal-notices), or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:  
[support@10xgenomics.com](mailto:support@10xgenomics.com)  
 10x Genomics  
 7068 Koll Center Parkway  
 Suite 401  
 Pleasanton, CA 94566 USA

