

DEMONSTRATED PROTOCOL

Isolation of Nuclei for Single Cell DNA Sequencing

Overview

This protocol outlines how to isolate, wash, and count single nuclei from frozen solid tissue samples for use with 10x Genomics Single Cell DNA protocol. Modifications to this protocol (e.g. homogenization, centrifugation speed/time, and filtration steps) may be required when working with different sample types. Refer to Appendix for a practice protocol for nuclei isolation from a cell suspension.

Consult Demonstrated Protocol Cell Preparation Guide (Document CG000053) for Tips & Best Practices when handling cells. Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.

Preparation - Buffers

Buffer	Composition	Final Conc.
Lysis Buffer (maintain at 4°C)	Tris-HCl	10 mM
	NaCl	10 mM
	MgCl ₂	3 mM
	Nonidet P40 Substitute	0.1%
	Prepare in nuclease-free water	
Nuclei Wash & Resuspension Buffer (maintain at 4°C)	BSA	0.04%
	Prepare in PBS	
Nuclear Stain (maintain at 4°C)	Ethidium Homodimer-1	2 %
	Glycerol	30%
	Prepare in nuclease-free water	

Specific Reagents & Consumables

Vendor	Item	Part Number
Thermo	RNase-Free Disposable Pellet Pestles	12-141-364
Fisher	Ethidium Homodimer-1	E1169
Scientific	UltraPure BSA (50 mg/ml)	AM2616
Millipore-Sigma	Trizma Hydrochloride Solution, pH 7.4	T2194
	Sodium Chloride Solution, 5M	59222C
	Magnesium Chloride Solution, 1M	M1028
	Nonidet P40 Substitute	74385
	PBS with 10% Bovine Albumin*	SRE0036
	(*alternate to Thermo Fisher Scientific product-AM2616)	
Bel-Art	Flowmi Cell Strainer, 40 µm	H13680-0040
Corning-Cellgro	Phosphate-Buffered Saline 1X without Calcium & Magnesium	21-040-CV

Protocol

Cell Preparation & Sourcing: This Protocol was demonstrated using commercially-sourced fresh frozen normal and tumor tissue from human breast and lung.

Tissue Homogenization & Nuclei Isolation

- Thaw the tube containing the tissue sample on ice to enable tissue transfer. Minimize the time the tissue is maintained on ice.
- Transfer to a petri dish on ice and cut the tissue into ~2–4 mm³ pieces. **Immediately** transfer one piece to a 2-ml microcentrifuge tube on ice. Leftover tissue pieces may be preserved for subsequent analysis by snap freezing the individual pieces on dry ice, followed by storage in vapor-phase liquid nitrogen.
- Add 500 µl chilled Lysis Buffer to the tube containing the tissue piece.
- Using a Pellet Pestle, homogenize the tissue by gently moving the pestle straight up and down 15x. **DO NOT** twist the pestle.

Homogenization by Pestle



Homogenized Tissue



- Centrifuge at 300 rcf for 2-3 sec. Transfer the supernatant containing the nuclei to a 2-ml microcentrifuge tube. Avoid transferring large pieces of the tissue.
- Centrifuge at 850 rcf for 5 min at 4°C.
- Remove the supernatant without disturbing the nuclei pellet.
- Using a **regular-bore** pipette tip, add 1 ml Nuclei Wash & Resuspension Buffer to the nuclei pellet and gently pipette mix 15x.
- Repeat** steps f-h for a total of 2 washes. At the end of the 2 washes, the nuclei suspension may be passed through a 40 µm Flowmi Cell Strainer to remove aggregates.
- Determine the nuclei concentration using a Countess II FL Automated Cell Counter (see Appendix).

-
- k. Using a **regular-bore** pipette tip, concentrate or dilute the nuclei suspension to achieve a target nuclei concentration of **~800–1200 nuclei/μl** (see Appendix).

The final nuclei suspension may contain some small sub-cellular debris; this is sample dependent and will not impact assay performance.

Results

Nuclei were isolated from multiple tissue types and depending on the sample, typical yield was **50,000–2,000,000 nuclei/sample**.

Nuclei isolated from a breast tumor sample using this protocol, were used as input sample in the Single Cell DNA protocol for CNV detection. The results are described in the Application Note – Discover Genetic Heterogeneity in Cancer (Document LIT000024).

Appendix

Practice Protocol

For nuclei isolation from fresh and cryopreserved cultured cells suspension.

- a. Aliquot **1x10⁶ cells** into a 2-ml microcentrifuge tube.
- b. Centrifuge at **300 rcf** for **5 min** at **4°C**.
- c. Remove the supernatant without disturbing the pellet.
- d. Add **500 μl** of chilled Lysis Buffer, pipette mix, and place on ice for **5 min**.
- e. Follow Tissue Homogenization & Nuclei Isolation steps f-k.

Count Nuclei

This protocol provides instructions for counting nuclei using a fluorescent dye – Ethidium Homodimer-1 and the Countess II FL Automated Cell Counter (with RFP light cube) to enable accurate quantification even in the presence of sub-cellular debris. The optimal cell concentration for the Cell Counter is **1,000–4,000 cells/μl**. Refer to manufacturer's instructions for details on operations.

- a. Vortex Nuclear Stain, centrifuge briefly, and aliquot **10 μl** in each tube.
- b. Using a **regular-bore** pipette tip, gently mix the nuclei suspension. **Immediately** add **10 μl** nuclei suspension to **10 μl** Nuclear Stain aliquot. Gently pipette mix **10x**.
- c. Transfer **10 μl** stained nuclei to a Countess II Cell Counting Slide chamber.
- d. Insert the slide into the Countess II FL Cell Counter. Image the nuclei using the RFP setting for fluorescent illumination and filtering. Optimize focus and exposure settings. Confirm the absence of large clumps using the bright-field mode.

- e. Note the RFP-positive concentration. Multiply by dilution factor 2 (from step b) to determine nuclei concentration.

For accurate nuclei counts, sample and count two different areas of the nuclei suspension.

Adjust Nuclei Concentration

If cell debris and large clumps are present, pass the sample through a **40 μm** Flowmi Cell Strainer before counting. Passing through the strainer will reduce nuclei concentration by **~25%**.

- a. Determine nuclei concentration using a Countess II FL Automated Cell Counter.
- b. Based on the nuclei count, concentrate or dilute the sample to achieve the target nuclei concentration.

To Concentrate the Nuclei:

- Centrifuge at **850 rcf** for **5 min** at **4°C**.
- Remove the supernatant without disturbing the nuclei pellet.
- Using a **regular-bore** pipette tip, add an appropriate volume of Nuclei Wash & Resuspension Buffer and gently pipette mix **15x**.

To Dilute the Nuclei:

- Using a **regular-bore** pipette tip, add an appropriate volume of Nuclei Wash & Resuspension Buffer and gently pipette mix **15x**.

© 2018 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
10x Genomics
7068 Koll Center Parkway
Suite 401
Pleasanton, CA 94566 USA

