

## DEMONSTRATED PROTOCOL

# Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols

with Feature Barcoding technology

## Overview

Cell surface proteins can be labeled using a specific protein binding molecule, such as an antibody conjugated to a Feature Barcode oligonucleotide. This protocol provides guidance for antibody-oligonucleotide conjugation and outlines cell surface protein labeling for use with:

- Chromium Single Cell 3' Reagent Kits v3 User Guide with Feature Barcoding technology for Cell Surface Protein (CG000185)
- Chromium Single Cell V(D)J Reagent Kits User Guide with Feature Barcoding technology for Cell Surface Protein (CG000186)

## Additional Guidance

Consult Demonstrated Protocol Cell Preparation Guide (Document CG00053) for Tips & Best Practices on handling cells and Technical Note Guidelines on Accurate Target Cell Counts (Document CG000091) for determining accurate cell counts.

**Cells carry potentially hazardous pathogens. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.**

## Preparation – Buffers

Buffers	Composition
Maintain at 4°C	
Labeling Buffer	PBS + 1% BSA
Resuspension Buffer	PBS + 0.04% BSA
Dextran Sulfate Solution	1% w/v (10 mg/ml) Dextran Sulfate Sodium Salt in Nuclease-free Water

## Specific Reagents & Consumables

For Antibody-Oligonucleotide Conjugation		
Vendor	Item	Part Number
Expedeon	Thunder-Link PLUS Conjugation Kit	425-0300
IDT	Custom DNA Oligos (see Table 1)	-
-	100 µg Purified Azide-free Antibody (1 mg/ml)	-
For Cell Surface Protein Labeling		
Vendor	Item	Part Number
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TotalSeq Antibody-Oligonucleotide Conjugate*	-
MP Biomedicals	Dextran Sulfate Sodium Salt	101516
Thermo Fisher Scientific	Dextran Sulfate Sodium Salt (alternative to MP Biomedicals product)	AC441490050
	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
Millipore Sigma	Phosphate-Buffered Saline (PBS) with 10% Bovine Albumin (alternative to Thermo Fisher product)	SRE0036
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV

\*TotalSeq-B for Single Cell 3' v3 protocol with Feature Barcoding technology for Cell Surface Protein

\*TotalSeq-C for Single Cell V(D)J protocol with Feature Barcoding technology for Cell Surface Protein

## Protocol Overview

### Option A. Custom Conjugated Antibodies

Follow manufacturer's instructions  
(Thunder-Link PLUS) for conjugation and purification

### Option B. Preconjugated Antibodies

BioLegend TotalSeq-B (for Single Cell 3' v3)  
OR  
BioLegend TotalSeq-C (for Single Cell V(D)J)

#### 1. Label Cells



Prepare Antibody Mix  
Centrifuge (14,000 rcf, 10 min)  
Maintain Antibody Mix supernatant at 4°C

Centrifuge cells for labeling (400 rcf, 5 min)  
Remove supernatant

Resuspend pellet in Labeling Buffer (50 µl)

Add Human TruStain FcX (5 µl)  
Add Dextran Sulfate solution (2 µl)  
Gently pipette mix

Incubate at 4°C (10 min)

Add Prepared Antibody Mix supernatant

Add Labeling Buffer  
(for a total volume of 100 µl)  
Gently pipette mix

Incubate at 4°C (30 min)

#### 2. Wash Cells



See Appendix for Labeling Buffer volumes when  
working with 15-ml tubes

Wash 1: Add Labeling Buffer to cells (1.4 ml)

Centrifuge (400 rcf, 5 min)  
Remove supernatant

Wash 2: Add Labeling Buffer to pellet (1.5 ml)

Gently pipette mix

Centrifuge (400 rcf, 5 min)  
Remove supernatant

Repeat 1x  
(Wash 3)

Add Resuspension Buffer to pellet  
Gently pipette mix

Determine cell concentration & viability

Proceed to:

- Chromium Single Cell 3' Reagent Kits v3 User Guide with Feature Barcoding technology for Cell Surface Protein (CG000185) OR
- Chromium Single Cell V(D)J Reagent Kits User Guide with Feature Barcoding technology for Cell Surface Protein (CG000186)



## Antibody-Oligonucleotide Conjugation Guidance

Choose antibody based on cell surface protein/s being labeled and a Feature Barcode oligonucleotide compatible with specific 10x Genomics protocol (Table 1) for conjugation. Alternatively, use compatible pre-conjugated antibodies from BioLegend or other vendors for labeling cells.

### Specific Reagents for Conjugation

Vendor	Item	Part Number
Expedeon	Thunder-Link PLUS Conjugation Kit	425-0300
IDT	Custom DNA Oligos (see Table 1)	-
-	100 µg Purified Azide-free Antibody (1 mg/ml)	-
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV

#### Oligonucleotide:

Use ≥10 nmole HPLC-purified and lyophilized Feature Barcode oligonucleotide for conjugation. The oligonucleotide must contain an amine group at 5'-end (5' amine modified; IDT code /5AmMC12/). Resuspend lyophilized oligonucleotide in PBS or other compatible buffer (see Thunder-Link PLUS) at 100 µM, i.e. 10 nmole dissolved in 100 µl buffer.



DO NOT use Tris buffers for resuspension as they are not compatible with conjugation.

#### Conjugation:

Follow manufacturer's instructions (Thunder-Link PLUS Oligo Conjugation Kit from Expedeon) for antibody-oligonucleotide conjugation. Thunder-Link PLUS is compatible with many purified antibodies.

#### Antibody-Oligonucleotide Ratio:

This protocol was demonstrated using 1:3 antibody-oligonucleotide ratio for conjugation. Optimization may be needed depending on the antibodies used.

#### Conjugate Purification:


The antibody-oligonucleotide conjugate purification is recommended to remove any unbound oligonucleotides. Follow the Thunder-link PLUS conjugate purification protocol.

#### Verification of Conjugation:

Verify the conjugation by comparing the antibody-oligonucleotide conjugate with an unconjugated antibody control resolved on a non-reducing, gradient SDS-PAGE gel.

- Run a known volume and concentration of unconjugated antibody next to a known volume of antibody-oligonucleotide conjugate on SDS-PAGE gel.
- Estimate the conjugate concentration and calculate the degree of conjugation by comparing the respective band intensities.

Figure 1. Verification of conjugation on a 4-12% gradient SDS-PAGE gel under non-reducing conditions.



OPTIONAL Use a BCA or Bradford Protein Assay Kit to calculate the final antibody concentration.

Table 1. Feature Barcode Oligonucleotide Sequence for Antibody Conjugation.

10x Genomics Protocol	Feature Barcode Oligonucleotide Sequence			
Single Cell 3' v3 – Cell Surface Protein (CG000185)	/5AmMC12/GTGACTGGAGTTCCAGACGTGTGCTCTCCGACTCTNNNNNNNNNN- TruSeq Read 2	10 nt	NNNNNNNNNNNNNNNN- Feature Barcode (15 nt)	NNNNNNNNNNGCTTTAAGGCCGGTCTAGCAA 9 nt Capture Sequence 1
Single Cell V(D)J – Cell Surface Protein (CG000186)	/5AmMC12/CGGAGATGTGTATAAGAGACAGNNNNNNNNNN- Nextera partial Read 2	10 nt	NNNNNNNNNNNNNNNN- Feature Barcode (15 nt)	NNNNNNNNNNCCCATATAAGAAA 9 nt Capture Sequence

See Appendix for an illustrative overview of antibody-oligonucleotide conjugate capture by 10x Gel Bead primers. Consult Barcode Whitelist for Custom Feature Barcoding conjugates (Document CG000193), for more information.

## Cell Surface Protein Labeling Protocol

### Specific Reagents & Buffers

Vendor	Item	Part Number
BioLegend	Human TruStain FcX (Fc Receptor Blocking Solution)	422301
	TotalSeq Antibody-Oligonucleotide Conjugate*	-
MP Biomedicals	Dextran Sulfate Sodium Salt	101516
Thermo Fisher Scientific	Dextran Sulfate Sodium Salt (alternative to MP Biomedicals product)	AC441490050
	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
Millipore Sigma	Phosphate-Buffered Saline (PBS) with 10% Bovine Albumin (alternative to Thermo Fisher product)	SRE0036
Corning	Phosphate-Buffered Saline, 1X without Calcium and Magnesium	21-040-CV

Buffers	Composition
Maintain at 4°C	
Labeling Buffer	PBS + 1% BSA
Resuspension Buffer	PBS + 0.04% BSA
Dextran Sulfate Solution	1% w/v (10 mg/ml) Dextran Sulfate Sodium Salt in Nuclease-free Water

\*TotalSeq-B for Single Cell 3' v3 protocol with Feature Barcoding technology for Cell Surface Protein

\*TotalSeq-C for Single Cell V(D)J protocol with Feature Barcoding technology for Cell Surface Protein

### 1. Label Cells

**Prepare Antibody Mix:** Add appropriate/manufacturer's recommended amount of antibody-oligonucleotide conjugates to a 1.5-ml microcentrifuge tube. Centrifuge the mix at **14,000 rcf** for **10 min** at **room temperature**. Transfer the supernatant (containing Antibody Mix) to a new tube and maintain at **4°C**.

This protocol was demonstrated using  $0.2-1 \times 10^6$  PBMCs.

- Centrifuge cells at **400 rcf** for **5 min** at **4°C**. Use of swinging-bucket rotor is recommended for higher cell recovery.
- Remove the supernatant without disturbing the pellet.
- Resuspend cell pellet in **50 µl** Labeling Buffer.
- Add **5 µl** Human TruStain FcX.
- Add **2 µl** Dextran Sulfate Solution. Gently pipette mix.
- Incubate for **10 min** at **4°C**.
- Add the prepared Antibody Mix supernatant.


- Add Labeling Buffer to the cell suspension to bring the total volume to **100 µl**. Gently pipette mix 10x (pipette set to 90 µl).
- Incubate for **30 min** at **4°C**.

### 2. Wash Cells

To eliminate non-specific binding with comparable efficiency, wash steps may be performed either in 1.5-ml microcentrifuge tubes or 15-ml tubes using indicated Labeling Buffer volumes. If concerned about dislodging the pellet during supernatant removal, perform wash steps in 15-ml tubes (see Appendix).

Optimization of centrifugation speed/time may be needed based on cell type.

- Wash by adding **1.4 ml** Labeling Buffer (for a total volume of **1.5 ml**) to the cells from step 1i.
- Centrifuge at **400 rcf** for **5 min** at **4°C**. Larger or fragile cell types may require slower centrifugation speeds.
- Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet.




Leaving behind excess supernatant may cause non-specific binding, which may result in increased background reads during sequencing.

- Using a **wide-bore** pipette tip, resuspend the cell pellet in **1.5 ml** Labeling Buffer.
  - Centrifuge at **400 rcf** for **5 min** at **4°C**.
  - Remove the supernatant without touching the bottom of the tube to avoid dislodging the pellet.
  - Repeat **d - f**.
  - Based on starting cell concentration and assuming ~50% cell loss, add an appropriate volume of Resuspension Buffer to obtain a cell concentration of 700-1,200 cells/µl and gently pipette mix using a **regular-bore** pipette tip until a single cell suspension is achieved.
  - Determine cell concentration and viability using a Countess II Automated Cell Counter or a hemocytometer.
  - Proceed **immediately** to:
    - Chromium Single Cell 3' Reagent Kits v3 User Guide with Feature Barcoding technology for Cell Surface Protein (CG000186)
- OR
- Chromium Single Cell V(D)J Reagent Kits User Guide with Feature Barcoding technology for Cell Surface Protein (CG000186)


## Appendix

To eliminate non-specific binding with comparable efficiency, wash steps may be performed either in 1.5-ml microcentrifuge tubes or 15-ml tubes using indicated Labeling Buffer volumes. Non-specific binding contributes to increased background reads during sequencing. If concerned about dislodging the pellet during supernatant removal, perform wash steps in 15-ml tubes.

### Illustrative Overview of Wash Steps in 1.5-ml Tubes



### Illustrative Overview of Wash Steps in 15-ml Tubes




## Appendix

Antibody-oligonucleotide conjugate capture by protocol specific Gel Bead primers is illustrated below.




### Illustrative Overview of Antibody-Oligonucleotide Conjugate Capture

#### Single Cell 3' v3 – Cell Surface Protein (CG000185)



#### Single Cell V(D)J – Cell Surface Protein (CG000186)



© 2018 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: [www.10xgenomics.com/trademarks](http://www.10xgenomics.com/trademarks). 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: [www.10xgenomics.com/patents](http://www.10xgenomics.com/patents). The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at <http://www.10xgenomics.com/legal-notices>, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**Contact:**  
[support@10xgenomics.com](mailto:support@10xgenomics.com)  
 10x Genomics  
 7068 Koll Center Parkway  
 Suite 401  
 Pleasanton, CA 94566 USA

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10X GENOMICS STANDARD WARRANTY, AND 10X GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

