

TECHNICAL NOTE

Chromium Connect: Consistent Automated Single Cell Gene Expression Library Generation

Introduction

Chromium Connect automates the preparation of sequencing-ready, single cell libraries from input cell samples, minimizing technical variation in single cell gene expression data. Generation of Chromium Single Cell 3' Gene Expression libraries on the Chromium Connect instrument includes automated Gel Beads-in-emulsion (GEM) generation, barcoding, and library preparation from single cell suspensions, along with additional functionalities for library quantification and pooling. This Technical Note highlights the consistency and reproducibility of single cell gene expression data derived from libraries generated using the Chromium Connect automated workflow.

Chromium Connect Automated Workflow

The Chromium Connect automated workflow (Figure 1) generates Chromium Single Cell 3' Gene Expression libraries from up to 8 single cell samples per run, using Chromium Connect specific automation-compatible reagents and consumables with minimal user interaction. The libraries can be sequenced and the data analyzed and visualized using Cell Ranger and Loupe Cell Browser for assessing single cell gene expression.




Figure 1. Chromium Connect with integrated components for an automated workflow for single cell library generation.

Methods

Peripheral Blood Mononuclear Cell (PBMC) single cell suspension was prepared as described in Demonstrated Protocol for Fresh Frozen Human Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing (Document CG00039). Cells were drawn from the same pool and Chromium Single Cell 3' Gene Expression libraries were generated targeting 10,000 cells/sample for each of the 8 channels, using either Chromium Connect instruments (up to 4) for the automated workflow or the Chromium Next GEM Single Cell 3' Reagent Kits v3.1 protocol (Document CG000204) for the manual workflow.

Results


Highly consistent and reproducible data derived from Chromium Single Cell 3' Gene Expression libraries generated using the Chromium Connect automated workflow is shown in Table 1. A comparison of single cell gene expression data generated using the Chromium Connect automated workflow and the manual workflow is depicted in Table 2.

Table 1. Chromium Single Cell 3' Gene Expression libraries generated using the Chromium Connect automated workflow show highly consistent and reproducible data.

Automated Workflow–Data Consistency & Reproducibility


QC Traces: Channel-to-channel consistency

PBMC single cell suspension samples (1 replicate in each of the 8 channels in an instrument) were used to generate single cell gene expression libraries using the Chromium Connect automated workflow. Consistent cDNA QC and library QC traces were observed across all 8 channels of the chip used in the automated workflow, as shown in the representative Agilent BioAnalyzer traces. All sample traces showed comparable size distribution and concentration, indicating highly consistent performance.



Gene Expression Profile: Channel-to-channel consistency

Each of the libraries generated across the 8 channels of the Chromium Connect instrument showed highly consistent gene expression profiles. As shown below, the correlation value of gene expression profile of libraries derived from each channel relative to channel 1 was >0.99 .



Gene Expression Profile: Instrument-to-instrument consistency

The gene expression profiles of PBMCs used for generating single cell gene expression libraries across 4 different Chromium Connect instruments showed high degree of consistency and correlation.





Table 2. Comparison of data derived from Chromium Single Cell 3' Gene Expression libraries generated using either the Chromium Connect automated workflow or the manual workflow.

Single Cell Gene Expression Data—Automated vs. Manual Workflow

Minimal User Variability

The Chromium Connect automated hands-off workflow minimizes user variability. Differences in user handling likely contribute to a wider distribution of median genes and Unique Molecular Identifier (UMI) counts recovered per cell in 96 libraries prepared using the manual workflow compared to 96 libraries generated using the Chromium Connect automated workflow (over multiple days, targeting 500-10,000 cells/sample).



Consistent Cell Type Profiling

The gene expression profiles compared in libraries generated from 10,000 PBMCs using either the manual or the Chromium Connect automated workflow identifies all major PBMC cell types. The single cell barcode tSNE plots show the gene expression data analyzed and visualized using Cell Ranger and Loupe Cell Browser, respectively. Major cell populations in each cluster were identified via marker gene expression. Comparable cell type profile among the two sets of data is also shown in the graph below. This demonstrates that the automated workflow generates consistent gene expression data for cell type profiling.





Table 2 contd. Comparison of data derived from Chromium Single Cell 3' Gene Expression libraries generated using either the Chromium Connect automated workflow or the manual workflow.


Single Cell Gene Expression Data—Automated vs. Manual Workflow

Accurate Cell Type Marker Detection

To investigate possible changes in gene expression due to automated workflow, expression of major blood cell type markers were compared in data generated using either the automated or the manual workflows. All the key markers were detected in both datasets (see t-SNE plots below).



Comparison of cell type marker expression data generated from 10,000 PBMCs using 4 channels in one Chromium Connect run highly correlated with the data generated using the manual workflow (4 channels/run by 1 user) with a correlation value of >0.99, as shown in the plot below.



© 2019 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
10x Genomics
6230 Stoneridge Mall Road
Pleasanton, CA 94588 USA

