

10x Genomics®

Sample Preparation Demonstrated Protocol

Fresh Frozen Human-Mouse Cell Line
Mixtures for Single Cell RNA Sequencing



Notices

Manual Part Number

CG00014 Rev E

Legal Notices

© 2017 10x Genomics, Inc. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, Inc., is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. "10x", "10x Genomics", "Changing the Definition of Sequencing", "Chromium", "GemCode", "Loupe", "Long Ranger", "Cell Ranger" and "Supernova" are trademarks of 10x Genomics, Inc. All other trademarks are the property of their respective owners. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Product(s) in practicing the methods set forth herein has not been validated by 10x, and such non-validated use is NOT COVERED BY 10X STANDARD WARRANTY, AND 10X HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE.

Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics, Inc., terms and conditions of sale for the Chromium™ Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics, Inc that it currently or will at any time in the future offer or in any way support any application set forth herein.

Customer Information and Feedback

For technical information or advice, please contact our Customer Technical Support Division online at any time.

Email: support@10xgenomics.com

10x Genomics

7068 Koll Center Parkway

Suite 401


Pleasanton, CA 94566 USA

Table of Contents

Demonstrated Protocol	1
1. Overview	2
2. Getting Started	3
2.1. Tips & Safety	3
2.2. General Materials	4
2.3. Specific Cell Sourcing & Culturing	5
2.4. Preparation – Buffers & Media	5
3. Cryopreservation of Human-Mouse Cell Mixtures	5
3.1. Cryopreservation Protocol	5
4. Thawing & Washing Human-Mouse Cell Mixtures	7
4.1. Preparation	7
4.2. Thawing, Washing & Counting Cells	7


Demonstrated Protocol

Fresh Frozen Human-Mouse
Cell Line Mixtures for Single
Cell RNA Sequencing



1. Overview

The ability of the 10x Genomics® Single Cell Solutions to partition single cells in a heterogeneous population can be verified by profiling a mixture of cells from two different species. Ideally, all sequence reads from a single Gel Bead-in-EMulsion (GEM) will be unambiguously mapped to the transcriptome of only one of the two species. The fraction of GEMs containing a mixture of human and mouse transcripts is used to infer doublet rates (below, red circles).



10x Genomics routinely uses a 1:1 mixture of human and mouse cells to validate the technical performance of the 10x Genomics Single Cell Solutions. For convenience, cell mixtures are cryopreserved in large batches, thawed when needed, and then washed and counted prior to loading into the system without the need for further cell culturing.

This Demonstrated Protocol outlines (i) how to generate cryopreserved 1:1 mixtures of human and mouse cells and (ii) how to thaw, wash, and count cryopreserved 1:1 mixtures of human and mouse cells in preparation for use in 10x Genomics Single Cell Protocols.

2. Getting Started

2.1. Tips & Safety

Best practices for handling any cell line include using sterile technique, nuclease-free reagents and consumables, minimizing pipetting steps, and using wide-bore pipette tips where possible to minimize cell damage.

To determine accurate cell counts, best practices include sampling the cell suspension at least twice and at least two counts on each sample (*i.e.* a minimum of four counts in total, based on two independent draws from the cell suspension).

CRITICAL!

*293T/17 cells contain Adeno and SV-40 viral DNA sequences, which requires classification at **Biosafety Level 2** based on U.S. Public Health Service Guidelines. NIH/3T3 cells are classified at Biosafety Level 1. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage and disposal of biological materials.*

2.2. General Materials

Supplier	Description	Part Number (US)
-	Refrigerated Benchtop Centrifuge for 15 ml and 50 ml tubes Microcentrifuge for 2 ml LoBind tubes Heated Water Bath, 2l	- - -
Rainin	Tips LTS 1ML Filter RT-L1000FLR Tips LTS W-O 1ML Fltr RT-L1000WFLR Pipet-Lite LTS Pipette L-1000XLS+	17007954 17014297 17014382
Sigma-Aldrich	Phosphate-Buffered Saline (PBS) with 10% Bovine Albumin (alternative to Thermo Fisher product)	SRE0036
Miltenyi	MACS SmartStrainers, 30 µm	130-098-458
Thermo Fisher Sci	UltraPure™ Bovine Serum Albumin (BSA, 50 mg/ml) (alternative to Sigma-Aldrich product) Trypan Blue Stain (0.4%) for use with the Countess® Automated Cell Counter Countess® II Automated Cell Counter Countess® II Automated Cell Counting Chamber Slides Gibco 0.25% Trypsin-EDTA (1X), phenol red Nunc™ Biobanking and Cell Culture Cryogenic Tubes, 1.8 ml Dimethyl Sulfoxide (DMSO), for molecular biology	AM2616 T10282 AMQAX1000 C10228 25200056 368632 AC327182500
Eppendorf	DNA LoBind Tubes, 2.0 ml*	022431048
Integra	PIPETBOY acu 2	155018
VWR	Sterile Polypropylene Centrifuge Tubes with Flat Caps, 50 ml Sterile Polypropylene Centrifuge Tubes with Flat Caps, 15 ml 10 ml Serological Pipette Seradigm Premium Grade Fetal Bovine Serum (FBS) 75 cm ² Cell Culture Flask, Canted Neck	82018-050 21008-103 89130-898 97068-085 46610-078
ATCC	DMEM	30-2002
Bioction	CoolCell® FTS30 Cell Freezing Container	BSC-170

*No substitutions are allowed. Items have been validated by 10x Genomics® and are required for Single Cell workflow, training and system operations.

2.3. Specific Cell Sourcing & Culturing

- a) Obtain vials of cryopreserved human 293T/17 (CRL-11268) and murine NIH/3T3 (CRL-1658) cells from ATCC or a licensed local distributor. Culture the cells as recommended by the distributor.

NOTE

293T/17 cells grow to a higher density than NIH/3T3 cells, and NIH/3T3 cells will often be the limiting reagent in this part of the Protocol. Approximately 5 times more NIH/3T3 culture flasks are required to achieve balanced cell numbers (e.g., 4 confluent T75 flasks of 293T/17s will require ~20 nearly confluent T75 flasks of NIH/3T3 cells).

2.4. Preparation – Buffers & Media

- a) Prepare 1 ml chilled (4°C) cryopreservation medium (DMEM + 20% FBS + 10% DMSO) for every 10^6 cells to be cryopreserved.
- b) Chill a CoolCell® FTS30 or equivalent cell freezing container in a 4°C refrigerator.

3. Cryopreservation of Human-Mouse Cell Mixtures

3.1. Cryopreservation Protocol

- a) Remove medium from cultured 293T/17 and NIH/3T3 cells, rinse culture vessels, detach cells using trypsin, spin down cells, resuspend cell pellets in growth medium by pipetting up and down 10 times, and strain cells to remove clumps and ensure single cell suspensions. (Consult Sample Preparation Demonstrated Protocol *Single Cell Suspensions from Cultured Cell Lines for Single Cell RNA Sequencing* (CG00054) for more details on harvesting adherent cells).
- b) Using a Countess® II Automated Cell Counter, count the 293T/17 and NIH/3T3 single cell suspensions after cell straining and measure the total volume of each cell culture. If the concentration is too high to achieve an accurate cell count, further dilute the cells with growth medium, thoroughly mix the resulting cell suspension, and recount. An accurate count for each cell type is important.
- c) Based on total volumes (V) and concentrations (C), calculate total cell number ($N = C \times V$) for each cell type to determine which cell type is limiting.

Example:

$$V_{293T/17} = 10 \text{ ml} \quad C_{293T/17} = 3 \times 10^6 \text{ cells/ml} \quad N_{293T/17} = 3 \times 10^7 \text{ cells}$$

$$V_{NIH/3T3} = 6 \text{ ml} \quad C_{NIH/3T3} = 2 \times 10^6 \text{ cells/ml} \quad N_{NIH/3T3} = 1.2 \times 10^7 \text{ cells}$$

NIH/3T3 cells are limiting

DEMONSTRATED PROTOCOL Fresh Frozen Human-Mouse Cell Line Mixtures

- d) Determine how much of the non-limiting cell type should be added to the limiting cell type to achieve a 1:1 ratio, then prepare the 1:1 Human-Mouse mixture. Calculate the total cell number (293T/17 + NIH/3T3 cells).

Example continued:

Since NIH/3T3 cells are limiting

$$V_{293T/17} = \frac{N_{NIH/3T3}}{C_{293T/17}} = \frac{1.2 \times 10^7 \text{ cells}}{3 \times 10^6 \text{ cells/ml}} = 4 \text{ ml}$$

4 ml 293T/17 cells is added to 6 ml of NIH/3T3 cells to generate 10 ml of a 1:1 mixture

The $N_{293T/17+NIH/3T3}$ in this mixture is 2.4×10^7 cells

- e) Centrifuge the 1:1 mixture at **300 rcf** for **5 min**.
- f) Resuspend the pellet in an appropriate volume of cryopreservation medium to achieve an overall cell concentration of 1×10^6 cells/ml in the cryovials.

Example continued:

After centrifugation, cells resuspended in 24 ml chilled cryopreservation medium

- g) Dispense **1 ml** aliquots of the 1:1 mixture into cryovials.
- h) Remove the chilled CoolCell® FTS30 from the 4°C refrigerator and place the cryovials inside.
- i) Place the chilled CoolCell® FTS30 in the **-80°C** freezer for at least **4 h**, ensure that the bottom and top vents are not obstructed to allow for adequate air flow.
- j) After **4 h**, transfer the cryovials to either liquid nitrogen or vapor-phase nitrogen for long-term storage.

4. Thawing & Washing Human-Mouse Cell Mixtures

4.1. Preparation

- a) Warm a water bath to 37°C prior to commencing the thawing protocol.
- b) Prepare 1X PBS with 0.04% BSA (400 µg/ml) solution for the final resuspension.

4.2. Thawing, Washing & Counting Cells

- a) Remove cryovial(s) from liquid or vapor-phase nitrogen storage and rapidly thaw in the water bath at **37°C** for **2 – 3 min**. Remove from the water bath when a tiny ice crystal remains.
- b) Pipette mix the cells and transfer the entire volume to a 2 ml tube.
- c) Centrifuge cells at **150 rcf** for **3 min**.

CRITICAL!

Depending on the rotor type, the cell pellet forms on the side or on the bottom of the conical vial. Know the expected position of the pellet, especially when working with small or limited cells, as the pellet can be difficult to see.

- d) Using a disposable transfer pipette with fine tip, remove supernatant without disrupting cell pellet.
- e) Using a **wide-bore** pipette tip, add **1 ml** 1X PBS with 0.04% BSA to each tube and gently pipette mix 5 times and invert tubes to resuspend cell pellet. Pool the tubes if necessary.
- f) Centrifuge cells at **150 rcf** for **3 min**.
- g) Remove supernatant without disrupting cell pellet.
- h) Repeat steps e – g for a total of 2 washes.
- i) Using a **regular-bore** pipette tip, add **1 ml** 1X PBS with 0.04% BSA. Gently pipette mix 10 – 15 times or until cells are completely suspended.

Repeat

NOTE

Do not invert the tube in this step, as cells can stick to the sides of the tube, thereby changing the cell concentration.

- j) Determine the cell concentration using a Countess® II Automated Cell Counter. The target cell concentration is $\sim 7 \times 10^5$ cells/ml (700 cells/µl).
- k) Proceed with the 10x Genomics® Single Cell Protocol.